

Determination of per- and polyfluorinated alkyl substances (PFAS) in drinking water using automated solid-phase extraction and LC-MS/MS

Changling Qiu, Xin Zhang, Rahmat Ullah, Wei Chen, Yan Liu, Chris Shevlin Thermo Fisher Scientific, Sunnyvale, CA

INTRODUCTION

Per- and polyfluorinated alkyl substances (PFAS) are a group of man-made chemicals including perfluorooctanoic acid (PFOA), perfluorooctanesulfonic acid (PFOS), and GenX chemicals that have been manufactured and used in a variety of industries globally.^{1,2} These compounds have a wide range of commercial product applications including industrial polymers, stain repellents, surfactants, waterproofing products, packaging, and aqueous film forming foams used for firefighting. PFAS are highly soluble in water, chemically stable, persistent in the environment, and can accumulate in the human body over time, leading to adverse human health effects.³

In November 2018, the United States Environmental Protection Agency (U.S. EPA) published Method 537.1 "Determination of selected per- and polyfluorinated alkyl substances in drinking water by solid phase extraction and LC/MS/MS".⁴ The method uses an offline solid-phase extraction (SPE) with liquid chromatography tandem mass spectrometry (LC-MS/MS) to extract, enrich, and determine 18 PFAS in drinking water. Currently most testing laboratories perform the sample extraction manually using a vacuum manifold, which is labor-intensive, time consuming, and the flow rate through the cartridge is difficult to control. There is a high demand for automation of the SPE procedure.

In this poster, we discuss the development of an analytical method using an automated SPE system, AutoTrace 280, and LC-MS/MS for determination of eighteen PFAS following the guidelines provided by U.S. EPA Method 537.1. We have demonstrated that the AutoTrace 280 system provides reliable automated SPE for determination of PFAS in large-volume (20 mL–4 L) aqueous samples.

MATERIALS AND METHODS

Equipment

- Thermo Scientific™ Dionex™ AutoTrace 280 PFAS™ System (1)
- Thermo Scientific™ Vanquish™ Flex Duo UHPLC system, fitted with Thermo Scientific™ PFC free kit (2)
- Thermo Scientific™ TSQ Fortis™ triple quadrupole mass spectrometer (3)
- Organomation Associates™ 12 Position N-EVAP Nitrogen Evaporator (4)



Method Workflow




Figure 1. U.S. EPA Method 537.1 procedure workflow

Table 1. Cartridge Conditioning and Sample Loading (program this method in solid-phase extraction mode)

No.	Method (programmed)	User intervention/information
1	Process six samples using the following method steps	
2	Condition cartridge with 7.5 mL of MeOH into solvent waste	
3	Condition cartridge with 7.5 mL of MeOH into solvent waste	
4	Condition cartridge with 9.0 mL of water into aqueous waste	
5	Condition cartridge with 9.0 mL of water into aqueous waste	
6	Load 270.0 mL of sample onto cartridge	Sample bottle actually contains 250 mL of sample. The method is programmed to deliver 270 mL sample as it accounts for the delay volume in the system. Waste automatically goes to aqueous waste.
7	Pause and Alert operator, resume when CONTINUE is pressed	Add 7.5 mL reagent water into sample bottle, swirl over the inner walls to rinse out any residual sample. Make sure the sample weights are at the bottom of the sample bottle submerged into sample.
8	Load 17.5 mL of sample onto cartridge	The method is programmed to consider the delay volume
9	Pause and Alert operator, resume when CONTINUE is pressed	Add 7.5 mL water into sample bottle, swirl over the inner walls to rinse out any residual sample. Make sure the sample weights are at the bottom of the sample bottle submerged into sample.
10	Load 21.5 mL of sample onto cartridge	The method is programmed to consider the delay volume and to pull all the aqueous phase from the tubes.
11	Dry cartridge with gas for 10.0 minutes	
12	End	

Step	Flow rate (mL/min)	SPE parameters	Instrument parameters
Cond flow	10.0	Push delay 5 s	Max elution vol. 20.0 mL
Load flow	10.0	Air factor 1.0	Exhaust fan on Yes
Rinse flow	10.0	Autowash vol. 1.0 mL	Beeper on Yes
Elute Flow	1.0		
Cond air push	15.0		
Rinse air push	20.0		
Elute air push	5.0		

Table 2. Sample Elution (program this method in solid-phase elute mode)

No.	Method (programmed)	User intervention/information
1	Process six samples using the following method steps	
2	Manually rinse sample container with 14.0 mL to collect	First, elute with 4.0 mL MeOH. The method is programmed to consider the delay volume.
3	Pause and alert operator, resume when CONTINUE is pressed	Add 4.0 mL methanol into sample bottle, swirl over the inner walls to rinse out any residual sample. Make sure the sample weights are at the bottom of the sample bottle submerged into that 4 mL methanol.
4	Manually rinse sample container with 18.0 mL to collect	Second, elute with 4.0 mL MeOH. The method is programmed to consider the delay volume and push out any residual methanol.
5	End	

Step	Flow rate (mL/min)	SPE parameters	Instrument parameters
Cond flow	1.0	Push delay 5 s	Max elution vol. 20.0 mL
Load flow	1.0	Air factor 1.0	Exhaust fan on Yes
Rinse flow	1.0	Autowash vol. 1.0 mL	Beeper on Yes
Elute flow	1.0		
Cond air push	15.0		
Rinse air push	20.0		
Elute air push	5.0		

*Do not detach the cartridges during methods one and two

Table 3. LC conditions and MS parameters

Parameter	Value	Parameter	Value
Analytical column	Accucore RP-MS, 2.1 x 100 mm, 2.6 µm	Ion source type	H-ESI
Isolator column	Hypersil BDS C18, 2.1 x 50 mm, 5 µm. This column was installed prior to the autosampler to remove any contaminants from the mobile phase.	Polarity	Negative
Column temp.	45 °C	Negative ion	2500 V
Flow rate	0.5 mL/min	Sheath gas	50 arbitrary units
Injection volume	5 µL	Aux gas	10 arbitrary units
Autosampler temp.	6 °C	Sweep gas	1 arbitrary units
Solvent A	Water containing 0.1% acetic acid	Ion transfer tube temp.	325 °C
Solvent B	Methanol containing 0.1% acetic acid	Vaporizer temp.	300 °C
Solvent C	20 mM ammonium acetate in water	Q1 resolution (FWHM*)	0.7
	Time (min)	%B	%C
	0	30	5
	1	30	5
	14	95	5
	17	95	5
	18	30	5
	21	30	5
Gradient			

*FWHM: Full width at half maximum




Figure 2. LC-MS/MS chromatograms of PFAS at 4 µg/L standard solution

Table 4. Retention time, asymmetry factor, and internal standards for method PFAS

Peak No.	Peak Name	Retention Time (min)	Asymmetry Factor	IS # ref
1	PFBS	4.56	1.09	¹³ C ₇ -PFOS
2	PFHxA	6.56	1.01	¹³ C ₇ -PFOA
3	¹³ C ₇ -PFHxA	6.56	0.96	¹³ C ₇ -PFOA
4	HFPO-DA	7.16	0.84	¹³ C ₇ -PFOA
5	¹³ C ₇ -HFPO-DA	7.16	0.84	¹³ C ₇ -PFOA
6	PFHpA	8.37	1.01	¹³ C ₇ -PFOA
7	ADONA	8.57	1.12	¹³ C ₇ -PFOS
8	PFHxS	8.58	0.95	¹³ C ₇ -PFOA
9	PFOA	9.65	1.06	¹³ C ₇ -PFOA
10	¹³ C ₇ -PFOA	9.66	0.98	--
11	PFNA	10.66	0.99	¹³ C ₇ -PFOA
12	PFOS	10.70	1.03	¹³ C ₇ -PFOS
13	¹³ C ₇ -PFOS	10.70	1.04	--
14	9Cl-PF3ONS	11.16	1.16	¹³ C ₇ -PFOS
15	PFDA	11.50	1.03	¹³ C ₇ -PFOA
16	¹³ C ₇ -PFDA	11.50	0.95	¹³ C ₇ -PFOA
17	NMeFOSAA	11.96	1.08	--
18	d ₃ -NMeFOSAA	11.97	1.05	d ₃ -NMeFOSAA
19	PFUnA	12.19	1.00	¹³ C ₇ -PFOA
20	NEFOSAA	12.34	0.93	¹³ C ₇ -PFOA
21	d ₃ -NEFOSAA	12.35	1.10	d ₃ -NMeFOSAA
22	11Cl-PF3OUdS	12.50	1.05	¹³ C ₇ -PFOS
23	PFDoA	12.78	1.07	¹³ C ₇ -PFOA
24	PFTFA	13.27	1.01	¹³ C ₇ -PFOA
25	PFTA	13.70	0.94	¹³ C ₇ -PFOA

Table 5. Precision and accuracy (n=6) of PFAS in fortified drinking water

Analyte	Fortified conc. (ng/L)	Mean recovery (%)	RSD (%)	Fortified conc. (ng/L)	Mean recovery (%)	RSD (%)
PFBS	16.0	107	3.3	80.0	98.3	3.6
PFHxA	16.0	108	2.3	80.0	106	2.6
HFPO-DA	16.0	84.1	7.5	80.0	88.6	6.3
PFHpA	16.0	113	2.7	80.0	117	1.3
PFHxS	16.0	120	3.4	80.0	123	2.1
ADONA	16.0	117	2.5	80.0	121	1.1
PFOA	16.0	113	2.5	80.0	119	1.6
PFNA	16.0	114	2.9	80.0	118	2.1
PFOS	16.0	113	4.5	80.0	117	2.9
9Cl-PF3ONS	16.0	96.1	4.1	80.0	103	2.6
PFDA	16.0	105	3.2	80.0	111	2.1
PFUnA	16.0	96.8	5.0	80.0	103	3.1
NMeFOSAA	16.0	103	5.2	80.0	110	5.2
11Cl-PF3OUdS	16.0	88.5	5.5	80.0	97.1	4.8
NEFOSAA	16.0	100	9.9	80.0	104	2.3
PFDoA	16.0	89.8	4.4	80.0	97.3	3.4
PFTFA	16.0	89.6	3.8	80.0	95.8	3.7
PFTA	16.0	89.0	4.8	80.0	98.1	3.3

CONCLUSIONS

This application note reports a method that can be used for the extraction and determination of 18 PFAS in drinking water with a PFAS-safe AutoTrace 280 extraction system and LC-MS/MS. The modified AutoTrace 280 extraction system ensures inertness and prevents PFAS from leaching into sample during extraction, while at same time delivering consistent and reliable performance. Both sample path cleaning in SPE and separation method precaution for the LC system maintained a low system background, meeting the EPA method requirement. The calculated LCMRLs ranged from 0.20 to 3.5 ng/L and the MDLs ranged from 0.30 to 2.5 ng/L, which were below or comparable to those values reported in U.S. EPA Method 537.1. At both 16.0 ng/L and 80.0 ng/L fortified concentration levels, all the recoveries were within the acceptable range of 70–130%. The calculated RSDs were all less than 10%, suggesting good precision. Thermo Scientific LC-MS/MS with the automatic extraction AutoTrace 280 system demonstrated an efficient, reliable, and sensitive method to fulfill the requirements of U.S. EPA Method 537.1.

REFERENCES

1. U.S. EPA. Basic Information about Per- and Polyfluoroalkyl Substances (PFAS). <https://www.epa.gov/pfas/basic-information-about-and-polyfluoroalkyl-substances-pfas>
2. National Institute of Health, National Institute of Environmental Health Sciences, Polyfluoroalkyl Substances (PFAS). <https://www.niehs.nih.gov/health/topics/agents/pfc/index.cfm#footnote8>
3. U.S. EPA. Drinking Water Health Advisories for PFOA and PFOS. <https://www.epa.gov/ground-water-and-drinking-water/drinking-water-health-advisories-pfoa-and-pfos>
4. U.S. EPA. Method 537.1, Detections of Selected Per- and Polyfluorinated Alkyl Substances in Drinking Water by Solid Phase Extraction and Liquid Chromatography/Tandem Mass Spectrometry (LC/MS/MS). Version 1.0, November 2018 <https://www.epa.gov/water-research/epa-drinking-water-research-methods>

TRADEMARKS/LICENSING

© 2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.