




Glyphosate

## Analysis of Underivatized Glyphosate and Other Polar Contaminants by LC-MS/MS

Shun-Hsin Liang, Ph.D

Senior Scientist



Per- And Polyfluoroalkyl Substances

**RESTEK**<sup>®</sup>

Pure Chromatography

[www.restek.com](http://www.restek.com)

# Outline

- The Raptor Polar X Column
- Glyphosate and Polar Contaminants Analysis
- Novel Solution for Ultrashort-Chain PFAS Analysis
- Conclusions

# Polar X – A Novel Phase Chemistry



## Polar X




2.7  $\mu\text{m}$  superficially porous particle (SPP) silica

- A Single ligand capable of **HILIC** and **Ion Exchange** retention
- Proper retention for polar compounds

# Polar X – A Novel Phase Chemistry

## Tune Your Retention to What You Need


### HILIC Retention



● Acidic, Neutral, and Basic Polar Analytes  
— Silica Surface




### Ion Exchange Retention



# Polar X – A Novel Phase Chemistry

Balanced Retention of Non-Polar, Polar, positively-Charged, and Negatively-Charged Compounds



# Outline

- The Raptor Polar X Column
- **Glyphosate and Polar Contaminants Analysis**
- Novel Solution for Ultrashort-Chain PFAS Analysis
- Conclusions





# Glyphosate and Polar Contaminants Analysis in Waters

- Glyphosate is a herbicide and crop desiccant, and main ingredient in the broadleaf weed killer Roundup.
- By volume, it is one of the most widely used herbicides in the world.
- In 2015 the WHO International Agency for Research on Cancer classified glyphosate as a probable carcinogen in humans


# Glyphosate and Polar Contaminants Analysis in Waters

## Analytical Challenges of Glyphosate Analysis

- Analysis of **glyphosate** in water often include its main metabolite, **AMPA**, and another commonly used herbicide **glufosinate**. The target detection limit is at low part-per-trillion (ppt) levels.
- They are small and very polar organophosphorus compounds and therefore pose a significant chromatographic challenge especially with conventional reversed-phase approaches.
- In addition to the chromatographic challenges, their phosphate functional groups makes them prone to metal chelation, which makes them susceptible to interact with any active metal in the sample flow path.
- Chelation will adversely affect peak shape and signal response regardless of the separation performance of the LC stationary phase.



Aminomethylphosphonic acid (AMPA)






# Glyphosate and Polar Contaminants Analysis in Waters

## New Hardware & System Passivation

- To ensure that Raptor Polar X Columns and Guard Column Cartridges are **Ready to Go** right out of the box, the columns have been treated with a passivation solution using a proprietary manufacturing process.
- However because most LCs have active metal in the flow path beyond the column, Restek developed a **LC Passivation Solution** to allow customers to quickly and easily passivate their LC instrument to enable them to achieve optimal analytical performance.

# Glyphosate and Polar Contaminants Analysis in Waters

## Simple System Passivation



Snap open the Ampule of Restek's LC Passivation Solution

Transfer to a Vial




Inject several times through your injector and tubing but not including your column. Just let the flow goes to waste.



MS/MS

Set your mobile phase flow rate to 0.4 mL/min and mobile phase composition to 50%. Divert your flow to waste so it does not enter your detector, and then make 5 to 10 full loop injections of the passivation solution. After the final full loop injection, allow the mobile phase to continue to flow for 1 minute to allow the passivation solution to flow completely through your LC flow path.

# Glyphosate and Polar Contaminants Analysis in Waters



Polar X: 2.7 $\mu$ m 30x2.1 mm

### Method Conditions:

|                |                                  |    |
|----------------|----------------------------------|----|
| MPA:           | 0.5% formic acid in water        |    |
| MPB:           | 0.5% formic acid in acetonitrile |    |
| Flow:          | 0.5 mL/min                       |    |
| Gradient:      | Time                             | %B |
|                | 0.0                              | 65 |
|                | 5.0                              | 10 |
|                | 6.51                             | 65 |
|                | 8.0                              | 65 |
| Oven Temp:     | 35 °C                            |    |
| Injection Vol: | 5 $\mu$ L                        |    |

100 ppb in water




Pure Chromatography

[www.restek.com](http://www.restek.com)

# Glyphosate and Polar Contaminants Analysis in Waters

**500  $\mu$ L injection of 20 ppt**



**Polar X: 2.7 $\mu$ m 30x2.1 mm**

| Method Conditions: |                                  |    |
|--------------------|----------------------------------|----|
| MPA:               | 0.5% formic acid in water        |    |
| MPB:               | 0.5% formic acid in acetonitrile |    |
| Flow:              | 0.5 mL/min                       |    |
| Gradient:          | Time                             | %B |
|                    | 0.0                              | 95 |
|                    | 1.0                              | 95 |
|                    | 3.0                              | 10 |
|                    | 9.0                              | 10 |
|                    | 10.5                             | 95 |
| Oven Temp:         | 30 °C                            |    |
| Injection Vol:     | 500 $\mu$ L                      |    |




Pure Chromatography

[www.restek.com](http://www.restek.com)

# Glyphosate and Polar Contaminants Analysis in Waters

Calibration Curves for Polar Anionic Herbicides in Water  
(Isotopically Labeled Internal Standards Used)




# Glyphosate and Polar Contaminants Analysis in Waters

## Glyphosate in Bottled Drinking Water

**Column** Raptor Polar X (cat.# 9311A32)  
**Dimensions:** 30 mm x 2.1 mm ID  
**Particle Size:** 2.7  $\mu\text{m}$   
**Pore Size:** 90  $\text{\AA}$   
**Temp.:** 30  $^{\circ}\text{C}$   
**Sample** Glyphosate (cat.# 32426)  
**Diluent:** Direct injection of bottled drinking water fortified at 100 ppt with glyphosate.  
**Inj. Vol.:** 500  $\mu\text{L}$   
**Mobile Phase**  
A: Water, 0.5% formic acid  
B: Acetonitrile, 0.5% formic acid

| Time (min) | Flow (mL/min) | %A | %B |
|------------|---------------|----|----|
| 0.00       | 0.5           | 5  | 95 |
| 1.0        | 0.5           | 5  | 95 |
| 3.0        | 0.5           | 90 | 10 |
| 9.0        | 0.5           | 90 | 10 |
| 9.01       | 0.5           | 5  | 95 |
| 10.5       | 0.5           | 5  | 95 |

**Detector** MS/MS  
**Ion Mode:** ESI-  
**Instrument** UHPLC




# Glyphosate and Polar Contaminants Analysis in Waters

Expanded Panel of 17 Polar Contaminants

1. AMPA
2. Bialaphos
3. Perchlorate
4. Glufosinate
5. MPPA
6. TFA
7. HEPA
8. Difluoroacetic acid
9. Chlorate
10. Ethephon
11. Glyphosate
12. Bromide
13. Bromate
14. N-acetyl AMPA
15. Fosetyl Aluminum
16. Phosphonic Acid
17. N-acetyl glufosinate

Polar X: 2.7µm 30x2.1 mm

| Method Conditions: |                                  |    |
|--------------------|----------------------------------|----|
| MPA:               | 0.5% formic acid in water        |    |
| MPB:               | 0.5% formic acid in acetonitrile |    |
| Flow:              | 0.5 mL/min                       |    |
| Gradient:          | Time                             | %B |
|                    | 0.0                              | 65 |
|                    | 5.0                              | 10 |
|                    | 11.5                             | 10 |
|                    | 11.51                            | 65 |
|                    | 13                               | 65 |
| Oven Temp:         | 35 °C                            |    |
| Injection Vol:     | 1 µL                             |    |



Pure Chromatography

www.restek.com

# Outline

- The Raptor Polar X Column
- Glyphosate and Polar Contaminants Analysis
- Novel Solution for Ultrashort-Chain PFAS Analysis
- Conclusions






## Novel Solution for Ultrashort-Chain PFAS Analysis

- Ubiquitous presence of ultrashort-chain PFAS in environmental waters has become a major concern in parallel to long-chain PFAS contamination.
- The highly polar ultrashort-chain PFAS poses a challenge to the current reversed-phase approach to PFAS chromatographic analysis.

# Novel Solution for Ultrashort-Chain PFAS Analysis



Detected in different water matrices (snow, rain, river, and tap water)

1. PFPrA is up to 45% of the total PFAS in rain and snow samples (USA, France, Japan)
2. C2-C3 accounted for >40% of the total PFAS in rain samples (Canada)
3. TFA is the most abundant PFAS in rainwater (Japan)
4. PFEtS and PFPrS are detected in aqueous film-forming foam (AFFF)

Analytical methods for investigation on the sources and environmental levels

# Novel Solution for Ultrashort-Chain PFAS Analysis

## Legacy PFAS

C8




Perfluorooctanoic acid (PFOA)




Perfluorooctanesulfonic acid (PFOS)

## PFOA and PFOS Alternatives

### Perfluoroalkyl ether carboxylic acids (PFECAs)



Hexafluoropropylene oxide dimer acid (HFPO-DA)




ammonium 4,8-dioxa-3H-perfluorononanoate (ADONA)

### Polyfluoroalkyl ether Sulfonates (PFESAs)




F53-B  
(9-chlorohexadecafluoro-3-oxanonane-1-sulfonate)  
(9Cl-PF3ONS)




F53-B  
(11-chloroeicosafluoro-3-oxanonane-1-sulfonate)  
(11Cl-PF3OUDs)

## UltraShort & Short-Chain PFAS




Trifluoroacetic acid (TFA)

C2




Perfluoroethane sulfonate (PFES)




Perfluoropropionic acid (PFPrA)

C3




Perfluoropropane sulfonate (PFPrS)




Perfluorobutanoic acid (PFBA)

C4




Perfluorobutane sulfonate (PFBS)




Perfluorohexanoic acid (PFHxA)

C6



Perfluorohexane sulfonate (PFHxS)


# Novel Solution for Ultrashort-Chain PFAS Analysis



**Polar X: 2.7 $\mu$ m 50x2.1 mm**

| LC Conditions : (Waters Acquity UPLC) |                                                   |
|---------------------------------------|---------------------------------------------------|
| Mobile Phase A                        | 10mM ammonium formate, 0.05% formic acid in water |
| Mobile Phase B                        | 0.05% formic acid in 60:40 acetonitrile:methanol  |

| Gradient     | Time (min) | %B |
|--------------|------------|----|
|              | 0.00       | 85 |
|              | 5.00       | 85 |
| Injection    | 10 $\mu$ L |    |
| Flow Rate    | 0.5 mL/min |    |
| Run Time     | 5 min      |    |
| Column Temp. | 40°C       |    |



400 ppt in 50:50 water:methanol

# Novel Solution for Ultrashort-Chain PFAS Analysis

*Thursday, August 20*

**Polyfluoroalkyl Substances (PFAS) in the Environment – Session 3**

**2:00 pm**

*“Novel Stationary Phase for Simultaneous Analysis of Ultrashort-Chain (C2, C3),  
Alternative, and Legacy PFAS”*



Pure Chromatography

[www.restek.com](http://www.restek.com)

# Comprehensive & Simultaneous Analysis of Different Classes of PFAS on the Raptor Polar X


*Perfluoroalkylcarboxylic Acid*

**Polar X: 2.7 $\mu$ m 50 x 2.1 mm**


**LC Conditions : (Waters Acquity UPLC)**

Mobile Phase A 10mM ammonium formate, 0.05% formic acid in water

Mobile Phase B 0.05% formic acid in 60:40 acetonitrile:methanol



*Polyfluoroalkyl Ether Sulfonic Acid*



*Perfluoroalkyl Ether Carboxylic Acid*

| Gradient     | Time (min)       | %B |
|--------------|------------------|----|
|              | 0.00             | 95 |
|              | 4.00             | 65 |
|              | 4.01             | 95 |
|              | 6.00             | 95 |
| Injection    | 10 $\mu$ L       |    |
| Diluent      | 50/50 Water/MeOH |    |
| Flow Rate    | 0.5 mL/min       |    |
| Run Time     | 6 min            |    |
| Column Temp. | 40°C             |    |




Pure Chromatography

www.restek.com

# Comprehensive & Simultaneous Analysis of Different Classes of PFAS on the Raptor Polar X

**Polar X: 2.7 $\mu$ m 50 x 2.1 mm**



**RESTEK**

Pure Chromatography

[www.restek.com](http://www.restek.com)

## Conclusions

- The Raptor Polar X is a unique stationary phase specifically designed for balanced retention and selectivity of a wide range of polar compounds.
- HILIC methods are simplified by rapid equilibration that increases productivity.
- Proper chromatographic retention allows the development of fast and simple LC-MS/MS method for high-throughput analysis of glyphosate and ultrashort-chain PFAS in waters.
- Sensitive glyphosate analysis in water can be achieved by direct large volume injection on the Polar X column.
- Ultrashort-chain and long-chain PFAS can be simultaneously analyzed on the Polar X column.



*Thanks for Your Attention*



Pure Chromatography

[www.restek.com](http://www.restek.com)