

# **Microplastics analysis – automatic analysis by Dynamic Particle Imaging Analysis.**

Ruth Marfil-Vega, PhD

# Acknowledgments

## Molecular Spectroscopy Team



Liang Zhao  
Product Specialist



Sudhir Dahal  
Product Manager

## Physical Testing Team



William Gabler  
Sr. Applications Scientist



Andrew Lim  
Product Specialist

# In today's presentation



1. Microplastics analysis – ASTM
2. Analysis by Infrared (IR) spectroscopy
3. Analysis by Dynamic Imaging Analysis (DIA)
4. IR spectroscopy vs Dynamic Particle Imaging
5. Conclusions
6. Q&A

ASTM

IR

DIA

IR vs DIA

Conclusions

Q&A

# ASTM

- ❑ 6 Work Items
- ❑ Focus on:
  - ❑ *Sample collection*
  - ❑ *Sample preparation*
  - ❑ *Sample analysis by*
    - ❑ IR
    - ❑ Raman
    - ❑ Pyr-GCMS
    - ❑ Dynamic Particle Imaging

## Focus of my presentation

- ❑ Current data gaps WK67565 (IR)
- ❑ Initial developments WK72349 (Dynamic Particle Imaging)

- ❑ WK67563 - New Practice for **Collection** of All Matrices of High and Low Turbidity Water Samples for the Identification and Quantification of Micro-plastic Particles and Fibers.
- ❑ WK67564 - New Practice for **Preparation of Water Samples** with High, Medium or Low Suspended Solids for the Identification and Quantification of Microplastic Particles and Fibers using **Raman Spectroscopy, IR Spectroscopy or Pyrolysis-Gas Chromatography/Mass Spectroscopy (GC/MS)**.
- ❑ WK67565 - Spectroscopic Identification and Quantification of Microplastic Particles and Fibers in all High and Low Turbidity Water Matrices including Municipal Wastewater Using **IR and Raman Spectroscopy**.
- ❑ WK67788 Identification of Microplastic Particles and Fibers in All Matrices of High and Low Turbidity Water Samples including Municipal Raw Wastewater using **Pyrolysis-GC/MS**.
- ❑ WK70831 The Production and Quantification of Microplastic Particles to be used in **Quality Control Reference Samples** in the Collection Practices, Preparation Practices and Identification Analysis Methods used on All Matrices of High and Low Turbidity Waters.
- ❑ WK72349 Determination of Microplastics Particle and Fiber Size, Distribution, Shape and Concentration in Waters with High to Low Suspended Solids Using a **Dynamic Image Particle Size and Shape Analyzer**.

ASTM

IR


DIA

IR vs DPI

Conclusions

Q&A

# Analysis of Microplastics



Variety of complimentary techniques available.

# How does IR work?

- ❑ An IR spectrometer measures the absorption of IR radiation (wavelengths  $\sim 1 \mu\text{m}$ -  $\sim 40 \mu\text{m}$ ) by molecules in various states. Modern FTIR spectrometers often use Fourier transformed interferograms to generate IR spectra.
- ❑ In an IR spectrum the energy of IR radiation is usually represented in unit of reciprocal centimeter ( $\text{cm}^{-1}$ ). Typically a Mid IR spectrum records from about  $400 \text{ cm}^{-1}$  to  $4000 \text{ cm}^{-1}$
- ❑ Peaks in an IR spectrum reflect specific vibrational/rotational modes in molecules, therefore, can be used for qualitative/quantitative analysis of chemical species.
- ❑ FTIR microscopes can be applied to investigate samples on the microscopic scale. There are three operational modes: transmission, reflection or attenuated total reflection (ATR).

ASTM

IR

DIA

IR vs DIA

Conclusions

Q&amp;A

# Experimental Conditions

| IR  |                                     |
|---|-------------------------------------|
| <b>Detector</b>                           | Liquid nitrogen cooled MCT detector |
| <b>Objective mirror</b>                   | 15x reflection objective mirror     |
| <b>Optical system mode</b>                | Reflection and ATR                  |
| <b>Resolution (cm<sup>-1</sup>)</b>       | 8                                   |
| <b>No. of Scans averaged per spectrum</b> | Variable, 64 (typical)              |
| <b>Intensity mode</b>                     | Absorbance                          |
| <b>Range (cm<sup>-1</sup>)</b>            | 700-4000                            |
| <b>Mirror speed (mm/s)</b>                | 9                                   |
| <b>Apodization</b>                        | HappGenzel                          |

| Sample ID | True diameter (μm) | Material         |
|-----------|--------------------|------------------|
| 10 μm     | 10.23              | Polystyrene      |
| 20 μm     | 20.10              | Polystyrene      |
| 50 μm     | 50.56              | Polymethacrylate |
| 100 μm    | 100                | Polystyrene      |

| Filter Type          | Description                                    |
|----------------------|--|
| Stainless steel mesh | 5 μm pore size                                 |
| Cellulose ester      | 25 mm diameter, 0-6 μm pore size               |
| Glass fiber          | 25 mm diameter, grade 934-AH, 1.5 μm pore size |

ASTM

IR


DIA

IR vs DIA

Conclusions

Q&A

# Results: Background and importance of filter material




- ❑ Analysis performed in reflection mode
- ❑ Spectra from cellulose and glass fiber filter presents absorbance across the measuring region
- ❑ Stainless steel is highly reflective and does not show absorbance
- ❑ Qualitative results demonstrates the negative impact of filter background on sensitivity

*Replication sample analysis: n=5*




# Results: Mapping and its value

100  $\mu\text{m}$



10  $\mu\text{m}$



- ❑ Analysis performed in reflection mode
- ❑ Qualitative initial screening:
  - ❑ *Identify areas of interest*
  - ❑ *General assessment of size distribution*
- ❑ Sizes of particles from 100 to 20  $\mu\text{m}$  were unequivocally determined
- ❑ Reflection from neighboring particles interfered in the analysis of 10  $\mu\text{m}$

*Replication sample analysis: n=1*

ASTM

IR


DIA

IR vs DIA

Conclusions

Q&A

## Results: ATR vs reflection mode



- ❑ ATR mode presents higher specificity than reflection mode
- ❑ Sizes of particles from 100 to 20 μm were unequivocally determined
- ❑ 10 μm particles were also identified
- ❑ Increase in sensitivity when using ATR mode was qualitatively demonstrated

*Replication sample analysis ATR: n=3*

ASTM

IR



DIA

IR vs DIA

Conclusions

Q&A


# Results: Quantitation of detection threshold



- ❑ Analysis performed in reflection mode (“worse” mode)
- ❑ Decreasing particle size mimicked by decreasing aperture size from 100x100 µm to 5x5 µm
- ❑ Results at 10x10 µm:
  - ❑  $S/N \geq 3$
  - ❑ *Relevant increase in S/N when increasing # of scans to 256*


*Replication sample analysis: n=4*

# How does DIA work?



# How does DIA work?


\* Measurement Parameters (13 Types)  
 Area circle equivalent diameter, perimeter equivalent diameter, maximum length, maximum perpendicular length, vertical Feret diameter, horizontal Feret diameter, particle perimeter, envelope perimeter, circularity, aspect ratio, horizontal bounding rectangle aspect ratio, particle area, average brightness



Number concentration  
(particles/mL)

Size distribution

Particle image



# Experimental Conditions

|                        |  |
|------------------------|--|
| Measurement Size Range | 5 - 100 $\mu\text{m}$  |
| Sample Volume          | 50 - 1000 $\mu\text{L}$  |
| Flow Cell Size         | 120 $\mu\text{m}$ deep x 1 mm wide   |
| Analysis Flow Rate     | 0.1 $\mu\text{L}/\text{min}$<br>(30 s – 10 min analysis)                                     |
| Measurements           | Particle size, distribution, count, concentration, shape, morphology, brightness (grayscale) |

| Sample ID             | True diameter ( $\mu\text{m}$ ) | Material         | Analysis time (frames) | Output filters (exclude <) |
|-----------------------|---------------------------------|------------------|------------------------|----------------------------|
| 10 $\mu\text{m}$      | 10.23                           | Polystyrene      | 5 min (2,100)          | <8 $\mu\text{m}$           |
| 20 $\mu\text{m}$      | 20.10                           | Polystyrene      | 5 min (2,100)          | <10 $\mu\text{m}$          |
| 50 $\mu\text{m}$      | 50.56                           | Polymethacrylate | 5 min (2,100)          | <10 $\mu\text{m}$          |
| 100 $\mu\text{m}$     | 100                             | Polystyrene      | 2 min (1,200)          | <10 $\mu\text{m}$          |
| NIST-20 $\mu\text{m}$ | 20.06                           | Polystyrene      | 5 min (1,800)          | <10 $\mu\text{m}$          |

ASTM

IR

**DIA**

IR vs DIA



Conclusions

Q&A




# Example Data Output

|  |                 |
|--|-----------------|
| File Name :                                | 06122020_50um_3 |
| Sample ID:                                 | 06122020_100_3  |
| Sample No.:                                |                 |
| Number of Frames: 2100                     |                 |
| Number of Particles: 2069                  |                 |
| Concentration(count/mL): 4899              |                 |
| Condition of Classification:<br>10.000 ≤ X |                 |
| True                                       |                 |
| Number of Particles: 2069 (100.000 %)      |                 |
| Concentration(count/mL): 4899.132          |                 |
| False                                      |                 |
| Number of Particles: 0 (0.000 %)           |                 |
| Concentration(count/mL): 0.000             |                 |
| Median Value(μm):                          | 50.646          |
| Modal Value(μm):                           | 51.400          |
| Mean Value(μm):                            | 47.814          |
| Standard Deviation:                        | 10.861          |
| Coefficient of Variation:                  | 0.227           |
| Value of Cumulative %                      |                 |
| 50.000 % (μm) :                            | 50.646          |
| Cumulative % of Interval                   |                 |
| 40.000- 60.000 μm (%) :                    | 90.875          |
| Comment                                    |                 |
| 500μL in 100mL stirred                     |                 |




Sample contamination



Particle agglomeration



## Results: Sizing particles



- ❑ Individual samples prepared with each microplastics standard
- ❑ Mean size calculated from all particles present in the sample
- ❑ Accurate and precise determination of particle size from 10 to 100 μm:
  - ❑ 10.6 μm ± 1.2%
  - ❑ 20.7 μm ± 0.2%
  - ❑ 49.8 μm ± 3.4%
  - ❑ 97.0 μm ± 5.2%

*Replication sample analysis: n=3*

ASTM

IR

DIA

IR vs DIA


Conclusions

Q&A



# Results: Counting particles


- ❑ Individual samples prepared with each microplastics standard
- ❑ Particles present in the sample counted
- ❑ Reproducible counting results from different size particles.
  - ❑  $25654 \pm 4.9\%$
  - ❑  $1865 \pm 11.3\%$
  - ❑  $2007 \pm 15.6\%$
  - ❑  $196 \pm 56.5\%$



*Replication sample analysis: n=3*



# Results: Impact of variable brightness conditions





NIST Traceable standard 20  $\mu\text{m}$ , acceptable concentration range: 2,430 – 3,630 particles/mL

- Light source intensity can be adjusted to improve brightness and contrast in dark or turbid solutions

High Brightness



Low Brightness



- Results (size and concentration) accurate, precise and reproducible under both brightness conditions.

*Replication sample analysis: n=4*

# Results: Analysis of simulated microplastics solution






Average diameter: 7.737  $\mu\text{m}$   
Concentration: 7141 count/mL  
(2585 particles observed)



# Results: Analysis of simulated microplastics solution

- Data can be processed based on particle brightness, in addition to particle size.



ASTM

IR

**DIA**

IR vs DIA

Conclusions

Q&A

# IR & DIA - applicability



|  | IR  | DIA  |
|--|---|--|
| Sample Matrix  | Filtered liquid sample on filter or dried sample on microscope slide  | Liquid sample (prep steps feasible)  |
| Minimum particle size (µm) with current measuring conditions | 10  | 10   |
| Information obtained   | <ul style="list-style-type: none"> <li>• Size, shape and agglomeration</li> <li>• Plastic type after entering unique spectral characteristics parameters in software</li> </ul> | <ul style="list-style-type: none"> <li>• Size, shape and agglomeration</li> <li>• Results can be filtered for flexible data analysis and stats analysis (ex. particle distribution)</li> </ul> |



# Conclusions

- ❑ Detection threshold for particle size by IR (in reflection and ATR modes) was qualitatively determined: 10  $\mu\text{m}$ .
- ❑ Stainless steel filters demonstrated to be the most suitable ones for IR analysis.
- ❑ Detection threshold for particle size by IR in reflection mode (less sensitive than ATR) was quantitatively determined: 10  $\mu\text{m}$ , with  $S/N \geq 3$ . Number of scans  $>124$  increased S/N.
- ❑ Initial demonstration of performance by DIA was completed:
  - ❑ *Accurate and precise determination of particle size from 10 to 100  $\mu\text{m}$ .*
  - ❑ *Reproducible counting results from different size particles.*
  - ❑ *Results (size and concentration) accurate, precise and reproducible under variable brightness conditions.*

ASTM

IR

DIA

IR vs DIA

Conclusions

Q&A



**For more information, contact me at:  
rmmarfilvega@shimadzu.com  
Phone number: 410-910-0884**



**Visit [www.OneLabOneEarth.com](http://www.OneLabOneEarth.com)**

ASTM

IR

DIA

IR vs DIA

Conclusions

**Q&A**