



Per- And Polyfluoroalkyl Substances

Novel Stationary Phase for Simultaneous Analysis of Ultrashort-Chain (C2, C3), Alternative, and Legacy PFAS



Shun-Hsin Liang, Ph.D
Senior Scientist, LC Applications
Shun-HsinLiang@Restek.com




Pure Chromatography

www.restek.com


Outline

- PFAS Panel
- Ultrashort-Chain PFAS
- Direct Injection Method Evaluation for Ultrashort-Chain/Alternative/Legacy PFAS (Polar X Column)
- Comprehensive PFAS Analysis (Polar X Column)
- Conclusions


PFAS (Per- and Polyfluoroalkyl Substances)



Perfluoroalkylcarboxylic Acid
(PFCA)




Perfluoroalkylsulfonic Acid
(PFSA)




X:2 Telomer Sulfonic Acid
(X:2 FTS)

Surfactant


Processing aid
Mist suppressant
Fire fighting foam
Cleaning Products




Perfluorooctanesulfonamide
(FOSA)



Perfluorooctanesulfonamido-
ethanol (FOSE)



Perfluorooctanesulfonamido-
acetic acid (FOSAA)



Telomer Carboxylic Acid
(X:2 FTA)

Raw Material

*Environmental Transformation
Product*

PFAS Testing Standards

Test Method	US EPA 537.1	ISO 25101:2009	DIN 38407-42	ASTM D7968-17a	ASTM D7979-17	US EPA 533	US EPA 8327	ISO 21675
Sample Matrix	Drinking Water	All water types	All water types	soil	All water types (- drinking water)	Drinking Water	Non-potable water	All water types
# of Analytes	18	24	2	21	21	25	24	30
Sample Prep	SPE	Direct injection	SPE	Direct injection	Direct injection	SPE	Direct injection	SPE
Sample Volume	250 mL	1000 mL	50 mL	2 g	5 mL	250 mL	5 mL	50 – 1000 mL*
Detection limit	Optional	Not shown	<ul style="list-style-type: none"> 0.01 ug/L 0.025 ng/L for treated waste water 	MDL (2.41 – 258.37 ng/kg)	MDL (0.7 – 106.8 ng/L)	LCMRL (1.7 – 20 ng/L)	<ul style="list-style-type: none"> MDL (0.7 – 4.6 ng/L) LLOQ is 10 ng/L 	LOQ: 0.2 ng/L (loose term)




Pure Chromatography


www.restek.com

PFOA and PFOS Alternatives

Perfluoroalkyl ether carboxylic acids (PFECAs)




Hexafluoropropylene oxide dimer acid
(HFPO-DA) GenX




ammonium 4,8-dioxa-3H-perfluorononanoate
(ADONA)

Polyfluoroalkyl ether sulfonates (PFESAs)




F-53B
(9-chlorohexadecafluoro-3-oxanonane-1-sulfonate)
(9Cl-PF3ONS)



F-53B
(11-chloroeicosafluoro-3-oxanonane-1-sulfonate)
(11Cl-PF3OUdS)

Updated 537.1 Method (Drinking Water)




Analyte	Peak # (Fig. 1)	RT (min)	IS# Ref
PFBS	1	7.62	2
PFHxA	2	10.42	1
HFPO-DA	4	11.38	1
PFHpA	6	13.40	1
PFHxS	7	13.58	2
ADONA	8	13.73	1
PFOA	9	15.85	1
PFOS	11	17.91	2
PFNA	13	17.92	1
9CI-PF3ONS	14	18.91	2
PFDA	15	19.69	1
NMeFOSAA	17	20.50	3
PFUnA	19	21.21	1
NEtFOSAA	20	21.26	3
11CI-PF3OUDs	22	21.84	2
PFDoA	23	22.52	1
PFTrDA	24	23.66	1
PFTA	25	24.64	1
¹³ C ₂ -PFHxA	3	10.42	1
¹³ C ₃ -HFPO-DA	5	11.40	1
¹³ C ₂ -PFDA	16	19.69	1
d ₅ -NMeFOSAA	21	21.24	3
¹³ C ₂ -PFOA-IS#1	10	15.85	-
¹³ C ₄ -PFOS-IS#2	12	17.91	-
d ₅ -NMeFOSAA-IS#3	18	20.49	-

Outline

- PFAS Panel
- Ultrashort-Chain PFAS
- Direct Injection Method Evaluation for Ultrashort-Chain/Alternative/Legacy PFAS (Polar X Column)
- Comprehensive PFAS Analysis (Polar X Column)
- Conclusions

Ultrashort-Chain PFAS



Detected in different water matrices (snow, rain, river, and tap water)

1. PFPrA is up to 45% of the total PFAS in rain and snow samples (USA, France, Japan)
2. C2-C3 accounted for >40% of the total PFAS in rain samples (Canada)
3. TFA is the most abundant PFAS in rainwater (Japan)
4. PFEtS and PFPrS are detected in aqueous film-forming foam (AFFF)


Analytical methods for investigation on the sources and environmental levels

Reversed-Phase LC/MS Method


Raptor C18
2.7 μ m, 100x3.0mm

LC Conditions : (Shimadzu Nexera X2)		
Mobile Phase A	5mM ammonium acetate in water	
Mobile Phase B	methanol	
Gradient	Time (min)	%B
	0.00	20
	7.00	95
	9.00	95
	9.01	20
	11.00	20
Injection	10 μ L	
Flow Rate	0.25 mL/min	
Run Time	11 min	
Column Temp.	40°C	

Spiked Reagent Water
(80 ppt)



Spiked Groundwater
(80 ppt)



RESTEK®

Pure Chromatography

www.restek.com

Outline

- PFAS Panel
- Ultrashort-Chain PFAS
- Direct Injection Method Evaluation for Ultrashort-Chain/Alternative/Legacy PFAS (Polar X Column)
- Comprehensive PFAS Analysis (Polar X Column)
- Conclusions

Direct Injection Method Evaluation



Polar X



- A Single ligand capable of **HILIC** and **Ion Exchange** retention
- Proper retention for polar compounds

RESTEK


Pure Chromatography

www.restek.com

Direct Injection Method Evaluation


Legacy PFAS

C8




PFOA and PFOS Alternatives

Perfluoroalkyl ether carboxylic acids (PFECAs)




Polyfluoroalkyl ether Sulfonates (PFESAs)




UltraShort & Short-Chain PFAS

C2




Perfluoroethane sulfonate (PFEtS)

C3




Perfluoropropane sulfonate (PFPrS)

C4




Perfluorobutane sulfonate (PFBS)

C6



Perfluorohexane sulfonate (PFHxS)

Direct Injection Method Evaluation



Polar X: 2.7 μ m 50x2.1 mm

LC Conditions : (Waters Acquity UPLC)

Mobile Phase A	10mM ammonium formate, 0.05% formic acid in water
Mobile Phase B	0.05% formic acid in 60:40 acetonitrile:methanol

Gradient	Time (min)	%B
	0.00	85
	5.00	85
Injection	10 μ L	
Flow Rate	0.5 mL/min	
Run Time	8 min	
Column Temp.	40°C	

400 ppt in 50:50 water:methanol

RESTEK®

Pure Chromatography

www.restek.com

Direct Injection Method Evaluation

*Precision & Accuracy
of Fortified Water Samples:*

1. Tap water
2. River water (Chicago)
3. Groundwater
4. POTW water (Effluent)



Direct Injection Method Evaluation

Analyte	Precursor Ion	Product Ion	IS for Quantification
TFA	113.03	69.01	$^{13}\text{C}_2\text{-PFHxA}$
PFPrA	163.03	119.01	$^{13}\text{C}_2\text{-PFHxA}$
PFBA	212.97	168.97	$^{13}\text{C}_2\text{-PFHxA}$
PFHxA	312.97	268.90	$^{13}\text{C}_2\text{-PFHxA}$
PFOA	412.90	368.91	$^{13}\text{C}_2\text{-PFOA}$
HFPO-DA	284.97	168.92	$^{13}\text{C}_2\text{-PFOA}$
ADONA	376.90	250.93	$^{13}\text{C}_2\text{-PFOA}$
PFEtS	198.87	79.92	$^{13}\text{C}_3\text{-PFBS}$
PFPrS	248.97	79.98	$^{13}\text{C}_3\text{-PFBS}$
PFBS	298.97	79.97	$^{13}\text{C}_3\text{-PFBS}$
PFHxS	398.90	79.97	$^{13}\text{C}_3\text{-PFBS}$
PFOS	498.84	79.97	$^{13}\text{C}_4\text{-PFOS}$
9CI-PF3ONS	530.78	350.85	$^{13}\text{C}_4\text{-PFOS}$
11CI-PF3OUDS	630.78	450.80	$^{13}\text{C}_4\text{-PFOS}$
$^{13}\text{C}_2\text{-PFHxA}$	314.97	269.93	-
$^{13}\text{C}_2\text{-PFOA}$	414.90	369.87	-
$^{13}\text{C}_3\text{-PFBS}$	301.90	79.97	-
$^{13}\text{C}_4\text{-PFOS}$	502.84	79.97	-

Direct Injection Method Evaluation

Polar X Column
2.7 μ m, 50x2.1mm

Sample Preparation:



(polypropylene vial)

- 250 μ L water sample or standard
- +
250 μ L methanol
- +
5 μ L internal standard

(10 ng/mL $^{13}C_2$ -PFHxA, $^{13}C_2$ -PFOA, $^{13}C_3$ -PFBS, $^{13}C_4$ -PFOS in methanol)

Direct Injection Method Evaluation

Quadratic regression (1/x weighted)


Calibration Range:

10 – 800 ng/L

20 – 800 ng/L (for TFA)

r^2 value > 0.996

deviations <20%



Pure Chromatography

www.restek.com


Direct Injection Method Evaluation

Analytes in Blank Water Samples

Samples	Detected Concentration (ng/L)														
	TFA	PFPRA	PFBA	PFHxA	PFOA	HFPO-DA	ADONA	PFEtS	PFPrS	PFBS	PFHxS	PFOS	9Cl-PF3ONS	11Cl-PF3OUDs	
Tap Water	164.2	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
River Water	193.3	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND	ND
Groundwater	1425.2	ND	ND	ND	5.4	ND	ND	ND	ND	6.7	3.9	ND	ND	ND	ND
POTW Water	352.8	9.6	15.3	93.5	20.4	ND	ND	ND	ND	6.8	6.7	9.6	ND	ND	ND

Direct Injection Method Evaluation

PFAS in Blank POTW Water



RESTEK®

Pure Chromatography



www.restek.com

Direct Injection Method Evaluation

Accuracy

% Recovery:



70 – 130 %



Precision

% RSD:

<20 %



TFA Contamination in Reagent Solvent




Outline


- PFAS Panel
- Ultrashort-Chain PFAS
- Direct Injection Method Evaluation for Ultrashort-Chain/Alternative/Legacy PFAS (Polar X Column)
- Comprehensive PFAS Analysis (Polar X Column)
- Conclusions

Comprehensive & Simultaneous Analysis of Different Classes of PFAS on the Raptor Polar X

Perfluoroalkylcarboxylic Acid



Polar X: 2.7 μ m 50 x 2.1 mm



LC Conditions : (Waters Acquity UPLC)

Mobile Phase A	10mM ammonium formate, 0.05% formic acid in water	
Mobile Phase B	0.05% formic acid in 60:40 acetonitrile:methanol	

Gradient	Time (min)	%B
	0.00	95
	4.00	65
	4.01	95
	6.00	95
Injection	10 μ L	
Diluent	50/50 Water/MeOH	
Flow Rate	0.5 mL/min	
Run Time	6 min	
Column Temp.	40°C	


RESTEK

Pure Chromatography

www.restek.com

Comprehensive & Simultaneous Analysis of Different Classes of PFAS on the Raptor Polar X

Polar X: **2.7 μ m 50 x 2.1 mm**



RESTEK

Pure Chromatography

www.restek.com

Conclusions

- Unique stationary phase provides proper chromatographic retention of small, polar ultrashort-chain PFAS.
- Fast and simple isocratic LC-MS/MS method allows high-throughput PFAS analysis in potable and non-potable waters.
- This method is suitable for labs interested in adding ultrashort-chain compounds to an existing PFAS assay.
- Multi-classes of PFAS can be simultaneously analyzed on the Polar X column.

Acknowledgement

LC R&D:

Xiaoning Lu
Vernon C. Bartlett
Connor Flannery
Ahren Green
Terry S. Reid

EPA (Region 5):

Larry Zintek



Pure Chromatography

www.restek.com

Thanks for Your Attention



Pure Chromatography

www.restek.com