

# **Supercritical Fluid Extraction: A Solution for the Extraction of PFAS in Environmental Samples**

Ruth Marfil-Vega, PhD

# Acknowledgments

Shimadzu Scientific Instruments  
Innovation Center



Yuka Fujito

William Hedgepeth

Pennsylvania Department of Environmental  
Protection  
Bureau of Laboratories



Department of Environmental  
Protection

# In today's presentation

1. Supercritical Fluid Extraction
2. Method development for PFAS in fish
3. Results – method performance
4. Results – samples
5. Conclusions
6. Q&A



SFE

Method

Performance


Samples

Conclusions

Q&A

# Supercritical Fluids

Supercritical CO<sub>2</sub> is a fluid state of carbon dioxide where it is held at or above its critical temperature (31.1 °C) and critical pressure (73.8 bar)



- Low viscosity of the mobile phase
- Applicable for Extraction (SFE) and LC separations (SFC)
- Superior to LC for chiral separations
- Green technique as less organic solvents are used

# Supercritical Fluid Configurations

Supercritical Fluid Extraction(SFE)  
pretreatment system  
: **Off-line SFE**



Supercritical Fluid Chromatograph (SFC)  
: **SFC(-MS)**




**On-line SFE-SFC(-MS)**



# Supercritical Fluid Extraction (SFE)

## What is it?

- ❑ SFE is a process that separates a component(s) from a matrix.
- ❑ Typically CO<sub>2</sub> is used to extract components from a solid matrix.
- ❑ Cosolvents like methanol and ethanol may be needed to help extract more polar components



# Benefits

## 1) **Automation:**

- reduces analyst errors and improves reproducibility

## 2) **Green technique:**

- reduces organic solvent use and total volume of waste generated


## 3) **Selectivity:**

- pressures and temperatures can be varied to allow step-wise extraction

## 4) **Speed:**

- supercritical fluid has a faster diffusion into matrix than liquids, reducing extraction times


# SFE vs conventional extraction



Automated sample preparation reduces analyst time and effort



Sorbent for dehydrating samples with high water content



Up 48 samples can be automatically extracted and analyzed by using Rack Changer



# Off-line SFE for “legacy” contaminants



## Off-line supercritical fluid extraction/gas chromatography-mass spectrometry analysis of pesticides in fish

William Hedgepeth<sup>1</sup>, Tairo Ogura<sup>1</sup>, Riki Kitano<sup>1</sup>, Jeff Dahl<sup>1</sup>, June Black<sup>2</sup>

<sup>1</sup>Shimadzu Scientific Instruments, Inc., Columbia MD, <sup>2</sup>Pennsylvania Department of Environmental Protection, Harrisburg PA

### Analytical Conditions

#### Extraction Conditions

Vessel Temp: 50 °C  
 System Pressure: 30.0 MPa  
 Flow rate: 1.0 mL/min 100% CO<sub>2</sub>  
 Static Extraction time: 25 min  
 Dynamic Extraction Time: 30 min  
 Trap Column: Shimadzu C18, 4.6 x 50 mm, 5µm  
 Column Temp: trap 20 °C, elute 50 °C  
 Column Rinse: Hexane, 2 mL/min



#### Analytical Conditions

GC-MS: GCMS-TQ8040 (Shimadzu)  
 Sampler: AOC-20i/s (Shimadzu)  
 Column: SH-Rxi-5MS (30m x 0.25mm I.D., df=0.25µm, Shimadzu)  
 IF Temp.: 290 °C IS Temp.: 230 °C Event Time: 0.3sec  
 Inj. Temp.: 275 °C  
 Flow Control: Linear Velocity, 43.5cm/sec  
 Inj. Mode: Splitless (High Press. Inj., 250kPa, 1.5min)  
 Oven Temp.: 50 °C (0.5min), 28 °C /min to 265 °C,  
 3 °C /min to 285 °C, 25 °C /min to 330 °C (1min)  
 Inj. Volume: 1 µL

Pesticide	Walleye “post-spiked” extract, ppb	Walleye “pre-spiked” extract with corn oil <sup>(*)</sup> , ppb	%Recovery
2,4,5,6-tetrachloro-m-Xylene	48.39	45.28	93.6
<i>alpha</i> -BHC	49.83	40.07	80.4
<i>gamma</i> -BHC	51.46	40.46	78.6
Chlordene	52.2	44.35	85
<i>Heptachlor</i>	47.63	39.07	82
Aldrin	53.21	41.77	78.5
Heptachlor epoxide	54.12	44.14	81.6
<i>trans</i> -Chlordane	49.83	43.21	86.7
<i>cis</i> -Chlordane	52.5	44.99	85.7
<i>trans</i> -Nonachlor	52.47	43.71	83.3
4,4'-DDE	51.25	45.98	89.7
<i>o,p'</i> -DDD	52.44	42.71	81.5
Dieldrin	53.53	44.09	82.4
Endrin	53.36	44.36	83.1
4,4'-DDD	61.15	41.64	68.1
<i>cis</i> -Nonachlor	53.65	41.34	77
4,4'-DDT	44.41	40.88	92.1
<i>Methoxychlor</i>	46.1	45.14	97.9
Mirex	50.74	38.79	76.4
Decachlorobiphenyl	41.74	41.09	98.5

(\*) oil added to improve recoveries from low fat samples

SFE

Method

Performance


Samples

Conclusions


Q&A

# SFE Instrument configuration


## Trapping method



## Direct collection method



# Experimental workflow



Add ground sample



# Instrument conditions

Nexera UC offline SFE	Value
Trap column	None
Mobile phase	A: CO <sub>2</sub> B: MeOH
Modifier concentration	20%
Flow rate	5 mL/min.
Time program	46 min (combination of 3 static and dynamic cycles)
Vessel temperature	60 °C
BPR pressure	20 MPa

LCMS-8050	Value
Column	Shim-pack GIST C18 2.7 um 100 x 2.1 mm XR-ODSII 3 x 75 mm for delay column
Mobile phase	A: 10 mM ammonium acetate in H <sub>2</sub> O B: MeOH
Flow rate	0.5 mL/min.
Gradient program	0 min: 20 %B; 9 min: 90 %B; 11 min: 90 %B; 11.5 min: 20 %B; stop 15 min
Oven temperature	35 °C
Injection volume	1 µL

SFE

Method

Performance


Samples

Conclusions

Q&A



# Method development

## Trapping vs direct collection




✓ Acceptable recoveries (>95%) for all target compounds with direct collection (spike: 50 ng on Hydroprotect)

Direct collection



# Method Development

## Modifier and additive concentration



### Modifiers:

- 0% methanol
- 10% methanol
- 20% methanol
- 40% methanol

### Additives:

- none
- 0.1% formic acid
- 10 mM ammonium acetate
- 5% H<sub>2</sub>O

(Spike: 50 ng on Hydroprotect)

- ✓ Best recoveries (>95%) for all targeted compounds with 20% methanol as modifier.
- ✓ Recoveries were within same range for the additives tested → no additive was selected.

SFE

Method

Performance

Samples

Conclusions

Q&amp;A

# Method Performance


## Calibration – Matrix matched

	Lowest Cal point (LOQ) ng/g	Highest Cal point ng/g	Linearity (R <sup>2</sup> )
PFBS	0.5	100	0.9999
PFHxA	0.5	100	0.9995
HFPO-DA	1.0	100	0.9997
PFHpA	1.0	100	0.9996
PFHxS	0.5	100	0.9999
ADONA	0.5	100	0.9997
PFOA	0.5	100	0.9997
PFNA	0.5	100	0.9997
PFOS	2	100	0.9999
9CI-PF3ONS	1.0	100	0.9995
PFDA	0.5	100	0.9998
N-MeFOSAA	2.0	100	0.9994
N-ETFOSAA	1.0	100	0.9999
PFUnA	1.0	100	0.9997
11CI-PF3OUdS	0.5	100	0.9999
PFDoA	1.0	100	0.9996
PFTriA	2.0	100	0.9997
PFTreA	1.0	100	0.9995

- ✓ Good linearity was observed in the 0.50 to 100 ng/g range with matrix matched calibration. PFAS spiked on 0.5 g of fish and 1 g of Hydroprotect. Unspiked fish analyzed as blank.

# Method Performance

## *Extraction efficiency and reproducibility*




- ✓ Good extraction efficiency (94-116%) and reproducibility were obtained.
  - %RSDs at 100 ng/g:  $\leq 5\%$  (except N-MeFOSA & N-EtFOSA:  $<10\%$ )
  - %RSD at 20 ng/g:  $\leq 10\%$  (except PFTriA: 12%)
  - %RSD at 2 ng/g: from 2% to 27% (except N-MeFOSA: 44%)



# Method Performance

## Matrix Effect



- ✓ No severe matrix effect observed. Recoveries for all target compounds except N-MeFOSAA were within 70-130%

*Spike: 50 ng on 1 g of Hydroprotect without or with 0.5 g of fish. Unspiked fish analyzed as blank.*

SFE

Method

Performance

Samples


Conclusions

Q&amp;A


# Sample analysis

- ✓ PFAS compounds were detected in Walleye and Large Mouth Bass, but not in the farm raised Trout sample.


(a) Wild caught Large Mouth Bass



(b) Wild caught Walleye



(c) Farm raised Trout



	Walleye (ng/g)	Large Mouth Bass (ng/g)	Trout (ng/g)
PFBS	1.00	1.62	n.d.
PFHxA	n.d.	n.d.	n.d.
HFPO-DA	n.d.	n.d.	n.d.
PFHxS	n.d.	n.d.	n.d.
PFHpA	n.d.	n.d.	n.d.
ADONA	n.d.	n.d.	n.d.
PFOA	1.00	1.43	n.d.
PFNA	2.37	1.13	n.d.
PFOS	51.65	77.34	n.d.
9Cl-PF3ONs	0.98	2.68	n.d.
PFDA	6.68	10.52	n.d.
N-MeFOSAA	n.d.	n.d.	n.d.
N-EtFOSAA	n.d.	n.d.	n.d.
PFUnA	5.65	14.24	n.d.
11Cl-PF3OUds	0.68	2.97	n.d.
PFDaA	2.79	4.48	n.d.
PFTriA	4.12	7.28	n.d.
PFTreA	1.39	2.33	n.d.

# Conclusions

- The suitability of SFE as a sample preparation technique for PFAS analysis was demonstrated.
  - *SFE provided excellent results for recovery, linearity, and reproducibility.*
  - *Current extraction time could be shorten (experiments to confirm this are pending).*
- Wild caught fish samples contained several target PFAS above the limit of quantification.
- SFE can be automated, resulting in increased lab productivity.



Contact me at:  
rmmarfilvega@shimadzu.com  
Phone number: 410-910-0884



Visit [www.OneLabOneEarth.com](http://www.OneLabOneEarth.com)

SFE

Method

Performance

Samples

Conclusions

Q&A