

Automated Sample Preparation for Determination of Perfluoroalkyl and Polyfluoroalkyl Substances (PFAS)



Angelika Köpf, LCTech GmbH

August 2020



LCTech GmbH

History and Development

SOLUTIONS BY **LC**Tech

- 1998 Foundation LCTech GmbH
- 2000 First own developed system from LCTech
 - TACS totally automated clean-up system for pesticide sample preparation
 - ...
- 2011 Start of a new period in LCTech - FREESTYLE
 - Robotic system for automated sample clean-up SPE, EVA, GPC...
 - 2016 FREESTYLE XANA – Automated sample clean-up for water samples
 - 2019 FREESTYLE SPE (XANA)-PFAS – Dedicated material for PFAS application



Perfluoroalkyl and Polyfluoroalkyl Substances (PFAS)

Recent Activities Europe, EFSA 2018



- Following Stockholm Convention POP Regulation (EG) No 850/2004 for persistent organic pollutants, EU No 757/2010 for update of above, e.g. that restricts PFOS worldwide EU published in 05/2018 the restriction for production/distribution, trading and import of PFOA/PFC with deadline for execution of 04.07.2020 (EU 2017/1000)
- Restriction will also be applied for all other matrix e.g. workman and other textiles with deadline of 04.07.2023



Perfluoroalkyl and Polyfluoroalkyl Substances (PFAS)

Recent Activities Europe, EFSA 2018



- New work for risk assessment PFAS
- First PFOS/PFOA (published in December 2018)
- Other PFAS: under development
- Based on human data, instead of animals

- EFSA recommendation after meeting of the CWG PFAS 16.10.2018;
EURL/NRL-Workshop on 30 - 31 Oct 2018

- In comparison to the Health Based Guidance Values of EFSA of 2008:
 - a factor of about 80 lower for PFOS (13 ng/kg bw/wk)
 - a factor of about 1750 lower for PFOA (6 ng/kg bw/wk)



Perfluoroalkyl and Polyfluoroalkyl Substances (PFAS)

Work EFSA Mentioned in 2018 to Look into



- Mixture approach compared to TEQ like Dioxins?
- Extend to look at a certain number of PFAS? Which ones?
- Lower LOQs in case possible, especially for food regularly consumed
- Looking for hot spots:
 - Animals foraging outside
 - Also private owners of chicken farms ...
- Bioassay to estimate contribution unknown PFAS by comparison with LC/MS



Perfluoroalkyl and Polyfluoroalkyl Substances (PFAS)

Recent Activities in the USA



- In November 2018, the United States Environmental Protection Agency (U.S. EPA) update a Review of Method US EPA 537.1:
"Determination of Selected Per- and Polyfluorinated Alkyl Substances in Drinking Water by Solid Phase Extraction and Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS)"
- The method uses an offline solid-phase extraction (SPE) with (LC-MS/MS) to extract, enrich, and determine 18 PFAS in drinking water.



Perfluoroalkyl and Polyfluoroalkyl Substances (PFAS)

Recent Activities in the USA



<u>Analyte^a</u>	<u>Acronym</u>
Hexafluoropropylene oxide dimer acid	HFPO-DA
N-ethyl perfluorooctanesulfonamidoacetic acid	NEtFOSAA
N-methyl perfluorooctanesulfonamidoacetic acid	NMeFOSAA
Perfluorobutanesulfonic acid	PFBS
Perfluorodecanoic acid	PFDA
Perfluorododecanoic acid	PFDoA
Perfluoroheptanoic acid	PFHpA
Perfluorohexanesulfonic acid	PFHxS
Perfluorohexanoic acid	PFHxA
Perfluorononanoic acid	PFNA
Perfluorooctanesulfonic acid	PFOS
Perfluorooctanoic acid	PFOA
Perfluorotetradecanoic acid	PFTA
Perfluorotridecanoic acid	PFTTrDA
Perfluoroundecanoic acid	PFUnA
11-chloroeicosafluoro-3-oxaundecane-1-sulfonic acid	11Cl-PF3OUdS
9-chlorohexadecafluoro-3-oxanone-1-sulfonic acid	9Cl-PF3ONS
4,8-dioxa-3H-perfluorononanoic acid	ADONA

Perfluoroalkyl and Polyfluoroalkyl Substances (PFAS)

Summary of the Method



- 250 mL water sample with added surrogates is loaded onto an offline SPE cartridge containing polystyrene-divinylbenzene (SDVB)
- Rinsing of bottle with MeOH, this rinsate applied for elution
- Concentration to dryness
- Adjusted to a 1-mL volume with 96:4 % (vol/vol) methanol:water
- addition of internal standards
- 10- μ L injection into an LC- MS/MS equipped with a C18 column



Perfluoroalkyl and Polyfluoroalkyl Substances (PFAS)

QC-Requirements: Initial QC to prove Analytical Setup



- 9.2.2. Initial Demonstration of Low System Background
- 9.2.3. Initial Demonstration of Precision (IDP) of 4-7 replicates
- 9.2.4. Initial Demonstration of Accuracy (IDA)
- 9.2.5. Initial Demonstration of Peak Asymmetry Factor
- 9.2.6. Minimum Reporting Level (MRL) Confirmation
- 9.2.7. Calibration Confirmation



Perfluoroalkyl and Polyfluoroalkyl Substances (PFAS)

QC-Requirements: Ongoing Requirements for Labs



- 9.3.1. Laboratory Reagent Blank (LRB)
- 9.3.2. Continuing Calibration Check (CCC)
- 9.3.3. Laboratory Fortified Blank (LFB)
- 9.3.4. Internal Standards (IS)
- 9.3.5. Surrogate Recovery



Perfluoroalkyl and Polyfluoroalkyl Substances (PFAS)

SPE Clean-up – Challenges in Daily Use



- Heavy workload
- Standardization
- Reproducibility

high demand for automation of the SPE procedure

- Easy handling
- Intuitive programming
- Plannable deadlines
- Meeting low detection limits
- it needs to be ensured that there is No blind values for PFAS Analysis brought in by equipment



Perfluoroalkyl and Polyfluoroalkyl Substances (PFAS)

SPE – Solution in Daily Use


SOLUTIONS BY **LC** Tech

- FREESTYLE SPE PFAS or
- FREESTYLE XANA PFAS

→ Minimised fluoroplastic components

→ Constantly low PFAS background

→ robust automation for 24/7 operation




FREESTYLE

SOLUTIONS BY **LC** *Tech*



XANA for Automated SPE in PFAS Analysis

Automation based on FREESTYLE BASIC



XANA for Automated SPE in PFAS Analysis

Easy Handling – Your Labware as Usual

SOLUTIONS BY **LC** *Tech*




Just hook the rack into the toothed plate.



XANA for Automated SPE in PFAS Analysis

Application of SPE Cartridges



Flexibility
just by adapters

XANA for Automated SPE in PFAS Analysis

Easy Handling – Simple Operation

SOLUTIONS BY **LC** Tech

- For samples > 100mL FREESTYLE XANA is applied
- 24 positions for 250 up to 1 L sample bottles
- Rack is easy to move and load
- Definition of *Analysis Batch* in US EPA 537.1 – a set of samples that is analysed on the same instrument during a 24-hour period, including no more than 20 Field Samples



XANA for Automated SPE in PFAS Analysis

Easy Handling – Simple Operation

SOLUTIONS BY **LC** *Tech*

- Sample bottles with special designed cap for one-hand operation.
- Angular position of the suction capillary to maximize the sample load.
- Suction capillary can be height-adjusted



XANA for Automated SPE in PFAS Analysis

Automated Rinsing of Sample Bottles

SOLUTIONS BY **LC** *Tech*

- Subsequently to sample loading Rinsing of the sample bottle and quantitative sample transfer by separate rinsing capillary is possible with up to 8 different solvents
- Fully automated
- Rinsing solution can also be applied after drying of columns e.g. to directly use this as eluting solvent



XANA for Automated SPE in PFAS Analysis

Simple Operation – Intuitive Software

Generate any SPE Method:

- Move next step from the basic tools with **drag & drop** to second column.
- The default button remains in the basic tools and can be chosen again.
- The collected steps in the second column build the method .
- Each button there can be individually changed and adopted by using the „bars“.



XANA for Automated SPE in PFAS Analysis

Simple Operation – Intuitive Software



Software allows the transfer of nearly any manual SPE method straight to the automation

The screenshot displays the XANA software interface for configuring an SPE step. The main window is titled "SPE Steps: Step number: 2" and "ID: 530". The "Load" step is selected, with a "Sample Volume" of 10.0 mL. The "Transfer Type" is "Transfer Sample-Aliquot over sample loop" and the "Input Vial Type" is "TypeI@16". The "Suction Speed" is 5 mL/min, "Dispensing Speed" is 2 mL/min, "Waiting time after dosage" is 0 s, and "Waiting Time after Step" is 5 s. The "Dispense" option is set to "into Waste". The "Calculated time" for this step is 2,2 min. The "SPE Steps: 4" section shows a "Total: 118.6 min".

Basic Tools: Conditioning, Emptying, Load, Washing, Drying, Eluting

SPE Steps: Conditioning, Load, Washing, Drying

Parkposition:

Memo:

Input vials					
Step N°	Owner	Count	Type	Basic type	ID
2	Load	1	TypeI@16	LAD	530

Result vials					
Step N°	Owner	Count	Type	Basic type	ID
4	Drying	1	TypeI@60	DRY	

XANA for Automated SPE in PFAS Analysis

Throughput- External Working Station

SOLUTIONS BY **LC** *Tech*

- polypropylene transfer tube system acc. to material chapter of the method
- Pumps 3 samples in parallel with flow rates: 1 – 30 mL/min.
- Conditioning, Washing, Rinsing and Drying of 3 columns in parallel with Up to 8 solvents for conditioning, rinsing and washing
- Detection of empty bottles, positions that are not taken aren't processed!



XANA for Automated SPE in PFAS Analysis

Workstation ensures parallel Loading

SOLUTIONS BY **LC** *Tech*

- Loading with positive pressure
- Tolerates back pressure up to 4 bar
- Batches of 3 columns processed at a time. 6 positions for loading, drying or eluting of columns (parallel loading of 3 samples and drying of 3 samples)



XANA for Automated SPE in PFAS Analysis

Planable Results – Meeting Timelines



Non-stop policy – this is an unique error handling

The screenshot shows the SPE Planner software interface. The main window is titled 'SPE Column LCTech_3ml' and 'ID: 603'. The 'Washing' step is selected, showing a volume of 4 ml, 2 repetitions, and a total volume of 12 ml. The calculated time for this step is 5.7 min. The 'SPE Steps: 6' section shows a total time of 24.6 min. The 'Input vials' and 'Result vials' tables are also visible.

Step N°.	Owner	Count	Type	Basic type	ID
3	User	1	Type1@00	LAD	603

Step N°.	Owner	Count	Type	Basic type
3	User	1	Type1@00	LAD

Pressure Monitoring

Once this monitoring alerts a blocked cartridge, this one will be taken out of the sequence (noted in report) and the rest of the sequence processed

XANA for Automated SPE in PFAS Analysis

Study with FREESTYLE XANA for US EPA 537.1



- Samples
 - Laboratory Fortified Blanks (LFB) – purified water spiked with targets and surrogates with 4 ng/L (3 replicates)
- Standards:
 - Internal standard mix (not required for sample spiking) (Wellington PN EPA-5371S)
 - Surrogate mix (Wellington PN EPA 537SS-R1)
 - PFAS native mix (Wellington PN EPA-537PDS-R1)
- SPE Clean-up: Bond Elut LMS 500 mg, 6 mL (part number 12255021)
- Clean-up performed at LC Tech GmbH Germany on FREESTYLE XANA
- Analysis performed at Agilent Inc. Wilmington on LC-MS/MS

Automated SPE applied in acc. to US EPA 537.1

Application of US EPA 537.1 - labware

Collection of eluates in conical polypropylene tubes acc. to list of material in method

SOLUTIONS BY **LC** *Tech*



Automated SPE applied in acc. to US EPA 537.1



Application of US EPA 537.1 - labware



From US EPA Note: Samples and extracts should not come in contact with any glass containers or pipettes as these analytes can potentially adsorb to glass surfaces -

Material:

250mL polypropylene bottles



Automated SPE applied in acc. to US EPA 537.1

Application of US EPA 537.1

SOLUTIONS BY **LC** *Tech*

- Adapter for 6 mL Cartridge
- Applied here:
 - Bond Elut LMS 500 mg, 6 mL (part number 12255021)
 - EPA methods 537/537.1 were validated using this cartridge



Automated SPE applied in acc. to US EPA 537.1

Analytical Set-up

SOLUTIONS BY **LC***Tech*

- LC/MS/MS: HPLC Instrument: Agilent 1290 Infinity II
- Delay column: Agilent InfinityLab PFC Delay Column (the purpose of this column is to separate potential PFAS contamination in the solvents, tubing, seals, etc. from the injected analytical peaks. This is commonly used for all PFAS analysis by LC).
- Analytical column: Agilent Zorbax RRHD EclipsePlus C18, 2.1 x 100 mm, 1.8 μm (part number 959758-902).
- Mass Spectrometer: Agilent 6470 Triple Quadrupole LC/MS/MS with Agilent Jet Stream Ion Source




Automated SPE applied in acc. to US EPA 537.1

Demonstration of Low Background



All analytes to be below 1/3 the lowest Standard; Proposed MRL of 0.6 ng/L




Automated SPE applied in acc. to US EPA 537.1

Additional Info Demonstration of Low Background



Despite the values are according to requirement, check for the cause of the contamination by evaporating all solvents offline.

Evaporation Experiment of the sum of all solvents in the protocol without FREESTYLE involved




Automated SPE applied in acc. to US EPA 537.1

Initial Demonstration of Precision



Analyse 4 to 7 replicate LFB fortified near the midrange calibration whereof the Relative Standard Deviation RSD within must be <20%




Automated SPE applied in acc. to US EPA 537.1

Initial Demonstration of Accuracy IDA




Calculate average recovery for replicates used in IDP whereof the average recovery to be $\pm 30\%$ of true value



Automated SPE applied in acc. to US EPA 537.1

Initial Calibration



Recoveries must be within 70 to 130% for all targets and surrogates with the lowest standard within 50 – 150%


At least 5 calibration points - here
 All compounds calibrated with 1st order curve, forced zero intercept, and 1/(concentration) weighting.

Automated SPE applied in acc. to US EPA 537.1

Internal Standard IS



Compare continuously internal standard peak area to average response in calibration standards. Peak areas must be within $\pm 50\%$ of the average areas in the calibration standards.





Automated SPE applied in acc. to US EPA 537.1

Surrogate Standard Recoveries SUR



Compare Surrogate Standard Recoveries which are added to CCCs, LRBs, LFBs, LFSMs, LFSMDs, FD, and FRB prior to extraction. Recoveries must be within 70 % to 130 %



Automated SPE applied in acc. to US EPA 537.1

The Automated Solution



- FREESTYLE XANA as automated sample delivery system in combination with material listed in chapter 6 and SPE cartridge Bond Elut is a working system that is applicable for US EPA quality objectives for water analysis on PFAS following US EPA 537.1.
- Automation of SPE helps to standardize, increase reproducibility, and meeting timelines thus ensuring to avoid blind values for PFAS Analysis
- FREESTYLE XANA means full automation without any manual interaction from loading an extract to the eluted sample of up to 24 samples in one batch



Special Thanks to



Sebastian Wierer and Dr. Uwe Aulwurm

LCTech GmbH Obertaufkirchen, Germany

wierer@lctech.de, aulwurm@lctech.de



Agilent Technologies

Matt Giardina and Tarun Anumol

Agilent Technologies, Inc., DE Wilmington, USA

matthew_giardina@agilent.com, tarun.anumol@agilent.com



Thank you for Your Attention!

SOLUTIONS BY **LC***Tech*



LC Tech GmbH
Daimlerstraße 4
84419 Obertaufkirchen
Germany

info@LCTech.de

