


Analysis of Persistent Organic Pollutants (POPs) by high-resolution GC/MS

Tarun Anumol, Matthew Curtis, Paul Contreras and Sofia Nieto
Agilent Technologies Inc.

Why do you need a GC/Q-TOF? Water Treatment plant samples in NY




Good reproducibility

Statistically no relevant features to account for spike in toxicity on day 2

Samples re-run on GC/MS with suitable extraction

7250 GC/Q-TOF

Resolving power: >25,000 at m/z 272
Mass accuracy: <2 ppm
IDL <60 fg OFN (EI)
Acquisition rate: 1-50 Hz



Standard/Low Energy EI Source




CI Source

Reproducible Spectral Performance


High Resolution and Mass Accuracy

Simultaneous High Resolution and Wide Dynamic Range





Traditional EI vs Low Energy EI

hexachlorobutadiene




Chemical Ionization vs Low Energy EI

2,4-Dinitrotoluene




Workflow Strategy




GC/Q-TOF Accurate Mass Library

GC/Q-TOF Pesticides & Environmental Pollutants
PCDL – now with **1000+ compounds**:

- High Resolution Spectra
- Expert curation
- Better compound alignment with EPA SVOC targets and Agilent GC/Q-TOF MRM database



Talk on Screening & Identification by GC-QTOF




PAHs


Profiling Environmental Contaminants in Water Using GC-Q/TOF + NEMC 2019, available online




Amines




Organophosphates



Phthalates



Nitroanilines



Chloronitrobenzenes

Compound: Sulfotep

Mass: 322.02274 Formula: C8H20O5P2S2 Save

Notes:

- Toxicology drug; Pesticide; Veterinary drug; Parasympathomim
- 同義語: 殺虫剤; 殺菌剤
- CAS RN: 135-05-1; Appendix A: Extremely Hazardous Chemicals
- Chinese National Food Safety Standard: Maximum Residue Limits for Pesticides in food (GB 2763-2014)
- JPL

Structure: MOL Text

Chemical structure diagram showing a bisphosphate ester with methyl groups and sulfur atoms.

Persistent Organic Pollutants (POPs) in Stockholm Convention

- The Stockholm Convention on POPs (2001) is a global treaty to protect **human health** from chemicals that remain in the environment and are **persistent, bioaccumulative** and transportable across the globe.
- **Include:**
 - Industrial chemicals ex. PCBs, hexachlorobenzene
 - Pesticides ex. Aldrin, DDT, endrin, toxaphene etc.
 - Pharmaceuticals
 - Solvents
 - By-products ex. **Dioxins & Furans**
- Initially 'dirty dozen' but new POPs include Perfluoro alkyl substances (**PFAS**), short chain chlorinated parafins (**SCCPs**), pesticides, Deca-BDE etc.
- Typically hydrophobic, therefore GC/MS analysis is a superior choice




Analysis of Dioxins by GC/Q-TOF



Analysis of Dioxins in US done mainly by EPA 1613B


Tetra- through Octo-Chlorinated Dioxins & Furans by HRGC/HRMS



PCDD and PCDF Nomenclature and Isomers

Dioxins

Polychlorinated dibenzodioxin (PCDD)




2,3,7,8 – TCDD

(2,3,7,8-Tetrachlorinated dibenzodioxin)

Furans

Polychlorinated dibenzofurans (PCDF)



2,3,4,7,8 – PCDF

(2,3,4,7,8-Pentachlorinated dibenzofurans)

Chlorine Atoms	PCDD Isomers	PCDF Isomers
1	2	4
2	10	16
3	14	28
4	22	38
5	14	28
6	10	16
7	2	4
8	1	1
Total	75	135



75 total PCDD isomers

135 total PCDF isomers

17 toxic isomer


Tetrachlorinated dibenzodioxins (TCDD)

Separation of Isomers



$x + y = 4$
22 isomers

1 toxic



Method Detection Limit (MDL) for 2,3,7,8-TCDD


Samples concentration of 10,000X (500 mL to 50 μ L)

MDL for 2,3,7,8-TCDD is 4.4 pg/L in EPA 1613B

Sample Concentration	Injected amount	15 eV, 1 μ A			70 eV, 1 μ A			70 eV, 10 μ A		
Spiked equivalent		RSD %	MDL (pg/L)	IDL (fg)	RSD %	MDL (pg/L)	IDL	RSD %	MDL (pg/L)	IDL
1 pg/L	10 fg	43%	1.6	16	84%	3.1	31	15%	0.6	6
2.5 pg/L	25 fg	17%	1.6	16	24%	2.3	23	9%	0.8	8
10 pg/L	100 fg	7%	63	63				5%	1.8	18


Dioxins concentration in Herring (6.4% fat)

Comparison between GC/Q-TOF & Magnetic Sector



Dioxins in Sediment reference



Comparison between GC/Q-TOF & Magnetic Sector




Short Chain Chlorinated Paraffins (SCCPs)

- SCCPs are emerging persistent organic pollutants, bioaccumulative, toxic and persistent in the environment
- It is a complex mixture of n-alkanes with different degree of chlorination
- Used as flame retardants in plastics and other materials, as well as in few other applications such as metal processing
- Represent a substantial challenge due to their self-interference as well as interference with other components of complex industrial matrices

Challenges for the Analysis of the SCCPs: Chromatographic Separation



Challenges for the Analysis of the SCCPs: Excessive Fragmentation



- Spectral fragmentation under normal conditions for 70 eV acquisitions
- The number below the m/z is the relative intensity

Automatic Molecular Formula Generation

Method Editor: Generate Formulas

Generate Formulas from Spectrum Peaks

Allowed Species Limits Charge State Fragment Formulas Scoring

Charge carrier to be assumed if not known

Positive ions: Negative ions:

-electron +electron

+H -H

+Na +Cl

+K +Br

+NH4 +HCOO

+C2H5 +CH3COO

+C3H5 +CF3COO

MS ion electron state: allow both even and odd



Group hits with same formula (but different charge carriers)

Elements and limits


Element	Minimum	Maximum
C	1	60
H	0	120
O	0	6
N	0	4

Molecular Formula Generation uses accurate m/z of the molecular ion, accurate m/z of the isotopes, isotope spacing and isotope intensity for confirmation.

Examples of SCCP Spectra in NCI





EIC Overlay of SCCP mixture



Quantitation Based on Pure Congeners

Sample				C12Cl6...	C12Cl6 Results					Qualifier (339.9905) Resu ^		
Name	Type	Level	Exp. Conc.	RT	Resp.	Final Conc.	Mass Accuracy	S/N	MZ	Mass Accuracy	S/N	
S8_StdCong_50ppb-1	Cal	3	50.0000	10.9...	78333	50.0713	0.9431	8669.78	339.9912	2.0387	1424	
S8_StdCong_50ppb-2	Cal	3	50.0000	10.9...	96498	62.4729	0.6964	8225.01	339.9908	0.7971	3467	
S9_StdCong_500ppb-1	Cal	4	500.0000	10.9...	758521	514.4659	-4.3322	20306...	339.9892	-3.9533	5015	
S9_StdCong_500ppb-2	Cal	4	500.0000	10.9...	714430	484.3629	-1.6121	20703...	339.9902	-0.8431	4520	
S1_SCCPmix_51p5%C...	Sample			10.9...	298472	200.3700	-1.1223	8.29	339.9906	0.3375	12	
S1_SCCPmix_51p5%C...	Sample			10.9...	314725	211.4666	-1.3633	9.60	339.9906	0.3773	13	



Challenges of SCCP Quantitation

Method Table

Time Segment: < <All> > Compound: < C11C15 > Reset Table View

Qualifier			
MZ	Rel. Resp.	Uncertainty	Int. Pams.
294.0295	100.0	20.0	

Qualifier			
Name	TS	Scan	Type
C11C15	1	Scan	Target

Qualifier			
MZ	Rel. Resp.	Uncertainty	Int. Pams.
292.0139	225.0	20.0	

RT	Int.	Int.
10.300	Spectrum Summation	...

Sample Information

Compound Information

- EIC(293.0217) Scan S2_SCCPmix_55%CI-1.D

Counts x10³

Acquisition Time (min)

293.0217, 292.0139
Ratio = 264.5 (117.6%)

Relative Abundanc. x10²

Acquisition Time (min)

Quant Results: Comparison of 55.5% and 63% SCCP Mixtures


Congener	RT range, min	Concentration, ppb			%		
		51.50%	55.50%	63%	51.50%	55.50%	63%
C ₁₀ Cl ₄	8.8-9.1	115.5	193.7	23	2.3	3.9	0.5
C ₁₀ Cl ₅	9-10.3	106.1	135.3	84.3	2.1	2.7	1.7
C ₁₀ Cl ₆	9.6-10.8	5.9	15.3	41.7	0.1	0.3	0.8
C ₁₀ Cl ₇	10.1-11.2	0.9	6.7	51.6	0.02	0.1	1
C ₁₀ Cl ₈	10-11.3	2.5	4.4	38	0.05	0.1	0.8
C ₁₁ Cl ₄	9.2-10	189.2	96.2	36.6	3.8	1.9	0.7
C ₁₁ Cl ₅	9.5-10.5	364.6	340.7	89.4	7.3	6.8	1.8
C ₁₁ Cl ₆	10-10.8	342	614.5	330.3	6.8	12.3	6.6
C ₁₁ Cl ₇	10.5-11.7	70.4	353.2	825.9	1.4	7.1	16.5
C ₁₁ Cl ₈	11-12.5	3.3	25.4	210.6	0.1	0.5	4.2
C ₁₂ Cl ₄	9.4-10.5	290.7	129.8	11.1	5.8	2.6	0.2
C ₁₂ Cl ₅	10-11.2	351.3	253.7	31.3	7.0	5.1	0.6
C ₁₂ Cl ₆	10.3-11.5	205.9	240.2	46.8	4.1	4.8	0.9
C ₁₂ Cl ₇	10.9-12.1	331.9	733.3	763.7	6.6	14.7	15.3
C ₁₂ Cl ₈	11.4-12.6	9.5	49.3	167.3	0.2	1	3.3
C ₁₃ Cl ₅	10.1-11.3	218.8	126.5	12.3	4.4	2.5	0.2
C ₁₃ Cl ₆	10.8-11.8	200.9	161.9	26.1	4	3.2	0.5
C ₁₃ Cl ₇	11.4-12.5	642.3	865.9	497.4	12.8	17.3	9.9
C ₁₃ Cl ₈	11.9-13	84.9	287.8	628.2	1.7	5.8	12.6
Total				70.7	92.7	78.2	

Low Energy EI Analysis of SCCP

- To improve sensitivity of detection and accuracy of quantitation for SCCP congeners with low chlorine content, the low energy EI approach was used.
- Traditional 70 eV EI results in a high degree of fragmentation of SCCP molecules, and does not provide enough unique ion clusters for individual identification.
- Multiple low electron energy settings were evaluated to determine the optimal value.
- The optimum combination of spectral tilt and signal response was achieved with an electron energy set at 22 eV.


Low Energy EI Analysis of SCCP

- Low energy EI data indicated a higher degree of fragmentation of the SCCP molecules compared to negative CI



Low Energy EI Analysis of SCCP

- However, this technique allowed more sensitive detection of the SCCP species with low chlorine atom number (such as C10Cl4)



Perfluoroalkyl Substances (PFAS)

Analysis of Fluorinated Alkyl Alcohols in Biosolids



Standards & Samples

Standards at 100 ng/ μ L in MTBE

Biosolid extracts prepared in EtOAc

Positive CI, 20% Methanol

MMI: 2 μ L cold, splitless injection

65 °C (0.01 min), 300 °C/min to 250 °C

Two min post-column backflush


Oven Ramp:

- 60 °C for 1 min
- 3 °C/min to 75°C for 0 min
- 20 °C/min to 210 °C for 0 min

Acronym	Type	r.t.	Exact Mass + H
4:2 FTOH	Target	4.078	265.0269
5:1 FTOH	Target	5.221	301.0081
6:2 FTOH	Target	5.773	365.0206
6:1 FTOH	Target	6.188	351.0049
7:2 sFTOH	Target	6.237	415.0174
7:1 FTOH	Target	6.981	401.0017
8:2 MFTOH	ISTD	7.449	469.0334
8:2 FTOH	Target	7.471	465.0142
8:1 FTOH	Target	7.598	450.9985
9:1 FTOH	Target	8.098	500.9953
10:2 FTOH	Target	8.470	565.0078
10:1 FTOH	Target	8.523	550.9921
11:1 FTOH	Target	8.886	600.9889
d7-MeFOSE	ISTD	11.868	565.0466
MeFOSE	Target	11.889	558.0026
d9-EtFOSE	ISTD	11.897	581.0748
EtFOSE	Target	11.928	572.0183


Chromatogram

17 FTOHs & FOSEs




Structure Elucidation

11:1 FTOH



Structural Elucidation: Why use a HRAM?

Nominal loss = 38 m/z; F₂ or something else?



Acronym	Observed Base Peak m/z	Molecular ion -F ₂ m/z	Δppm	Molecular ion -H ₂ O, -HF m/z	Δppm
4:2 FTOH	227.0104	227.0301	86.7726	227.0102	-0.8810

Loss Elucidation

Why use a 'Q'-TOF

Mass accuracy facilitated in structural elucidation of fragment ions

- Nominal mass loss of m/z 38 for n:2 FTOHs (n=4,6,8,10...)
 - Loss of F_2 ruled out due to mass accuracy
 - Loss of H_2O and HF confirmed
 - **MS/MS confirmed H_2O and HF Losses**
- Loss of HF for n:1 FTOH
- Loss of H_2O for perfluorosulfonamidoethanols

Identifying PFAS with loss elucidation

Mass Accuracy

Acronym	Formula	Loss	Exact Mass-Loss	Observed Mass 2	Δ ppm 2
4:2 FTOH	C ₆ H ₅ F ₉ O	-H ₂ O, -HF	227.0102	227.0105	-1.3215
6:2 FTOH	C ₈ H ₅ F ₁₃ O	-H ₂ O, -HF	327.0038	327.0040	-0.6116
8:2 FTOH	C ₁₀ H ₅ F ₁₇ O	-H ₂ O, -HF	426.9972	426.9975	-0.7026
10:2 FTOH	C ₁₂ H ₅ F ₂₁ O	-H ₂ O, -HF	526.9910	526.9915	-0.9488
7:2 sFTOH	C ₉ H ₅ F ₁₅ O	-H ₂ O, -HF	377.0006	377.0011	-1.3263
5:1 FTOH	C ₆ H ₃ F ₁₁ O	-HF	281.0017	281.0019	-0.7117
6:1 FTOH	C ₇ H ₃ F ₁₃ O	-HF	330.9988	330.9990	-0.6042
7:1 FTOH	C ₈ H ₃ F ₁₅ O	-HF	380.9954	380.9956	-0.5249
8:1 FTOH	C ₉ H ₃ F ₁₇ O	-HF	430.9923	430.9927	-0.9281
9:1 FTOH	C ₁₀ H ₃ F ₁₉ O	-HF	480.9894	480.9896	-0.4158
10:1 FTOH	C ₁₁ H ₃ F ₂₁ O	-HF	530.9860	530.9864	-0.7533
11:1 FTOH	C ₁₂ H ₃ F ₂₃ O	-HF	580.9834	580.9835	-0.1721
MeFOSE	C ₁₁ H ₈ F ₁₇ N O ₃ S	-H ₂ O	539.9920	539.9943	-4.2593
EtFOSE	C ₁₂ H ₁₀ F ₁₇ NO ₃ S	-H ₂ O	554.0077	554.0100	-4.1516

Summary

- The GC/Q-TOF is a complementary tool for environmental monitoring that is capable of unique identification of regulated and emerging environmental contaminants.
- The 7250 GC/Q-TOF is a versatile instrument that can analyze various classes of POPs including Dioxins, SCCPs, PCBs and PFAs on the same instrument

Acknowledgements

Dr. Peter Haglund – Umea University

Dr. Anthony Macherone – Agilent Technologies Inc.

Dr. Pierre Dumas - Institut National de Santé Publique du Québec (INSPQ)

Dr. Shoji Nakayama - National Institute for Environmental Studies, Japan