

A New SPE Sorbent for Improved EPH Fractionation




**Brian Kinsella
Applications Chemist
UCT**

NEMC 2020

EPH Fractionation



- Extractable Petroleum Hydrocarbons (EPHs)
- Assess risk of petroleum products in environment
 - Water, Sediment, Soil
- Separate aliphatic & aromatic hydrocarbons using SPE
 - Quantify separately
 - Improved toxicological profile / risk assessment




EPH Fractionation




Approved Methods:


- No fixed Federal / EPA method
- Most labs follow MA DEP method
- States can have different requirements
 - MA DEP Method:
 - 40–140% Recovery for Analytes and Surrogates
 - <5% Breakthrough for Naphthalene and 2-Methylnaphthalene
 - NJ DEP Method:
 - Based off MADEP EPH Method
 - Requires monitoring of additional analytes
 - TCEQ Method:
 - Uses Pentane instead of Hexane
 - 60–140% Recovery for Analytes
 - 70–130% Recovery for Surrogates
 - <10–20% Crossover



Naphthalene



2-Methylnaphthalene



Pyrene

EPH Fractionation



Sample Extraction:

- Soil/Sediment/Solids – prepared according to EPA Method 3540 (Soxhlet extraction) or EPA Method 3546 (microwave extraction)
- Water samples – prepared according to EPA Method 8015D (SPE) or EPA Method 3510 (LLE)
- Sample extract should be in *n*-Hexane for EPH fractionation
 - May necessitate solvent exchange

Analysis:

- Typically GC-FID
 - Non-selective
 - Extraction cartridges should be very clean (no significant background)

Current EPH Methods




- Solid-Phase Extraction with silica gel
 - Weak electrostatic interaction of aromatic ring and -OH group

SPE Cartridge:

- Special heat-treated silica
- UCT cartridges are capped at both ends to ensure silica integrity
- No moisture - deactivates sorbent!

UCT Products:

- XRSIHT15M25 (5g/25mL)
- XRSIHT13M15 (3g/15mL)




SPE Procedure (5g/25mL Cartridge)




1. Condition

- 30 mL Hexane:





2. Load Sample (1 mL Hexane)



● Aliphatic

● Aromatic




Note: Some procedures use DCM & Hexane


SPE Procedure (5g/25mL Cartridge)




3. Elute Aliphatic Fraction



● Aliphatic



● Aromatic



4. Elute Aromatic Fraction



5. Concentrate



Problems Encountered



- Deactivation of silica gel (hygroscopic!)
 - Inconsistent results
 - Low recoveries
 - Breakthrough of aromatic fraction into aliphatic fraction
 - Naphthalene and 2-Methylnaphthalene (most polar)
- Tedious
 - Do not let sorbent go dry!
 - Use stopcocks to control solvent level in SPE cartridge
- Capacity (polar aromatics)
- Optimization of Hexane volume for each batch of Si-OH
- Storage - unused cartridges must be stored in desiccator




New Sorbent



➤ Silver ions functionalized onto a solid support

- Improved stability
- Excellent reproducibility



➤ High Selectivity

- Strong retention of aromatic & unsaturated compounds
- Forms a charge-transfer complex with silver ions
- No breakthrough of aromatic fraction into the aliphatic fraction
 - Naphthalene and 2-Methylnaphthalene



- **High Capacity** – of aromatic hydrocarbons (25 mg/g)
- **Simplified Procedure** – simple, fast protocol
 - No optimization of solvent volume for each batch of cartridges
- **More Environmentally Friendly** – less solvent, no DCM
- **Moisture Stability** – no deactivation of sorbent
- **Lot-To-Lot Reproducibility**
- **Three Formats Available** – 0.5g, 1g, 2g (6mL cartridge)
- **Amenable to Automation** – improved reproducibility & increased sample throughput

Compounds Evaluated



	Aliphatics	Aromatics
1	1-Chlorooctadecane (surrogate)	o-Terphenyl (surrogate)
2	(C9) n-Nonane	2-Bromonaphthalene (fractionation surrogate)
3	(C10) n-Decane	2-Fluorobiphenyl (fractionation surrogate)
4	(C12) n-Dodecane	1,2,3-Trimethylbenzene
5	(C14) n-Tetradecane	2-Methylnaphthalene
6	(C16) n-Hexadecane	Acenaphthene
7	(C18) n-Octadecane	Acenaphthylene
8	(C19) n-Nonadecane	Anthracene
9	(C20) n-Eicosan	Benzo[a]anthracene
10	(C21) n-Heneicosane	Benzo[a]pyrene
11	(C22) n-Docosane	Benzo[b]fluoranthene
12	(C24) n-Tetracosane	Benzo[ghi]perylene
13	(C26) n-Hexacosane	Benzo[k]fluoranthene
14	(C28) n-Octacosane	Chrysene
15	(C30) n-Triacontane	Dibenz[a,h]anthracene
16	(C32) n-Dotriacontane	Fluoranthene
17	(C34) n-Tetratriacontane	Fluorene
18	(C36) n-Hexatriacontane	Indeno[1,2,3-cd]pyrene
19	(C38) n-Octatriacontane	Naphthalene
20		Phenanthrene
21		Pyrene

GC Conditions Used



GC-MS Conditions:

Parameter	Conditions
GC-MS system	Thermo Trace 1300 GC & ISQ MS
GC column	Restek Rxi®-5sil MS, 30m x 0.25mm, 0.25µm with Integra-Guard
GC liner	4 mm split liner with deactivated glass wool
Injection	1 µL split (1:100) at 300°C
Carrier gas	Ultra-high purity Helium at a constant flow of 1.2 mL/min
Oven temp. program	Initial temperature at 50°C, hold for 3 min; ramp at 10°C/min to 320°C, hold for 10 min
Temperatures	Transfer line = 275°C; Ion source = 275°C
Full scan range	35 - 600 amu


GC-FID Conditions:

Parameter	Conditions
GC-FID system	HP5890 Series II
GC column	Restek Rrx®-5 MS, 30m x 0.32mm, 0.5µm
GC liner	4 mm split/splitless liner with deactivated glass wool
Injection	1 µL splitless at 280°C
Carrier gas	Ultra-high purity Helium at a constant flow of 1.3 mL/min
Oven temp. program	Initial temperature at 55°C, hold for 0.3 min; ramp at 10°C/min to 320°C, hold for 10 min
FID temp.	320°C
Make-up gas	Ultra-high purity Hydrogen at a constant flow of 30 mL/min
Air flow	350 mL/min


Chromatogram - Aliphatic Fraction (Hexane)



GC-MS:




GC-FID:




Chromatogram - Aromatic Fraction (Acetone)



GC-MS:



GC-FID:



Fusion Ag⁺ Procedure (1g/6mL Cartridge)




1. Condition


- 4 mL Acetone
- 4 mL Hexane




2. Load Sample (1 mL Hexane)



● Aliphatic



● Aromatic




- No need to close stopcock in between steps
- Simply perform each step until dripping stops

Fusion Ag⁺ Procedure




3. Elute Aliphatic Fraction

4 mL
Hexane



4. Elute Aromatic Fraction

4 mL
Acetone



5. Concentrate
or
Direct Analysis

Fusion Ag⁺ Procedure





		0.5g, 6mL	1g, 6mL	2g, 6mL
Condition	1. Acetone	2 mL	4 mL	6 mL
	2. Hexane	2 mL	4 mL	6 mL
Elute	1. Hexane	2 mL	4 mL	6 mL
	2. Acetone	2 mL	4 mL	6 mL

Complete Separation - Aliphatics



- 200 µg/mL Aliphatics Standard + Surrogates

Hexane (Aliphatic) Fraction



Acetone (Aromatic) Fraction



*Included in aliphatic standard

Complete Separation - Aromatics



- 200 µg/mL Aromatics Standard + Surrogates

Hexane (Aliphatic) Fraction



Acetone (Aromatic) Fraction



Recovery - Aliphatics




	Sample 1	Sample 2	Sample 3	Sample 4	Mean Rec (%)	RSD (%)
1-Chlorooctadecane (surr.)	92	92	95	95	93	1.69
(C9) n-Nonane	96	97	99	98	98	1.37
(C10) n-Decane	96	97	99	98	98	1.31
(C12) n-Dodecane	98	98	100	100	99	1.21
(C14) n-Tetradecane	96	97	99	98	98	1.27
(C16) n-Hexadecane	95	96	98	98	97	1.56
(C18) n-Octadecane	91	92	94	94	93	1.67
(C19) n-Nonadecane	93	93	96	96	94	1.64
(C20) n-Eicosan	111	112	119	119	115	3.71
(C22) n-Docosane	91	92	94	94	93	1.74
(C24) n-Tetracosane	87	88	90	90	89	1.76
(C26) n-Hexacosane	89	89	92	92	90	1.80
(C28) n-Octacosane	87	88	90	90	89	1.84
(C30) n-Triacontane	82	83	85	85	84	1.92
(C36) n-Hexatriacontane	80	80	83	84	82	2.53

Recovery - Aromatics




	Sample 1	Sample 2	Sample 3	Sample 4	Mean Rec (%)	RSD (%)
o-Terphenyl (surr.)	92	98	111	102	101	8.05
2-Bromonaphthalene (frac. surr.)	92	98	112	102	101	8.27
2-Fluorobiphenyl (frac. surr.)	93	98	112	103	101	8.07
2-Methylnaphthalene	92	97	111	101	100	8.15
Acenaphthene	83	87	98	91	90	7.24
Acenaphthylene	92	98	112	103	101	8.21
Anthracene	93	99	113	103	102	8.32
Benz[a]anthracene	91	97	111	102	100	8.34
Benzo[a]pyrene	91	94	106	99	97	6.98
Benzo[b]fluoranthene	91	97	112	103	101	8.74
Benzo[k]fluoranthene	90	96	110	100	99	8.48
Chrysene	91	97	111	102	100	8.57
Dibenz [a,h] anthracene	85	93	108	99	96	10.07
Dibenzo[ghi]perylene	86	94	108	99	97	9.83
Fluoranthene	90	96	109	100	99	8.13
Fluorene	92	97	111	102	100	8.22
Indeno[1,2,3-cd]pyrene	84	90	103	95	93	8.56
Naphthalene	108	115	136	124	121	9.89
Phenanthrene	91	97	111	101	100	8.20
Pyrene	91	97	112	102	101	8.52

Batch-to-Batch Consistency - Aliphatics



n=3 per batch

Batch-to-Batch Consistency - Aromatics




n=3 per batch

Time Savings vs. HT Silica




Fractionation Step		Extraction Time (min)		
		Fusion-Ag ⁺ 0.5g, 6mL	Fusion-Ag ⁺ 1g, 6mL	Heat Treated Silica 2g, 6mL
Conditioning	Acetone	0.8	1.8	
	DCM (+5 min hold)			5.8
	Hexane	0.8	1.7	1.9
Sample Loading	Hexane	0.3	0.5	0.5
Elute Aliphatics	Hexane	0.75	1.7	1.4
Elute Aromatics	Acetone	0.8	1.5	
	DCM			2.8
Total extraction time		3.5	7.2	12.4




Solvent Savings vs. HT Silica



Fractionation Step		Solvent Consumption (mL)		
		Fusion-Ag+ 0.5g, 6mL	Fusion-Ag+ 1g, 6mL	Heat Treated Silica 2g, 6mL
Conditioning	Acetone	2	4	
	DCM			3
	Hexane	2	4	6
Sample Loading	Hexane	1	1	1
Elute Aliphatics	Hexane	2	4	4
Elute Aromatics	Acetone	2	4	
	DCM			6
Total solvent consumption		9	17	20
Total chlorinated solvent used		0	0	9



Summary



Simplified Procedure

High Selectivity

Less Solvent (no DCM)

High Capacity

Improved Productivity

Complete Fractionation

Excellent Analytical Performance

Lot-to-Lot Consistency

A decorative horizontal bar at the top of the slide, transitioning from dark blue on the left to dark red on the right.

Questions?

A decorative horizontal bar at the bottom of the slide, transitioning from dark blue on the left to dark red on the right.