

Comprehensive Microcystins Analysis and Identification Using Innovative Sample Preparation and Novel LC/MS Techniques

Tarun Anumol, Ph.D.
Director, Global Environment & Food Markets
Agilent Technologies Inc.

Contact: tarun.anumol@agilent.com



Introduction

Cyanobacteria are naturally occurring in most lakes and rivers

In warm temperatures and presence of nutrients, they can form harmful algal blooms (HABs) and release cyanotoxins

Cyanotoxins can be classified into 3 categories:

- Cyclic peptides, including **microcystins (MCs)** and nodularins
- Alkaloids
- Lipopolysaccharides

MCs are the most commonly occurring in water and potent hepatotoxins. (100s of variations)

Red tide, the toxic algae bloom that kills wildlife, returns to southwest Florida

By **Doug Stanglin** USA TODAY
Published 12:20 p.m. ET Nov. 13, 2019 | Updated 2:56 p.m. ET Nov. 13, 2019



'Red tide' toxic algae bloom kills sea life and costs Florida millions
Florida is watching the approach of a red tide invasion to its beaches which costs the tourist and fishing industry millions of dollars in losses. USA TODAY

Harmful algae bloom affecting water supply: 400,000 in Toledo, Ohio await test results

POSTED 9:37 PM, AUGUST 3, 2014, BY CNN WIRE SERVICE

Toxic algae bloom now stretches 650 miles along Ohio river

Microcystins

Global Guidelines

Organization	Compound	Drinking Water (ug/L)
World Health Organization	LR	1.0
EU Drinking Water Directive (proposed)	LR	1.0
Canada	LR*	1.5
USEPA	total	8.0
Ohio		0.3 µg/L bottle-fed infants and preschool age children 1.6 µg/L schoolage children and adults
Oregon		0.3 µg/L age 5 and younger 1.6 µg/L age 6 and older
Minnesota		0.1 ug/L

* “maximum acceptable concentrations (MAC); protective of total microcystins”


EPA Method 544: DETERMINATION OF MICROCYSTINS AND NODULARIN IN DRINKING WATER BY SOLID PHASE EXTRACTION AND LIQUID CHROMATOGRAPHY/TANDEM MASS SPECTROMETRY (LC/MS/MS)

General Microcystins Structure

- 7 amino acids in cyclic structure
- Most variations differ in R² & R⁴ location


Naming convention uses amino acid abbreviations for these positions

Superscript number indicates amino acid order within molecule



Microcystin	R ²	R ⁴	Formula	Neutral Mass
LR	Leucine	Arginine	C ₄₉ H ₇₄ N ₁₀ O ₁₂	994.5488
Desmethyl LR	Leucine	Arginine	C ₄₈ H ₇₂ N ₁₀ O ₁₂	980.5331
RR	Arginine	Arginine	C ₄₉ H ₇₅ N ₁₃ O ₁₂	1037.5658
YR	Tyrosine	Arginine	C ₅₂ H ₇₂ N ₁₀ O ₁₃	1044.5280
LA	Leucine	Alanine	C ₄₆ H ₆₇ N ₇ O ₁₂	909.4848
LW	Leucine	Phenylalanine	C ₅₄ H ₇₂ N ₈ O ₁₂	1024.5270
LF	Leucine	Tryptophan	C ₅₂ H ₇₁ N ₇ O ₁₂	985.5161
HtyR	Homotyrosine	Arginine	C ₅₃ H ₇₄ N ₁₀ O ₁₃	1058.5437


Microcystin-LR Structure



C₄₉H₇₄N₁₀O₁₂
Mass = 994.5488

Position	Abbreviation	Amino Acid
R ¹	Ala ¹	Alanine
R ²	Leu ² (L)	Leucine
R ³	MeAsp ³	Methylaspartic acid
R ⁴	Arg ⁴ (R)	Arginine
R ⁵	Adda ⁵	3-amino-9-methoxy-2,6,6-trimethyl-10-phenyldeca-4(E),6(E)-dienoic acid
R ⁶	Glu ⁶	Glutamic acid
R ⁷	Mdha ⁷	N-methyldehydroalanine

Optimizing Microcystin Analysis: Considerations



EPA Method 544

Analysis of Microcystins & Nodularin by Offline SPE + LC/MS/MS

Collection and preservation: 0.5 L

- TRIZMA, 2-chloroacetamide, ascorbic acid, EDTA

Add surrogate (C_2D_5 -LR)

Filter; rinse bottle with 10% MeOH

- Soak filter in 80% MeOH; freeze 1 – 16 hours
- Rinse filter and combine with filtered water

Offline SPE to extract MC's

- Elute with 10 mL of 90% MeOH
- * Evaporate to dryness
- * Reconstitute with 1 mL of 90% MeOH

* “the laboratory is permitted to modify the evaporation technique, separation technique, LC column, mobile phase composition, LC conditions and MS and MS/MS conditions”

Sample Preparation for this study

✓ Collection and preservation: 0.5 L (well water + tap water)

- TRIZMA, 2-chloroacetamide, ascorbic acid, EDTA

✓ Add surrogate (C_2D_5 -LR)

Filter; rinse bottle with 10% MeOH

- Soak filter in 80% MeOH; freeze 1 – 16 hours
- Rinse filter and combine with filtered water

✓ Offline SPE to extract MC's

- Elute with 10 mL of 90% MeOH
- *Evaporate to dryness **2-3 mL**
- *Reconstitute ~~with 1 mL of 90% MeOH~~ **to 5.0 mL with MeOH**

* “the laboratory is permitted to modify the evaporation technique, separation technique, LC column, mobile phase composition, LC conditions and MS and MS/MS conditions”

MRM Parameters^a



**1260 Infinity II
LC + Ultivo
LC/MS/MS**

Microcystin	Precursor ^b (m/z)	Product Ion (m/z)	Collision Energy (V)
LR	995.6	135.2	80
		213.2	80
Desmethyl-LR	981.5	135.2	80
		213.2	80
RR	520.0	135.2	30
		213.2	40
YR	1045.5	135.2	80
		213.2	70
LA	910.5	135.2	70
		213.2	70
LY	1002.5	135.2	80
		213.2	70
LW	1025.5	135.2	80
		213.2	60
LF	986.5	135.2	70
		213.2	50
HtyR	1059.5	135.2	80
		213.2	70

^a The fragmentor and cell acceleration voltages were 150 V and 2 V, respectively, for all transitions; minor differences on Ultivo LC/MS


^b All precursors are singly charged, except RR which is doubly charged.



**1260 Infinity II LC +
6470 MS/MS**


Chromatography

EPA 544 Method Flexibility (Sec. 1.6): “Analytes must be adequately resolved chromatographically in order to permit the mass spectrometer to dwell on a minimum number of compounds eluting within a retention time window.”




Agilent Poroshell
SB-C18
2.1 x 50 mm 2.7
µm LC column

Low Calibrators (ng/mL in vial) 10-uL injections




MC-LR vs Drinking Water Guideline Levels 0.5 L SPE; reconstitute in 5 mL (100x conc)




Guidelines for World Health Organization: Microcystin-LR = 1 ug/L in water
(1 ug/L * 0.5 L) in 5 mL final volume = 0.1 ug/mL = 100 ng/mL in vial

Ion Suppression Evaluation

C₂D₅-LR surrogate in solvent blank, calibrators, and water samples



Quantitative Spike Recovery Results Tap water spikes @ 0.02 ug/L (n=4)



Offline SPE vs Online SPE

Sample collection



Filtration (for surface & WW analysis)



Surrogate addition



Extraction (SPE)



Evaporation



LC-MS Analysis

Benefits of Online SPE

Lower sample collection & transport

More reproducibility due to automation

Much lower labor cost

Faster time to results

Sensitivity

Lower solvent usage

Lower internal std usage


Workflow issues

Risk of contamination between samples



Inter-operator consistency

Reproducibility & Accuracy

Online SPE Setup & Configuration



Solvent selection valve for up to 3 solvents (i.e. Wash, Prepare, Elute solvents)



Online SPE & Direct Injection Conditions

Online SPE

Parameter	Value										
Analytical column	Agilent InfinityLab Poroshell 120 EC-C18, 3.0 × 100 mm, 2.7 μm (p/n 695975-302)										
Column temperature	45 °C										
Trapping columns	2 × PLRP-S cartridges, 2.1 × 12.5 mm (p/n 5982-1271); 6 mL screw top vials (p/n 9301-1377), screw caps (p/n 9301-1379), preslit septa (p/n 5188-2758)										
Injection volume	1,800 μL										
Draw and eject speed	1,000 μL/min; 4,000 μL/min										
Needle wash	MeOH, 5 seconds										
Mobile phases	A) 0.1 % FA + 1 mM NH ₄ F in H ₂ O B) 0.1 % FA in ACN										
Stop time	14 minutes										
Post time	3 minutes										
Flow rate	0.4 mL/min										
Gradient	<table border="1"> <thead> <tr> <th>Time (min)</th> <th>%B</th> </tr> </thead> <tbody> <tr> <td>0</td> <td>30</td> </tr> <tr> <td>4</td> <td>30</td> </tr> <tr> <td>12</td> <td>95</td> </tr> <tr> <td>14</td> <td>95</td> </tr> </tbody> </table>	Time (min)	%B	0	30	4	30	12	95	14	95
Time (min)	%B										
0	30										
4	30										
12	95										
14	95										





Direct Injection

Parameter	Value								
Analytical column	Agilent InfinityLab Poroshell 120 EC-C18, 3.0 × 100 mm, 2.7 μm (p/n 695975-302)								
Column temperature	45 °C								
Injection volume	10 μL								
Draw and eject speed	100 μL/min; 400 μL/min								
Needle wash	Methanol (MeOH), 5 seconds								
Mobile phases	A) 0.1 % formic acid (FA) + 1 mM ammonium fluoride (NH ₄ F) in water (H ₂ O) B) 0.1 % FA in acetonitrile (ACN)								
Stop time	10 minutes								
Post time	3 minutes								
Flow rate	0.4 mL/min								
Gradient	<table border="1"> <thead> <tr> <th>Time (min)</th> <th>%B</th> </tr> </thead> <tbody> <tr> <td>0</td> <td>30</td> </tr> <tr> <td>8</td> <td>95</td> </tr> <tr> <td>10</td> <td>95</td> </tr> </tbody> </table>	Time (min)	%B	0	30	8	95	10	95
Time (min)	%B								
0	30								
8	95								
10	95								




Method Performance for Online SPE & Direct Injection

MCs	Direct injection		Online SPE	
	LOD (ppt)	LOQ (ppt)	LOD (ppt)	LOQ (ppt)
RR	10	50	0.05	0.2
HtyR	50	100	0.2	1
LR	100	500	0.5	2
HilR	50	100	0.2	1
LA	50	100	0.2	1
LF	50	100	0.2	1



Spike Recovery Experiments for Online SPE & Dir. Injection



Avg. RSD (%)


Offline SPE: 8.7%

Online SPE: 5.1%


Dir. Injection: 4.3%

Offline SPE vs Online SPE time-savings Assuming a batch of 20 samples

Offline Extraction




Online Extraction




Analysis done in less than half the time!

Suspect peaks can show up on LC/MS/MS that can be misidentified

MC-YR identification in lake water sample



RT = 2.24 min
1045.5 -> 135.2
1045.5 -> 213.2
Qual Ratio = 24%




RT = 1.91 min
1045.5 -> 135.2
1045.5 -> 213.2
Qual Ratio = 2%

- Is this MC-YR?
- Is this another MC?
- Can it be something else with same precursor mass of 1045.5?
- How do we find out?






Calibration Curves for Quantification by LC-Q/TOF

0.1 – 25 µg/L





Cpd	R ²	Cal 1 Accuracy
YR	0.9993	118 %
LR	0.9997	110 %
RR	0.9992	119 %




Is it a Microcystin?

Fragmentation – m/z 135



135.0804 = $C_9H_{11}O^+$






135.1168 = $C_{10}H_{15}^+$

- The characteristic ion for microcystins, often used for quantitation, is nominal m/z 135
- There are actually two product ions at m/z 135 that are formed in microcystins from the Adda⁵ group
- The presence in the lake sample indicates this is an unknown microcystin

Is it a Microcystin?

Fragmentation – m/z 213






- m/z 213.0870 corresponds to the sequence Glu⁶ + Mdha⁷, as shown in the spectrum for YR (top)
- This ion is absent in the unknown spectra in lake

The Unknown is NOT MC-YR

What else is different in the Unknown Spectra?

Fragment at m/z 199




- Desmethylation of Mdha⁷ would result in a loss of CH₂ (-14.0157), producing an ion at m/z 199.0713
- The lake spectrum is consistent with a desmethylated microcystin from Mdha⁷, and the compound is not YR


What is generating such fragments from 1045.5 m/z?

How can a methyl group be added elsewhere to a MC-YR like backbone?

Amino Acid substitutions at the R2 & R4 positions result in most MC analogues



Is the unknown a demethylated HtyR?



Microcystin	R ²	R ⁴	Formula	Neutral Mass
LR	Leucine	Arginine	C ₄₉ H ₇₄ N ₁₀ O ₁₂	994.5488
Desmethyl LR	Leucine	Arginine	C ₄₈ H ₇₂ N ₁₀ O ₁₂	980.5331
RR	Arginine	Arginine	C ₄₉ H ₇₅ N ₁₃ O ₁₂	1037.5658
YR	Tyrosine	Arginine	C ₅₂ H ₇₂ N ₁₀ O ₁₃	1044.5280
LA	Leucine	Alanine	C ₄₆ H ₆₇ N ₇ O ₁₂	909.4848
LW	Leucine	Phenylalanine	C ₅₄ H ₇₂ N ₈ O ₁₂	1024.5270
LF	Leucine	Tryptophan	C ₅₂ H ₇₁ N ₇ O ₁₂	985.5161
HtyR	Homotyrosine	Arginine	C ₅₃ H ₇₄ N ₁₀ O ₁₃	1058.5437

Is the unknown demethylated HtyR?


Look for fragments that include R2 (Hty)

Compare ions from authentic standard of HtyR (top) to the unknown (bottom), to show the presence of Hty in the R² position and desmethylation in R⁷ position

- Fragment Ion : R⁷ – Ala¹ – R² – MeAsp³ – Arg⁴

- HtyR (from 1059.5 m/z)
- R⁷ = Mdha; R² = Hty
- Calc m/z : 617.3042
- Measured: 617.3059 (2.8 ppm mass error)


- dm-HtyR (1045.5 m/z)
- R⁷ = Dha; R² = Hty
- Calc m/z : 603.2885
- Measured: 603.2884 (0.2 ppm mass error)



Is the unknown demethylated HtyR?

Look for fragments that include R2 (Hty)

No Demethylation occurring in R³ and no contribution of R⁷ so fragments should be same



- Fragment *Ion* : Ala¹ – R² – MeAsp³ – Arg⁴

- HtyR: R² = Hty
- Calc m/z : 534.2671
- Measured: 534.2692 (2.3 ppm mass error)


- dm-HtyR: R² = Hty
- Calc m/z : 534.2671
- Measured: 534.2683 (4.0 ppm mass error)



Confirmation that unknown is MC-DesMe-HtyR

Additional unknowns detected – Lake A

7 additional new MCs found in lake sample without commercially available standards



RT	Name	Formula	Mass	m/z	Area	Score	Diff (Tgt, ppm)	Ions	Flags (Tgt)	Flag Severity (Tgt)
4.66	MC-YR	C52 H72 N10 O13	1044.5273	1045.5346	1853141	98.48	-0.65	5	multiple IDs	Information
4.69	MC-LR	C49 H74 N10 O12	994.5487	995.5549	36788	90.41	-0.09	3	multiple IDs	Information
4.74	MC-HtyR-DAsP3-Dha7	C51 H70 N10 O13	1030.5117	1031.5185	422076	96.79	-0.65	8	multiple IDs	Information
4.76	MC-HtyR	C53 H74 N10 O13	1058.5434	1059.5501	94974	98.43	-0.26	4		Pass
4.84	MC-DesMe-LR	C48 H72 N10 O12	980.533	981.5405	653497	99.12	-0.07	5	multiple IDs	Information
4.97	MC-HphR-Dha7	C52 H72 N10 O12	1028.533	1029.54	548460	99.56	-0.15	8	multiple IDs	Information
4.99	MC-LR-DAsp3-Dha7	C47 H70 N10 O12	966.516	967.5236	116668	94.3	-1.5	4		Pass
4.99	MC-LR	C49 H74 N10 O12	994.5482	995.5559	32185	92.69	-0.62	4	multiple IDs	Information

- A library of 100+ microcystins was built
- MS domain confirmation
- Auto searching possible using isotope pattern

Conclusions

LC/TQ

- Excellent tool for the targeted analysis of microcystins when analytical standards are available
- Quantification down to sub $\mu\text{g/L}$ easily possible

Online SPE

- Reduces labor and significantly increases throughput with ng/L quantification possible.

LC/QTOF

- Complimentary tool that can confirm suspect compounds using accurate mass of the molecular ion adducts, as well as MS/MS fragments.
- Quantification of all MCs upto $0.1 \mu\text{g/L}$ possible

Databases

- Can be compiled to include the chemical formula for additional microcystins based on reported analogues in the literature

Acknowledgements & Resources

Ralph Hindle

Kathy Hunt

Chang Jiang



Analysis of microcystins and nodularin in drinking water using an Agilent Ultivo triple quadrupole LC/MS



Figure 1. Agilent Ultivo Integrated into LC Stack.

Abstract

Abstract text is present but illegible in the image.

Authors
Tarun Anumol¹,

5991-9087EN



Quantitation of Microcystins in Water by Direct Injection and Online SPE LC/MS/MS Systems

Authors

Chang Jiang, Pei-bin Hu
Agilent Technologies Inc.
Chengdu, China
Tarun Anumol
Agilent Technologies Inc.,
Wilmington, DE

Abstract

This Application Note compares direct injection and online SPE methods using liquid chromatography (LC) coupled with tandem mass spectrometry (MS/MS) to analyze microcystins (MCs) in water. These include MC-RR, MC-HyR, MC-LR, MC-HiR, MC-LA, and MC-LF. The online SPE LC/MS/MS system eliminates time-consuming and laborious offline SPE extraction steps. Both the direct injection and online SPE method achieved limits of detection (LODs) of low or sub-ng/L levels, which are much lower than the provisional guidelines for drinking water expressed by the World Health Organization (WHO).

5994-0007EN



Identification of Unknown Microcystins in Alberta Lake Water

Application Note

Environmental

Abstract

Detection, characterization, and tentative identification of very low levels of unknown microcystins in lake water are possible in the absence of analytical standards using a combination of triple quadrupole LC/MS and LC/Q-TOF analysis and a Personal Compound Database (PCD) compiled from the World Health Organization (WHO) list of microcystins.

Authors

Ralph Hindle
Vogon Laboratory Services Ltd.
Cochrane, Alberta
Canada

Xu Zhang and David Kinniburgh
Alberta Centre for Toxicology
Alberta Health & Wellness
University of Calgary
Calgary, Alberta
Canada

5991-4444EN

Contact: tarun.anumol@agilent.com



Agilent

Trusted Answers