



Comprehensive and Fast Multi-Elemental Analysis Using Inductively Coupled Plasma Optical Emission Spectroscopy (ICP-OES)

Presented by: Mike Mourgas
Applications Scientist
Trace Elemental Analysis

A Complete Portfolio of Elemental Analysis Solutions

Features include:

- Single or multi-element analysis
- Low-cost, high-throughput systems
- Space-saving design
- Low gas consumption



Flame AAS



Graphite Furnace AAS



ICP-OES



Single Quadrupole ICP-MS



Triple Quadrupole ICP-MS



High Resolution ICP-MS

Measure elemental impurities:

- Any sample matrix
- Measurement at concentrations from parts per trillion to percent

Introducing the New Thermo Scientific iCAP PRO Series ICP-OES



Perform Your Environmental Analysis Like a Pro with Thermo Scientific™ iCAP™ PRO Series ICP-OES

Thermo Scientific iCAP PRO Series ICP-OES: Solutions for Every Environmental Sample



iCAP PRO ICP-OES Entry level ICP-OES

- Maximum ease of use for aqueous samples



iCAP PRO X ICP-OES Routine ICP-OES

- Optimized settings for aqueous and organic samples



iCAP PRO XP ICP-OES Workhorse ICP-OES

- Full flexibility for all samples
- Enhanced UV sensitivity



iCAP PRO XPS ICP-OES High Speed ICP-OES

- Best speed and performance for high-throughput laboratories

Four New Models Optimized for Performance, Flexibility and Speed

Thermo Scientific iCAP PRO Series ICP-OES – Core Technologies



Optical System

- New high-energy, compact Echelle polychromator
- Full wavelength range in one exposure
- Resolution < 7 pm at 200 nm



Vertical Torch & POP Interface

- New vertical torch and inner torch box designed for enhanced robustness
- Axial and radial POP nozzles
- New POP window



CID Detector

- New CID 821
- Continuous wavelength coverage
- Immunity from blooming
- Non-destructive readout
- High sensitivity
- Full Frame imaging






Sample Introduction System

- Clip-in components
- Lock-in-place torch
- Large peristaltic pump on all models
- Drain sensor



Vertical Torch and POP Interface Designed for Enhanced Robustness



Vertical Torch & Inner Torch Box

- Vertical torch for all models
- Adjustable radial viewing height for radial and duo instruments
- Removable inner torch box – Designed for corrosion resistance, easy to remove and clean when needed


Optimized Exhaust Flow

- Proprietary exhaust flow design through torch box ensures:
 - Optimal cooling of torch
 - Stability of the plasma
 - Minimal deposition of the torch and center tube

POP interface and window


- Durable ceramic cones – Purge gas exits from optical system through cones to remove interferences that can compromise sensitivity
- New POP window – Protects fore optics from dust, dirt and contamination from the plasma, easy to remove and clean when needed

Optical System and New CID Detector for Enhanced Speed, Sensitivity and Robustness



Improved speed and sensitivity

- **Simultaneous** measurement of the whole spectrum in one acquisition!
- **Full wavelength coverage** – Range of 167-852 nm using iFR Mode or 167-240 nm with eUV mode
- **High quality optics** – Maximum light throughput to maintain speed



Reduced interferences

- **Greater order separation** – Reduced order overlap and interferences (Pb in 10,000 ppm Cu)
- **Resolution** – <7pm at 200 nm
- **Improved stability** – Optimized air flows to prevent sample deposition



Robust optics

- **Mirror protection** – Beam blocker to prevent extreme UV exposure
- **Sealed optics** – Moving parts outside the optical tank
- **Efficient purge** – Just 3.4 L/min

iCAP PRO Series ICP-OES Detector

New Camera – Thermo Scientific CID821

Large chip allows simultaneous readout
Array size of 2048 x 2048 pixels




Improved speed

Pre-exposure not required
Ultra-fast signal readout, 30-40% faster than previous generation detector

Reduced interferences


Maximum order separation (e.g., Cu on Pb)
Anti-blooming advantage over CCD detectors
Off-axis window eliminates optical scatter



iCAP 7000 Series
ICP-OES
CID86


iCAP PRO Series
ICP-OES CID821

Features for Enhanced Ease of Use



New Features

- **New Torch Box Door** – Sliding door for easy access to the torch interface
- **New LED instrument status panel** – Allows user to know status of the instrument or analysis from across the laboratory



Small Footprint and Easy Installation

- **Smallest footprint on the market** – Dimension: 24.2”x27.2”x36.7” (LxWxH), optimizes laboratory bench space
- **Easy access** – Improved connection area, can be pushed close against the wall



Sample Introduction System

- **Sample introduction kits** – Five application specific kits
- **Clip-in sample introduction system** – easy to assemble/disassemble
- **Ceramic torch** – for enhanced robustness when analyzing high matrix sample continuously over long periods

Improved Speed – Reduced Warm-Up Time

Reduced warm-up time


- **Efficiency** – Helps users to increase their lab efficiency
- **Speed** – Fast warm-up with system stability and quick purge
- **Software** – Software stabilization enables fast measurements **from standby**

intelligent Full Range - iFR

- **Standard purge** – Fast warm-up from system standby
- **Quick start up** – Just five minutes from standby to first measurement

enhanced UV - eUV

- **Standard purge for UV** – No special requirements from standby (gas on)
- **Enhanced purge from power off** – Enhanced UV, fast warm-up from power off (no power, no gas)



Ready to go and use
much faster than
iCAP 7000 Series
ICP-OES


Increase your
laboratory
efficiency

Applies to the Thermo Scientific™ iCAP™
PRO XP ICP-OES and the Thermo
Scientific™ iCAP™ PRO XPS ICP-OES

Thermo Scientific™ Qtegra™ Intelligent Scientific Data Solution™ (ISDS) Software

• Benefits of the Qtegra ISDS Software

- Installed on over > 6000 instruments
- Intuitive, streamline workflow platform
- Continues to support existing liquid sampling accessories
- A range of new features added for ease of use
- 21 CFR Part 11 compliance tool set
- Same software platform as the Thermo Scientific ICP-MS instruments for easy cross-training and operation between techniques



Method Validation – A Must for Confirming Method Optimization and Data Quality

- **Method Validation**

- Required to ensure that the method optimized results in accurate, precise and reproducible results.

- **Method validation achieved through:**

- MDL or IDL determinations – DLs determined for the instrument and each sample method
- Calibration curve – correlation coefficient (R^2) of at least 0.995, the closer to 1.0 the better
- Linear Dynamic Range (LDR) determination
- Precision tests – measured by %RSD between 3 sample replicates (short term precision)
- Accuracy tests – matrix spiking of samples and analysis of a Certified Reference Material (CRM)
- Repeat above tests over several days (long term precision and accuracy)
- Quality Control Protocol - built in the Qtegra ISDS Software
 - For routine, continuous monitoring of the analysis
 - Includes automatic actions that can be selected upon failure of a QC standard



Application Data for the iCAP PRO Series ICP-OES

- Analysis of water samples according to USEPA Method 200.7 using the iCAP PRO XPR ICP-OES
- 25% NaCl measurements using the iCAP PRO XP Duo vs. Radial

Analysis of Water Samples by US EPA Method 200.7

Analysis of water samples by USEPA Method 200.7 using the iCAP PRO XPS ICP-OES



- This note describes the analysis of water samples by Method 200.7
- **Challenge:**
 - Compliance with comprehensive method quality control protocol

National Primary Drinking Water Regulations					
Contaminant	MCL (mg/L)	MCLG (mg/L)	Contaminant	MCL (mg/L)	MCLG (mg/L)
Antimony	0.006	0.006	Copper	1.3	1.3
Arsenic	0.01	0	Lead	0.015	0
Barium	2.0	2.0	Mercury	0.002	0.002
Beryllium	0.004	0.004	Selenium	0.05	0.05
Cadmium	0.005	0.005	Thallium	0.002	0.002
Chromium	0.1	0.1	Uranium	0.03	0

National Secondary Drinking Water Regulations	
Contaminant	MCL (mg/L)
Aluminum	0.05 – 0.2
Copper	1
Iron	0.3
Manganese	0.05
Silver	0.1
Zinc	5

Analysis of Water Samples by US EPA Method 200.7

Instrumentation

- iCAP PRO XPS Duo ICP-OES
- Teledyne CETAC ASX-560 autosampler

Samples and Standards

- All calibration and QC standards prepared using 1000 mg/L stock standards in 1.5% HNO₃
- Internal Standard: 5 mg/L Yttrium added online
- Samples - Drinking water, trench water and well water

Analysis

- MDL study
 - Reagent blank spiked at low concentrations
- LDR study
- Interference Study
 - Single element SIC solutions
- Run sequence



Calibration
QCS
IPC
Check Blank
10 Samples
IPC
Check Blank



**152 total samples,
including QC and
calibration
standards**

Parameter	Setting
Pump Tubing	Sample: Tygon® orange/white Drain: Tygon® white/white
Pump Speed	45 rpm
Spray Chamber	Glass Cyclonic
Nebulizer Gas Flow	Glass Concentric
Coolant Gas Flow	12 L/min
Auxiliary Gas Flow	0.5 L/min
Nebulizer Gas Flow	0.5 L/min
Center Tube	2 mm
RF Power	1150 W
Repeats	3
Radial view height	10 mm
Exposure Time	Axial View – 10 seconds Radial View – 10 seconds

*Analysis using iFR mode for all 31 elements

Analysis of Water Samples by US EPA Method 200.7 – Results

MDL and LDR Results

- **MDL**
 - All results shown are below the levels for the National Primary and Secondary Drinking Water standards.
- **LDR**
 - LDR results are more than sufficient for the analysis of typical water samples.

*UMCR 3 – Third Unregulated Contaminant Water Rule

Analyte	Wavelength (nm)	Plasma View	LDR (mg/L)	MDL (µg/L)	Level of Interest (µg/L)
Ag	328.608	Axial	>10	0.84	100
Al	308.215	Radial	>1000	21	50 – 200
As	193.759	Axial	>100	2.1	10
Ba	455.403	Axial	>2	0.47	2000
Be	234.861	Axial	>10	0.08	4
Cd	226.502	Axial	>10	0.25	5
Co	228.616	Axial	>10	0.75	1*
Cr	284.325	Axial	>10	0.29	100 / 0.2*
Cu	224.700	Axial	>10	0.51	1300
Fe	258.940	Radial	>1000	3.7	300
Mn	257.610	Axial	>10	0.06	50
Mo	203.844	Axial	>10	0.90	1*
Pb	220.353	Axial	>100	3.2	15
Sb	206.833	Axial	>100	3.3	5
Se	196.090	Axial	>100	4.8	50
Tl	190.856	Axial	>10	0.50	0.2*
Zn	213.856	Axial	>2	0.02	5000

Analysis of Water Samples by US EPA Method 200.7 – Results

Sample Results


All spike recoveries for the different water samples were within the required recovery range of 85 – 115%.

Analyte	Drinking Water			Trench Water			Well Water		
	Unspiked (mg/L)	Spiked (mg/L)	Recovery (%)	Unspiked (mg/L)	Spiked (mg/L)	Recovery (%)	Unspiked (mg/L)	Spiked (mg/L)	Recovery (%)
As	<MQL	0.200	100.0	<MQL	0.203	101.5	<MQL	0.203	101.5
Ca	40.52	42.43	95.5	49.92	57.39	99.6	46.20	53.66	99.5
Cd	<MQL	0.199	99.5	<MQL	0.196	98.0	0.01	0.199	99.0
Cu	0.024	0.319	98.3	<MQL	0.291	97.0	0.007	0.296	96.3
Fe	0.045	0.239	97.0	1.360	8.701	97.9	27.40	34.82	98.9
Hg	<MQL	0.196	98.0	<MQL	0.196	98.1	<MQL	0.197	98.5
K	2.747	7.795	101.0	12.56	15.31	110.0	1.401	4.116	108.6
Mg	4.271	11.60	97.7	7.863	14.95	9.45	6.953	14.02	94.2
Na	14.24	19.67	108.6	145.31	170.8	102.0	92.85	118.3	101.8
P	0.015	1.644	108.6	0.102	1.730	108.5	1.185	2.742	103.8
Pb	<MQL	0.197	98.5	<MQL	0.192	96.0	0.077	0.266	94.5
Sb	<MQL	0.200	100.0	<MQL	0.195	97.5	<MQL	0.197	98.5
Se	<MQL	0.193	96.5	<MQL	0.193	96.5	<MQL	0.197	98.5
SO4	40.43	43.82	113.2	77.69	93.13	103.1	1.295	16.80	103.5
Tl	<MQL	0.198	99.0	<MQL	0.281	93.7	<MQL	0.283	94.3
Zn	0.0009	0.22	109.6	0.0013	0.22	109.4	0.282	0.48	99.0

Analysis of Water Samples by US EPA Method 200.7

Sample Results

- Stability demonstrated by the analysis of the Instrument Performance Check (IPC) Solution
- 9 IPC solutions analyzed over the run, all results within the 90 – 110% recovery as shown in the graph



Analysis of Water Samples by US EPA Method 200.7 – Summary


- Compliance with the comprehensive Method 200.7 QC protocol for the determination of 31 elements in various water samples was demonstrated using the iCAP PRO XPS Duo ICP-OES.
- The productivity tools of Qtegra ISDS Software combined with the speed of the iCAP PRO XPS ICP-OES using the iFR analysis mode resulted in an analysis time of **1 minute and 58 seconds per sample, or 30 samples per hour.**
- Method 200.7 QC protocol was easily implemented in the LabBook using the built-in QC features of the Qtegra ISDS Software.
- The optimized vertical torch combined with the high-resolution optics minimized physical and spectral interferences.




Analysis of 25% NaCl to Demonstrate High Matrix Tolerance

Challenges:

- Devitrification of quartz glass ICP torch
- Clogging of center tube due to salt deposition
- Clogging of nebulizer
- High RF power, high nebulizer gas flow, etc., required for high matrix samples
- Poor data stability or failure over long runs



(A) Devitrification and (B) premature failure of quartz glass torch



Poor stability/failure in analysis in less than 2 hrs. with normal set-up


Analysis of 25% NaCl – Sample Introduction and Instrumentation

- **Sample introduction system**

- Burgener Mira Mist nebulizer
- Baffled spray chamber
- Sheath gas adaptor with additional gas supply
- Ceramic D-torch


- **iCAP PRO XP ICP-OES**

- Dedicated Radial – traditional choice for high salt solutions
- Duo – New Vertical Torch design on the iCAP PRO Series ICP-OES



A. Ceramic torch outer tube

B. Torch with quartz outer tube for displaying inner ceramic tube



Sheath gas adaptor

Analysis of 25% NaCl - Instrument Parameters

- Both configurations of the iCAP PRO XP ICP-OES were used for the analysis to compare results
- Additional gas is argon for the Sheath Gas Adaptor


Parameter	Setting for iCAP PRO XP ICP-OES Radial	Setting for iCAP PRO XP ICP-OES Duo
Pump Tubing	Sample Tygon® orange/ white	Sample Tygon® orange/ white
	Drain Tygon® white/white	Drain Tygon® white/white
Spray Chamber	Baffled cyclonic	Baffled cyclonic
Nebulizer	Burgener Mira Mist	Burgener Mira Mist
Center Tube	2.0 mm (ceramic)	2.0 mm (ceramic)
Torch	Ceramic D-Torch Radial	Ceramic D-Torch Duo
Pump Speed	45 rpm	45 rpm
Flush Pump Speed	100 rpm	100 rpm
Pump Stabilization Time	10 s	10 s
Wash Time	30 s	30 s
Nebulizer Gas Flow	0.55 L·min ⁻¹	0.55 L·min ⁻¹
Auxiliary Gas Flow	0.5 L·min ⁻¹	1.5 L·min ⁻¹
Coolant Gas Flow	12.0 L L·min ⁻¹	12.0 L·min ⁻¹
Additional Gas	0.15 L L·min ⁻¹	0.15 L·min ⁻¹
RF Power	1400 W	1350 W
Radial Viewing Height	11 mm	11 mm
Exposure Time	iFR 10 s	Radial iFR 10 s, Axial iFR 10 s

Results – Sensitivity and Long-Term Stability Experiment

Method Detection Limits (MDL) obtained from different configurations of the iCAP PRO XP ICP-OES

Element and wavelength (nm)	Radial only MDL ($\mu\text{g}\cdot\text{L}^{-1}$)	Duo – Radial view MDL ($\mu\text{g}\cdot\text{L}^{-1}$)	Duo – Axial view MDL ($\mu\text{g}\cdot\text{L}^{-1}$)
Al 167.079	2.33	7.46	4.91
Ba 455.403	1.17	1.68	0.75
Co 228.616	7.79	4.78	3.37
Cr 205.560	2.73	3.94	1.58
Cu 324.754	7.82	5.65	1.93
Fe 259.940	9.8	13.06	6.58
Mn 257.610	3.88	8.15	2.31
Ni 221.647	4.71	5.81	2.58
Sr 407.771	7.7	8.48	4.41
Zn 213.856	2.54	2.85	0.72

Long-term (18 hours) experiment on the iCAP PRO XP ICP OES, Radial only instrument




Analytes – Al, Ba, Co, Cr, Cu, Fe, Mn, Ni, Sr, Zn

--- ±10%

Results – Long-Term Stability Experiments on the Duo Instrument

Long-term (~10 hours)
stability experiment on
the iCAP PRO XP
ICP- OES Duo



Analysis of 25% NaCl – Summary

- Increased operation time with minimal downtime
- Reduced maintenance requirements and improved robustness – cost effective. No ICP torch failure, no cleaning of torch and cones necessary
- No problem analyzing high matrix samples like 25% NaCl solutions, i.e., seawater-like samples over uninterrupted long time periods (at least 10 hours)
 - Radial only ✓
 - Duo ✓
- Dilution steps for high matrix samples eliminated – vertical ceramic torch and sheath gas adaptor easy to install
- Better detection limits and enhanced sensitivity in Axial mode

Conclusion – Benefits for Environmental Analysis



Speed Sensitivity

- More samples in less time means productivity
- Savings in operational costs
- Meet DL requirements for environmental analysis

Ease of Use

- Easy cross-training
- Less learning curve
- Convenience in operation

Data Quality

- Confidence in reported results
- Correct results the first analysis
- Less sample and QC failures

Robustness

- Stability throughout analytical run
- More uptime
- Less maintenance