


# Critical Review and Screening of Laboratory Supplies for PFAS Analysis in Water Samples

Landon Wiest, Shun-Hsin Liang, Mike Chang; Restek Corporation


## Introduction

The ubiquitous nature of PFAS in the environment makes ensuring a contaminant-free workflow (Figure 1) essential. In this presentation, we demonstrate that Resprep S-DVB SPE cartridges and related sample preparation products are consistently free of background interferences. In addition, a PFAS delay column effectively removes any contamination that may be present in the instrument. Using the materials and procedure presented here, EPA Method 537.1 requirements for cleanliness, accuracy, and precision were reliably met.


**Figure 1: Sample preparation following EPA Method 537.1**



**Figure 2: Reagent Blanks (LRB) from the Entire Extraction Workflow with Multiple Lots of S-DVB SPE**



**Figure 3: Laboratory Fortified Blank (LFB) at Midrange (40ppt)**



**Table I: Precision and Accuracy Results for Method 537.1 PFAS Analysis (n = 4)**

Analyte	%RSD*	Mean Recovery**
Perfluorobutanesulfonic acid (PFBS)	11.9%	91.1%
Perfluorohexanoic acid (PFHxA)	7.96%	99.4%
Hexafluoropropylene oxide dimer acid (HFPO-DA)	6.34%	94.4%
Perfluoroheptanoic acid (PFHpA)	4.19%	92.7%
Perfluorohexanesulfonic acid (PFHxS)	11.9%	89.4%
4,8-Dioxa-3H-perfluorononanoic acid (ADONA)	5.18%	96.6%
Perfluorooctanoic acid (PFOA)	5.21%	91.6%
Perfluorononanoic acid (PFNA)	6.79%	97.2%
Perfluorooctanesulfonic acid (PFOS)	6.78%	87.8%
9-Chlorohexadecafluoro-3-oxanone-1-sulfonic acid (9Cl-PF3ONS)	8.59%	85.1%
Perfluorodecanoic acid (PFDA)	6.96%	93.6%
N-methyl perfluorooctanesulfonamidoacetic acid (N-MeFOSAA)	10.1%	82.8%
N-ethyl perfluorooctanesulfonamidoacetic acid (N-EtFOSAA)	16.5%	106%
Perfluoroundecanoic acid (PFUnA)	2.30%	97.5%
11-Chloroeicosafuoro-3-oxaundecane-1-sulfonic acid (11Cl-PF3OUdS)	5.47%	87.6%
Perfluorododecanoic acid (PFDoA)	5.73%	99.0%
Perfluorotridecanoic acid (PFTrDA)	12.7%	89.1%
Perfluorotetradecanoic acid (PFTA)	8.90%	89.7%

## Methods and Results

- Laboratory reagent blanks were prepared for three different lots of Resprep S-DVB SPE cartridges
- Entire workflow, including manifolds, sample vials, solvents, reservoirs and connectors were free of PFAS contaminants tested in EPA Method 537.1 (Figure 2) via LC-MS/MS on a Shimadzu 8045.
- The use of a PFAS delay column provided excellent chromatography, by retaining any PFAS contaminants in the LC
- Raptor C18 analytical LC column gives excellent selectivity for PFAS compounds (Figure 3).
- Recoveries and %RSDs meet all acceptance criteria for EPA Method 537.1 (Table I)

## Conclusions

- Eliminating PFAS contaminants from workflows is difficult.
- The above conditions and products provide an entire sample preparation workflow and HPLC method that is free of PFAS contaminants.
- The entire workflow meets the criteria for drinking water analysis as set forth by the EPA 537.1 Method.

**PATENTS & TRADEMARKS**  
 Restek patents and trademarks are the property of Restek Corporation. (See [www.restek.com/Patents-Trademarks](http://www.restek.com/Patents-Trademarks) for full list.) Other trademarks in Restek literature or on its website are the property of their respective owners. Restek registered trademarks are registered in the U.S. and may also be registered in other countries.