

US EPA Method 524.2: Software Tools for Successful Method Development and Analysis of Purgeable Organic Compounds in Drinking Water

Anastasia A. Andrianova, Ph.D.
GC/MS Applications Scientist

Bruce D. Quimby, Ph.D.
Senior GC/MS Applications Scientist



Volatile Organic Compound (VOC) Testing

EPA method 8260

Surface and ground water as well as solid sample matrices

EPA method 624

Wastewater samples

The procedures are very similar to 8260, however, the reported compound lists are somewhat different

EPA method 524.2

Drinking water

The conditions used in this method enables a lower detection limit to be reached. This provides results that meet drinking water regulatory limits or MCLs

Purge and Trap vs. Headspace Analysis for VOCs

EPA 524.2 requires the use of purge and trap

METHOD 524.2

MEASUREMENT OF PURGEABLE ORGANIC COMPOUNDS IN WATER BY CAPILLARY COLUMN GAS CHROMATOGRAPHY/MASS SPECTROMETRY

1. SCOPE AND APPLICATION

1.1 This is a general purpose method for the identification and simultaneous measurement of purgeable volatile organic compounds in surface water, ground water, and drinking water in any stage of treatment (1,2). The method is applicable to a wide range of organic compounds, including the four trihalomethane disinfection by-products, that have sufficiently high volatility and low water solubility to be removed from water samples with purge and trap procedures. The following compounds can be determined by this method.

6.2 PURGE AND TRAP SYSTEM - The purge and trap system consists of three separate pieces of equipment: purging device, trap, and desorber. Systems are commercially available from several sources that meet all of the following specifications.

US EPA Method 524.2 VOCs in Water

EPA method 524.2 provides procedures and requirements for the quantitation of VOCs in surface water, ground water, and drinking water by GC/MS

Four specific groups of VOCs

- Trihalomethanes (THMs)
- Gases
- Aromatics
- Products of dehydrohalogenation

Dichlorodifluoromethane	Methyl acrylate	Toluene	1,4-Dichlorobut-2-ene
Chloromethane	Bromochloromethane	trans-1,3-Dichloropropene	Propylbenzene
Chloroethene	Methylacrylonitrile	Ethyl methacrylate	2-Chlorotoluene
Bromomethane	Tetrahydrofuran	1,1,2-Trichloroethane	Mesitylene (1,3,5-Trimethylbenzene)
Ethyl chloride	Trichloromethane	Tetrachloroethylene	tert-Butylbenzene
Trichloromonofluoromethane	1,1,1-Trichloroethane	1,3-Dichloropropane	1,2,4-Trimethylbenzene
Ethyl ether	1-Chlorobutane	2-Hexanone	1-Methylpropyl benzene
1,1-Dichloroethene	Carbon Tetrachloride	Dibromochloromethane	1,3-Dichlorobenzene
Acetone	1,1-Dichloropropene	1,2-Dibromoethane	p-Cymene (4-Isopropyltoluene)
Iodomethane	Benzene	Chlorobenzene	1,4-Dichlorobenzene
Carbon disulfide	1,2-Dichloroethane	1,1,1,2-Tetrachloroethane	1,2-Dichlorobenzene-d4 (SURR)
Allyl chloride	Fluorobenzene (ISTD)	Ethylbenzene	1,2-Dichlorobenzene
Methylene chloride	Trichloroethylene	m+p-Xylene	n-Butylbenzene
Acrylonitrile	1,2-Dichloropropane	o-Xylene	Hexachloroethane
trans-1,2-Dichloroethylene	Dibromomethane	Styrene	1,2-Dibromo-3-chloropropane
Methyl tert-butyl ether	Methyl methacrylate	Tribromomethane	Nitrobenzene
1,1-Dichloroethane	Bromodichloromethane	Isopropylbenzene	1,2,4-Trichlorobenzene
2,2-Dichloropropane	2-Nitropropane	p-Bromofluorobenzene (SURR)	1,1,2,3,4,4-Hexachlorobuta-1,3-diene
cis-1,2-Dichloroethylene	cis-1,3-Dichloropropene	Bromobenzene	Naphthalene
2-Butanone	2,2-Dimethoxybutane	1,1,2,2-Tetrachloroethane	1,2,3-Trichlorobenzene
Propanenitrile	Methyl Isobutyl Ketone (MIBK)	1,2,3-Trichloropropane	

80 target compounds

Parameters for GC, MSD, and P&T Necessary for Success

Hardware:

- GC/MS system coupled with a purge and trap (P&T) concentrator, and a liquid vial autosampler

MS Tuning:

- BFB Autotune


Software tools to streamline:

- Quantitation method development
- Data processing
- Reporting

US EPA Method 524.2 Instrument Configuration

Consumables:

- 1.0 mm inlet liner (5190-4047)
- DB-624UI, 20m x 0.18 mm x 1 μ m column (121-1324UI)
- 6 mm diameter stainless steel drawout lens (G3163-20530)



Tune Evaluation: Meeting BFB Tune Check Criteria

EPA 524.2 requires that the spectrometer must produce a mass spectrum that meets all criteria in Table 3 when 25 ng or less of 4-bromofluorobenzene (BFB) is introduced into the GC.



TABLE 3. ION ABUNDANCE CRITERIA FOR 4-BROMOFLUOROBENZENE (BFB)

Mass (M/z)	Relative Abundance Criteria
50	15 to 40% of mass 95
75	30 to 80% of mass 95
95	Base Peak, 100% Relative Abundance
96	5 to 9% of mass 95
173	< 2% of mass 174
174	> 50% of mass 95
175	5 to 9% of mass 174
176	> 95% but < 101% of mass 174
177	5 to 9% of mass 176

Tune Evaluation Report

Data Path: D:\MassHunter\GCMS\1\data\8860\Purge&Trap Cal_try_2_01-24-2019\blank_ISTD 5ppb.D
 Acq on: 1/24/2019 5:50:34 PM
 Operator:
 Sample: blank_ISTD 5 ppb
 Inst Name: Opportunity 8860
 ALS Vial: 1
 Method:

BFB




Target Mass	Rel. To Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Pass/Fail
50	95	15	40	17.6	1428	Pass
75	95	30	60	49.6	4029	Pass
95	95	100	100	100.0	8125	Pass
96	95	5	9	7.3	591	Pass
173	174	0	2	0.6	41	Pass
174	95	50	100	88.2	7164	Pass
175	174	5	9	7.9	569	Pass
176	174	95	101	99.8	7151	Pass
177	176	5	9	6.8	486	Pass

Tune Evaluation: Meeting BFB Tune Check Criteria





MassHunter Quantitative Analysis in Environmental Quant Mode (EnviroQuant)



Tune Evaluation Report

Data Path: D:\MassHunter\GCMS\1\data\8860\Purge&Trap Cal_try_2_01-24-2019\blank_ISTD 5ppb.D
 Acq on: 1/24/2019 5:50:34 PM
 Operator:
 Sample: blank_ISTD 5 ppb
 Inst Name: Opportunity 8860
 ALS Vial: 1
 Method:

BFB





Target Mass	Rel. To Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Pass/Fail
50	95	15	40	17.6	1428	Pass
75	95	30	60	49.6	4029	Pass
95	95	100	100	100.0	8125	Pass
96	95	5	9	7.3	591	Pass
173	174	0	2	0.6	41	Pass
174	95	50	100	88.2	7164	Pass
175	174	5	9	7.9	569	Pass
176	174	95	101	99.8	7151	Pass
177	176	5	9	6.8	486	Pass

For more information: Agilent application note 5991-0029EN

Tune Evaluation: Meeting BFB Tune Check Criteria



MassHunter Quantitative Analysis in Environmental Quant Mode (EnviroQuant)



Apps



TuneEvaluationUI.exe

A stand-alone Tune Evaluation app (installed automatically with MassHunter Quantitative Analysis)

Tune Evaluation Report

Data Path: D:\MassHunter\GCMS\1\data\8860\Purge&Trap Cal_try_2_01-24-2019\blank_ISTD 5ppb.D
 Acq on: 1/24/2019 5:50:34 PM
 Operator:
 Sample: blank_ISTD 5 ppb
 Inst Name: Opportunity 8860
 ALS Vial: 1
 Method:

BFB




Target Mass	Rel. To Mass	Lower Limit%	Upper Limit%	Rel. Abn%	Raw Abn	Pass/Fail
50	95	15	40	17.6	1428	Pass
75	95	30	60	49.6	4029	Pass
95	95	100	100	100.0	8125	Pass
96	95	5	9	7.3	591	Pass
173	174	0	2	0.6	41	Pass
174	95	50	100	88.2	7164	Pass
175	174	5	9	7.9	569	Pass
176	174	95	101	99.8	7151	Pass
177	176	5	9	6.8	486	Pass

For more information: Agilent application note 5991-0029EN

Total Ion Chromatogram of Method 524.2 50 µg/L Standard, ISTD and Surrogates (5 µg/L)


80 target compounds
2 surrogates
1 ISTD





GC column DB-624 UI
PN 121-1324UI
20 m x 0.18 mm x 1 µm

Setting up EPA 524.2 Analysis with MassHunter




1. A VOC standard was analyzed and the deconvoluted spectra were searched against the NIST spectral library for compound identification using MassHunter Unknowns Analysis



Setting up EPA 524.2 Analysis with MassHunter

1. A VOC standard was analyzed and the deconvoluted spectra were searched against the NIST spectral library for compound identification using MassHunter Unknowns Analysis 
2. A user library was created using MassHunter Library Editor. It included the identified compounds with their retention times and deconvoluted spectra 


Setting up EPA 524.2 Analysis with MassHunter

1. A VOC standard was analyzed and the deconvoluted spectra were searched against the NIST spectral library for compound identification using MassHunter Unknowns Analysis 
2. A user library was created using MassHunter Library Editor. It included the identified compounds with their retention times and deconvoluted spectra 
3. A quantitation method was created from the user spectral library using the compound information and retention times included in the library using MassHunter Quantitative Analysis 
4. Real-world samples were analyzed with the quantitation method. Spectral deconvolution allowed for confirmation of compound identifications.

The concentration of the VOCs in the real-world sample was determined against the initial calibration

Step 1: Analyzing a Standard

Analyze a calibration standard with Unknowns Analysis against the NIST library



Step 1: Analyzing a Standard

Analyze a calibration standard with Unknowns Analysis against the NIST library

Samples

Sample Name	File Name	Components	HR
50 ppb_ISTD 5	50 ppb_ISTD 5ppb.D	603	

Components

File Name	Component RT	Compound Name	Match Factor	Best Hit
50 ppb_ISTD...	1.2055	Dichlorodifluorometh...	96.2	✓
50 ppb_ISTD...	1.3624	Chloromethane	97.5	✓
50 ppb_ISTD...	1.4621	Ethene, chloro-	95.5	✓
50 ppb_ISTD...	1.7522	Methane, bromo-	96.3	✓
50 ppb_ISTD...	1.8575	Ethyl Chloride	95.8	✓
50 ppb_ISTD...	2.1205	Trichloromonofluoro...	96.5	✓
50 ppb_ISTD...	2.4653	Ethyl ether	97.4	✓
50 ppb_ISTD...	2.6828	Ethene, 1,1-dichloro-	96.4	✓
50 ppb_ISTD...	2.7934	Acetone	93.0	✓
50 ppb_ISTD...	2.7934	Hydrogen isocyanate	92.7	✓
50 ppb_ISTD...	2.8410	Methane, iodo-	99.2	✓
50 ppb_ISTD...	2.9071	Carbon disulfide	97.1	✓
50 ppb_ISTD...	3.1450	Allyl chloride	96.2	✓
50 ppb_ISTD...	3.3099	Methylene chloride	97.1	✓
50 ppb_ISTD...	3.3105	Methane-d, trichloro-	70.2	✓
50 ppb_ISTD...	3.6961	2-Propenentile	90.2	✓
50 ppb_ISTD...	3.7090	Ethylene, 1,2-dichloro...	93.5	✓
50 ppb_ISTD...	3.7091	Ethylene, 1,2-dichloro...	97.3	✓
50 ppb_ISTD...	3.7866	2-Ethylhexanal ethyl...	85.9	✓
50 ppb_ISTD...	3.7866	Propane, 2-methoxy...	96.7	✓
50 ppb_ISTD...	4.3701	Ethane, 1,1-dichloro-	92.3	✓
50 ppb_ISTD...	5.1907	Propane, 2,2-dichloro-	89.0	✓
50 ppb_ISTD...	5.2017	Ethylene, 1,2-dichloro...	91.6	✓
50 ppb_ISTD...	5.2019	Ethane, 1,1-dichloro-	95.6	✓
50 ppb_ISTD...	5.2027	Cyclobutane, 1,1-dic...	65.3	✓
50 ppb_ISTD...	5.2960	2-Butanone	90.8	✓
50 ppb_ISTD...	5.2963	...	97.8	✓

Method

Libraries: O:\MassHunter\Library\NIST17.L

Search criteria: Pre-search type: Normal

Forward-Reverse Search: Pure Weight Factor: 0.7

Match factor: Use RT Match

RT penalty function: Trapezoidal, Multiplicative, Additive

RT mismatch penalty: Max RT penalty: 20

RT calibration file: [Empty field]

Library search against NIST


Mass Spectrum

Acquisition Time (min): 11,000 to 17,000

Mass-to-Charge (m/z): 150 to 600

Step 2: Creating a Spectral Library

Analyze a calibration standard with Unknowns Analysis against the NIST library



Step 2: Creating a Spectral Library

Analyze a calibration standard with Unknowns Analysis against the NIST library

The screenshot displays the Agilent MassHunter Unknowns Analysis interface. The 'Samples' table shows 603 components and 105 hits. The 'Components' table lists various compounds with their retention times and areas. A context menu is open over the table, with 'Export...' highlighted. An 'Export Component Table' dialog box is also visible, with 'All components/hits' selected and 'Automatically name compounds for non-hit components' checked. Three spectra are shown: a 'Deconvoluted spectrum' (TIC scan), a 'NIST spectrum' (mass spectrum), and an 'Apex spectrum' (mass spectrum).

Sample Name	File Name	Components	Hits
50 ppb_ISTD 5	50 ppb_ISTD 5ppb.D	603	105

File Name	Component RT	Compound Name	Component Area
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	452380.8
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	408459.8
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	419776.6
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	385388.3
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	260548.9
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	858659.5
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	455981.9
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	848692.4
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	134271.0
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	103826.8
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	449762.6
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	816115.0
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	744666.1
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	714108.8
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	343251.6
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	294670.3
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	824938.9
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	845869.1
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	1010895.0
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	1133587.9
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	608656.9
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	834972.8
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	857980.3
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	874059.6
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	362350.2
50 ppb_ISTD 5	1.2005	Disubstituted benzene...	227610.1
50 ppb_ISTD 5	3.3099	Methylene chloride	...
50 ppb_ISTD 5	3.3105	Methane, trichloro-	...

Export Component Table

Export from:

- All components/hits
- Selected components/hits

Export to:

Library

Automatically name compounds for non-hit components

Prefix: Unknown

- Add index
- Add retention time

OK Cancel



Deconvoluted spectrum

NIST spectrum

Apex spectrum

Step 2: Creating a Spectral Library

Export deconvoluted spectra to the Library Editor and review



Step 3: Creating Quantitative Analysis Method

New method from acquired scan data with library search

Agilent MassHunter Quantitative Analysis (for GCMS and LCMS) - Purge&Trap Cal_try 2_01-24-2019 - App Note_Corrected_Cal_try -2_ADDED more QUAL IONS_P-T_8860-01-24-2019.batch.bin

File Edit View Analyze Method Update Library Report Tools Help

New Method from Acquired Scan Data with Library Search...

Ethene, 1,1-dichloro- ISTD: Benzene, fluoro-

Ethene, 1,1-dichloro- Area	Qualifier (96.0) Re...	Qualifier (98.0) Re...	Qualifier (63.0) Re...	Benzene, fluoro- (ISTD)...	Qualifier (77.0) Re...
982	40.3		32.6	6.613	206229
					1.6
					1.6
					2.0
					1.7
					1.6
					1.8
					1.8

New Method from Acquired Scan Data with Library Search

Library Method:
\\andriano\Documents\2019 Volatiles_Purge&Trap\P&T UA lib method.m
Choose... New... Edit...

Sample Path:
Purge&Trap Cal_try 2_01-24-2019\50 ppb_ISTD 5ppb.D
Browse...

Number of Qualifiers to add:
3

OK Cancel

Step 3: Creating Quantitative Analysis Method

Define ISTD and surrogates

Agilent MassHunter Quantitative Analysis (for GCMS and LCMS) - [New Method]

File Edit View Analyze Method Update Library Report Tools Help

Analyze Batch Layout: [Icons] Restore Default Layout

Method Tasks: New / Open Method Workflow Method Setup Tasks

- Compound Setup
- Retention Time Setup
- ISTD Setup
- Concentration Setup
- Qualifier Setup
- Calibration Curve Setup
- Globals Setup

Save / Exit

- Validate
- Save
- Save As...
- Exit

Manual Setup Tasks

Outlier Setup Tasks

Advanced Tasks

Method Table

Time Segment: <All> Compound: Benzene, fluoro- Reset Table View

Quantifier	Name	RT	Scan	Type	MZ	CAS#	Match Factor
	1-Propene, 1,1-dichloro-	6.038	Scan	Target	75.0	563-58-6	100.0
	Benzene	6.274	Scan	Target	78.0	71-43-2	99.9
	Ethane, 1,2-dichloro-	6.308	Scan	Target	62.0	107-06-2	99.5
	Benzene, fluoro-	6.614	Scan	Target	96.0	462-06-6	100.0
	Trichloroethylene	7.023	Scan	Target	130.0	79-01-6	100.0
	Propane, 1,2-dichloro-	7.261	Scan	ISTD	63.0	78-87-5	100.0
	Methane, dibromo-	7.379	Scan	Surrogate	174.0	74-95-3	100.0
	Methyl methacrylate	7.453	Scan	Matrix Spike	100.0	80-62-6	100.0
	Methane, bromodichloro-	7.568	Scan	Target	83.0	75-27-4	99.2
	Propane, 2-nitro-	7.820	Scan	Target	43.0	79-46-9	99.5
	1-Propene, 1,3-dichloro-, cis-	8.044	Scan	Target	75.0	10061-01-5	100.0
	2,2-Dimethoxybutane	8.105	Scan	Target	89.0	3453-99-4	99.9
	Methyl Isobutyl Ketone	8.237	Scan	Target	58.0	108-10-1	99.7
	Toluene	8.380	Scan	Target	91.0	108-88-3	100.0
	1-Propene, 1,3-dichloro-, trans-	8.619	Scan	Target	75.0	10061-02-6	100.0
	Methacrylic acid, ethyl ester	8.748	Scan	Target	69.0	97-63-2	100.0
	Ethane, 1,1,2-trichloro-	8.797	Scan	Target	97.0	79-00-5	99.8
	Tetrachloroethylene	8.933	Scan	Target	164.0	127-18-4	99.9
	Propane, 1,3-dichloro-	8.961	Scan	Target	76.0	142-28-9	99.8
	2-Hexanone	9.082	Scan	Target	58.0	591-78-6	99.6
	Methane, dibromochloro-	9.181	Scan	Target	129.0	124-48-1	99.6
	Ethane, 1,2-dibromo-	9.284	Scan	Target	109.0	106-93-4	99.7
	Benzene, chloro-	9.787	Scan	Target	112.0	108-90-7	100.0
	Ethane, 1,1,1,2-tetrachloro-	9.875	Scan	Target	133.0	630-20-6	99.1
	Ethylbenzene	9.909	Scan	Target	91.0	100-41-4	100.0
	m-p-Xylene	10.028	Scan	Target	91.0	108-38-3	100.0
	o-Xylene	10.418	Scan	Target	91.0	95-47-6	100.0
	Styrene	10.431	Scan	Target	104.0	100-42-5	100.0
	Methane, tribromo-	10.600	Scan	Target	173.0	75-25-2	100.0
	Benzene, (1-methylethyl)-	10.791	Scan	Target	105.0	98-82-8	100.0
	p-Bromofluorobenzene	10.933	Scan	Target	174.0	460-00-4	100.0
	Benzene, bromo-	11.074	Scan	Target	158.0	108-86-1	100.0
	Ethane, 1,1,2,2-tetrachloro-	11.083	Scan	Target	83.0	79-34-5	100.0

Step 3: Creating Quantitative Analysis Method

Add target deconvolution (optionally)

Agilent MassHunter Quantitative Analysis (for GCMS and LCMS) - Method - <C:\Users\andriano\Documents\2019 Volatiles_Purge&Trap\Purge&Trap Cal_try 2_01-24-2019\>

File Edit View Analyze Method Update Library Report Tools Help

Analyze Batch Layout: Restore Default Layout

Method Tasks

- New / Open Method
- Workflow
 - Target Deconvolution Setup
 - Screening - GC
 - Screening - LC
- Method Setup Tasks
 - Compound Setup
 - Retention Time Setup
 - ISTD Setup
 - Concentration Setup
 - Qualifier Setup
 - Calibration Curve Setup
- Globals Setup
- Save / Exit
 - Validate
 - Save
 - Save As...
- Exit
- Manual Setup Tasks
- Outlier Setup Tasks
- Advanced Tasks

Method Table

Target Deconvolution Setup

Reference Library:
Setup Reference Library...

Library Method:
C:\Users\andriano\Documents\2019 Volatiles_Purge&Trap\P&T UA lib method.uamethod.xml
Edit... New... Choose...

Spectrum Setup:
 Deconvoluted scan as Spectrum Extraction Override
 Show reference spectrum
 Show override spectrum
 Show match scores

Outlier Setup:
Library Match Score Minimum: 60
Min. Percent Purity: 60
Alternative Peak Criteria: Deconvoluted Library Match Score

OK Cancel

Step 3: Creating Quantitative Analysis Method

Create calibration levels from the analyzed calibration standards

Sample			
Name	Data File	Type	Level
blank_ISTD 5 ppb	blank_ISTD 5ppb.D	Blank	
0-5 ppb_ISTD 5 ppb	0-5 ppb_ISTD 5ppb.D	Cal	0.5
0-25 ppb_ISTD 5 ppb	0-25 ppb_ISTD 5ppb.D	Cal	0.25
1 ppb_ISTD 5 ppb	1 ppb_ISTD 5ppb.D	Cal	1
5 ppb_ISTD 5 ppb	5 ppb_ISTD 5ppb.D	Cal	5
10 ppb_ISTD 5 ppb	10 ppb_ISTD 5ppb.D	Cal	10
25 ppb_ISTD 5 ppb	25 ppb_ISTD 5ppb.D	Cal	25
50 ppb_ISTD 5 ppb	50 ppb_ISTD 5ppb.D	Cal	50

Step 3: Creating Quantitative Analysis Method

Create calibration levels from the analyzed calibration standards

The screenshot shows the Agilent MassHunter software interface. The 'Method' menu is open, and the 'Create Levels from Calibration Samples' option is highlighted with a yellow box. A yellow arrow points to the 'Method' menu. The 'Method Table' window is visible, showing a table of calibration standards.

Name	Data
blank_ISTD 5 ppb	blank_ISTD 5ppb
0-5 ppb_ISTD 5 ppb	0-5 ppb_ISTD 5ppb
0-25 ppb_ISTD 5 ppb	0-25 ppb_ISTD 5ppb
1 ppb_ISTD 5 ppb	1 ppb_ISTD 5ppb
5 ppb_ISTD 5 ppb	5 ppb_ISTD 5ppb
10 ppb_ISTD 5 ppb	10 ppb_ISTD 5ppb
25 ppb_ISTD 5 ppb	25 ppb_ISTD 5ppb
50 ppb_ISTD 5 ppb	50 ppb_ISTD 5ppb

Sample					
Name	Data File	Type	Level	Acq. Method File	Ac
1 ppb_ISTD 5 p...	1 ppb_ISTD 5pp...	Cal	1	BFB_AA_RTL_u...	1/2

Quantifier					
Name	TS	Scan	Type	CF	
Dichlorodifluoro...	1	Scan	Target	Linear	
Chloromethane	1	Scan	Target	Linear	
Ethene, chloro-	1	Scan	Target	Linear	
Methane, bromo-	1	Scan	Target	Linear	
Ethyl Chloride	1	Scan	Target	Linear	
Trichloromonoflu...	1	Scan	Target	Linear	
Ethyl ether	1	Scan	Target	Linear	
Ethene, 1,1-dichl...	1	Scan	Target	Linear	
Acetone	1	Scan	Target	Linear	
Methane, iodo-	1	Scan	Target	Quadratic	
Carbon disulfide	1	Scan	Target	Linear	
Allyl chloride	1	Scan	Target	Linear	
Methylene chlori...	1	Scan	Target	Linear	
Acrylonitrile	1	Scan	Target	Linear	
Ethylene, 1,2-dic...	1	Scan	Target	Linear	
Methyl tert-butyl...	1	Scan	Target	Linear	
Ethane, 1,1-dichl...	1	Scan	Target	Linear	
Propane, 2,2-dic...	1	Scan	Target	Linear	
Ethylene, 1,2-dic...	1	Scan	Target	Linear	
2-Butanone	1	Scan	Target	Linear	
Propanenitrile	1	Scan	Target	Linear	

Step 3: Creating Quantitative Analysis Method

Create calibration levels from the analyzed calibration standards


Quantifier						
Name	TS	Scan	Type	MZ	RT	
Ethene, 1,1-dichl...	1	Scan	Target	61.0	2.683	

Qualifier		
MZ	Rel. Resp.	Uncertainty
96.0	66.3	100.0
98.0	41.3	100.0
63.0	32.6	100.0

Calibration		
Level	Conc.	Response
0.25	0.2500	2296
0.5	0.5000	5760
1	1.0000	11233
5	5.0000	54171
10	10.0000	103157
25	25.0000	262043
50	50.0000	518953

The method is ready to analyze calibration standards and samples

ICAL Review in MassHunter Quantitative Analysis 10.1



Initial Calibration (ICAL)

0.25–50 µg/L for 68 compounds (85%)

0.50–50 µg/L for 74 compounds (93%)

79 compounds out of 80 met the EPA criteria of less than 20% RSD with linear calibration

A quadratic fit with R^2 of 0.9963 was used for iodomethane

The %RSDs for the internal standard and surrogate compounds introduced by the AQUATek LVA were <5% RSD

Table 4. ICAL for Method 524.2 From 0.25 to 50 µg/L.

Compound	RT (min)	0.25 µg/L	0.5 µg/L	1 µg/L	5 µg/L	10 µg/L	25 µg/L	50 µg/L	Avg RRF	%RSD
		RRF	RRF	RRF	RRF	RRF	RRF	RRF		
Fluorobenzene (ISTD)	6.613	ISTD	ISTD	ISTD	ISTD	ISTD	ISTD	ISTD	1.000	0.6
Dichlorodifluoromethane	1.207	0.087	0.150	0.178	0.136	0.157	0.153	0.153	0.145	19.6
Chloromethane	1.362	0.148	0.178	0.143	0.159	0.170	0.162	0.161	0.160	7.5
Chloroethene	1.462	0.244	0.182	0.212	0.180	0.191	0.183	0.184	0.196	12.0
Bromomethane	1.752				0.149	0.145	0.119	0.117	0.132	12.8
Ethyl Chloride	1.858	0.066	0.118	0.104	0.099	0.115	0.113	0.112	0.104	17.2
Trichloromonofluoromethane	2.120	0.299	0.353	0.342	0.337	0.334	0.337	0.336	0.334	5.0
Ethyl Ether	2.466	0.062	0.096	0.115	0.099	0.105	0.105	0.106	0.098	17.4
1,1-Dichloroethene	2.683	0.223	0.280	0.274	0.262	0.252	0.255	0.250	0.257	7.3
Acetone	2.793				0.021	0.022	0.020	0.020	0.021	5.1
Iodomethane	2.841	0.186	0.160	0.177	0.171	0.209	0.275	0.297	0.211	0.9963*
Carbon Disulfide	2.907	0.603	0.515	0.516	0.477	0.489	0.483	0.489	0.510	8.5
Allyl Chloride	3.145	0.074	0.089	0.071	0.091	0.099	0.089	0.089	0.086	11.8
Methylene Chloride	3.310	0.225	0.182	0.178	0.179	0.174	0.171	0.171	0.183	10.5

Method Detection Limits (MDLs)

<0.10 µg/L for 58 compounds (73%)

<0.15 µg/L for 66 compounds (83%)

8 trials at 0.25 µg/L for 66 compounds

For the compounds with higher reporting limits – 8 trials at 0.5 and 1 µg/L

$$\text{MDL} = s \times t_{(n-1, 1-\alpha=99)} = s \times 2.998$$

- 1.2 Method detection limits (MDLs) (3) are compound, instrument and especially matrix dependent and vary from approximately 0.02 to 1.6 µg/L.

Table 5. Calculated MDLs for VOCs.

MDLs in pg level


Compound	RT (min)	Calculated Concentration in the Sample (µg/L)										Average concentration (µg/L)	SD	MDL
		Spike (µg/L)	Samp. 1	Samp. 2	Samp. 3	Samp. 4	Samp. 5	Samp. 6	Samp. 7	Samp. 8				
Fluorobenzene (ISTD)	6.613	5.00	5.00	5.00	5.00	5.00	5.00	5.00	5.00	5.00	5.00	5.00	N/A	N/A
Dichlorodifluoromethane	1.207	0.25	0.23	0.24	0.24	0.24	0.26	0.22	0.24	0.23	0.24	0.24	0.010	0.031
Chloromethane	1.362	0.25	0.26	0.26	0.25	0.24	0.23	0.29	0.25	0.26	0.26	0.26	0.017	0.052
Chloroethene	1.462	0.25	0.17	0.17	0.14	0.14	0.18	0.15	0.16	0.19	0.16	0.16	0.017	0.051
Bromomethane	1.752	1.00	0.88	1.03	0.90	0.91	0.78	0.74	0.75	0.97	0.87	0.106	0.316	
Ethyl Chloride	1.858	0.25	0.36	0.34	0.29	0.35	0.29	0.38	0.37	0.37	0.34	0.035	0.103	
Trichloromonofluoromethane	2.120	0.25	0.25	0.24	0.18	0.23	0.18	0.26	0.20	0.21	0.22	0.030	0.091	
Ethyl Ether	2.466	0.25	0.28	0.21	0.29	0.27	0.33	0.23	0.27	0.30	0.27	0.038	0.114	
1,1-Dichloroethene	2.683	0.25	0.24	0.24	0.27	0.23	0.24	0.21	0.24	0.25	0.24	0.019	0.057	
Acetone	2.793	1.00	1.04	1.21	1.51	1.33	1.37	1.25	1.03	1.08	1.23	0.173	0.518	

MDLs for Gases

$$\text{MDL} = s \times t_{(n-1, 1-\alpha=99)} = s \times 2.998$$


Table 5. Calculated MDLs for VOCs.

	Compound	RT (min)	Average concentration (µg/L)	SD	MDL
1	Dichlorodifluoromethane	1.207	0.24	0.010	0.031
2	Chloromethane	1.362	0.26	0.017	0.052
3	Chloroethene	1.462	0.16	0.017	0.051
4	Bromomethane	1.752	0.87	0.106	0.316
5	Ethyl Chloride	1.858	0.34	0.035	0.103
6	Trichloromonofluoromethane	2.120	0.22	0.030	0.091




US EPA Method 524.2: Real World Samples

Tested tap water sample from Eastern Pennsylvania: MassHunter Quantitative Analysis




US EPA Method 524.2: Real World Samples

Tested tap water sample from Eastern Pennsylvania: MassHunter Quantitative Analysis




US EPA Method 524.2: Real World Samples

Tested tap water sample from Eastern Pennsylvania: MassHunter Unknowns Analysis







VOCs found:




Components						
Component RT	Compound Name	Match Factor	Best Hit	Formula	Library RT	Delta RT
2.1257	Trichloromonofluoromethane	71.8	<input checked="" type="checkbox"/>	CCl3F	2.1205	-0.0052
5.6403	Trichloromethane	95.2	<input checked="" type="checkbox"/>	CHCl3	5.6396	-0.0007
6.6152	Benzene, fluoro-	97.5	<input checked="" type="checkbox"/>	C6H5F	6.6141	-0.0011
7.5693	Methane, bromochloro-	95.3	<input checked="" type="checkbox"/>	CHBrCl2	7.5681	-0.0012
8.9355	Tetrachloroethylene	91.9	<input checked="" type="checkbox"/>	C2Cl4	8.9325	-0.0029
9.1814	Methane, dibromochloro-	94.3	<input checked="" type="checkbox"/>	CHBr2Cl	9.1810	-0.0004
10.6010	Methane, tribromo-	75.5	<input checked="" type="checkbox"/>	CHBr3	10.5996	-0.0014
10.9340	p-Bromofluorobenzene	98.3	<input checked="" type="checkbox"/>	C6H4BrF	10.9330	-0.0010
12.4527	1,2-Dichlorobenzene-D4	87.0	<input checked="" type="checkbox"/>	C6D4Cl2	12.4523	-0.0004



Trichloromonofluoromethane Found in a Tap Water Sample from Eastern Pennsylvania: Review in Unknowns Analysis





Trichloromonofluoromethane Found in a Tap Water Sample from Eastern Pennsylvania: Review in Unknowns Analysis



Library spectrum

Components

Component RT	Compound Name	Match Factor	Best Hit	Formula	Library RT	Delta RT
2.1257	Trichloromonofluoromethane	71.8	<input checked="" type="checkbox"/>	CCl3F	2.1205	-0.0052



US EPA Method 524.2: Real World Samples

VOCs found in a tap water sample from four sources in Pennsylvania

Compound	RT (min)	Concentration ($\mu\text{g/L}$)			
		Southern Pennsylvania	Eastern Pennsylvania	Southeastern Pennsylvania	City of Philadelphia
Trichloromonofluoromethane	2.120		0.30		
Trichloromethane (Chloroform)	5.640	1.05	7.15	12.56	14.06
Bromodichloromethane	7.568		5.15	4.81	5.77
Toluene	8.380	0.29			
Tetrachloroethylene	8.933		0.36		
Dibromochloromethane	9.181		4.49	1.03	1.44
Tribromomethane	10.600		1.26		

Summary

The software tools offer simple and rapid means to streamline quantitation method development, data processing, and reporting for the EPA 524.2 Method via:

- Built-in software tools for BFB tune check evaluation
- Creating user library of VOCs with deconvoluted spectra and retention times
- Creating a quantitation method from the user spectral library
- Analysis of the real-world samples using spectral deconvolution

ICAL was performed over the range of 0.25 to 50 $\mu\text{g/L}$

Several VOCs were identified and quantified in **real world tap water samples** at concentrations varying over the range of 0.3 to 14.1 $\mu\text{g/L}$, in many cases much lower than current EPA 524.2 MCLs

EPA Method 524.2 with Agilent 8860/5977B

Application note 5994-0833EN

- Demonstrated applicability of an Agilent 8860/5977B GC/MSD system coupled with a Teledyne Tekmar Lumin purge and trap (P&T) concentrator, and an AQUATek LVA (liquid vial autosampler)
- ICAL range 0.25–50 µg/L
- MDLs in ppt levels
- Deconvoluted spectra for compound identification
- Several VOCs were identified and quantified in real-world drinking water samples

Application Note
Environmental



US EPA Method 524.2: Successful Measurement of Purgeable Organic Compounds in Drinking Water by Agilent 8860/5977B GC/MSD

Authors

Bruce D. Quimby and
Anastasia A. Andrianova
Agilent Technologies, Inc.

Abstract

An Agilent 8860/5977B GC/MSD system coupled with a Teledyne Tekmar Lumin purge and trap (P&T) concentrator and an AQUATek liquid vial autosampler (LVA) was successfully used for the analysis of volatile organic compounds (VOCs) to the requirements of United States Environmental Protection Agency (US EPA) method 524.2. The analysis of VOCs in water following this or similar methods are widely used around the world as part of insuring the safety of potable water supplies.

Thank you!