

# High Throughput Analysis of Drinking Water by ICP-MS According to EPA Method 200.8

Siqi Sun, PhD  
Analytik Jena US LLC



## Method: EPA 200.8

- Drinking water is the world's most important resource
- Quality control for drinking water is regulated by international and national regulations and norms
- The U.S. Environmental Protection Agency (EPA) has released the Method 200.8, which specifies criteria for the determination of trace elements in waters and wastes by ICP-MS
- The method applies to the determination of the 21 listed elements in ground water, surface waters, and drinking water.

### Analyte

---

|            |      |
|------------|------|
| Aluminum   | (Al) |
| Antimony   | (Sb) |
| Arsenic    | (As) |
| Barium     | (Ba) |
| Beryllium  | (Be) |
| Cadmium    | (Cd) |
| Chromium   | (Cr) |
| Cobalt     | (Co) |
| Copper     | (Cu) |
| Lead       | (Pb) |
| Manganese  | (Mn) |
| Mercury    | (Hg) |
| Molybdenum | (Mo) |
| Nickel     | (Ni) |
| Selenium   | (Se) |
| Silver     | (Ag) |
| Thallium   | (Tl) |
| Thorium    | (Th) |
| Uranium    | (U)  |
| Vanadium   | (V)  |
| Zinc       | (Zn) |

---


## Challenge: **Speed vs Precision**

---

- For labs with high sample throughput, the more samples analyzed per hour the lower the cost per sample
  - Accuracy and precision often come at the cost of speed
  - How can we speed up our analyses without sacrificing the required precision?
- **With Higher Sensitivity**
- Higher sensitivity translates to more counts/second
  - More counts/second means fewer seconds to achieve the required counts for precise analysis
- With the PlasmaQuant MS, >82 samples per hour was achieved with the method required precision

## Sample Introduction

- The PlasmaQuant® MS was equipped with:
  - ▶ oneFAST, ESI introduction system
  - ▶ ASX 560 (240 sample capacity)
  - ▶ standard aqueous sample nebulizer, spray chamber, and torch
  
- Non-measurement steps (sample uptake, washout, autosampler movement) are minimized using the ESI oneFAST system
  - Provides a significant boost to sample throughput



## Experiment Parameters


| Parameter            | Specification                |
|----------------------|------------------------------|
| Plasma Gas Flow      | 9 L/min                      |
| Auxiliary Gas Flow   | 1.45 L/min                   |
| Nebulizer Gas Flow   | 1.01 L/min                   |
| Spray Chamber Temp   | 3°C                          |
| RF Power             | 1450 W                       |
| Sampling Depth       | 5.0 mm                       |
| Dwell Time           | 20 ms (50 ms for Be, As, Se) |
| Scans per Replicate  | 7 (peak hopping, 1 pt/peak)  |
| No. of Replicates    | 6                            |
| Stabilization Delay  | 19 s                         |
| Sample Load Time     | 7 s                          |
| iCRC Gas Flow        | He, 120 mL/min               |
| Detector Attenuation | None                         |

■ Total argon usage is <12 L/min

■ Dwell time is increased for ions with high ionization energies or low analyte isotope abundance

■ iCRC technology moves interference management at the skimmer cone and ahead of the ion optics

➤ Increases quadrupole cleanliness and lifetime




## The Instrument: PlasmaQuant MS

- High throughput analysis was carried out with PQ MS
- Discrete Dynode Electron Multiplier (DDEM) as the unique all digital detector
  - In pure 'pulse-counting' mode
  - 10 orders of linear dynamic range
  - No separate analog electronics
  - No regular cross-calibrations
- At high concentration, signal can be automatically attenuated

| Isotope                               | Sensitivity of PQ MS Mcps/ppm | Sensitivity of PQ MS Elite Mcps/ppm |
|---------------------------------------|-------------------------------|-------------------------------------|
| <sup>7</sup> Li                       | 65                            | 250                                 |
| <sup>9</sup> Be                       | 20                            | 50                                  |
| <sup>59</sup> Co                      | 400                           | 1000                                |
| <sup>115</sup> In                     | 500                           | 1500                                |
| <sup>205</sup> Tl                     | 200                           | 750                                 |
| <sup>232</sup> Th                     | 300                           | 1000                                |
| <sup>238</sup> U                      | 300                           | 1000                                |
|                                       |                               |                                     |
| <b>CeO<sup>+</sup>/Ce<sup>+</sup></b> | <b>&lt; 2%</b>                | <b>&lt; 2%</b>                      |
| <b>Ce<sup>2+</sup>/Ce<sup>+</sup></b> | <b>&lt; 2%</b>                | <b>&lt; 2%</b>                      |
|                                       |                               |                                     |
| <b>Background</b>                     | <b>&lt; 0.5 cps</b>           | <b>&lt; 1 cps</b>                   |

## The Instrument: PlasmaQuant MS

- High throughput analysis was carried out with PQ MS
- Discrete Dynode Electron Multiplier (DDEM) as the unique all digital detector
  - In pure 'pulse-counting' mode
  - 10 orders of linear dynamic range
  - No separate analog electronics
  - No regular cross-calibrations
- At high concentration, signal can be automatically attenuated



## Samples preparation


- Samples are tap water from Jena, Germany as well as two certified reference materials
- Only high purity reagents were used in the sample preparation
  - 1:2 dilution in deionized water
  - Samples were acidified to contain 1% nitric acid
- Calibration solutions were prepared from
  - multi-element standard
  - single-element standards for Ag, Sb, Hg






## Accuracy

- Certified reference materials were measured in order to verify the accuracy of the method
- Results were within 91 - 103 % of the specified value.
- Additionally, two lab-fortified matrices (+1 ppb and +10 ppb) were measured.
- The recovery rates of the LFM's were between 91 - 104 %.




# Precision



- The relative standard deviation (RSD) was used to assess measurement precision
- An average RSD of 1.5% was achieved
- Precision and accuracy of the certified reference materials and lab-fortified matrices are shown
- RSD is <2 % for low concentrations, <1 % for higher concentrations
- The required RSD <5 % for a 10 or 100 ppb tuning solution standard (depending on sensitivity) was achieved for all elements verifying instrument stability

## Speed

- On the PlasmaQuant MS, 60 samples per hour can be measured with very high precision (mean RSD <1.5 %)
- Results surpass requirements of the U.S. EPA method 200.8, meeting even stricter requirements and more demanding regulations.
- If that level of precision is not the main priority, the spectrometer's unmatched sensitivity allows to further boost sample throughput to **>80 samples per hour** with still competitive precision (mean RSD 2.2 %)



## Conclusion

- Unmatched sensitivity of the PlasmaQuant MS enables high sample throughput without sacrificing precision
- The requirements of the U.S. EPA 200.8 regulation are surpassed by far and the performance is maintained for hours proving the instrument's robustness
- Consuming only 50 % argon volume compared to a conventional ICP-MS, and requiring considerably less consumables and maintenance.
- Minimal running costs and high sample throughput gives the lowest costs per sample, making the PlasmaQuant MS the ideal solution for laboratories who need to routinely measure a large number of samples.

