

On-line sample concentration coupled with LC-MS/MS Orbitrap HRMS for the EU Water Framework Directive: Target quantitation and non-targeted screening of PFAS and other priority environmental contaminants in water

John Quick¹, Ed George², Neville Llewellyn², Olaf Scheibner², Cristina Jacob², and Martin

Jech²

¹ALS - Coventry UK; ²ThermoFisherScientific



RIGHT SOLUTIONS · RIGHT PARTNER



Outline

- Introduction to ALS
- Background to UK Chemical Investigation Programme (CIP)
- Past experience with Exactive–Plus for CIP (phase 2)
- Development of a new method for CIP (phase 3) with Q Exactive Focus and Exploris 120
- Non-Target Screening
- Conclusions & future plans



ALS Global locations



Right Solutions • Right Partner

ALS Life Sciences in UK

- Food and Pharma (1000 staff)
 - Chatteris HUB laboratory
 - Carlisle
 - Mirfield
 - Ely
 - Shrewbury
 - Rotherham
 - Trowbridge
 - Sittingbourne
 - Newton Abbot
 - Clonmel (Ireland)
- Environmental Water (500 staff)
 - Coventry HUB laboratory
 - Wakefield
 - Otterbourne
 - Glasgow
 - Dublin (Ireland sales office)
- Environmental Land (250 staff)
 - Chester/Hawarden HUB laboratory
 - Aberdeen



ALS



ALS - Background

- ALS Environmental
 - Busy, high throughput commercial environment
 - Efficient, robust methods required
- 2015 Set of low level methods developed to meet requirements of phase 2 of the UK Chemical Investigation Program (CIP)
 - Ultra-trace level work low pg/L levels required for some compounds
 - Investment in modern instrumentation including Exactive-Plus
- 2016/2018 CIP methods extended to fully cover WFD suite and matrices
 - Accredited to 17025
 - 2019 Successful completion of CIP2 project
 - 2020 Start of phase 3 of CIP
 - New compounds added
 - Larger numbers of more "challenging" matrix types
 - Purchase of Thermo QE Focus

UK Chemical Investigation Program



- The UKWIR Chemicals Investigation Programme (CIP) is the UK water Industry's response to current and emerging legislation on trace substances in the water environment. It brings together the water and wastewater companies in England and Wales with the various regulators in a collaborative programme.
- Phase 1 2010 to 2014
 - £25 million, 46 chemicals, effluent screening, process sampling, sewerage catchment sampling.
- Phase 2 2015 to 2020
 - £140 million, 74 chemicals, 600 STW sites, over 3 million determinations
 - Pharmaceuticals, pesticides, personal care products, EDCs, flame retardants & industrial chemicals.
 - Focus on effluent impact assessment.
- Phase 3 2020 to 2022
 - Builds on CIP2, more new and emerging compounds, trend monitoring, microplastics and AMR also included.
 - Split into 14 work packages (CHEM1 to CHEM14)
 - Start delayed until September 2020 due to COVID.

CIP 2 (2015) METHOD - EXACTIVE PLUS/EQUAN

- Analytes
 - PFOS, PFOA, DEHP, Triclosan
- Matrices
 - Surface waters, sewage effluents
- Method- RSH EQuan 850 Autosampler for On-line Sample Preparation
 - Load 5mls of sample onto Hypersil Gold AQ column (20 x 2.1mm)
 - Trapping column to reduce PFOA background
 - Backflush onto Phenylhexyl analytical column (100 x 2.1mm)
 - Mobile phase Methanol/Water with formic acid and ammonium formate
 - Gradient elution, 12 minute runtime
 - Exactive Plus Simultaneous pos/neg acquisition @ 70,000 resolution
- Performance
 - All targets for uncertainty and LODs were met.
 - Sub-ppt LODs for PFOS/PFOA obtained
 - Method ran robustly throughout the duration of the project



What is TriPlus RSH EQuan 850 Autosampler?



Automated high-throughput LC-MS solution for the analysis of contaminants in

environmental water, drinking water and beverages at low ng/L levels

- Online sample pre-concentration
 - $\,\circ\,$ 2 LC pumps: Loading and Eluting
 - $\,\circ\,$ 2 injection valves and 1 switching valve
 - o 2 LC columns
- High injection volumes
 - o 1 20 mL
- Standard injection volumes
 - $o \ 1{-}100 \ \mu L$



Find out more at thermofisher.com/EQuan850

Right Solutions • Right Partner

8

Offline vs. Online Sample preparation



2 days Ittration Condition Load sample Wash Euton Evaporation to dryness Reconstitution LCMS/MS (LCMS/MS Ittration Ittration Ittration Ittration Ittration Ittration Ittration Final Association to a sample Preparation Interfering compounds Ittration Ittration Ittration Ittration Analysis time is reduced from hours to cninutes Ittration Ittration Ittration

Conventional Method: Offline Sample Preparation

9 Proprietary & Confidential

Online Sample Pre-Concentration: How it Works?



EQuan Valve Standard Valve Standard Valve EQuan Valve Step 3: Elute sample from trap /aste Waste column Waste Analytical colun In from In from Loading Pump (000)Loading Pump (Quaternary) C (Quaternary) Large sample loop Large sample loop Waste Waste Step 2: Load sample onto trap In from Eluting In from Eluting column Pump (binary) Pump (binary) Switching Valve Switching Valve EQuan Valve Standard Valve **EQuan Valve Standard Valve** Step 1: Inject and load the sample into large sample loop Waste Waste up to 20 mL Analytical co In from In from Loading Pump Loading Pump 0 (000)0 (0000) (Quaternary) (Quaternary) Large sample loop Large sample loop Waste Waste Step 4: Elute sample from In from Eluting In from Eluting analytical column Pump (binary) Pump (binary) Switching Valve Switching Valve

EQuan Injection Mode: Red Flow Path

10 10

CIP2 – PFOS on EXACTIVE–PLUS



Sensitivity at 0.2ng/L



PFOS in Sewage Effluent – 1.8ng/L

R	r: 0.00 100 _∃	-12.01 SM: 7G	6.40	NL: 1.26E6				
	95			498.9252-498.9352 F: FTMS - p ESIFul ms				
	85			[100.00-1500.00] MS Stanley_FE_Blank_04				
	80	PFOS at 1.8ng/L in Final Effluent						
	75							
	65							
-	60		1					
- Provide	55							
A A A	45		()					
o la C	40							
	35							
	25							
	20							
	15							
	5							
	여루	1 2 3 4 5		12				
1	Time (min)							

Calibration curve – 0.2 – 100.0ng/L



Performance Statistics

ANALYTE	PFOS	PFOA		
UNITS	ng/L	ng/L		
Method LOD	0.0676	0.0539		
Method MRL	0.09	0.09		
EFFLUENT				
Spike Recovery	96.9%	94.8%		
Spike RSD	4.15%	2.38%		
Uncertainty	11.41%	10.00%		
RIVERINE				
Spike Recovery	94.4%	95.9%		
Spike RSD	5.53%	5.17%		
Uncertainty	16.65%	14.47%		



CIP – Phase 3

- CHEM 14 work package new and emerging contaminants
- Matrices surface waters, sewage effluents & sewage influents
- LC-MS amenable compounds shown below

		Required I	LOD (ng/L)
Name	CAS	STW Effluent	STW Influent
Ranitidine	66357-35-5	100	1000
Benzotriazole	95-14-7	100	1000
Imidacloprid	138291-41-3	1	10
Bisphenol S	80-09-1	10	100
Climbazole	38083-17-9	10	100
Perfluorooctanoic acid	335-67-1	0.1	1
Perfluorooctanesulfonic acid	1763-23-1	0.1	1
Perfluoropentanoic acid	2706-90-3	1	10
Perfluorobutanesulfonic acid	375-73-5	1	10
Perfluorohexanoic acid	307-24-4	1	10
GEN-X	62037-80-3	1	10
Perfluoroheptanoic acid	375-85-9	1	10
Perfluorohexanesulfonic acid	355-46-4	1	10
6:2 Fluorotelomer Sulfonamide Alkylbetaine (6:2 FTAB)	34455-29-3	1	10
6:2 Fluorotelomer Sulfonate (6:2 FTS)	27619-97-2	1	10
Fipronil	120068-37-3	10	100
Triclocarban	101-20-2	100	1000

Q EXACTIVE FOCUS



- Justification
 - Development of extended method for CIP phase 3 required
 - Options to improve sensitivity and selectivity over Exactive Plus targeted SIM, PRM
 - Demonstrated reliability of Exactive Plus during CIP2
 - Possibility for developing a wide ranging environmental screening method.
- Purchase
 - QE Focus purchased January 2020
 - Method for the analysis of the LC-MS amenable compounds for the CHEM14 work package of CIP 3 developed in conjunction with Thermo.

Q EXACTIVE FOCUS









- Meet LOD and uncertainty of measurement targets as defined in the CIP technical specification.
- Single method if at all possible
- Robust operation across different matrix types
- On-line sample preparation (EQUAN)

QE Focus Method Summary

- Analytes
 - 17 compounds as listed includes PFAS compounds, pesticides industrial chemicals and pharmaceuticals.
 - Isotope labelled surrogates available for all except 1 (6:2FTAB)
- Matrices
 - Surface waters, sewage effluents and influents
- Sample preparation
 - On-line trace enrichment using EQUAN
 - Centrifugation to remove particulates
 - Sample volume 5ml
 - Loading column Thermo Scientific Hypersil Gold aQ pre-concentration column (20x2.1, 12µm)





QE Focus Method Summary – continued

HPLC

- Column Thermo Scientific Acclaim RSLC PolarAdvantage 2.1x150mm. Particle size 2.2um.
- Mobile Phase Water, 5mM Amm.Form, 1% Formic Acid / Methanol, 5mM Amm.Form, 1% Formic Acid
- Run Gradient, 20 minutes
- MS
 - Full Scan MS
 - Pos/neg switching
 - Scan range 110 600 m/z
 - Resolution 70,000 (QE Focus), 120,000 (Exploris 120)
 - AGC target 1e6
 - t–SIM
 - Isolation window 1.0 m/z
 - Resolution 70,000 (QE Focus), 120,000 (Exploris 120)
 - AGC target 2e5
- Quantification
 - Full MS used for all compounds apart from 6:2FTAB (t-SIM)





Sensitivity and Linearity





Sensitivity and Linearity





Right Solutions • **Right Partner**

19



Method Validation Results

Compound	LOD, ug/L				MRL	Low Std		High Std	
	Treated Sewage	Surface Water	Ground Water	Crude Sewage	ug/L	%RSD	%Bias	%RSD	%Bias
Ranitidine	0.00890	0.00230	0.000804	0.0225	0.1	3.10	-2.89	2.87	-0.63
Benzotriazole	Benzotriazole 0.944 0.0140		0.00399	0.361	0.1	5.43	-1.24	2.26	-0.59
Imidacloprid	0.00205	0.000462	0.000293	0.00644	0.001	4.57	-2.81	1.53	-0.43
Bisphenol S	0.00426	0.00243	0.00284	0.0233	0.01	4.52	-3.87	2.30	-0.57
PFPeA	0.00191	0.000623	0.000332	0.00467	0.001	4.11	-2.35	2.13	-0.22
Climbazole	0.00297	0.00170	0.00149	0.0110	0.1	3.91	-1.39	1.98	-0.29
PFBS	0.000426	0.0000973	0.000179	0.00176	0.001	2.37	-2.52	2.21	-0.48
PFHxA	0.00126	0.000292	0.000192	0.00159	0.001	4.24	-2.99	1.32	-0.70
GenX	0.00129	0.000334	0.000147	0.00816	0.001	4.07	-2.81	1.13	-1.87
PFHpA	0.000591	0.000230	0.000158	0.00137	0.001	2.37	-2.81	2.21	-0.35
PFHxS	0.000424	0.0000788	0.000107	0.00112	0.001	2.37	-3.63	1.94	-0.94
6:2FTAB	0.00535	0.000507	0.000995	0.0221	0.001	8.26	0.76	8.76	-4.30
6:2FTS	0.00249	0.000326	0.000176	0.00176	0.001	2.21	-1.73	2.21	1.39
PFOA	OA 0.00108 0.000232		0.0000668	0.000594	0.00009	2.37	-4.04	2.01	-0.38
PFOS	0.00830	0.000151	0.0000233	0.0200	0.00009	10.01	-3.95	4.45	-0.77
Fipronil	0.000939	0.000216	0.000234	0.00227	0.01	6.47	-3.14	2.03	0.28
Triclocarban	0.00132	0.00110	0.00115	0.00736	0.1	4.83	-1.99	1.63	0.30



Method Validation Results – continued

Compound	Treated Sewage 10%		Treated Sewage 60%		Surface Water 60%		Crude Sewage 60%	
	%RSD	%Rec	%RSD	%Rec	%RSD	%Rec	%RSD	%Rec
Ranitidine	3.99	94.7	4.91	101	3.98	98.5	6.93	99.3
Benzotriazole	10.61*	106.0*	9.61*	95.7 [*]	4.26	103.1	5.49	104.4
Imidacloprid	3.14	96.9	1.61	101	1.64	99.6	6.95	98.8
Bisphenol S	3.09	97.1	3.03	99.8	2.61	99.3	6.18	97.7
PFPeA	6.46	100	3.43	102	2.18	102	6.62	99.7
Climbazole	2.84	96.5	3.73	101	2.35	101	7.37	98.7
PFBS	2.93	96.8	2.24	99.2	1.34	99.1	6.60	98.4
PFHxA	2.33	97.1	3.33	101	2.28	100	6.38	97.7
GenX	4.78	95.9	3.04	99.9	1.58	101	6.92	98.5
PFHpA	2.29	95.7	2.42	99.2	1.49	99.1	7.52	97.8
PFHxS	3.46	94.8	2.68	96.8	1.08	97.3	6.73	96.1
6:2FTAB	18.19	93.9	9.73	98.2	9.70	128	12.50	109
6:2FTS	2.60	95.3	2.32	99.1	1.69	101	6.95	94.9
PFOA	4.54	84.0	5.70	105	2.53	102	2.48	102
PFOS	*	*	5.38	118	2.92	101	11.27	118
Fipronil	2.55	95.5	2.79	100	1.77	101	5.83	99.2
Triclocarban	5.28	96.2	3.93	98.5	2.49	98.9	7.25	98.2



Real samples

6:2 FTAB in sewage effluent - 11ng/L



6:2 FTS in sewage effluent - 10ng/L





Real samples

Fipronil in sewage effluent - 57ng/L C19617379_3_01 (Method Settings) 09/05/20 00:32:52 RT: 12.83 - 20.83 SM: 7G RT: 16.83 AA: 3734433 BP: 436.9284 NL: 4.78E5 Base Peak m/z= 434.922-434.9336+ 436.9259-436.9303 F: FTMS -p ESI Full ms [110.0000-600.0000] MS ICIS 100₃ 90-80 C19617379_3_01 Relative Abundance 70-20-10-19.38 19.76 0.0000 434,9324 17.65 17.99 0.0000 436.9285 0-20 18 19 14 15 16 17 13 Time (min)

Bisphenol S in sewage effluent - 750ng/L



Non Target Screening



Food for thought.....

- Typically an Orbitrap acquisition is Full Scan with a range of 100-1000
- We can routinely acquire to 1 ppm mass accuracy
- If we have a targeted list of 500 compounds obtained from Extracted lon chromatograms – the this leaves 99.95% of data unused.....and this is just at one snapshot of the retention time axis.

So what can we learn from the rest of the data and how do we use it?



Non Target Screening



What else can we get out of our Orbitrap data?

 Identification of Unknowns and emerging compounds?



Food and Drink
 Authenticity

 Extractables and Leachables?





 Identification of Pesticide Metabolites?



• Identify Markers for Environmental states or catchment change?



Non Target Screening – *Peak Picking*

m/z Features to Compounds to Answers

For the questions raised above, the answers lie in our Full Scan MS1 and related MS2 data.

Close examination of the TIC time-axis will show a highly complex picture of 'm/z features' with some <u>exhibiting chromatographic behaviour.</u>

Are m/z features with the same chromatographic profile, at the same retention time, related? For the same Rt – are they the same m/z feature in other related samples?



XIC overlay of tryptophan from Compound Discoverer software showing associated ion species. Several ions detected in positive polarity at the same retention time are assembled to reflect one compound within one sample injection of rat plasma extract. Associated ions include [M+H]r, [M+K]r, [M+H-NH];r. [2M+H]r, [2M+N]r. *Peak Picking*: If we can 'assemble' these features into compounds both within the sample and those common to related samples;

- and if we can assign elemental composition, we can begin to move forward -



Non Target Screening – *Peak Annotation*





Thermo Scientific[™] Compound Discoverer[™] software solves this challenge.

Using:

- Sophisticated feature assembly algorithm.
- Comprehensive integrated set of libraries, databases: mzCloud, ChemSpider Custom mass lists Compound Class Scoring
- Tools to assist Structural Identification: FISh, mzCloud Slimilarity, mzLogic Fine Isotope Trace, Mass Defect.
- Integrated Statistical analysis tools: Differential analysis, PCA, PLS-DA, ANOVA, Volcano, box-whisker, and trend-line plots
- · Template-based and Customisable node-based workflows

.....The software streamlines the identification of unknowns; the determination of real differences between samples; and the elucidation of degradation and biological pathways.

Non Target Screening – *Peak Annotation*







Non Target Screening

mzCloud ID Search Results in Compound Discoverer



Non Target Screening





• Results are always directly linked to raw data.

Non Target Screening - Illicit Drugs



- Compound Discoverer used to detect parent compounds as well as related metabolites, to investigate different metabolism pathways that shed light on the original source of illicit drugs detected in wastewater effluent.
- Here we look at effluent data collected from six UK sewage plants, where cocaine and related metabolites were detected. We will focus on Benzoylecgonine, a dominant metabolite in urinary excretion.





Non Target Screening



Example of Cocaine detection showing mzCloud mirror plot



Non Target Screening – Illicit Drugs



Example of Benzoylecgonine detection showing mzCloud mirror plot



Non Target Screening – Illicit Drugs



- For the six effluent samples, the ratio of Cocaine with respect to Benzolylecgonine was plotted below.
- The ratio in WTPPX clearly shows cocaine dominating, which suggests that the drugs were not from an excretionary route but rather likely from direct disposal.



Conclusions / Further Work



- A method has been developed using HRMS orbitrap with on-line trace enrichment for the low level analysis of a suite of pollutants spanning multiple compound classes in environmental waters using a sample volume of just 5mL.
- Low limits of detection in many cases <1ng/L were obtained together with excellent precision and bias statistics.
- Software tools and libraries are critical aspects in looking at non-targeted workflows as shown in the previous examples
- Next steps will focus on the development of a wide ranging, low level screening method covering a large number of pesticides, pharmaceuticals, PFAS compounds and other chemicals of environmental concern.
- Some new features of the Equan 850 will help simplify and automate the method even further

Future work- Equan850 Exploris120







Equan 850 Direct and On-line SPE and Automatic Calibration Preparation



Right Solutions • **Right Partner**

37

Acknowledgements



- Sam Towers ALS Environmental
- David Evans ALS Environmental
- Neville Llewellyn Thermo Fisher

Any Questions?