

# Targeted quantitation workflow with Orbitrap GC-MS: PAHs and PCBs in soil extracts

**Suresh Seethapathy, Ph.D.**

Senior Applications Scientist  
Thermo Fisher Scientific  
Centers of Excellence – GC/GC-MS

11-Aug 2021

 The world leader in serving science



# Agenda

1 PAHs and PCBs: Importance, challenges

---

2 Sample preparation

---

3 Orbitrap Exploris GC: Schematic and hardware

---

4 Instrument method, tuning, and calibration

---

5 Results and discussion


---

6 Summary

---

# PCBs and PAHs: Importance

- Polychlorinated biphenyls (PCBs)
  - Coolants and lubricants (i.e. transformer oils, capacitors etc.,)
- Polyaromatic hydrocarbons (PAHs)
  - Combustion products (i.e. diesel soot, stack emissions, bushfires)
- Substituted PAHs (oxyPAHs, methylPAHs & polyaromatic (N,S,O) heterocycles (PASHs, PAOHs, PANHs)
  - Combustion products and indirectly from PAHs




1= PCB (general structure), 2=Benzo(a)pyrene (PAH),  
 3=Benzantrone (oxyPAH), 4= 2-methylnaphthalene (methylPAH), 5=Carbazole (PANH), 6=Dibenzofuran (PAOH), 7=Dibenzothiophene (PASH)

# PCBs and PAHs: Challenges and method considerations

- Gas Chromatography (GC)
  - Mass Spectrometer (MS) - EI/CI
- Multiple chromatographic methods
- Time-consuming sample preparation
  - Soxhlet extraction (slow [~48 hours] consuming large amounts of solvents)
- Chromatographic separation of critical pairs
  - i.e. Benzo(b)fluoranthene/ benzo(k)fluoranthene
- Long run times >40 mins

# Sample prep (QuEChERS extraction and dSPE clean-up)



Two hours sample preparation time up to 20X faster vs. Soxhlet extraction

# Thermo Scientific Orbitrap Exploris GC series schematic

ExtractaBrite  
EI/CI ION  
SOURCE

Robust design, vent free source  
and column change  
EI/veV/CI

Advanced Quadrupole  
Technology (AQT)

Isolation widths down to  
0.4 Da improves  
sensitivity and selectivity

Ion Routing  
Multipole

Facilitates higher  
energy collisional  
dissociation (HCD)  
fragmentation


C-Trap

AdvancedActive Beam  
Guide (AABG)  
with axial gradient

Prevents neutrals from entering  
the quadrupole and improves  
robustness

High-Field Orbitrap  
Mass Analyzer

Resolution up to 240,000 at  
 $m/z$  200 and acquisition rates up to 40 Hz



# Thermo Scientific Orbitrap Exploris GC Mass Spectrometer: Schematic/hardware

ExtractaBrite  
EI/CI ION  
SOURCE

Advanced Active Beam  
Guide (AABG)  
with axial gradient


Advanced Quadrupole  
Technology (AQT)

Isolation widths down to  
0.4 Da improves  
sensitivity and selectivity

C-Trap

Ion Routing  
Multipole

High-Field Orbitrap  
Mass Analyzer



# Thermo Scientific Orbitrap Exploris GC Mass Spectrometer: Schematic/hardware



ExtractaBrite  
EI/CI ION  
SOURCE

Advanced Active Beam  
Guide (AABG)  
with axial gradient

Advanced Quadrupole  
Technology (AQT)

C-Trap

Ion Routing  
Multipole



High-Field Orbitrap  
Mass Analyzer




# Instrument method

- Column: **Thermo Scientific™ TraceGOLD™ TG-5 Sil MS (30m × 0.25mm, 0.25µm)**
- Inlet and mode: SSL, Splitless
- Flowrate: 1.2 mL/min
- Inject volume: **1 µL**
- Thermo Scientific™ Orbitrap Exploris™ GC MS
- Transfer line temperature: 320 °C
- Source temperature: 350 °C
- Electron energy (eV): 70
- Acquisition mode: Full scan (FS)
- Mass range (*m/z*): 50-550
- Mass resolution : FS 60,000(FWHM @ *m/z* 200, scan speed 7.4 Hz)
- Lock mass (*m/z*: 207.03235)

| No | Retention time [min] | Rate [°C/min] | Target value [°C] | Hold time [min] |
|----|----------------------|---------------|-------------------|-----------------|
| 1  | 0.000                | Run           |                   |                 |
| 2  | 1.000                | 0.00          | 40.0              | 1.00            |
| 3  | 9.750                | 28.00         | 285.0             | 0.00            |
| 4  | 16.417               | 3.00          | 305.0             | 0.00            |
| 5  | 22.917               | 30.00         | 350.0             | 5.00            |
| 6  |                      | New Row       |                   |                 |
| 7  | 22.917               | StopRun       |                   |                 |

# Tuning and calibration

Tune and calibrate in under 5 minutes and be confident system is at peak performance for all users



# Thermo Scientific™ Chromeleon™ CDS Software

**Multivendor**  
Instrument control

**Scalability**  
Workstation to enterprise & cloud

**Compliance & Data Integrity**

**Integration**

**MS**

**Productivity**  
Up to 1/3 faster\*





**Ease of Use**  
Workflows

**MS**

Built for **the Lab.**  
Built for **IT.**  
Why compromise?  
**Chromeleon 7.3 CDS**


# Selectivity

50 pg/uL calibration standard in n-hexane




Retention time (minutes)

# Sensitivity MDLs




# Linearity

| Compound                    | Compound type | R <sup>2</sup> | AVCF % RSD | Compound                | Compound type | R <sup>2</sup> | AVCF % RSD  |
|-----------------------------|---------------|----------------|------------|-------------------------|---------------|----------------|-------------|
| Naphthalene                 | PAH           | 0.9999         | 1.6        | 1,2-Dimethylnapthalene  | methylPAH     | 0.9993         | 4.5         |
| Acenaphthylene              | PAH           | 0.9987         | 5.4        | 1,8-Dimethylnapthalene  | methylPAH     | 0.9998         | 2.6         |
| Acenaphthene                | PAH           | 0.9995         | 4.0        | PCB-28                  | PCB           | 0.9997         | 2.5         |
| Biphenyl                    | PAH           | 0.9998         | 2.6        | PCB-52                  | PCB           | 0.9991         | 2.8         |
| Fluorene                    | PAH           | 0.9981         | 9.0        | PCB-101                 | PCB           | 0.9998         | 3.2         |
| Phenanthrene                | PAH           | 0.9995         | 3.8        | PCB-118                 | PCB           | 0.9998         | 3.7         |
| Anthracene                  | PAH           | 0.9981         | 4.3        | PCB-153                 | PCB           | 0.9998         | 1.6         |
| Fluoranthene                | PAH           | 0.9998         | 3.0        | PCB-138                 | PCB           | 0.9991         | 2.8         |
| Pyrene                      | PAH           | 0.9997         | 3.2        | PCB-180                 | PCB           | 0.9997         | 4.3         |
| Benzo[a]anthracene          | PAH           | 0.9999         | 1.7        | Benzo(b)thiophene       | PASH          | 0.9998         | 3.2         |
| Chrysene                    | PAH           | 0.9997         | 3.1        | Dibenzothiophene        | PASH          | 0.9988         | 3.7         |
| Benzo[b]fluoranthene        | PAH           | 0.9998         | 2.6        | 1&2-Ethylnapthalene     | ethylPAH      | 0.9996         | 3.7         |
| Benzo[k]fluoranthene        | PAH           | 0.9994         | 4.5        | Quinoline               | PANH          | 0.9988         | 4.0         |
| Benzo[a]pyrene              | PAH           | 0.9987         | 5.4        | 1-Indanone              | PAOH          | 0.9993         | 4.7         |
| Indeno[1,2,3-cd]pyrene      | PAH           | 0.9964         | 9.3        | Dibenzofuran            | PAOH          | 0.9993         | 5.3         |
| Dibenzo[a,h]anthracene      | PAH           | 0.9978         | 7.3        | Carbazole               | PAOH          | 0.9980         | 4.7         |
| Benzo[ghi]perylene          | PAH           | 0.9989         | 5.1        | 9, 10-Anthraquinone     | PAOH          | 0.9951         | 12.9        |
| 1-Methylnapthalene          | methylPAH     | 1.0000         | 1.1        | 2-Methylanthraquinone   | PAOH          | 0.9981         | 6.5         |
| 2-Methylnapthalene          | methylPAH     | 0.9999         | 1.8        | 9-Fluorenone            | oxyPAH        | 0.9997         | 5.0         |
| 2, 7-Dimethylnapthalene     | methylPAH     | 0.9999         | 1.5        | Benzanthrone            | oxyPAH        | 0.9985         | 6.0         |
| 1,3 +1,6-Dimethylnapthalene | methylPAH     | 0.9999         | 2.0        | 5, 12-Napthacenequinone | oxyPAH        | 0.9963         | 9.6         |
| 2,3-Dimethylnapthalene      | methylPAH     | 0.9999         | 1.8        |                         | <b>Min</b>    | <b>0.9951</b>  | <b>1.1</b>  |
|                             |               |                |            |                         | <b>Max</b>    | <b>1.0000</b>  | <b>12.9</b> |
|                             |               |                |            |                         | <b>Mean</b>   | <b>0.9991</b>  | <b>4.2</b>  |



oxyPAH


R<sup>2</sup>=0.9997  
% RSD=5.0




PAH

R<sup>2</sup>=0.9989  
% RSD=4.1

# Quantification of PAHs and PCBs in soil




|                            |                                      |
|----------------------------|--------------------------------------|
| Amount found in sample     | 0.6 $\mu\text{g}/\text{Kg}$          |
| Ion ratio % deviation      | 0.7 %                                |
| Measured mass ( $m/z$ )    | 255.96074                            |
| Theoretical mass ( $m/z$ ) | 255.96078                            |
| Chemical formula           | $\text{C}_{12}\text{H}_7\text{Cl}_3$ |
| Mass error (ppm)           | 0.2                                  |



|                            |                                   |
|----------------------------|-----------------------------------|
| Amount found in sample     | 0.8 $\mu\text{g}/\text{Kg}$       |
| Ion ratio % deviation      | 4.8%                              |
| Measured mass ( $m/z$ )    | 180.05698                         |
| Theoretical mass ( $m/z$ ) | 180.05697                         |
| Chemical formula           | $\text{C}_{13}\text{H}_8\text{O}$ |
| Mass error (ppm)           | 0.1                               |


# Robustness: Peak area repeatability




QuEChERS Soil Extract spiked at 10 ng/Kg




# Robustness: Peak area repeatability




# Robustness: QC standard ion ratio stability




# Robustness: QC standard mass accuracy stability




# Thermo Scientific Compound Discoverer 3.2 Software deconvoluted soil EI spectrum



# NIST search results



# Screening for additional contaminants (EI deconvolution and NIST matches)



- **Orbitrap Exploris GC mass spectrometer:** Best mass accuracy and high resolution leads to easy, reliable, and compliant data generation.
- **Easy full-scan data acquisition:** Retrospective analysis and comprehensive targeted and non-targeted workflow of environmental sample analysis with Thermo Scientific™ Compound Discoverer™ software.
- **Comprehensive method consolidation:** PAHs and PCBs analysis in soil in < 20 min. Up to 20x, realized using modified QuEChERS method compared to tradition Soxhlet extraction methods, saving cost and time.
- **Sensitivity:** MDL values calculated for 45 native compounds ranging from 118 to 475 fg OC (corresponding to 0.1–0.5 µg/kg in sample). LOQs ranged from 0.5 to 5.0 µg/kg in soil as determined.
- **Linear dynamic range:** Better than 6 orders of magnitude.
- **Robustness:** Stable peak area, ion ratios, and mass accuracies were achieved for QuEChERS soil extracts and QC standards, studied over a three week period.
- Overall, an excellent solution for both the chemists running the system and lab managers keeping track of the return on investment.

# Thank you

# Questions?

