

# Streamlining the Sample Preparation workflow for POP Analysis

Angelika Köpf, LCTech GmbH




LCTech GmbH

Experts in Sample Preparation since 1998

SOLUTIONS BY **LC**Tech




From Sample to Chromatography



## Streamlining Sample prep Workflow in POP Analysis

### Typical Sample Preparation Workflow

SOLUTIONS BY **LC***Tech*



## Streamlining Sample prep Workflow in POP Analysis

### Common Challenges in sample prep of POP Analysis



- Reproducible and efficient Extraction, to avoid downtime
- Some analytes should not be evaporated to dryness as they get lost
- Tedious cleanup based on SPE, which is error prone from column preparation to proper fractionation
- Every transfer step can cause errors
- High cost on consumables and solvents
- Avoid blind values and Cross Contamination



# Streamlining Sample prep Workflow in POP Analysis


## Solution by LCTech – this years Example Dioxin/PCB



## Streamlining Sample prep Workflow in POP Analysis

### X-TRACTION for Extraction

SOLUTIONS BY **LC** *Tech*



Tray for 5 solvents

Touchscreen

Core of the system,  
heating bar  
compartment with cell

Benchtop unit – no hood  
Pre-programmed – easy to use  
Operation by touchscreen (reports)

Easy to install

1 unit – 1 sample – up to 6 units

## Streamlining Sample prep Workflow in POP Analysis

### Expandable Design

SOLUTIONS BY **LC** Tech



**Throughput** with runtime 15 min + 3 min change = 18 min = 3,3 per unit / hour  
30 per unit / 9 hour workday  
180 per 6 units / 9 hour workday

## Streamlining Sample prep Workflow in POP Analysis

### Unique LCTech Solution – What makes it that Special

---



- **Easiest, quick handling** of cells, no tools needed, intuitive and safe
- **Low cost** investment due to entry level single station (main system) and expandable design (up to five add-on systems)
- **Fast and efficient protocols**
- **Green chemistry!** Low solvent consumption in process and rinsing
- **Low cost** per sample due to **less consumables** (only one filter) and **low wear and tear**

The word "UNIQUE" in a white, sans-serif font, centered within a dark blue, rounded rectangular banner that has a slight 3D effect with a white shadow on the right side.

UNIQUE





## Streamlining Sample prep Workflow in POP Analysis

### Quick and save handling

SOLUTIONS BY **LC** *Tech*



- ✓ Time saving with unique easy closing by magnet
  - ✓ No tools needed
  - ✓ No risk to damage screw
- ✓ Besides the glass fibre filter all parts are reusable for many times

## Streamlining Sample prep Workflow in POP Analysis

### Efficient Extraction using special concept


---

SOLUTIONS BY **LC** *Tech*

- Fill in sample



- No need to overfill with expensive Hydromatrix
- No need to fill with solvent
- Universal cell size



## Streamlining Sample prep Workflow in POP Analysis

### Working Principle of X-TRACTION

---

SOLUTIONS BY **LC** Tech

- The extraction system works with an air gap in the cell. This allows for an extraction, without generating the high pressure, that other units pose at the sample with the same efficiency  
→ **LPFE** - **L**ow **P**ressure **F**luid **E**xtraction.
- Working in **low pressure ranges**, leads to decreased wear and tear of instrument parts, higher longevity and an easy and safe handling.



## Streamlining Sample prep Workflow in POP Analysis

### Efficient process by Design supporting principle

---

SOLUTIONS BY **LC** *Tech*

- Solvent loading



- Solvent heating



# Streamlining Sample prep Workflow in POP Analysis

## Solution by LCTech – this years Example Dioxin/PCB



## Streamlining Sample prep Workflow in POP Analysis

### Controlled and unattended Evaporation

---

SOLUTIONS BY **LC***Tech*

- Vacuum Centrifuge ideally for parallel processing, space saving
- Ideal Flexibility to insert containers of Extract or Fraction
- D-EVA (Dioxin-EVAporation) consists of:
  - Centrifuge esp. for Dioxins & PCB analysis with special software and Default settings
  - Cold trap
  - Rotor
  - Special LCTech sensor



## Streamlining Sample prep Workflow in POP Analysis

### Controlled and unattended Evaporation

---

- Challenge:
  - To stop evaporation at defined low end volume needed for dioxin analysis
  - To avoid the impact of memory effect and “burning” of the material.
- ➔ **Special LCTech sensor and software** developed for automatic STOP
- ➔ Energy supply via light that can be stopped



## Streamlining Sample prep Workflow in POP Analysis

### Default methods

---

SOLUTIONS BY **LC** *Tech*

- Parallel processing of up to 26 toluene samples (dioxins)
- Process time: **37 min**
- Final volumes of about **30 – 100  $\mu$ L or nearly dryness**





# Streamlining Sample prep Workflow in POP Analysis


## Solution by LCTech – this years Example Dioxin/PCB



## Streamlining Sample prep Workflow in POP Analysis

### Cleanup unit DEXTech


SOLUTIONS BY **LC***Tech*



- Bench top unit, all easy accessible
- Easy setup, no fume hood
- Ready to use programm
  
- High quality columns
- No Cross Contamination
- Highest degree of automation
- Attractive cost for ready to use columns

# Streamlining Sample prep Workflow in POP Analysis


## Alox Plus Method – conventional or DCM Free



	Universal	SMART
Process Time (includes loading approx. 3 min and conditioning approx. 15 min)	50 min	34 min
Total Solvent Consumption (incl. conditioning and fractions)	230 mL	118 mL
	Fraction 1	Fraction 2
Fraction Volume	24 mL	10 mL

# Streamlining Sample prep Workflow in POP Analysis

## Results





## Streamlining Sample prep Workflow in POP Analysis

### BCR – 536 – PCB in Harbour Sediment

Results after entire sample preparation and analysis workflow

native [µg/kg]	mean n=4 [µg/kg]	cert. Value [µg/kg]	recovery [%]
PCB-#28	47.3	44.0	108
PCB-#52	42.9	38.0	113
PCB-#101	42.9	44.0	98
PCB-#118	23.0	27.5	84
PCB-#105	3.0	3.5	87
PCB-#153	53.3	50.0	107
PCB-#138	25.7	27.0	95
PCB-#156	2.7	3.0	90
PCB-#180	21.4	22.4	96
PCB-#128	5.3	5.4	98
PCB-#170	18.1	17.2	105
PCB-#163	13.0	13.4	97



LC Tech Laboratory, Obertaufkirchen

## Streamlining Sample prep Workflow in POP Analysis

### BCR – 677 - PCDD/F in sewage sludge



Results after entire sample preparation and analysis workflow

native [pg/g]	mean n=5 [ng/kg]	cert. Value [ng/kg]	recovery [%]
2,3,7,8-TCDF	41	45	91
1,2,3,7,8-PeCDF	23	24,8	93
2,3,4,7,8-PeCDF	16	16,9	95
1,2,3,4,7,8-HxCDF	13	14,5	90
1,2,3,6,7,8-HxCDF	5,7	6,1	93
2,3,4,6,7,8-HxCDF	6,4	5,6	114
1,2,3,7,8,9-HxCDF	1,1	0,8	131
1,2,3,4,6,7,8-HpCDF	59,1	61,6	96
1,2,3,4,7,8,9-HpCDF	4,7	6,3	74
1,2,3,4,6,7,8,9-OCDF	158	177	89
2,3,7,8-TCDD	1,6	1,5	103
1,2,3,7,8-PeCDD	4	4,1	97
1,2,3,4,7,8-HxCDD	nd	nd	n.d
1,2,3,6,7,8-HxCDD	239	235	102
1,2,3,7,8,9-HxCDD	73,6	79	93
1,2,3,4,6,7,8-HpCDD	3221	3500	92
1,2,3,4,6,7,8,9-OCDD	12921	12700	102



LC Tech Laboratory, Obertaufkirchen

## Streamlining Sample prep Workflow in POP Analysis

### Excellent Precision for PCDD/F



	Extraction 1	Extraction 2	Extraction 3	Extraction 4	Mean	Deviation	RSD %
	Fresh weight ng/kg	Fresh weight ng/kg	Fresh weight ng/kg	Fresh weight ng/kg	Fresh weight ng/kg		
<b>Polychlorinated Dibenzofurane</b>							
2,3,7,8-TeCDF	0.4292	0.4263	0.3877	0.3862	<b>0.4074</b>	<b>0.0204</b>	<b>5.0</b>
1,2,3,7,8-PeCDF	0.2477	0.2428	0.2462	0.2615	<b>0.2496</b>	<b>0.0071</b>	<b>2.9</b>
2,3,4,7,8-PeCDF	0.4973	0.5043	0.4909	0.4892	<b>0.4954</b>	<b>0.0059</b>	<b>1.2</b>
1,2,3,4,7,8-HxCDF	0.3864	0.3791	0.3631	0.3696	<b>0.3746</b>	<b>0.0089</b>	<b>2.4</b>
1,2,3,6,7,8-HxCDF	0.2518	0.256	0.2447	0.247	<b>0.2499</b>	<b>0.0044</b>	<b>1.7</b>
2,3,4,6,7,8-HxCDF	0.2585	0.2687	0.2752	0.2619	<b>0.2661</b>	<b>0.0064</b>	<b>2.4</b>
1,2,3,7,8,9-HxCDF	0.1844	0.1819	0.1881	0.1832	<b>0.1844</b>	<b>0.0023</b>	<b>1.3</b>
1,2,3,4,6,7,8-HpCDF	0.3229	0.3448	0.3284	0.3236	<b>0.3299</b>	<b>0.0088</b>	<b>2.7</b>
1,2,3,4,7,8,9-HpCDF	0.2241	0.2323	0.2236	0.2165	<b>0.2241</b>	<b>0.0056</b>	<b>2.5</b>
OCDF Octachlordibenzofuran	0.9901	1.043	0.9885	1.001	<b>1.0057</b>	<b>0.0221</b>	<b>2.2</b>
<b>Polychlorinated Dibenzodioxin</b>							
2,3,7,8-TeCDD	0.2163	0.1917	0.1949	0.1634	<b>0.1916</b>	<b>0.0188</b>	<b>9.8</b>
1,2,3,7,8-PeCDD	0.247	0.2377	0.2255	0.2275	<b>0.2344</b>	<b>0.0086</b>	<b>3.7</b>
1,2,3,4,7,8-HxCDD	0.1831	0.1759	0.19	0.2084	<b>0.1894</b>	<b>0.0121</b>	<b>6.4</b>
1,2,3,6,7,8-HxCDD	0.2067	0.2032	0.2019	0.1904	<b>0.2006</b>	<b>0.0061</b>	<b>3.1</b>
1,2,3,7,8,9-HxCDD	0.2422	0.2373	0.2526	0.2351	<b>0.2418</b>	<b>0.0067</b>	<b>2.8</b>
1,2,3,4,6,7,8-HpCDD	0.2567	0.2704	0.2598	0.2401	<b>0.2568</b>	<b>0.0109</b>	<b>4.2</b>
OCDD Octachlordibenzodioxin	1.472	1.409	1.389	1.353	<b>1.4058</b>	<b>0.0432</b>	<b>3.1</b>

Quality Control Feed Material FMP

Data provided by Dr. Bernsmann, Chemical and Veterinary Analytical Institute, Münsterland-Emscher-Lippe (CVUA-MEL), Muenster, Germany



# Streamlining Sample prep Workflow in POP Analysis

## Excellent Precision for PCB


Data provided by Dr. Bernsmann, Chemical and Veterinary Analytical Institute,  
Münsterland-Emscher-Lippe (CVUA-MEL), Muenster, Germany

Quality Control Feed Material FMP SOLUTIONS BY 

	Extraction 1	Extraction 2	Extraction 3	Extraction 4	Mean	STABWN	RSD %
	Fresh weight   ng/kg	Fresh weight   ng/kg	Fresh weight   ng/kg	Fresh weight   ng/kg	Fresh weight   ng/kg		
<b>non-ortho PCB</b>							
PCB 77	74.58	77.17	77.39	73.8	75.74	1.57	2.1
PCB 81	6.957	7.275	6.708	6.934	6.969	0.2019	2.9
PCB 126	4.561	4.628	4.575	4.43	4.549	0.0728	1.6
PCB 169	2.542	2.546	2.516	2.405	2.5023	0.0573	2.3
<b>mono-ortho PCB</b>							
PCB 105	370.4	378.8	370.9	368.1	372.1	4.0376	1.1
PCB 114	33.57	33.36	33.13	31.39	32.9	0.8643	2.6
PCB 118	886.1	934.4	905.6	899.5	906.4	17.64	1.9
PCB 123	23.05	23.69	23.32	23.00	23.3	0.2739	1.2
PCB 156	493.1	503.8	496.7	498.8	498.1	3.8710	0.8
PCB 157	53.14	53.23	52.49	51.88	52.7	0.5455	1.0
PCB 167	229.7	224.5	224.4	223.6	225.6	2.4213	1.1
PCB 189	117.2	122.3	118.3	116.9	118.7	2.1568	1.8
<b>ndl-PCB</b>	µg/kg	µg/kg	µg/kg	µg/kg	µg/kg	µg/kg	
PCB 28	0.7508	0.7508	0.7764	0.7697	0.7619	0.0114	1.5
PCB 52	1.705	1.705	1.719	1.713	1.7105	0.0059	0.3
PCB 101	2.498	2.498	2.53	2.477	2.5008	0.0189	0.8
PCB 138	3.771	3.771	3.742	3.674	3.7395	0.0396	1.1
PCB 153	4.603	4.603	4.653	4.611	4.6175	0.0208	0.4
PCB 180	4.437	4.437	4.483	4.361	4.4295	0.0438	1.0

# Streamlining Sample prep Workflow in POP Analysis

## Exceptional Precision PCDD/F




Quality Control Feed Material FMP

Data provided by Dr. Bernsmann, Chemical and Veterinary Analytical Institute, Münsterland-Emscher-Lippe (CVUA-MEL), Muenster, Germany

# Streamlining Sample prep Workflow in POP Analysis

## Exceptional Precision for PCB



Quality Control Feed Material FMP

Data provided by Dr. Bernsmann, Chemical and Veterinary Analytical Institute, Münsterland-Emscher-Lippe (CVUA-MEL), Muenster, Germany

## Streamlining Sample prep Workflow in POP Analysis

### Customer Carry Over Test



PCDD/Fs	Carry Over (%)
2378-TetraCDF	0.000999476
2378-TetraCDD	0.000933642
12378-PentaCDF	0.000894312
23478-PentaCDF	0.001000214
12378-PentaCDD	0.000927239
123478-HexaCDF	0.000897438
123678-HexaCDF	0.000799541
234678-HexaCDF	0.001123207
123478-HexaCDD	0.000950550
123678-HexaCDD	0.001042591
123789-HexaCDD	0.001110944
123789-HexaCDF	0.001051582
1234678-HeptaCDF	0.001760539
1234678-HeptaCDD	0.003740431
1234789-HeptaCDF	0.001014591
OctaCDD	0.008554827
OctaCDF	0.002093968

PCB	Carry Over (%)
PCB-28 **	0.037920007
PCB-52 **	0.057359327
PCB-101 **	0.042526266
PCB-81	0.001353596
PCB-77	0.005647445
PCB-123	0.000957830
PCB-118	0.032228256
PCB-114	0.000778819
PCB-153 **	0.051951924
PCB-105	0.016001291
PCB-138 **	0.023677053
PCB-126	0
PCB-167	0.003419919
PCB-156	0.004653067
PCB-157	0.003789838
PCB-180 **	0.017443246
PCB-169	0.004602942
PCB-189	0

For the test two certified reference materials (soil) have been used

(BCR - 529 Industrial soil and SQCo68 - 50 g)

Each time the high contaminated soil sample was performed including with the DEXTech Plus system followed with a blank sample.

## Streamlining Sample prep Workflow in POP Analysis

### Separate Cross contamination test in D-EVA

---

SOLUTIONS BY **LC** *Tech*

- Centrifuge glasses in 90° and 180° orientation
- Standard 4 x highest calibration level
- No detectable signal /  
➔ cross-contamination by GC-HRMS (Thermo DFS)



## Streamlining Sample prep Workflow in POP Analysis

### Advantages in this automated processing

---




- Complete workflow from sample to injection, whilst still having control over each chapter
- Unmatched easy and safe handling of X-TRACTION, quick, low cost and solvent saving
- Cleanup on DEXTech brings great performance with safe operation, Time saving, solvent saving and cost saving flexibility of protocols and high quality Columns
- D-EVA saves the sensitive fractions with unique safety stop
- Design and material chosen ensure no blind values and no Cross Contamination
- All in all Robust LCTech engineering for reliable operation and low maintenance



## Streamlining Sample prep Workflow in POP Analysis

Any Questions?

SOLUTIONS BY **LC***Tech*



[koepf@LCTech.de](mailto:koepf@LCTech.de)



Thank you for your attention!

---

SOLUTIONS BY **LC** *Tech*



LCTech GmbH  
Daimlerstraße 4  
84419 Obertaufkirchen  
Deutschland

info@LCTech.de  
www.LCTech.de  
www.LCTech-online.com





## X-TRACTION

# EPA Compliance

---



*Text of EPA 3534A (SW 846):*

*2.2 The extraction cell containing the sample is heated to the extraction temperature (see Sec. 11.8), **pressurized with the appropriate solvent system**, and extracted for 5-10 minutes (or as recommended by the instrument manufacturer). Multiple extractions are 3545A - 3 Revision 1 February 2007 recommended for some groups of analytes. The solvent systems used for this procedure vary with the analytes of interest and are described in Sec. 7.7.*

*6.1.2 Other system designs may be employed, provided that adequate performance can be demonstrated for the analytes and matrices of interest*

*11.9.4 Optimize the conditions, as needed, according to the manufacturer's instructions. **In general, the pressure is not a critical parameter, as the purpose of pressurizing the extraction cell is to prevent the solvent from boiling at the extraction temperature and to ensure that the solvent remains in intimate contact with the sample.***

Proof on Performance data , that X-TRACTION is as efficient

