



# Validation of Controls for Wastewater Surveillance for SARS CoV-2

Brian Swalla Ph.D., IDEXX Laboratories, Water Division


**IDEXX**

# IDEXX has invested heavily in developing materials, protocols, and data to support wastewater surveillance

## IDEXX Reagents




## IDEXX Protocols



## IDEXX Validation Data

- 100+ samples
- 15+ different wastewater treatment facilities
  - Geographically diverse
  - Different processes
- Testing over several months
- Multiple analysts
- Different material & reagent lots
- Varied lab equipment


# IDEXX has experience implementing wastewater surveillance around the world with various types of laboratories



## Customer Types

- Public Health Labs
- National Public Health Organizations
- Independent Laboratories
- Major contract laboratories
- Academic Research Labs
- Universities (on campus testing)
- Wastewater Utilities


# Validation data allows labs with robust quality procedures to produce results that can be used for public health decisions



**Because of the urgency of the pandemic, validation has often been limited in scope, which can lead to unintended consequences**

# IDEXX's experience with wastewater surveillance has demonstrated the importance of validation

- Matrix Recovery
- Fecal Normalization
- Internal Control



## Matrix Recovery Control

- Importance of material choice

## Human Fecal Control

- Importance of probe design

## Internal Control

- Ability to assess inhibition



# A matrix recovery control is critical for assessing yield loss during concentration

Matrix Recovery


Fecal Normalization

Internal Control

## Sources of variation in concentration

- Relatedness of spike to SARS-CoV-2
- Inherent variability of wastewater matrix
- Variability of wastewater inputs (non-Human sources)
- Run-to-run variation

Sample preparation      **Sample concentration**      RNA Extraction      RNA Measurement



Spike known amount of exogenous material

Measure remaining material to calculate yield

# IDEXX testing has demonstrated that certain materials perform differently in wastewater

Matrix Recovery

Fecal Normalization

Internal Control

Collection Date	Site	Sample	BRSV Recovery %		
			Replicate		Average
			A	B	
Day 1	A	260	69.3%	49.1%	59.2%
Day 3		259	34.5%	38.7%	36.6%
Day 1	B	275	65.7%	83.5%	74.6%
Day 2	C	266	15.9%	15.1%	15.5%
Day 4		262	22.7%	20.2%	21.5%
Day 6		280	15.4%	14.6%	15.0%
Day 2	D	274	33.6%	36.6%	35.1%
Day 3		247	12.6%	13.4%	13.0%
Day 5		283	30.9%	28.7%	29.8%
Day 4	E	264	32.2%	31.4%	31.8%
Day 5	F	268	17.2%	15.1%	16.2%
Day 6	G	281	25.5%	26.5%	26.0%
Average Recovery:					31.2%



Higher variation between sites

- IDEXX data show BRSV can be used reliably as a matrix recovery control
- IDEXX also tested Bovine Coronavirus (BCoV)
- BCoV showed lower stability in wastewater and more variation in results
- However, BCoV is still widely used

Recovery of Bovine Respiratory Syncytial Virus (BRSV) from wastewater using PEG-based concentration


# A human fecal control can account for differences in relative human waste input over time

Matrix Recovery

**Fecal Normalization**


Internal Control

## Without human fecal control



**Poor public health decisions**

## With human fecal control



**Appropriate public health decisions**




# IDEXX testing demonstrated the importance of the original probe design

Matrix Recovery

Fecal Normalization


Internal Control

sample 1



MGB = blue

sample 2



no MGB = red

- The original PMMoV probe design included a “minor groove binder,” or MGB
- Some publications reported results for wastewater not using MGB
- IDEXX testing demonstrated that MGB is important for probe performance
- Methods need to be examined critically to evaluate risks


# Internal controls must be validated to show no impact on quantification

Matrix Recovery

Fecal Normalization

Internal Control

- An internal control (IC) can indicate validity for each sample result, based on:
  - Inhibition potential
  - Successful purification and amplification
- The IC reaction **must** be multiplexed with the SARS-CoV-2 reaction
- There is a potential the IC reaction could affect quantification
- Validation is the **only** way to assess this risk



No impact of Multiplex Internal Control Reaction on SARS-CoV-2 Detection or Quantification

# IDEXX Laboratories, Inc.



- Headquarters in Westbrook, Maine
- 700,000+ sq. ft. manufacturing and offices
- USDA licensed; ISO 9001, 14001, 17025 certified
- Offices in 17 countries, serving customers in 120 countries
- A S&P 500 and Nasdaq 100 company

- Worldwide market leader
- Specialist in water microbiology
- Accepted or approved in 50+ countries; over 100 different approvals/acceptances by more than 20 entities globally
- **Helping protect water quality for an estimated 2.5 billion people every day**



# IDEXX's experience in validated wastewater testing enabled rapid validation of wastewater surveillance methods

## >20 years of wastewater expertise

- More than 1,500 testing systems used in wastewater facilities around the globe
- Methods approved by or included in ~20 regulatory bodies and standards



## >10 years of real-time PCR expertise

- Creator of the RealPCR tests, with ~20 kits commercialized globally
- The OPTI SARS-CoV-2 RT-PCR Test Kit has received FDA emergency use authorization
- >1 million PCR tests performed each year





# Questions

[Brian-Swalla@idexx.com](mailto:Brian-Swalla@idexx.com)