

Understanding the Need for Mass Balance Measurements in Microplastics and the Three GC/MS Techniques Used to Get Them

- A Short Overview -

Kurt Thaxton
www.gerstel.com

Today's Goal: Clarification!

- ▷ What is Mass Balance, and why care?
- ▷ How is it measured? (short answer: by some type of GC-MS)
- ▷ What are the three ways of measuring Mass Balance, and how are they different?

For what I want to measure related to Microplastics, what way is the best for me?

No “Catch-All” Analysis for MP’s

- ▷ Particle number
- ▷ Particle size
- ▷ Morphology
- ▷ State of degradation/oxidation
- ▷ Polymer ID(s)
- ▷ Mass Balance

Mass Balance (aka *Mass Fraction*) is the relative composition of the MP sample (such as 35% polyethylene, 40% polypropylene, etc.)

Additives don’t contribute much to the overall mass, so they are not counted in mass balance (but can be used for other purposes)

Polymer ID, Mass Balance, and Additives together can potentially identify the source of the MP

How is Mass Balance Determined?

- ▶ Thermal Extraction GC-MS (TE GC-MS)
- ▶ Pyrolysis GC-MS (Py GC-MS)
- ▶ Thermal Extraction/Desorption GC-MS (TED GC-MS)

Note

- ▶ GERSTEL makes instruments for all three ways of determining mass balance by GC-MS
- ▶ Each way has advantages and disadvantages
- ▶ We do not endorse one over the other; you should choose a technique based on your needs


Method One : Thermal Extraction (TE)

- ▷ Heats a microplastic sample in an empty 3.5" L x ¼" OD tube to ~300 °C in the carrier-gas stream.
- ▷ Is a "special case" of Thermal Desorption GCMS (TD GC-MS)
- ▷ Extracts adsorbed/absorbed materials from the surface of the solid sample
- ▷ Results, if quantitated, are expressed as mass/mass (i.e. ng/mg)
- ▷ The sample is not consumed and only materials from the surface are released; most of the sample remains afterwards
- ▷ Sample size is typically 10-50 mg



TE and Microplastics

- ▶ Instruments used for TD typically have a max temp of 350 to 450 °C. When a TD is used for TE, this is not high enough to decompose the MP polymer and ID/quantify it directly
- ▶ For some polymers, certain additives, such as antioxidants in PET, are always present at a known concentration
- ▶ The non-polymer species can serve as ‘marker compounds’ for the polymer itself, and quantitation of the polymer is possible by quantifying the non-polymer marker



TE and Microplastics


- ▷ TE temperatures are high enough to also ‘extract’ out some other useful species, such as anti-oxidants, UV blockers, etc., but only if they are compatible with the TE and GCMS temperatures used
- ▷ The sample is NOT consumed in TE; the majority of it is left behind in the tube, and only ng to μg quantities are actually injected to the GCMS
- ▷ The tube allows for relatively large samples, including whole filtrate samples should they fit inside the tube

TE analysis of MP's can be done on a Pyrolysis- or TED GC-MS system by heating the sample to only lower temps (such as 320 °C)



Thermal Extraction vs Pyrolysis: Lake Erie Microplastic Particle

GERSTEL



- ▷ Thermal Extraction, 320 °C, 10 min
 - Some SVOC's that may or may not be related to the polymer
 - C₉ alkene is weak, but is a decomp product of one type of PP (but not all)
 - No library match to commercial libraries to any polymer
- ▷ Smart-Ramped Pyrolysis 300 - 800 °C
 - Polymer decomposition pattern for polypropylene clearly visible
 - Library match to PP: 93%
 - Multiple peaks can be used as markers

Thermal Extraction GC-MS: Good, but...

- ▷ Larger, representative sample sizes are possible (whole filter media)
- ▷ Additive information, although not used for mass balance, is readily accessible
- ▷ *Marker Compounds need to be SVOC additives or minor decomposition products that are unique to the polymer in the MP particle. ID and quantification is indirect, inferred.*
- ▷ The lower temperatures used in TD/TE instruments makes carry-over and contamination from difficult matrices more likely, for larger samples in particular



GERSTEL TD 3.5+


The utility for all common classes of common MP polymers (~10 or more) has yet to be demonstrated.

Method Two: PYROLYSIS and Microplastics

GERSTEL



How does pyrolysis GC-MS work?

- ▷ Samples are placed in pyrolysis tubes
- ▷ Can be particles, punches of filter media, or cryo-milled sediment
- ▷ Pyrolysis occurs between 600 and 1000 °C
- ▷ Non-volatile sample material is decomposed
- ▷ The sample is fully consumed in the process
- ▷ 100% of the decomposition products are introduced to the GC/MS
 - A trap can optionally be used for trapping (focusing) and/or splitting of the sample
- ▷ Material(s) in the sample are identified by their pyrolysis fragments





The resulting pyrograms show not what the sample *is*, but what it *became* when heated.

Pyrolysis Examples of Common Plastics




Sample weight:
~ 100 µg



Classic Standard Pulsed Pyrolysis

- ▷ Fast heating, fast pyrolysis
- ▷ Required if going direct to column
- ▷ Same temperature used for ALL polymers in the sample
 - The pulse temp is a compromise
 - Secondary reactions due to overheating add complexity
- ▷ Requires method development
 - Works best for known samples
- ▷ Not suitable for complex unknown samples (can't optimize temp)




A newer alternative: Smart-Ramped Pyrolysis (SRP)

GERSTEL

- ▷ Pyrolysis using a temperature ramp
- ▷ Trapping the pyrolysates is necessary before introduction to the column
- ▷ Polymers are pyrolyzed without overheating and secondary reactions are eliminated
- ▷ One combined GC/MS run follows


Improvements over Pulsed Pyrolysis

- ▷ Better signal and less noise
- ▷ Don't have to choose a single pyro temp
- ▷ Even unknown samples can be analyzed
- ▷ Polymer and additive information in one run are possible



Lake Erie River Sediment: Smart-Ramped Pyrolysis



- ▷ Several polymers present:
 - Methyl Methacrylate / PMMA
 - Styrene (PS)
 - Phenol (Polyphenol)
 - Phthalic Anhydride (Polyester)
- ▷ Monomers can be used as “marker compounds”
- ▷ Markers can be used for ID and for quantitation.
- ▷ Note the presence of some additives as well (Bis(2-ethylhexyl)phthalate)



Fractionated Pyrolysis: Primary MP's in Facial Scrub

GERSTEL

- ▷ Direct SRP Pyrogram up to 600 °C
 - ▷ Very complex but the use of marker compounds allows the MP polymers to be 'pulled out' (more later) →
- ▷ Fractionated Pyrolysis: TE followed by Pyrolysis
 - ▷ A: 100 °C (VOCs)
 - ▷ B: 300 °C (SVOCs, additives)
 - ▷ C: 600 °C (Well defined pyrogram of MP polymer)
- ▷ Fractionated pyrolysis is an easy way to get polymer and additive information from MP's
- ▷ Interpretation is simplified, but data in 2-3 files




PYRO and Microplastics

GERSTEL

Publications

- ▶ Microplastics (MPs) by Pyrolysis GC-MS (AppNote 212)
 - MPs in Filtered Great Lakes Water and Sediment
 - MPs in Body Care Product (Facial Scrub)
- ▶ Microplastics from fish stomach content, Texas Gulf Coast ...
 - E. Hendrickson et al., Mar. Pollut. Bull., 2018, 137, 91-95
- ▶ Microplastics from Western Lake Superior ...
 - E. Hendrickson et al., Environ. Sci. Technol., 2018, 52, 1787-1796
- ▶ Simultaneous Determination of Plastic Particle Identity and Adsorbed Organic Compounds by TD-Pyrolysis GC-MS
 - Molecules 2020, 25, 4985; doi:10.3390/molecules25214985



Pyrolysis GC-MS: Good, but....

GERSTEL

- ▷ Samples are small in mass (0.1 - 1 mg, ideally)
- ▷ Samples small in size (few square mm²)
- ▷ For filters, several punches must be taken
- ▷ Pyrolyzing the whole filter risks overloading the GCMS (pyrolyzers are connected directly to the GCMS)



GERSTEL Pyrolyzer

Thus many smaller samples from a single filtrate need to be run to get a representative sample, or sample must be cryomilled

Method Three: Automated TED-GC-MS

GERSTEL



Thermal Extraction and Pyrolysis of Larger Samples,
Offline in a TGA, followed by Mass Balance
Determination via Evolved Gas Analysis – TD GC-MS

Example: Water Sample on a Filter Crucible


GERSTEL

Large enough for a representative sample

No punching or cutting of paper filters

No other filter media needed

mesh size 5-6 μm




Example: Water Sample on a Filter Crucible

GERSTEL

Analysis of Microplastics in drinking Water and beverages with limited matrix load

- ▷ Filtration set-up
- ▷ For aqueous samples with low matrix load
- ▷ When filtration is complete, the crucible placed directly in TED system for analysis



TED-GC-MS




1st step: Thermal Extraction in TGA



Environmental sample:
Water Filtrate or
Sediment, Soil, Filter Paper
= Matrix + Microplastics

NOTE THE SIZE

Sample size up to 100 mg



Thermogravimetric Analysis (TGA) System

Collection of pyrolysis products on Polydimethylsiloxane (PDMS) using a GERSTEL Twister or similar









Off-line to GC-MS:


Autosampler moves the Twister to a TD GC-MS for analysis

TED GC-MS

2nd Step: analysis of the pyrolysis products




Eight example polymer pyrolysis products and example marker compounds



Polymer	LOD in μg
PE	2,2
PP	0,14
PS	0,08
PET	0,24
PA6	0,24
PA 6.6	3,4
PMMA	0,12
SBR	0,06

Example Calibration for PE samples





E. Dümichen, et al. Water Research 85 (2015) 451

- ▷ Good correlation between spiked mass of polyethylene and the peak area of typical PE marker substances
- ▷ Note the use of four unique marker compounds to make data more reliable
- ▷ Marker compounds are created during pyrolysis, and don't exist "naturally" in the sample
- ▷ Unlikely to have false positives

Application example: street drainage

Styrene-Butadiene rubber (SBR)
from tire wear \approx 4-10 % by mass




P. Eisentraut, E. Dümichen, A. S. Ruhl, M. Jekel, M. Albrecht, M. Gehde, U. Braun, *Environ. Sci. Technol. Lett.*, 2018, 5, 10 608-613.

Polymer particle mixture with Py-GC/MS and TED-GC/MS




Marker Compounds from
ChromIdent Software



*20190619_002_Polymermix_fest_SRP800_01

	UnionMatches	ForwardMatches	Intersection M.	ReverseMatches	MarkerP...	Ambigu...	Unidentif...
> LDPE	uSI 17.3	fSI 42.1	iSI 100	rSI 41	7	10	79
> PA	uSI 14.9	fSI 39.5	iSI 100	rSI 36.4	6	10	81
> PET	uSI 20.1	fSI 39.5	iSI 100	rSI 51.9	15	0	81
> PP	uSI 14.1	fSI 44.5	iSI 100	rSI 31.7	17	2	77
> PS	uSI 14	fSI 35.4	iSI 100	rSI 39.7	6	6	83
> PVC	uSI 6.4	fSI 17.7	iSI 100	rSI 36.1	1	2	93
> SBR	uSI 11.6	fSI 22.8	iSI 100	rSI 47.7	0	5	91



*20191017_004_Polymermix_STD02_Tube02

	UnionMatches	ForwardMatches	Intersection M.	ReverseMatches	MarkerP...	Ambigu...	Unidentif...
> LDPE	uSI 21.8	fSI 31	iSI 100	rSI 70.7	39	12	469
> PA	uSI 8.6	fSI 20.6	iSI 100	rSI 42.9	10	11	497
> PET	uSI 13.1	fSI 21.1	iSI 100	rSI 62	22	1	496
> PP	uSI 20.8	fSI 35.4	iSI 100	rSI 58.6	62	2	454
> PS	uSI 9	fSI 19.6	iSI 100	rSI 43.6	16	4	499
> PVC	uSI 6.4	fSI 11.6	iSI 100	rSI 55.2	5	2	512
> SBR	uSI 6.6	fSI 11.6	iSI 100	rSI 53.5	4	3	512

70.7% match with larger sample;
39 marker peaks found



Contents lists available at ScienceDirect

Journal of Chromatography A

journal homepage: www.elsevier.com/locate/chroma



Automated thermal extraction-desorption gas chromatography mass spectrometry: A multifunctional tool for comprehensive characterization of polymers and their degradation products

E. Duemichen^{a,*}, P. Eisentraut^a, M. Celina^b, U. Braun^a


^a Bundesanstalt für Materialforschung und -prüfung (BAM), Unter den Eichen 87, 12205 Berlin, Germany
^b Sandia National Laboratories, Organic Materials Science Dept., 1853, Albuquerque, NM, 87185-1411, USA



Contents lists available at ScienceDirect

Water Research

journal homepage: www.elsevier.com/locate/watres



Analysis of polyethylene microplastics in environmental samples, using a thermal decomposition method



Erik Dümichen^a, Anne-Kathrin Barthel^a, Ulrike Braun^{a,*}, Claus G. Bannick^b, Kathrin Brand^{b,c}, Martin Jekel^c, Rainer Senz^d

^a BAM Federal Institute for Material Research and Testing, Unter den Eichen 87, 12205 Berlin, Germany
^b UBA Umweltbundesamt, Wörlitzer Platz 1, 06844 Dessau-Roßlau, Germany
^c Technical University of Berlin, Water Urban Area, Strasse des 17. Juni, 10623 Berlin, Germany
^d Beuth University of Applied Sciences, Luxemburger Straße 10, 13353 Berlin, Germany



Contents lists available at ScienceDirect

Environmental Pollution

journal homepage: www.elsevier.com/locate/envpol



Comparison of different methods for MP detection: What can we learn from them, and why asking the right question before measurements matters?[☆]

Anna M. Elert[†], Roland Becker, Erik Duemichen, Paul Eisentraut, Jana Falkenhagen, Heinz Sturm, Ulrike Braun

[†] Federal Institute for Material Research and Testing (BAM), Unter den Eichen 87, 12205 Berlin, Germany

GERSTEL



Contents lists available at ScienceDirect

Chemosphere

journal homepage: www.elsevier.com/locate/chemosphere



Fast identification of microplastics in complex environmental samples by a thermal degradation method



Erik Dümichen^{a,*}, Paul Eisentraut^a, Claus Gerhard Bannick^b, Anne-Kathrin Barthel^b, Rainer Senz^c, Ulrike Braun^a

^a BAM Bundesamt für Materialforschung und -prüfung, Unter den Eichen 87, 12205 Berlin, Germany
^b UBA Umweltbundesamt, Corrensplatz 1, 14195 Berlin, Germany
^c BHT Beuth University of Applied Sciences, Luxemburger Straße 10, 13353 Berlin, Germany



Cite This: Environ. Sci. Technol. Lett. XXXX, XXX, XXX–XXX

pubs.acs.org/journal/estlc

Letter

Two Birds with One Stone—Fast and Simultaneous Analysis of Microplastics: Microparticles Derived from Thermoplastics and Tire Wear

Paul Eisentraut,[†] Erik Dümichen,[†] Aki Sebastian Ruhl,^{‡,§} Martin Jekel,[‡] Mirko Albrecht,[§] Michael Gehde,[§] and Ulrike Braun^{†,*}

[†] Bundesanstalt für Materialforschung und -prüfung, Unter den Eichen 87, 12205 Berlin, Germany
[‡] Technische Universität Berlin, Straße des 17. Juni 135, 10623 Berlin, Germany
[§] Technische Universität Chemnitz, Reichenhainer Straße 70, 09126 Chemnitz, Germany

Supporting Information

ABSTRACT: Analysis of microplastic particles in environmental samples needs sophisticated techniques and is time intensive due to sample preparation and detection. Alternatives to the most common (micro-) spectroscopic

Microplastic Analysis using TED-GC-MS



The logo for GERSTEL, featuring the word "GERSTEL" in a bold, black, sans-serif font. The text is enclosed within a dark red rectangular border that has a slight 3D effect, with the top and bottom lines being thicker than the sides.

So what does it all mean?

The Short Answers

GERSTEL

- ▶ All three techniques can determine the mass balance of MP samples on filter media; all three can determine the ID of polymers AND additives
- ▶ Thermal Extraction GC-MS analyzes larger representative samples but is indirect: polymers must have unique marker additives and samples should be relatively clean (e.g., PET MP's in drinking water); MP analysis range is limited / undetermined
- ▶ Pyrolysis GC-MS offers direct MP ID and analysis of a wide range of polymers, can also do additives if needed, and can handle a wider range of matrices; sample sizes are smaller & multiple reps required for representative sampling
- ▶ TED GC-MS can do a wide range of polymers, additives runs, larger representative samples, and being off-line to the GC-MS makes it the most robust; more investment (two instruments and training) are needed



Thermal Desorber used as a Thermal Extractor



Pyrolyzer



Thermal Extraction/Desorption (TED-GCMS) System

Instrumentation Comparison

GERSTEL

	Thermal Extraction (~315 °C)	Pyrolysis (600 to 1000 °C)	TED (600 to 1000 °C)
Sample Intro	On-line to GCMS	On-line to GCMS	Off-line from GCMS
Process	Thermally extract to trap, then to column (through valve / transfer line or direct to column)	Pyrolyze direct to GCMS (direct to column, or optional trapping step before column)	Extraction/Pyrolysis in TGA; trap on PDMS and transfer to TD-GCMS (PDMS trap is desorbed in TD; re-trapping before column recommended)
Bake-out	Up to 450 °C, N ₂ or He	Up to 1000 °C, N ₂ or He	Up to 1000 °C, in N ₂ , He, O ₂
Sample Size	Typically 10-50 mg	0.1 - 10 mg (0.1 – 1 mg typically; all goes to inlet)	Typically 10-50 mg
Representative Sample	One Run	Multiple Runs	One Run
Type of Data	Mass Spectral	Mass Spectral	MS and TGA both

Analytical Comparison

GERSTEL

	Thermal Extraction (~310 °C)	Pyrolysis (600 to 1000 °C)	TED (600 to 1000 °C)
Sample Prep (filter media)	Dry, solvent rinse of filter, dry again, analyze	Dry, punch out correct sized samples, analyze (multiple punches very likely needed)	Dry, punch out correct sized samples, analyze (one, larger punch is typically sufficient)
Polymer ID	Yes, Indirect / Inferred (marker additives)	Direct (un-zipped monomer or targeted degradant markers)	Direct (un-zipped monomer or targeted degradant markers)
Range of Polymers	Limited to uniqueness of additives (typ. PET)	Eight or more common MP polymers (PET, PE, PP, PS, ...)	Eight or more common MP polymers (PET, PE, PP, PS, ...)
Mass Balance	Yes, Indirect / Inferred quant through additives	Yes, direct quant through monomers / degradants	Yes, direct quant through monomers / degradants
Additives <small>(note: not needed for mass balance)</small>	Yes, Direct, one step	Yes, Direct, one or two steps (300 °C and then 600 to 1000 °C), if two steps needed	Yes, Direct, one or two steps (300 °C and then 600 to 1000 °C), if two steps needed

Cost and Complexity

GERSTEL

	Thermal Extraction (~310 'C)	Pyrolysis (600 to 1000 'C)	TED (600 to 1000 'C)
Instruments	Three (TD, GC, and MS)	Three (PYRO, GC, and MS)	Four (TGA, TD, GC, and MS)
Cost	\$\$	\$\$	\$\$\$
Support	One Provider	One Provider	Two Providers (TGA and TD-GCMS)
Whole Filtrate in One Run	Possible	No	Possible
Carry-over Risk (sludge, tissue, ...)	Highest (TE only = lowest temp technique)	High (PYRO only less risky) or Higher if TE step used	High (or Moderate with O ₂ cleaning step)

Thanks to:

GERSTEL

BAM

Bundesanstalt für
Materialforschung
und -prüfung

GERSTEL

Ulrike Braun
Erik Dümichen
Paul Eisentraut
et al.

Eike Kleine-Benne
Huan-Xiao Zhou,
Jackie Whitecavage,
John R. Stuff