

# A New Rapid, Simple, and Efficient Extraction Method of PFAS from Soil


Alicia D. Stell, Ph.D.  
CEM Corporation

# PFAS: Per- and Polyfluoroalkyl Substances




# Sample Preparation is the Bottleneck

## Time Spent on Typical Analysis




R.E. Majors, LC/GC Magazine

# Typical PFAS Extraction from Soil




# EDGE Extraction




# Standard versus EDGE Extraction

## Standard



## EDGE



# EDGE Advantages

- Automation
- In-cell Cleanup
- Simple
- Multi-matrix Multi-residue Method
- PFAS free



# PFAS Free Instrumentation

- PEEK and Polypropylene Tubing
- Side Enclosure
- Centrifuge Tube Collection
- Nitrogen Option
- PFAS free consumables





# EDGE Sample Prep



Layered in Q-Cup (sample holder)

# EDGE Sample Prep



Layered in Q-Cup (sample holder)

Q-Disc

S1 (C9+G1+C9)

# EDGE Sample Prep



Layered in Q-Cup (sample holder)

Sample	5 g: Soil
Q-Disc	S1 (C9+G1+C9)

# EDGE Sample Prep






Layered in Q-Cup (sample holder)	
Spike	Wellington Laboratories: PFAC30PAR and MPFAC-24ES
Sample	5 g: Soil
Q-Disc	S1 (C9+G1+C9)

# EDGE Method

## Edit Method - PFAS project

Settings	Name	PFAS project
Cycles	Q-Discs® Type	S1 >
Wash	Collection	Combined >

### Notes

  9:14 09/01/2020 

# EDGE Method

Edit Method - PFAS project									
Settings	Cycle	Solvent	Top Volume (mL)	Bottom Volume (mL)	Temp (°C)	Hold Time	Rinse Solvent	Rinse Volume (mL)	Vial
Cycles	1	80:20 MeOH:Water + 0.3% NH4OH >	10	0	65	4:00 >	80:20 MeOH:Water + 0.3% NH4OH >	0	1
Wash	2	80:20 MeOH:Water + 0.3% NH4OH	10	0	65	4:00	80:20 MeOH:Water + 0.3% NH4OH	0	1
[Large blue redacted area]									



9:15  
09/01/2020



# EDGE Method





**Cycle 1**

Solvent	80:20 MeOH:Water + 0.3% NH <sub>4</sub> OH >	
Volume (mL)	Top <input type="text" value="10"/>	Bottom <input type="text" value="0"/>
Rinse Solvent	80:20 MeOH:Water + 0.3% NH <sub>4</sub> OH >	
Volume (mL)	<input type="text" value="0"/>	
Temperature (°C)	<input type="text" value="65"/>	
Hold Time	4:00 >	
Bubbling	<input type="checkbox"/>	




OK Cancel

9:15  
/01/2020

# EDGE Method

Edit Method - PFAS project									
Settings	Wash	Solvent	Volume (mL)	Hold Time	Temperature (°C)				
Cycles	1	Methanol >	30	0:30 >	65				
Wash	2	Methanol	30	0:30	65				
	3	80:20 MeOH:Water + 0.3% NH4OH	10	--:--	----				
[Redacted area]									

9:15  
09/01/2020





# Minimal Post work



Dilute to known volume




Neutralize

# Instrument Conditions

## Acquisition Parameters

HPLC System	ExionLC™ HPLC System
MS/MS System	5500+ QTRAP
Ion Source	Turbo V™
Injection Volume	1 µL
Analytical Column	Phenomenex. Gemini C18 3 µm, 3 x 100 mm
Delay Column	Phenomenex. Gemini C18 5 µm, 3 x 50 mm
LC Flow Rate	500 µL/min
Mobile Phases	Water & MeOH (both with 10 mM ammonium acetate)
Source & MS Parameters	TEM = 600 C GS1= 60, GS2 = 60 ISV = -4500 CUR = 35, CAD = 8




### Flow program


Flow program  Simple

	Time	Flow	A.Con	B.Con	B.Curve
1		0.5000	90.0	10.0	0
2	0.50	0.5000	90.0	10.0	0
3	1.00	0.5000	45.0	55.0	0
4	5.50	0.5000	1.0	99.0	0
5	7.00	0.5000	1.0	99.0	0
6	7.10	0.5000	90.0	10.0	0
7	10.00	0.5000	90.0	10.0	0

# Example Chromatogram – 2 ng/mL standard



# Example Chromatogram – 2 ng/mL standard



# Soil Data: High Spike 20 ppb

PFAS	% Recovery	% RSD (n=3)
1H, 1H, 2H, 2H-perfluorodecane sulfonic acid (8:2 FTS)	91	9.7
1H, 1H, 2H, 2H-perfluorooctane sulfonic acid (6:2 FTS)	90	11
1H,1H,2H,2H-perfluorohexane sulfonic acid (4:2 FTS)	101	5.8
N-ethylperfluoro-1-octanesulfonamidoacetic acid (EtFOSAA)	95	12
N-methylperfluoro-1-octanesulfonamidoacetic acid (MeFOSAA)	87	12
perfluoro-1-butanesulfonic acid (PFBS)	97	13
perfluoro-1-decanesulfonic acid (PFDS)	91	14
perfluoro-1-heptanesulfonic acid (PFHpS)	96	12
perfluoro-1-nonanesulfonic acid (PFNS)	91	16
perfluoro-1-octanesulfonamide (PFOSA)	85	20
perfluoro-1-pentanesulfonic acid (PFPeS)	88	11
perfluorohexanesulfonic acid (PFHxS)	97	11

# Soil Data: High Spike 20 ppb

PFAS	% Recovery	RSD (n=3)
perfluoro-n-butanoic acid (PFBA)	87	12
perfluoro-n-decanoic acid (PFDA)	89	15
perfluoro-n-dodecanoic acid (PFDoA)	87	10
perfluoro-n-heptanoic acid (PFHpA)	86	10
perfluoro-n-hexanoic acid (PFHxA)	88	11
perfluoro-n-nonanoic acid (PFNA)	88	10
perfluoro-n-octanoic acid (PFOA)	91	10
perfluoro-n-pentanoic acid (PFPeA)	86	9.5
perfluoro-n-tetradecanoic acid (PFTeDA)	91	11
perfluoro-n-tridecanoic acid (PFTrDA)	87	11
perfluoro-n-uoecanoic acid (PFUdA)	87	12
perfluorooctanesulfonic acid (PFOS)	98	8.9


# No Carryover – Sample Immediately After High Spike

## 2 ng/mL std


## High Spike

## Blank


L1 std - HFPO-DA (Standard) 285.0 / 169.0...04-CEM samples.wiff), (sample Index: 3)  
Area: 1.943e5, Height: 8.332e4, RT: 4.41 min




Sample 12 - HFPO-DA (Quality Control) 28...4-CEM samples.wiff), (sample Index: 20)  
Area: 1.466e6, Height: 6.350e5, RT: 4.37 min




Sample 13 - HFPO-DA (Unknown) 285.0 / 1...4-CEM samples.wiff), (sample Index: 21)  
Area: 7.439e2, Height: 2.750e2, RT: 4.37 min




L1 std - PFHxS (Standard) 399.0 / 80.0 - (...4-04-CEM samples.wiff), (sample Index: 3)  
Area: 2.304e5, Height: 7.853e4, RT: 4.70 min




Sample 12 - PFHxS (Quality Control) 399.0...04-CEM samples.wiff), (sample Index: 20)  
Area: 1.688e6, Height: 6.155e5, RT: 4.67 min




Sample 13 - PFHxS (Unknown) 399.0 / 80.0...4-CEM samples.wiff), (sample Index: 21)  
Area: 1.267e3, Height: 3.684e2, RT: 4.66 min




L1 std - PFOA (Standard) 413.0 / 369.0 - (...4-04-CEM samples.wiff), (sample Index: 3)  
Area: 6.766e5, Height: 3.003e5, RT: 5.12 min



Sample 12 - PFOA (Quality Control) 413.0...04-CEM samples.wiff), (sample Index: 20)  
Area: 5.445e6, Height: 2.360e6, RT: 5.09 min



Sample 13 - PFOA (Unknown) 413.0 / 369.0...4-CEM samples.wiff), (sample Index: 21)  
Area: 4.427e3, Height: 1.716e3, RT: 5.09 min




# ERA Soil CRM Data

PFAS	% Recovery	RSD (n=2)
PFBA	88	13
DONA	96	18
4:2 FTS	106	15
HFPO-DA	86	13
PFDoA	92	7.5
PFHpS	99	19
PFHxS	108	20
PFOSA	95	24
PFOS	105	15
PFUnA	98	16



# Conclusions

- EDGE instrumentation is PFAS free
- No contamination from sample preparation or consumables
- Good recoveries and RSD values for high spiked and CRM soil samples
- No carryover from high spike sample
- EDGE is a rapid, simple, and efficient technology for the extraction of PFAS from soil



# We Are Where You Are



- 60 associates in the United States and Canada
- 24-hour telephone support

# Extraction Team



Alicia Stell, PhD  
Lead R&D Scientist



Candice Cashman, PhD  
Senior Scientist



Brittany Fessler  
Product Specialist



Benedict Liu  
Applications Scientist

[Molecular.support@cem.com](mailto:Molecular.support@cem.com)



# Questions

[Molecular.Support@cem.com](mailto:Molecular.Support@cem.com)

CEM