


# A Rapid Method for the Analysis of Air Toxics Based on US EPA TO-15

**Terry Jeffers**

Senior Application Scientist

August 9, 2021

 The world leader in serving science



## What are HAPs?

- HAPs are pollutants that are known or suspected to cause serious health effects or adverse environmental effects.
- Examples include:



Benzene  
Found in gasoline



Dichloromethane  
Used as a solvent and paint stripper

## Background

- US EPA Method TO-15 is established method for hazardous air pollutants (HAPs).
- Utilises canister sampling method and thermal desorption gas chromatography mass spectrometry (TD-GC-MS).
- Ingress of water can negatively affect results.
- Standard run times limit sample throughput.
- Other challenges include the wide range of sample compound temperatures and varying temperature and humidity of samples.

# Objectives

- Demonstrate an alternative chromatographic approach to speed up analysis compared to US EPA TO-15.
- Demonstrate the combined use of an innovative trap-based water removal device as part of a robust TD-GC-MS system.
- Demonstrate both typical and 0.5 – 50 ppb Calibration data for US EPA TO-15

# Equipment

## Sample introduction




Markes™ CIA Advantage™ canister sampler coupled to a Markes™ Kori-xr™ water removal device and a Markes™ Unity-xr™ thermal desorber

## Sample analysis



Thermo Scientific™ ISQ™ 7000 mass spectrometer (MS) coupled to Thermo Scientific™ TRACE 1310 gas chromatograph (GC)

# Sample flow path



# Sample Preparation Method conditions

## Markes™ CIA Advantage-xr™

Canister sampling volume	Up to 1000 mL
Water removal	Kori-xr™
IS loop fill	1 min

## Markes™ Unity-xr™

Trap type	Cold trap, "TO-15/TO-17 Air toxics"
Trap temperature	-30 °C to 300 °C
Flow path	120 °C
Loop equilibration	0.1 min
IS loop injection	1.0 min at 50 mL/min
Sample flow	50 mL/min



# GC-MS Method conditions

## Thermo Scientific™ TRACE 1310 GC

Inlet temperature	260 °C
Injection mode	Splitless, 3 min
Split flow	20 mL/min
Carrier gas flow	He, 0.8 mL/min
GC oven temperature program	Initial 35 °C for 3 min, Ramp 14 °C/min to 100 °C, Ramp 20 °C/min to 230 °C, hold for 4 min
Total run time	18 min

## Thermo Scientific™ ISQ 7000 mass spectrometer

Transfer line temperature	230 °C
Mode/range	Full scan, 35-260 amu
Ion source	Thermo Scientific™ ExtractaBrite™
Ion source temperature	310 °C
Ionisation mode	El @ 70 eV
Solvent delay	1.09 min
Dwell/scan time	0.15 s
Emission current	25 µA

## Column

TG-VMS	20 m x 0.18 mm x 1.0 µm Film
--------	------------------------------





# Thermo Scientific ISQ 7000 MS - Never Vent

What is NeverVent technology?




Increase GC-MS Uptime with the vacuum probe interlock **(VPI)**

Through the VPI, no need to vent mass spec system for extracting the wireless ExtractaBrite ion source



Extends the capability of the Vacuum Probe Interlock **(VPI)** design with the newly introduced source plug, **V-Lock**

## Standard Preparation



- Bromochloromethane, 1,4-Difluorobenzene, Chlorobenzene-d5 used as Internal Standards
- 4-bromofluorbenzene – Tune Mix


# Result – Sample to Sample Cycle Time




Sample Preparation



Sample Separation & Detection




Sample Results



25 minutes – Sample to Sample

# Result - Chromatography





Separation of 50 ppb standard achieved in under 13 minutes

# Result – BFB Tune Check

TABLE 4. 4-BFB KEY IONS AND ION ABUNDANCE CRITERIA

Mass	Ion Abundance Criteria
50	15 to 40% of mass 95
75	30 to 60% of mass 95
95	Base Peak, 100% Relative Abundance
96	5 to 9% of mass 95
173	<2% of mass 174
174	>50% of mass 95
175	5 to 9% of mass 174
176	>95% but < 101% of mass 174
177	5 to 9% of mass 176


> 20 Scans across Peak



No.	Name	Eval. Result	Operator	Ref. Value 1	Ref. Value 2	Result
1	m/z 50 - 15 to 40% of m/z 95	20.4	between	15	40	Passed
2	m/z 75 - 30 to 60% of m/z 95	53.8	between	30	60	Passed
3	Base Peak (m/z 95)	95.1	=	95		Passed
4	m/z 96 - 5 to 9% of m/z 95	6.4	between	5	9	Passed
5	m/z 173 - Less than 2% of m/z 174	0.6	<	2		Passed
6	m/z 174 - Greater than 50% of m/z 95	81.4	>	50		Passed
7	m/z 175 - 5 to 9% of m/z 174	7.5	between	5	9	Passed
8	m/z 176 - 95 to 101% of m/z 174	96.3	between	95	101	Passed
9	m/z 177 - 5 to 10% of m/z 176	6.8	between	5	10	Passed
<b>Overall Result:</b>						<b>Passed</b>

Tune Check – Top 3 scans across apex of peak

# Result – Internal Standard Chromatography



## Internal Standards and Tuning Component

# Result – Internal Standard Stability


Injection Name	Area Counts	Area Counts	Area Counts
	MS Quantitation Bromochloromethane	MS Quantitation 1,4-Difluorobenzene	MS Quantitation Chlorobenzene-d5
Blank-1	48866	236955	215257
ICAL 0-5PPbv	37589	200913	192575
ICAL 1-00PPbv	40130	199230	194095
ICAL 2-00PPbv	44729	201277	189526
ICAL 3-00PPbv	46003	206774	184438
ICAL 4-00PPbv	47968	210992	189314
ICAL 5-00PPbv	47224	205001	181838
ICAL 8-00PPbv	49643	218012	187905
ICAL 20-00PPbv	48914	218198	190552
ICAL40-00PPbv	48219	212215	182758
ICAL 50-00PPbv	47111	216549	177721
rinse 1	45987	207411	177078
rinse-2	46395	207897	178401
rinse-3	45286	207740	178508
MDL -1	34067	182574	173666
MDL -2	33658	177160	174686
MDL -3	33827	176371	174661
MDL -4	34406	176554	177583
MDL -5	35442	184933	181716
MDL -6	35523	182928	182437
MDL -7	35703	185743	182678
MDL -8	35878	190579	187870
Validation-1	47590	218333	194671
Validation-2	48807	215986	198117
Validation-3	49122	224642	199111
Validation-4	51248	231040	209124
Validation-5	50179	226201	207141
Validation-6	49717	224492	204090
Validation-7	52020	231826	213275
<b>Average</b>	<b>43836</b>	<b>206156</b>	<b>188993</b>
<b>Standard Deviation</b>	<b>6368</b>	<b>17851</b>	<b>11978</b>
<b>% RSD</b>	<b>15%</b>	<b>9%</b>	<b>6%</b>

## % RSD

- **Bromochloromethane = 15 %**
- **1,4-Difluorobenzene = 9 %**
- **Chlorobenzene-d5 = 6 %**

## Internal Standard Area Count Stability

# Results - Chromatography




Peak number	HAP	Retention time (min)	Quantitation ion (m/z)
1	Dichlorodifluoromethane	1.33	85
2	1,1-Difluoroethane	1.35	65
3	1,2-Dichlorofluoromethane	1.42	135
4	Chloromethane	1.48	50
5	Vinyl chloride	1.53	62
6	1,3-Butadiene	1.53	54
7	Bromomethane	1.78	94
8	Chloroethane	1.88	64
9	Trichlorofluoromethane	1.98	101
10	1,1-Dichloroethene	2.42	61
11	Carbon disulfide	2.42	76
12	Freon 113	2.47	101
13	Isopropyl Alcohol	2.93	45
14	Methylene chloride	2.98	49
15	Acetone	3.04	43
16	Hexane	3.25	57
17	MTBE	3.30	73
18	tert-butanol	3.41	59
19	Diisopropyl ether	3.70	45
20	ETBE	4.07	59
21	cis-1,2-Dichloroethene	4.31	61
22	Bromochloromethane	4.50	128
23	Chloroform	4.59	83
24	Carbon Tetrachloride	4.69	119
25	Tetrahydrofuran	4.73	42
26	1,1,1-Trichloroethane	4.76	97
27	2-Butanone	4.89	43
28	Benzene	5.12	78
29	TAME	5.26	73
30	1,2-dichloro-ethane	5.30	62
31	Trichloroethylene	5.66	130
32	1,4-Difluorobenzene	5.71	114

EIC for a 50-ppb standard from 1.23 min to 6 min



# Results - Chromatography




Peak number	HAP	Retention time (min)	Quantitation ion (m/z)
33	1,2-Dichloropropane	6.13	63
34	Bromodichloromethane	6.20	83
35	1,4-Dioxane	6.39	88
36	Cis-1,3-Dichloropropene	6.79	75
37	Toluene	7.00	91
38	Tetrachloroethylene	7.34	166
39	4-Methyl-2-pentanone (MIBK)	7.38	43
40	Trans-1,3-Dichloropropene	7.39	75
41	1,1,2-trichloroethane	7.54	97
42	Dibromochloromethane	7.69	129
43	1,1-Dibromoethane	7.89	107
44	2-Hexanone	8.15	43
45	Chlorobenzene-d5	8.37	117
46	Chlorobenzene	8.38	112
47	Ethylbenzene	8.42	91
48	1,1,1,2-Tetrachloroethane	8.44	133
49	M,p-Xylene	8.56	91
50	O-Xylene	8.92	91
51	Bromoform	8.97	173
52	Styrene	8.97	104
53	4-Bromofluorobenzene	9.42	95
54	1,1,1,2-Tetrachloroethane	9.61	83
55	4-Ethyltoluene	9.64	105
56	1,3,5-Trimethylbenzene	9.72	105
57	1,2,4-Trimethylbenzene	10.04	105
58	1,3-Dichlorobenzene	10.28	146
59	1,4-Dichlorobenzene	10.36	146
60	1,2-Dichlorobenzene	10.70	146
61	Hexachlorobutadiene	11.89	225
62	1,2,3-Trichlorobenzene	11.90	180
63	Naphthalene	12.17	128

EIC for a 50-ppb standard from 6 min to 13 min


## Result – Method Validation

- Multi point Calibration
  - 10 calibration points
  - 0.5 – 50 ppbv
  
- Minimum Detection Limit Study
  - 0.5 ppbv standard
  
- Initial Demonstration of Capability (IDC)
  - 8.0 ppbv standard


# Results - Calibration



# Results – Sensitivity




# Results – Replicate precision



# Chromeleon CDS Benefits

- Data Review – at a glance



thermoscientific

APPLICATION NOTE

10729

## A rapid method for the analysis of air toxics based on US EPA TO-15

Authors: Terry Jeffers<sup>1</sup>, David Lee<sup>2</sup> and Dwain Cardona<sup>3</sup>

<sup>1</sup>Thermo Fisher Scientific, West Palm Beach, FL, USA

<sup>2</sup>Thermo Fisher Scientific, Runcorn, UK

<sup>3</sup>Thermo Fisher Scientific, Austin, TX, USA

Keywords: US EPA, VOCs, TO-15 volatiles, air toxics, gas chromatography, single quadrupole mass spectrometry, selected ion monitoring, thermal desorption gas chromatography mass spectrometry, TD-GC-MS, canister, environmental lab, air analysis, ambient air monitoring

### Goal

The following application demonstrates an alternative chromatographic approach to US EPA Method TO-15 for the rapid determination of toxic organic compounds in ambient air using a combined TD-GC-MS solution from Markes International and Thermo Fisher Scientific.

Method linearity, RRF variation, method detection limit (MDL), and precision were assessed to evaluate method performance.

### Introduction

US EPA Method TO-15 is an established method used for the measurement of hazardous air pollutants (HAPs) also known as air toxics. This subset of volatile organic compounds (VOCs) is collected and analyzed using



canister sampling methods and thermal desorption gas chromatography mass spectrometry (TD-GC-MS). Despite the popularity of canister sampling, maintaining system integrity and obtaining consistent results can be difficult for routine air analysis labs. Ingress of water to the analytical instrument can negatively impact analyte response and repeatability, as well as reduce the lifetime of the column and detector. In addition to this, traditional canister pre-concentration technologies are challenged by the wide ranges of sample compound concentrations and by the varying temperatures and humidity at sampling locations.

In this application note, we demonstrate the combined use of an innovative trap-based water removal device, alternative column choice, and a robust TD-GC-MS configuration to overcome the challenges experienced in the analysis of volatile air toxics in accordance with US EPA Method TO-15 and provide advantages in laboratory productivity.

ThermoFisher  
SCIENTIFIC

## Conclusions

- Alternative column dimensions reduces analysis times, increasing sample throughput.
- Effective chromatographic separation and accurate quantitation of 60 HAPs in under 13 minutes is achieved.
- Sample to Sample run time is 25 minutes
- Excellent linearity was obtained over a concentration range of 0.5 to 50 ppb, comfortably meeting the requirements of <30% RRF RSDs.
- Replicate precision for all target compounds was below the requirement of <25% from n=7 replicates
- Sensitivity assessment demonstrated that lower method detection limits were achieved than required by the published method



# Thank you

