

A Novel EI Source Optimized for Use with Hydrogen Carrier Gas in GC/MS and GC/MS/MS

NEMC 2022

Operational issues impacting the
environmental laboratory industry

Angela Smith Henry
GC,MS Supplies Application Chemist
Anastasia Andrianova
MSD Applications Chemist
Bruce Quimby
MSD Applications Chemist
Amanda McQuay
CSD R&D Applications Chemist



Helium

Helium remains the best carrier gas to use for GC/MS analysis, but...

- Supply shortages are becoming common
- Costs are rising
- Questions about sustainability are being raised



Reasons to use Hydrogen Carrier Gas


- Readily available (H_2 already used for FID and other detectors)
- Sustainable
- Lower Cost
- Cleans source during use
- Available on demand (H_2 generator) or by cylinder
- Faster analysis
- Lower temperature separation possible
- Move to “more efficient” columns
 - 30 m x 0.25 mm x 0.25 μm \rightarrow 20m x 0.18 mm x 0.18 μm

But....





Source-Induced Problems with Hydrogen Carrier Gas: Nitrobenzene Conversion

Expected mass spectrum



What really happens with H₂




Conversion of nitrobenzene to aniline with H₂ carrier gas



It's Not Just Nitrobenzene...

Compound Class	Example Compound(s)	H ₂ Carrier Gas + Conventional GC/MS Source
Nitro compounds	Nitrobenzene, fenitrothion, ethalfluralin	Unacceptable
Heavily chlorinated compounds	DDT, endrin, heptachlor, BHC compounds, pentachlorophenol	Unacceptable
Polycyclic aromatic hydrocarbons (PAHs)	Benzo(b)fluoranthene, benzo[g,h,i]perylene,	Neutral
Alkanes >C24	Tetratriacontane, hexadecane, tetracontane	Neutral
Pesticides	Deltamethrin, fipronil, permethrin, captan	Unacceptable
Fragrance/flavor compounds	Musk ketone, musk ambrette, linalool	Unacceptable
Volatile organic compounds	1,4-dioxane, trichloromethane, bromodichloromethane	Neutral



Permethrin




Phenoxytoluene



Musk ambrette




Deltamethrin



(m-phenoxyphenyl)-acetonitrile



Agilent HydroInert Source for Hydrogen Carrier Gas on GC/MS




- Allows for the use of Hydrogen Carrier Gas with better supply and reduced costs
- Faster, shorter separations
- Reduces loss of sensitivity and spectral anomalies
- Reduced source cleaning and maintenance


Back to Nitrobenzene and H₂ Carrier Gas



Existing extractor source




HydroInert source







HydroInert source prevents hydrogenation

H₂ Cleavage Reactions in the Source: Captan






Conventional source + H₂

HydroInert + H₂




Less hydrogenation in HydroInert source and better NIST match scores

H₂ Cleavage Reactions in the Source: Deltamethrin



- (m-phenoxyphenyl)-acetonitrile identified in some runs with LMS 74.4
- Deltamethrin identified in most runs with an LMS range of 70-77



Deltamethrin identified in **all** runs with HydroInert (LMS 90+)

We Retain Our Mass Spectra, but What About Calibration Curves?

Can we retain our calibration range?


- 0.1 – 100 µg/mL

Helium + conventional source + 9 mm extraction lens
Nitrobenzene
Average RF %RSD: 6.33% RSD


H₂ + HydroInert source + 9 mm lens
Nitrobenzene
Average RF %RSD: 6.59% RSD

Yes, we can retain calibration range, too
...with consumable and method parameter
changes.


Nitrobenzene - 13 Levels, 13 Levels Used, 13 Points, 13 Points Used, 0 QCs



Nitrobenzene - 13 Levels, 13 Levels Used, 13 Points, 13 Points Used, 0 QCs




What Happens with Alkanes and H₂ Carrier Gas in GC/MS?



- Improved peak shape (reduced tailing) with HydroInert source



- Shorter chain alkanes do well with H₂
- Longer chain alkanes have issues with tailing

Do I Have Similar Responses Between “Normal” Helium Data and H₂ + HydroInert Source Data?



- Using the same column, injection parameters, and tune (atune.u)
- Faster constant flow rate with H₂
- Benefits of H₂
 - Narrower peaks
 - Increased resolution
- Retained sensitivity

Comparing Helium and HydroInert + H₂ Results





Compound	Response Factor (RF) from EPA Method 8270E	RF He GC/MS	RF H ₂ HydroInert GC/MS	RF H ₂ HydroInert GC/MS/MS
Acenaphthene	0.9	1.3	1.1	0.2
Acetophenone	0.01	1.2	0.4	1.0
Benzo(b)fluoranthene	0.7	1.4	1.2	1.2
Bis(2-chloroethoxy)methane	0.3	0.4	0.3	0.7
Di-n-butyl phthalate	0.01	1.3	0.8	0.8
Hexachlorocyclopentadiene	0.05	0.3	0.2	0.03
2,4-Dinitrophenol	0.01	0.2	0.1	0.02
4-Nitrophenol	0.01	0.2	0.14	0.05
2,4,6-Trichlorophenol	0.2	0.3	0.2	0.2
62 compounds		3% low RF	11% low	19% low



Better resolution with H₂ and the majority of RFs match to EPA 8270 RF guidelines

Sensitivity Comparison of Sources with H₂ Carrier Gas

Carbon tetrachloride





Bromoform





What Happens to Our Compound Classes?

Compound Class	Example Compound(s)	H ₂ Carrier Gas + Conventional GC/MS Source	H ₂ Carrier Gas + HydroInert Source
Nitro compounds	Nitrobenzene, fenitrothion, ethalfluralin	Unacceptable	Differentiating
Heavily chlorinated compounds	DDT, endrin, heptachlor, BHC compounds, pentachlorophenol	Unacceptable	Differentiating
Polycyclic aromatic hydrocarbons (PAHs)	Benzo(b)fluoranthene, benzo[g,h,i]perylene,	Neutral	Neutral
Alkanes >C24	Tetratriacontane, hexadecane, tetracontane	Neutral	Neutral
Pesticides	Deltamethrin, fipronil, permethrin, captan	Unacceptable	Differentiating
Fragrance/flavor compounds	Musk ketone, musk ambrette, linalool	Unacceptable	Differentiating
Volatile organic compounds	1,4-dioxane, tichloromethane, bromodichloromethane	Neutral	Differentiating




Musk ambrette



Headspace GC/MS: VOCs with H₂ Carrier and HydroInert

8697 HS/8890/5977C VOCs with H₂ Carrier and HydroInert



Pulsed split 20:1

Agilent J&W DB-624 UI 20 m x 0.18 mm x 1.0 μm

0.95 mL/min

5 g Na₂SO₄ added to each vial

Cal from 0.25 ppb to 25 ppb Scan mode

Cal from 0.05 ppb to 25 ppb SIM mode

VOCs: Library Match Scores vs NIST20

	RT (min)	LMS NIST20	
Dichlorodifluoromethane	0.507	92	
Chloromethane	0.616	97	
Ethene, chloro-	0.698	91	
Methane, bromo-	0.891	96	
Ethyl Chloride	0.945	92	
Trichloromonofluoromethane	1.068	96	
Ethyl ether	1.199	97	
Ethene, 1,1-dichloro-	1.289	98	
Acetone	1.316	87 *	■
Methane, iodo-	1.351	99	
Carbon disulfide	1.380	95	
Allyl chloride	1.433	97	
Methylene chloride	1.479	97	
2-Propenenitrile	1.574	90	
Ethylene, 1,2-dichloro-, (E)-	1.587	99	
Propane, 2-methoxy-2-methyl-	1.592	98	
Ethane, 1,1-dichloro-	1.747	97	
Ethylene, 1,2-dichloro-, (Z)-	1.968	95	
Propane, 2,2-dichloro-	1.970	80 **	■
2-Butanone	1.975	72 *	■
Propanenitrile	1.994	67 *	■
2-Propenoic acid, methyl ester	2.010	97	
Methylacrylonitrile	2.053	95	
Methane, bromochloro-	2.060	97	
Trichloromethane	2.088	98	
Tetrahydrofuran	2.091	96	
Ethane, 1,1,1-trichloro-	2.169	98	
Butane, 1-chloro-	2.206	97	

1-Propene, 1,1-dichloro-	2.233	96	
Carbon Tetrachloride	2.237	96	
Benzene	2.317	94	
Benzene, fluoro-	2.425	98	
Trichloroethylene	2.578	99	
Propane, 1,2-dichloro-	2.673	98	
Methyl methacrylate	2.714	98	
Methane, dibromo-	2.723	98	
Methane, bromodichloro-	2.786	98	
Propane, 2-nitro-	2.884	93	
Chloromethyl cyanide	2.888	63 *	■
1-Propene, 1,3-dichloro-, (Z)-	2.986	98	
Toluene	3.146	99	
1-Propene, 1,3-dichloro-, (E)-	3.240	98	
Methacrylic acid, ethyl ester	3.284	98	
Ethane, 1,1,2-trichloro-	3.329	98	
Tetrachloroethylene	3.411	91	
Propane, 1,3-dichloro-	3.412	90	
Methane, dibromochloro-	3.524	98	
Ethane, 1,2-dibromo-	3.586	99	
Benzene, chloro-	3.836	99	
Ethane, 1,1,1,2-tetrachloro-	3.876	96	
Ethylbenzene	3.892	98	
Benzene, 1,3-dimethyl-	3.953	99	
o-Xylene	4.165	89	
Styrene	4.169	96	
Methane, tribromo-	4.266	99	
Benzene, (1-methylethyl)-	4.364	98	
p-Bromofluorobenzene	4.446	97	
Ethane, 1,1,2,2-tetrachloro-	4.521	97	



Benzene, bromo-	4.529	97	
Propane, 1,2,3-trichloro-	4.547	84	
2-Butene, 1,4-dichloro-, (E)-	4.554	65 **	■
Benzene, propyl-	4.592	98	
Benzene, 1-chloro-2-methyl-	4.638	98	
Mesitylene	4.691	91	
Benzene, 1-chloro-4-methyl-	4.698	95	
Benzene, tert-butyl-	4.875	97	
Ethane, pentachloro-	4.881	86	
Benzene, 1,2,4-trimethyl-	4.902	98	
Benzene, (1-methylpropyl)-	5.001	98	
Benzene, 1,3-dichloro-	5.059	99	
p-Cymene	5.085	97	
Benzene, 1,4-dichloro-	5.110	99	
1,2-Dichlorobenzene-D4	5.314	78 **	■
Benzene, n-butyl-	5.321	96	
Benzene, 1,2-dichloro-	5.324	92	
Ethane, hexachloro-	5.475	97	
Propane, 1,2-dibromo-3-chloro-	5.775	98	
Benzene, nitro-	5.896	94	←
Benzene, 1,2,4-trichloro-	6.269	99	
1,3-Butadiene, 1,1,2,3,4,4-hexachloro-	6.379	91	
Naphthalene	6.412	99	
Benzene, 1,2,3-trichloro-	6.557	99	

Average LMS = 94

* Low-responding compound


** Spectral interference from overlapping peaks

Iodomethane: Good Linearity from 0.05 ppb – 25 ppb





Methane, iodo- - 8 Levels, 8 Levels Used, 8 Points, 8 Points Used, 0 QCs

$y = 0.416384 * x + 4.598609E-004$
 $R^2 = 0.99922165$
Type: Linear, Origin: Ignore, Weight: 1/x
Avg. RF RSD = 3.268834



Nitrobenzene: Linearity from 0.25 ppb – 25 ppb



Being a polar analyte, nitrobenzene is a low responder in static headspace analysis.

It does, however, exhibit good linearity.

MDLs: SIM HS/GC/MSD VOAs with H₂ Carrier and HydroInert

Name	RT Min min	Conc for MDL	MDL µg/L
Dichlorodifluoromethane	0.508	0.10	0.011
Chloromethane	0.615	0.10	0.022
Ethene, chloro-	0.698	0.05	0.008
Methane, bromo-	0.891	0.10	0.029
Ethyl Chloride	0.945	0.05	0.010
Trichloromonofluoromethane	1.067	0.05	0.008
Ethyl ether	1.198	0.05	0.017
Ethene, 1,1-dichloro-	1.288	0.05	0.006
Acetone	1.317	[cont]	
Methane, iodo-	1.350	0.05	0.006
Carbon disulfide	1.379	0.05	0.003
Allyl chloride	1.432	0.05	0.014
Methylene chloride	1.478	0.05	0.007
Acrylonitrile	1.572	[0.25]	
Ethylene, 1,2-dichloro-, trans-	1.586	0.05	0.007
Methyl tert-butyl ether	1.592	0.05	0.003
Ethane, 1,1-dichloro-	1.745	0.05	0.003
Ethylene, 1,2-dichloro-, cis-	1.966	0.05	0.007
Propane, 2,2-dichloro-	1.969	0.10	0.017
Propanenitrile	1.993	[0.25]	
2-Propenoic acid, methyl ester	2.008	0.10	0.029
Methylacrylonitrile	2.052	0.10	0.032
Methane, bromochloro-	2.059	0.10	0.019
Trichloromethane	2.086	0.05	0.011
Tetrahydrofuran	2.090	0.05	0.030
Ethane, 1,1,1-trichloro-	2.168	0.05	0.007
Butane, 1-chloro-	2.205	0.05	0.007

1-Propene, 1,1-dichloro-	2.231	0.05	0.007
Carbon Tetrachloride	2.235	0.05	0.015
Benzene	2.315	0.05	0.004
Ethane, 1,2-dichloro-	2.316	0.05	0.005
Trichloroethylene	2.577	0.05	0.006
Propane, 1,2-dichloro-	2.671	0.05	0.011
Methyl methacrylate	2.713	0.05	0.033
Methane, dibromo-	2.722	0.05	0.009
Methane, bromodichloro-	2.785	0.05	0.011
Propane, 2-nitro-	2.883	0.10	0.041
Chloromethyl cyanide	2.887	[0.25]	
1-Propene, 1,3-dichloro-, cis-	2.985	0.05	0.003
Toluene	3.145	0.05	0.003
1-Propene, 1,3-dichloro-, trans-	3.239	0.05	0.005
Methacrylic acid, ethyl ester	3.283	0.05	0.008
Ethane, 1,1,2-trichloro-	3.328	0.05	0.034
Tetrachloroethylene	3.410	0.05	0.005
Propane, 1,3-dichloro-	3.412	0.05	0.007
Methane, dibromochloro-	3.524	0.05	0.008
Ethane, 1,2-dibromo-	3.585	0.05	0.005
Benzene, chloro-	3.835	0.05	0.002
Ethane, 1,1,1,2-tetrachloro-	3.875	0.05	0.007
Ethylbenzene	3.892	0.05	0.005
m-Xylene	3.953	0.05	0.001
o-Xylene	4.164	0.05	0.004
Styrene	4.169	0.05	0.005
Methane, tribromo-	4.266	0.05	0.003
Benzene, (1-methylethyl)-	4.364	0.05	0.004
p-Bromofluorobenzene	4.446	ISTD	
Ethane, 1,1,2,2-tetrachloro-	4.521	0.05	0.006

Benzene, bromo-	4.530	0.05	0.003
Propane, 1,2,3-trichloro-	4.548	0.05	0.024
2-Butene, 1,4-dichloro-, trans-	4.555	[0.25]	
Benzene, propyl-	4.592	0.05	0.008
Benzene, 1-chloro-2-methyl-	4.638	0.05	0.006
Mesitylene	4.692	0.05	0.008
Benzene, tert-butyl-	4.876	0.05	0.004
Ethane, pentachloro-	4.881	0.05	0.009
Benzene, 1,2,4-trimethyl-	4.903	0.05	0.007
Benzene, (1-methylpropyl)-	5.001	0.05	0.004
Benzene, 1,3-dichloro-	5.060	0.05	0.003
p-Cymene	5.086	0.05	0.009
Benzene, 1,4-dichloro-	5.110	0.05	0.004
1,2-Dichlorobenzene-D4	5.313	Surr	
Benzene, n-butyl-	5.322	0.05	0.012
Benzene, 1,2-dichloro-	5.325	0.05	0.003
Ethane, hexachloro-	5.476	0.05	0.008
Propane, 1,2-dibromo-3-chloro-	5.775	0.05	0.017
Benzene, nitro-	5.896	[0.25]	
Benzene, 1,2,4-trichloro-	6.270	0.05	0.007
1,3-Butadiene, 1,1,2,3,4,4-hexachloro-	6.380	0.05	0.006
Naphthalene	6.413	0.05	0.003
Benzene, 1,2,3-trichloro-	6.558	0.05	0.006

Most MDLs are in the ppt range

Summary

Developed a source for H₂ carrier gas with:

- Decreased in-source reactions with H₂
- Improved S/N
- Improved sensitivity

Gained:

- The ability to retain existing quantitative methods, libraries, and MRM transitions or SIM ions.

Still must consider:

- Time for method development
- What are the best column dimensions for H₂ and the analysis?
- Adjust inlet/injection parameters to avoid overload



Acknowledgements

Thanks to:

Bill Mock at Pace Analytical National Center for Testing Innovation Laboratory

Agustin Pierri at Weck Labs

Marc Centner at ALS Australia

Laura Miles and team at Markes International

Thank you!

Any Questions?

