

# Optimized Workflow of POPs Analysis in Environmental and Food Matrices Using Semi-Automated Cleanup and GC-MS/MS

Ruud Addink (FMS) and Dale Walker (Agilent Technologies)



# Overview

- Introduction
- Automated Hyphenated System Configurations for different Matrices and Extractions
- Semi Automated Systems
- GC/MS/MS Data



# Challenges of POPs Sample Prep

- **Labor intensive, prone to error**
- **Compliance with regulatory procedures and accreditation (lengthy method validation)**
- **Strict QA/QC requirements**
- **Sample matrix complexity**
- **Native background and interferences (can be orders of magnitude higher than analytes)**
- **Pico/femto-gram analyses require ultra pure extract and excellent instrument sensitivity**

# Automated Sample Prep

- **Advantages of Automated Sample Prep**

- Rapid Turn Around Time: 30 to 45 Minutes for 6 Samples
- Cleaner Background Interferences: Closed Loop System
- Quality Results: Certified Pre-packaged Columns
- Green Technology: Lower solvent and power use
- QA/QC & Accreditation Requirements: Easier to Manage
- Computerized Method: Instrumentation based prep

# Sample Processing Workflow

- Analysis of various matrices for PCDD/Fs and PCBs using extraction, cleanup and concentration.
- Soxhlet extraction (typically up to 16-24 h).
- Concentration step
- Preparative multi-column chromatography involving various solvents and steps.
- Can include acid-base-neutral silica, pure acidified silica, alumina, Florisil and carbon columns. Use of H<sub>2</sub>SO<sub>4</sub> acid mixed with silica; NaOH mixed with silica Neutral Silica
- Fractionation and Concentration

# Automation

- Advantages of automated sample prep are:
- Reduced time:
  - Automated Pressurized Liquid Extraction (PLE) takes 60 min start-to-finish
  - Manual Soxhlet up to 24 h.
- Reduced cost: less labor involved, shorter turnover time per sample, less electricity use for PLE than Soxhlet.
- Reduced volume: less solvent used.





# Automating the Workflow

- Automated Pressurized Liquid Extraction (PLE) for sample extraction is fast (60 min), efficient (120 °C, 1500 psi), green (less power), reliable (long track record).
- Users' choice of fully automated or semi automated clean up system: fast (30-45 min), low solvent usage



# Automated Cleanup

- Fully automated sample load and elution.
- Load Sample Extracts in hexane directly onto the system with no Manual Pretreatment
- Easy to perform QC sample simultaneously with a Real sample.
  - 2 samples per module
- Different column configuration: silica-carbon-alumina.
- Uses no DCM, only Hexane and Toluene.
- Total Clean Up time 30-45 min.
- Low volumes




# Columns


- Silica - PCB/PBDE-free Acid, Base and Neutral silica gel column (mini, classical, classical plus, high capacity, XL).
- Alumina – PCB/PBDE-free basic alumina column.
- Carbon – PCB/PBDE-free carbon/celite column.
- Packed in disposable Teflon tubes; individually sealed in Mylar packaging; production in clean room environment.

# Column Kits


Classical  
Kit




Classical Plus  
Kit




Mini Kit



High-Capacity  
Kit



XL Capacity  
Kit




# SuperVap Concentrator



# SuperVap Evaporation

- System pre-heated to 45-60 °C.
- Samples evaporated at stable T under 8 psi nitrogen.
- 1 mL extract vial transferred to GC vial (can have direct-to-vial feature).
- Recovery standards added (nonane/dodecane).
- Extract taken to 10 uL volume with a gentle stream of nitrogen at ambient temperature.

Glass Evaporation tube



GC vial




# Using the GC-MS/MS




Agilent 7010




# GC-MS/MS Precision




# GC-MS/MS Limit of Detection



# EP-110 13C Recoveries DD/F Matrices



# EP-110 13C PCBs Recoveries Matrices




# Semi-Automated Sample Cleanup

- Simple to run, no computerized instrumentation
- Fast: 30 - 45 min
- Closed loop system to give a clean background, low level detection
- Use certified pre-packaged columns
- Green technology, only vacuum pump uses power
- Economical column kits, choice of low fat and high fat column kits
- No electronics or mechanical equipment to fail




# Characteristics of Semi-Automated System (EZPrep)






## Stage 1: to waste



## Stage 2: collect




# Attributes

- **Closed loop system:**
  - **Eliminates background contaminants**
  - **No washing needed.**
  - **Capped solvent reservoirs**
- **Optimized for solvent reduction while obtaining highest possible recoveries**
- **Uses Hexane and Toluene, no Dichloromethane**
- **Easy sample loading on top of silica column via syringe vial**
- **Columns connect easy with SNAP connections**

# Procedure


- **Order of column assembly: Acid Silica-Carbon-Alumina**
- **Conditioning all columns with hexane (waste)**
- **Load sample**
- **Elute columns with hexane and transfer all target compounds to carbon and alumina (waste)**
- **Discard silica columns and remove carbon and alumina columns**
- **Elute carbon and alumina individually with toluene (reverse)**
- **Collect carbon fraction (PCDD/F and co-planary PCBs) and alumina fraction (mono- and di-ortho PCBs)**

# $^{13}\text{C}$ PCBs Recoveries EZPrep (%)






# $^{13}\text{C}$ PCDD/F Recoveries EZPrep (%)






# Comparison of Native Data with Automated System



# 13C PCBs recoveries in soil



10 g soil in  
toluene, n=6

# Conclusions (1)

- Combination of PLE, SuperVap and EP-110/EZPrep give flexibility to (semi-) automate extraction, concentration and cleanup for a wide variety of compounds.
- **PLE** (Pressurized Liquid Extraction)
  - Reduced time:
    - Automated Pressurized Liquid Extraction (PLE) takes 60 min start-to-finish
    - Manual Soxhlet up to 24 h.
    - Up to 8 samples in Parallel
  - Reduced cost: less labor involved, shorter turnover time per sample, less electricity use for PLE than Soxhlet.
  - Reduced volume: less solvent used.
- **EP110**
  - Uses No DCM
  - Low solvent use 100-250 mL
  - Total time from sample till data between 3-4.5 h
  - Clean up step time between 30 and 45 min
  - Modular 2 to 8 sample configurations in parallel
- **PLE and EP-110 can be configured in one system**

## Conclusions (2)

- **Ezprep Semi-Automated Sample Cleanup**
  - Is a low solvent and fastest clean up system for Dioxins and PCBs (30 - 45 min)
- High sample throughput → 18 samples/hour
  - 6 samples in parallel per station
  - 3 stations fit in one hood
- System gives excellent recoveries for PCDD/F, PCB and PBDEs comparable to automated systems
- Use of certified pre-packaged columns guarantees low native background

# Agilent 7010 MS/MS Dioxin Analysis

- **Provides a complete Dioxin/PCB workflow with FMS Sample Prep Systems**
- **Easy to Use and Implement**
- **Cost Effective Alternative to High-Res Mass Spec**
  - Lower instrumental purchase cost
  - Lower cost of real estate
  - Lower operational cost from lower training requirements and lower maintenance costs

Questions?  
[r.addink@fms-inc.com](mailto:r.addink@fms-inc.com)