

1,4-Dioxane in drinking water: Increasing the sensitivity for trace-level analysis

Environmental Measurement Symposium 2022

Rachael Szafnauer¹, Jan Peter Mayser¹, Lucy Hearn¹, Nathan Shafer², Rebecca Cole¹

¹Markes International Ltd, UK. ²Markes International, Inc., USA.



Markes International

World-leaders in thermal desorption



SCHAUBURG
Analytics

MARKES
international

SepSolve
Analytical




Sample extraction
& enrichment



Thermal
desorption




Sampling
technologies



GCxGC
separation



Time-of-flight mass
spectrometry



Data mining &
chemometrics

1,4-Dioxane: What is it and where can it be found?


As a solvent
stabilizer for
TCA and other
chlorinated
solvents




As a solvent in
lacquers,
paints, pharma,
etc.



As a
contaminant in
ethoxylated
products
(surfactants,
detergents)



In cosmetics as
a stabilizer



1,4-Dioxane: Should we be worried?

PERSISTENT POLLUTANTS

New York restricts 1,4-dioxane in cleaning and personal care products

State is first in US to limit level of this persistent pollutant in consumer goods

by Cheryl Hogue
December 1, 2019

1,4-Dioxane in Cosmetics: A Manufacturing Byproduct

[Share](#) [Tweet](#) [LinkedIn](#) [Email](#) [Print](#)

The FDA has received questions about 1,4-dioxane, a contaminant that may occur in trace amounts in certain cosmetics. The following information is from responses to those questions, scientific literature, and other public sources.

POLLUTION

1,4-Dioxane: Another forever chemical plagues drinking-water utilities

Highly miscible in water, the likely carcinogen is challenging to remove

by Cheryl Hogue
November 8, 2020 | A version of this story appeared in *Volume 98, Issue 43*



CAMPAIGN for SAFE COSMETICS

A PROJECT OF BREAST CANCER PREVENTION PARTNERS

[Regulations](#)

[Chemicals of Concern](#)

[Health & Science](#)

[Reports](#)

[Safer Cosmetics Companies](#)


[Home](#) / [Get the Facts](#) / [1,4-dioxane](#)

1,4-DIOXANE

1,4-Dioxane: An emerging contaminant of concern


- The USEPA acknowledges that people may be exposed to 1,4-dioxane via drinking water, as well as from ambient air and soil
- Toxic Substances Control Act (TSCA) manages the exposure from water
- 1,4-Dioxane is one of the first 10 high-priority chemical assessments the EPA conducted under 2016 revisions to TSCA
- Classified as “likely to be carcinogenic to humans” by all routes of exposure: US EPA, 2019

1,4 Dioxane




Why a focus on drinking water?

1,4-Dioxane Occurrence in 4864 Public Water Systems Included in UCMR3



1,4-Dioxane detected in 21% of public water systems but detection rates declined over time




Exposure not solely related to groundwater-based conceptual model for 1,4-dioxane releases


- 1,4-Dioxane has been detected at many public water sites > 0.35 $\mu\text{g/L}$ (red, reference concentration, RC) and even more sites > 0.07 $\mu\text{g/L}$ (blue, min. reporting level, MRL)

Why a focus on drinking water?

1,4-Dioxane Occurrence in 4864 Public Water Systems Included in UCMR3



1,4-Dioxane detected in 21% of public water systems but detection rates declined over time




Exposure not solely related to groundwater-based conceptual model for 1,4-dioxane releases


- 1,4-Dioxane has been detected at many public water sites $> 0.35 \mu\text{g/L}$ (red, reference concentration, RC) and even more sites $> 0.07 \mu\text{g/L}$ (blue, min. reporting level, MRL)
- **Since all routes of exposure are dangerous, why the huge focus on drinking water?**

Why a focus on drinking water?

1,4-Dioxane Occurrence in 4864 Public Water Systems Included in UCMR3



1,4-Dioxane detected in 21% of public water systems but detection rates declined over time




Exposure not solely related to groundwater-based conceptual model for 1,4-dioxane releases

- 1,4-Dioxane has been detected at many public water sites $> 0.35 \mu\text{g/L}$ (red, reference concentration, RC) and even more sites $> 0.07 \mu\text{g/L}$ (blue, min. reporting level, MRL)
- Estimated 1- to 3-day half-life due to photooxidation, short-lived in the atmosphere, while this is not the case when present in water

Guidelines in the US

- No enforceable federal drinking water standard for 1,4-dioxane
- Many states have started regulating¹ (right), but these vary as do the types of reporting required
- EPA has established non-enforceable screening levels for residential water use at 0.46 µg/L
- Intended to provide technical information to state agencies and public health officials
- US EPA Method 522 developed for determination of 1,4-dioxane in drinking water by SPE



A global issue: Guidelines in Germany

Regulation: Anfrage an Bayerischen Landtag 17/16517 26.06.2017

LOQ:
0.58 µg/L


- SPME method
HS-SPME-GC-MS

LOQ:
0.025 µg/L

- German water regulation
change (LOQ in surface water)

LOQ:
0.007 µg/L

- LWV has change to new
SPE GC-MS/MS method



German body for the Registration, Evaluation, Authorisation and Restriction of Chemicals (REACH)

Challenges with existing analytical techniques

Determination of 1,4-dioxane in drinking water

- A variety of techniques used
 - Headspace solid-phase microextraction (HS–SPME)
 - Purge-and-trap (P&T)
 - Solid phase extraction (SPE)
- HS–SPME and P&T
 - ✓ Use smaller sample volumes (10s of mLs) and can be automated
 - ✗ Can suffer with limited sensitivity and carryover issues due to sample foaming
- SPE
 - ✓ Can reach much lower detection limits thanks to improved preconcentration, with GC–MS in SIM or triple quad MS
 - ✗ Manual, lengthy extraction process using hazardous solvents such as dichloromethane (DCM)
 - ✗ Potential for matrix interferences




Need a method which is:


- Short
- Automated
- Reproducible

1,4-Dioxane: A challenging compound

- Polar, leading to high water solubility
- Similar boiling point to water, 101°C
- Poor liquid/liquid extraction
- Has a 50-100x lower RRF (Relative Response Factor) than other, similar volatility compounds
- Sensitivity levels are inherently lower




1,4-Dioxane





How do we tackle this?

Using focusing trap technology and multi-step enrichment (MSE)



- Can we modify the 'old' HS-SPME method with **cryogen-free sample enrichment**?
- Recent advances of SPME
 - SPME Arrow technology provides **increased phase capacity** to sample more of each analyte

	Sorption phase surface	Sorption phase volume
	62.8 mm ²	11.8 µL
	9.4 mm ²	0.6 µL

Automated sample extraction and enrichment



Headspace-trap
(& classical headspace)



SPME (Arrow)-trap
(& classical SPME)

Automated sample extraction and enrichment




Headspace-trap
(& classical headspace)



SPME (Arrow)-trap
(& classical SPME)




HiSorb high-capacity
sorptive extraction




Tube-based thermal
desorption (TD)

How does the trap work?




Multi-step enrichment (MSE®) SPME Arrow-trap



Increasing sensitivity: Experiments with SPME fiber

MSE–SPME–trap

- Using the 'old' HS-SPME GC-MS methodology
- Increased the extraction yield using enrichment
- However:
 - 6 enrichments meant that one sample took >2h
 - Not a commercially viable option for high-throughput lab
- ★ Increase phase capacity with SPME Arrow



1,4-Dioxane peak with no. of SPME fiber enrichments (from same vial) noted


MSE–HS–SPME Arrow–trap sampling

1,4-Dioxane in drinking water



- Calibration: 0.01 – 0.5 µg/L (10-500 ppt) in water (10 mL) saturated with NaCl (2 g)
- Headspace–SPME Arrow extractions using PDMS/CWR/DVB phase
- Sampling: 3 extractions (same vial) for 10 mins each
- Real-world sample: Tap water from the Markes International Bridgend site


Calibration: 1,4-Dioxane



- 0.01 to 0.5 µg/L 1,4-dioxane (10 – 500 ppt)
- Extracted ion chromatograms (88 ion) overlaid


Going beyond regulatory limits

- Excellent linearity
 - 6-point calibration ($n = 3$ per level) $R^2 > 0.99$
- Good reproducibility
 - Relative standard deviations (RSDs) 3 – 13%
 - US EPA Method 522 indicates < 20%
- Detection of lowest calibration point 10 ng/L below both regulations
- Sample and analysis time = <1 hour (overlap mode used to shorten this further)



Tap water sample

- Collected from Markes' UK site
- Little-to-no detection of 1,4-dioxane: < 10 ng/L calibration point
- Confirmed with replicates performed ($n = 10$)



Conclusions

MSE–SPME Arrow–trap analysis



- Low limits of quantification:
 - **46x lower** than the US EPA screening level for tap water (0.46 µg/L) and **35x lower** than the assessed cancer risk level (0.35 µg/L)²
 - **2.5x lower** than the more stringent German REACH regulation LOQ of 0.025 µg/L (25 ppt)
- Sample enrichment: *more* analyte extracted for detection enabling these lower limits to be reached with excellent peak shape and sensitivity
- Fully automated and solvent-free technique with potential for use in high-throughput laboratories
 - Time to result:
>24h manual SPE (no analysis time) → <1h for MSE–SPME Arrow–trap

Thank you for your attention



Contact Markes



enquiries@markes.com



UK: +44 (0)1443 230935

USA: +1 866-483-5684 (toll-free)

Germany: +49 (0)69 6681089-10

P.R. China: +86 21 5465 1216



www.markes.com
www.markes.com.cn



@MarkesInt



https://uk.linkedin.com/company/
markes-international

