

# Innovative and Sensitive Method for the Determination of 1,4-Dioxane in Drinking Water by HS-GC/MS.

Authors: Miles Snow, Scott Lacombe and Lee Marotta

PerkinElmer, Chromatography Centre of Excellence, Bolton, ON, Canada




## Overview

- Water quality is a global concern
- 1,4 Dioxane is highly soluble in water and thus transports through the aquifer with ease
- 1,4 Dioxane is toxic and known cancer agent
- Drinking limits are very low with some locations at 1ppb or lower
- “*1,4-Dioxane is a **likely human carcinogen** and has been found in groundwater at sites throughout the United States. The physical/chemical properties and behavior of 1,4-dioxane create challenges for its characterization and treatment. It is highly mobile and does not readily biodegrade in the environment.*”




# USA 1,4-Dioxane Regulations

## 1,4-DIOXANE REGULATIONS IN DRINKING WATER



The information is current as of June 21, 2022.

- Images from <https://www.jdsupra.com/legalnews/state-by-state-regulation-of-1-4-8538403/>



# Contract Lab Reported DL

- All regulated states are less than 10ug/L, thus 8260 and 8270 are challenged
- It is expected the limits will likely continue to drop.

General Reporting Limit Guidelines						
Matrix	Method	0.2ug/L	0.4ug/L	1ug/L	10ug/L	100ug/L
GW SW	8260					
GW SW	8270					
GW SW	8260 SIM					
GW SW	8270 SIM					
GW SW	522 MOD					
DW	522					
	HSTrap					

# Methods for the Chromatographic Analysis of 1,4 Dioxane

---

- USEPA 8260 Purge and Trap
- USEPA 8270 Liquid/liquid
- USEPA 522 with SPE
- Method proposed using HSTrap

# USEPA8260 Purge and Trap



Utilizes Purge and Trap with excellent recoveries with most VOC. 1,4 dioxane is high water solubility with a very large partition coefficient thus the need for heated purge. The heated purge temperature significantly increases the water load on the trap/column thus lead to interferences by water.

Advantages-- can be used on existing 8260 equipment

# Liquid-Liquid with USEPA8270

- If care is taken to prevent loss of 1,4 dioxane during sample extraction, good numbers can be achieved. Recommended to use 1,4 dioxane-d8 as a surrogate to correct for sample loss from the sample prep step. This is incorrectly called Internal Standard correction.
- Advantages-- can be used on existing 8270 semi volatile systems with ease



## SPE with USEPA 522

- EPA 522 require the sample be passed through a SPE cartridge for sample enrichment. The SPE cartridge is washed with an elution solvent which is run with a liquid Injection on an 8270 like system
- Advantage---great enhancement in signal with reduced background
- Disadvantage---significant labour required for sample prep and use of expensive isotopically labeled standard



# 1,4 Dioxane using HSTrap as modified USEPA 8260




Less sample prep than  
8270 and 522

Better DL than 8260 and  
existing modifications

High Precision


# HSTrap – How it works Standby



# HSTrap – How it works Pressurization/Trap Load

Trap Load


HS Trap Operation



9

# HSTrap – How it works Desorb

Splitless Desorb and Hold



# Standard Preparation

General Solution Dilution									
	PPB	PPB		uL		uL	uL		PPB
Analyte	Conc of Std	Desired Conc	Dilution factor	final volume		volume of std needed	actual used		actual final conc
1,4 Dioxane	2,000,000	20000	100	10000		100.00	100		20000.00
	20000	200	100	10000		100	100		200.00
	20000	200	100	10000		100	10		20.00
	20000	50	400	10000		25	100		200.00
	20000	20	1000	10000		10	10		20.00
	20000	10	2000	10000		5	5		10.00
	200	2	100	10000		100	100		2.00
	200	0.5	400	10000		25	25		0.50
	200	0.1	2000	10000		5	5		0.10
			#DIV/0!			#DIV/0!			#DIV/0!
			#DIV/0!			#DIV/0!			#DIV/0!
			#DIV/0!			#DIV/0!			#DIV/0!
			#DIV/0!			#DIV/0!			#DIV/0!
			#DIV/0!			#DIV/0!			#DIV/0!



Note: For this study, the standards are V/V so therefore ug/L values are 3% higher

# Instrument conditions



## Chromatography Conditions

GC Parameters		HeadSpace Parameters	
Instrument	PerkinElmer Clarus 690 GC	Instrument	PerkinElmer TurboMatrix 40 HTrap
Carrier Gas	Helium	Carrier Gas	Helium
Columns	Elite-624 SIL MS 30 m 0.18 mm 1.0 u	Oven temperature	80° C
Injector Parameters		Needle temperature	110° C
Carrier flow split	1.8 mL/minute(1.0 min) ballastic to 0.9 mL/min OFF	Transfer line	120° C
Type	PSS with 1mm liner	Trap low	40° C
Temperature	180° C	Trap high cycles	325° C
Oven Parameters		Pressurization time	1
Oven initial temperature	40° C	Decay time	2.0 minutes
Oven initial hold	3 minutes	Dry purge	5 minutes
Ramp	10C/min	Trap hold	4 minutes
Final tempurature	200° C	Desorb time	0.5 minutes
Final time	0 minutes	Thermostat time	25 minutes
Oven maximum	?	GC cycle time	20 minutes
Equilibration time	?	Outlet split	OFF
MS Detector Parameters		Shaker	ON
Mode	EI+	Column pressure	33.6 psi
Electon Energy	70 eV	Vial pressure	33.6 psi
Multiplier Voltage	1700 V	Desorb pressure	33.6 psi
Emission current	100 uA		
Transfer line	200° C		
Source temp	200° C		
Solvent delay	2 minutes		
SIFI SIR	58,64,88,96 da 0.2sec dwell 0.04 second ISD		
SIFI Scan	16-100 daltons 0.1 sec/scan 0.06 sec ISD		

# 0.5ppb 1,4 Dioxane and 1,4 Dioxane-d8 extracted ion from Scan in SIFI



# 0.5ppb 1,4 Dioxane and 1,4 Dioxane-d8 SIM ions from SIFI



# 1,4 Dioxane Calibration curve LOG

20ppb

Compound 1 name: 14Dioxane


Correlation coefficient: r = 0.999806, r'2 = 0.999611

Calibration curve: 1.00269 \* x + 3.80040

Response type: External Std, Area

Curve type: Linear, Origin: Exclude, Weighting: Null, Axis trans: Log

20m 0.18mm 624 Sil MS column\_5min DP 40, 05-Nov-2021 + 15:02:00



# %RSD with External Standard



Compound 1: 14Dioxane

#	Name	ID	Type	Std Conc	RT	Response	!					
108	dioxaneCA new salt 2	Analyte			6.255	149933.7	19.57					
109	dioxaneCA new salt 2	Analyte			6.255	162534.3	21.21		all 10			
110	dioxaneCA new salt 2	Standard		20	6.255	152417.2	19.89		average	153123		
111	dioxaneCA new salt 2	Standard		20	6.255	154108.6	20.11		stdev	4073.909		
112	dioxaneCA new salt 2	Analyte			6.255	153615.6	20.05		%RSD	2.660547		
113	dioxaneCA new salt 2	Analyte			6.255	151696.1	19.8					
114	dioxaneCA new salt 2	Analyte			6.274	153057.8	19.97					
115	dioxaneCA new salt 2	Analyte			6.274	147620.5	19.26					

Compound 2: 14Dioxane-d8

#	Name	ID	Type	Std Conc	RT	Response	!					
108	dioxaneCA new salt 2	Analyte			6.199	152423.8	19.41					
109	dioxaneCA new salt 2	Analyte			6.199	165174.1	21.04		all 10			
110	dioxaneCA new salt 2	Standard		20	6.199	155916.1	19.86		average	156431.1		
111	dioxaneCA new salt 2	Standard		20	6.199	158152.2	20.14		stdev	4061.379		
112	dioxaneCA new salt 2	Analyte			6.199	156876.9	19.98		%RSD	2.596273		
113	dioxaneCA new salt 2	Analyte			6.199	155899.9	19.86					
114	dioxaneCA new salt 2	Analyte			6.199	156552.8	19.94					
115	dioxaneCA new salt 2	Analyte			6.199	150453.1	19.16					



Excellent repeatability



## Recalculate using d-8 as Internal Standard

- As expected, the use of internal standard greatly improves the % RSD to less than 0.4%.

1,4 Dioxane	d8	
149933.656	152423.8	0.983663
162534.266	165174.1	0.984018
152417.156	155916.1	0.977559
154108.594	158152.2	0.974432
153615.609	156876.9	0.979211
151696.063	155899.9	0.973035
153057.797	156552.8	0.977676
147620.516	150453.1	0.981173
average	0.978846	
stdev	0.003738	
%RSD	0.381864	



## Conclusions

---

- Sample run times are faster than most traditional methods
- Less sample prep vs 8270/522
- External standard provides excellent results without extra cost associated with the use of an expensive d8 IS.
- The HSTrap method provides excellent detection limits as compared to the current environmental methods