Thermo Fisher S C I E N T I F I C

Reproducible trace analysis of PCBs in environmental matrices using triple quadrupole GC-MS/MS

Andy Fornadel, PhD Product Marketing Manager – Americas

4th August 2022

The world leader in serving science



Content



Introduction to analysis of PCBs

2 GC-MS/MS for targeted analysis of PCBs

3 Example of PCB analysis using GC-MS/MS



Introduction

Polychlorinated biphenyls (PCBs) are a group of industrial organic chemicals characterized by:

Thermo Fishei

- non-flammability
- chemical stability
- high boiling point
- electrical insulating properties
- persistence in environment

Many of these properties make PCBs historically useful in electrical insulators, hydraulic equipment, paints, rubbers, and other industrial applications

PCB congeners

Currently 209 known PCBs congeners that can be divided into two groups according to their structural and toxicological characteristics: Thermo Fi

- non-dioxin like PCBs (non-DL-PCB):
 - majority of the PCB congeners
 - lower degree of toxicity

- dioxin-like PCBs (DL-PCBs):
 - the 12 most toxic congeners
 - classified as POPs
 - regulated under the Stockholm Convention for POPs since 2001

PCB Regulation

 Because of toxicity/mutagenicity and persistence in the environment, particularly in soils and waters, PCBs are highly regulated in many countries Thermo Fi

- PCB clean-up and remediation can prove challenging due to their tenacity
- Although the goal for PCBs in water is zero in the US, the maximum contaminant level is set at 0.5 ppb.

GC-MS/MS for analysis of PCBs



- EPA developed Method 1668 using GC-HRMS
- Advances in GC-MS/MS allows for reliable application to PCB analysis.
 - Different extraction and clean-up processes
 - More sensitive/selective mass analyzers and column phases
 - Re-evaluation of performance metrics
- Important parameters to consider:
 - Working range, particularly on the low-end
 - Linearity, reproducibility, robustness
 - Applicability for routine analysis.

Traditional Challenges

Challenges



 Complicated sample preparation



 Sufficient chromatographic separation between PCBs needed to avoid isobaric interferences



Long run times for sufficient separations. Low throughput and high costs.

Analytical Parameters

AI/AS 1610 autosampler parameters				
Injection type	Standard			
Sample mode	Standard			
Fill strokes	10			
Sample depth	Bottom			
Injection mode	Fast			
Pre-injection delay time (s)	0			
Post-injection delay time (s)	0			
Pre-injection wash cycles	0			
Post-injection wash cycles	4			
Post-injection solvent wash volume (µL)	6.0			
Sample wash cycles	1			
Sample wash volume (µL)	1.0			
Injection volume (µL)	1.0			

.....

SSL parameters					
Injection temperature (°C)	280				
Liner	Thermo Scientifc [™] LinerGOLD [™] splitless/split liner single taper with wool (P/N 453A1925-UI)				
Inlet module and mode	SSL, splitless				
Split flow (mL/min)	75				
Splitless time (min)	1.2				
Septum purge flow (mL/min)	5, constant				
Carrier gas, flow (mL/min)	He, 1.2				

TRACE 1610 GC parameters					
Oven temperature program					
Temperature (°C)	90				
Hold time (min)	1.00				
Rate (°C/min)	25				
Temperature 2 (°C)	270				
Rate 2 (°C/min)	4				
Temperature 3 (°C)	330				
Hold time (min)	2				
GC run time (min)	25.20				
Column					
TRACE TR-PCB 8 MS	50 m, 0.25 mm, 0.25 μm (P/N 26AJ148P)				

TSQ 9610 mass spectrometer parameters					
Transfer line temperature (°C)	280				
lon source type and temperature (°C)	NeverVent AEI, 300				
lonization type	El				
Emission current (µA)	50				
Aquisition mode	timed-SRM				
Q1 and Q3 resolution	Normal (0.7 amu)				
Tuning parameters	AEI SmartTune				
Collision gas and pressure (psi)	Argon at 70				
Detector gain	X 7				

Sample preparation for environmental matrices



Aqueous samples (1L) containing less than 1% solid

Separatory funnel extraction

Clean-up with acid and/or base silica gel and Alumina chromatography

Concentration of extract to 10-50 μ L prior to GC injection

Solid, semi-solid, and multi-phase samples (10 g dry weight)

Soxhlet extraction

Clean-up with acid and/or base silica gel and Alumina chromatography

Concentration of extract to 10-50 μ L prior to GC injection

t-SRM for improved selectivity

Compound	Retention	Precursor Ion	Product Ion	Collision Energy	Compound	Retention	Precursor Ion	Product Ion	Collision Energy
	time (min)	(Da)	(Da)	(V)	DOD 440	time (min)	(Da)	(Da)	(V)
PCB-1	8.5	188.04	153.04	22	PCB-118	14.05	323.90	253.95	22
PCB-1	8.5	190.04	153.04	22	PCB-118	14.05	325.90	255.95	22
PCB-1L13C	8.5	200.08	165.10	22	PCB-118L13C	14.05	335.92	265.99	22
PCB-1L13C	8.5	202.08	165.10	22	PCB-118L13C	14.05	337.92	267.99	22
PCB-3	9.03	188.04	153.04	22	PCB-114	14.28	323.90	253.95	22
PCB-3	9.03	190.04	153.04	22	PCB-114	14.28	325.90	255.95	22
PCB-3L13C	9.03	200.08	165.10	22	PCB-114L13C	14.28	335.92	265.99	22
PCB-3L13C	9.03	202.08	165.10	22	PCB-114L13C	14.28	337.92	267.99	22
PCB-4	9.19	222.00	152.06	22	PCB-105	14.75	323.90	253.95	22
PCB-4	9.19	224.00	152.06	22	PCB-105	14.75	325.90	255.95	22
PCB-4L13C	9.19	234.04	164.10	22	PCB-105L13C	14.75	335.92	265.99	22
PCB-4L13C	9.19	236.04	164.10	22	PCB-105L13C	14.75	337.92	267.99	22
PCB-19	9.88	255.96	186.02	22	PCB-126	15.64	323.90	253.95	22
PCB-19	9.88	257.96	186.02	22	PCB-126	15.64	325.90	255.95	22
PCB-19L13C	9.88	268.00	198.02	22	PCB-126L13C	15.64	335.92	265.99	22
PCB-19L13C	9.88	270.00	198.02	22	PCB-126L13C	15.64	337.92	267.99	22
PCB-15	10.36	222.00	152.06	22	PCB-202	15.68	427.80	357.80	25
PCB-15	10.36	224.00	152.06	22	PCB-202	15.68	429.80	357.80	25
PCB-15L13C	10.36	234.04	164.10	22	PCB-202L13C	15.68	439.80	369.90	25
PCB-15L13C	10.36	236.04	164.10	22	PCB-202L13C	15.68	441.80	369.90	25
PCB-54	10.54	289.92	219.98	22	PCB-167	16	357.80	287.90	24
PCB-54	10.54	291 92	219.98	22	PCB-167	16	359.80	289.95	24
PCB-54I 13C	10.54	301.96	232.02	22	PCB-167L13C	16	369.90	299.51	24
PCB-54L13C	10.54	303.96	232.02	22	PCB-167L13C	16	371.90	301.95	24
PCB-104	11 34	323.90	253.95	22	PCB-156	16.6	357.80	287.90	24
PCB-104	11.34	325.00	255.95	22	PCB-156	16.6	359.80	289.95	24
PCB-104 13C	11.34	335.02	265.00	22	PCB-156I 13C	16.6	369.90	299.51	24
PCB-104L13C	11.34	337.02	267.99	22	PCB-156L13C	16.6	371.90	301.95	24
PCB-104L130	11.04	255.06	196.02	22	PCB-157	16.82	357.80	287.90	24
	11.04	257.06	196.02	22	PCB-157	16.82	359.80	289.95	24
PCD-37	11.04	257.50	100.02	22	PCB-157L12C	16.82	360.00	200.55	24
PCB-37L13C	11.04	200.00	190.02	22	PCB-157L13C	16.92	371.00	201.05	24
PCB-37L13C	11.04	270.00	198.02	22	DCB 160	17.02	257.90	297.00	24
PCB-155	12.2	357.80	287.90	24	PCD-109	17.00	357.60	207.90	24
PCB-155	12.2	359.80	289.95	24	PCB-109	17.00	309.60	209.90	24
PCB-155L13C	12.2	369.90	299.51	24	PCB-169L13C	17.00	369.90	299.51	24
PCB-155L13C	12.2	371.90	301.95	24	PCB-169L13C	17.00	371.90	301.95	24
PCB-101L13C	12.59	335.92	265.99	22	PCB-208	10.10	461.70	391.80	25
PCB-101L13C	12.59	337.92	267.99	22	PCB-208	18.18	463.70	393.80	25
PCB-111L13C	13.13	335.92	265.99	22	PCB-208L13C	18.18	473.80	403.80	25
PCB-111L13C	13.13	337.92	267.99	22	PCB-208L13C	18.18	475.80	405.80	25
PCB-77	13.42	289.92	219.98	22	PCB-189	18.92	391.80	321.90	25
PCB-77	13.42	291.92	219.98	22	PCB-189	18.92	393.80	323.90	25
PCB-77L13C	13.42	301.96	232.02	22	PCB-189L13C	18.92	403.80	333.90	25
PCB-77L13C	13.42	303.96	232.02	22	PCB-189L13C	18.92	405.80	335.90	25
PCB-81	13.69	289.92	219.98	22	PCB-205	20.04	427.80	357.80	25
PCB-81	13.69	291.92	219.98	22	PCB-205	20.04	429.80	357.80	25
PCB-81L13C	13.69	301.96	232.02	22	PCB-205L13C	20.04	439.80	369.90	25
PCB-81L13C	13.69	303.96	232.02	22	PCB-205L13C	20.04	441.80	369.90	25
PCB-188	13.8	391.80	321.90	25	PCB-206	20.6	461.70	391.80	25
PCB-188	13.8	393.80	323.90	25	PCB-206	20.6	463.70	393.80	25
PCB-188L13C	13.8	403.80	333.90	25	PCB-206L13C	20.6	473.80	403.80	25
PCB-188L13C	13.8	405.80	335.90	25	PCB-206L13C	20.6	475.80	405.80	25
PCB-123	13.97	323.90	253.95	22	PCB-209	21.03	495.70	425.80	25
PCB-123	13.97	325.90	255.95	22	PCB-209	21.03	497.70	427.80	25
PCB-123L13C	13.97	335.92	265.99	22	PCB-209L13C	21.03	507.70	437.80	25
PCB-123L13C	13.97	337.92	267.99	22	PCB-209L13C	21.03	509.70	439.80	25



Time [min]

Chromatographic resolution of critical pairs



Thermo Fisher

SCIEN



Linearity for accurate quantitative analyis

Compound	Retention Time (min)	Coefficient of determination (R ²)	AvCF %RSD
PCB-1	8.49	0.99993	1.5
PCB-3	9.03	0.99996	1.2
PCB-4	9.19	0.99987	2.1
PCB-19	9.88	0.99992	1.6
PCB-15	10.35	0.99996	1.2
PCB-54	10.53	0.99987	2.1
PCB-104	11.34	0.99987	2.1
PCB-37	11.83	0.99989	1.9
PCB-155	12.21	0.99168	17.3
PCB-77	13.40	0.99956	3.9
PCB-81	13.68	0.99925	5.0
PCB-188	13.80	0.99974	2.9
PCB-123	13.96	0.99967	3.3
PCB-118	14.04	0.99972	3.1
PCB-114	14.27	0.99962	3.6
PCB-105	14.74	0.99969	3.2
PCB-126	15.62	0.99990	1.9
PCB-202	15.67	0.99977	2.7
PCB-167	15.93	0.99707	10.0
PCB-156	16.62	0.99473	13.8
PCB-157	16.80	0.99437	14.3
PCB-169	17.84	0.99434	14.4
PCB-208	18.16	0.99991	1.8
PCB-189	18.90	0.99979	2.6
PCB-205	20.00	0.99960	3.6
PCB-206	20.57	0.99989	1.9
PCB-209	21.02	0.99978	2.7



Assessment of Sensitivity

LOQ = 0.05 ng/mL(i) ion ratios within \pm 15% the expected values (ii) absolute peak area RSD <15% (iii) the response factor (RF) RSD <15%

			I OD-5	T. 6	т.5	55
PCB-1	PCB-54	PCB-118	PCB-4	5.4	5.2	64
			PCB-19	2.8	2.9	62
Creen Grand Concern pp (BOO			PCB-15	2.9	2.3	64
			PCB-54	1.8	2.8	64
			PCB-104	7.3	6.3	95
			PCB-37	2.4	2.4	61
			PCB-155	13.5	6.5	124
			PCB-77	3.0	1.9	63
			PCB-81	3.1	2.9	63
			PCB-188	6.4	6.2	157
			PCB-123	7.5	7.5	102
PCB-157	PCB-205	PCB-209	PCB-118	4.1	3.5	101
	1 00 200	. 02 200	PCB-114	4.9	4.1	96
			PCB-105	4.8	5.7	95
			PCB-126	13.0	11.6	90
			PCB-202	5.4	6.6	64
			PCB-167	4.8	5.3	123
			PCB-156	10.5	8.2	123
	1. Sec. 1. Sec		PCB-157	9.1	5.4	122
			PCB-169	5.9	10.4	129
	1. Sec.		PCB-208	9.9	6.6	98
			PCB-189	11.2	10.9	154
			PCB-205	8.0	6.4	64
	Time [min]	 A state of the sta	PCB-206	8.8	6.2	97
			PCB-209	8.8	6.6	110
andy.fornadel@thermofisher.com 4-August-202	2	**********				

Thermo Fisher SCIENT

Calculated

IDL (fg OC)

Average

Measured

IR (n=10)

Absolute

Peak Area

%RSD (n=10)

2.2

Compound

PCB-1

DCB-3

....

RF

%RSD

2.4

Expected

IR

PL.B-I

g.0ę0

0e0

Reproducibility

	OC absolute neak area across the sequence	Compound	QC absolute peak area %RSD
	do absolute peak area across the sequence	PCB-1	6.9
12.0	0	PCB-3	7.0
		PCB-4	7.7
		PCB-19	7.2
10.0	0	PCB-15	8.0
		PCB-54	7.3
		PCB-104	7.6
		PCB-37	8.3
8.0		PCB-155	6.6
		PCB-77	8.5
\cap		PCB-81	8.7
S 6.0		PCB-188	7.2
₩		PCB-123	9.0
		PCB-118	8.7
		PCB-114	8.5
4.0		PCB-105	8.9
		PCB-126	8.5
		PCB-202	7.9
2.0		PCB-167	8.6
		PCB-156	9.0
		PCB-157	8.7
		PCB-169	8.7
0.0		PCB-208	8.3
		PCB-189	9.5
•	$\langle 0, \langle 0, \langle 0, \langle 0, \rangle \rangle \rangle \langle 0, \langle 0, \langle 0,$	PCB-205	10.4
		PCB-206	9.3
		PCB-209	9.3

........

.

. .

.....

•

.

Reproducibility



....

Conclusions

- The TRACE TR-PCB column ensured chromatographic separation of the target analytes in about 21 minutes with calculated resolution of the critical pair PCB-123 / PBDE-118 of 3%.
- The column thin film phase, high thermal stability, and low column bleed ensured elution of the high boiling point PCBs (e.g., PCB-209) with improved peak shapes.
- We demonstrated linearity over a concentration range of 0.10 to 2,000 ng/mL with coefficient of determination of R² >0.990 and AvCF %RSDs <20.
- IDLs ranged from 3 fg to 19 fg OC (corresponding to 0.15 pg/L to 0.95 pg/L in water samples and to 0.015 to 0.095 ng/kg in soil samples) and LOQ set at 0.05 ng/mL.

......

00001

Extended robustness demonstrated over 100 injections without need for maintenance or re-tuning.



7 andy.fornadel@thermofisher.com | 4-August-2022