

Strategies for Ultra Low-level Detection and Quantification of Short- and Long- Chain Per- and Polyfluoroalkyl Substances (PFAS) by Direct Injection LC-MS/MS


Ed George

Market Manager

LSMS Environmental, Food, and Beverage

August 4, 2022

NEMC 2022 Crystal City, VA



 The world leader in serving science

Executive Summary

- Part 1: A single LC-MS/MS method was developed to measure various PFAS classes, including perfluorinated ether acids in surface by direct injection
 - A Thermo Scientific TSQ Altis Plus MS and Vanquish Flex UHPLC system provided excellent quantitative performance for measuring PFAS in surface water samples via 25 µL direct injection LC-MS/MS down to low single digit ng/L concentrations.
 - Nearly all (40 of 43) targeted PFAS compounds had LODs at or below 1 ng/L for neat PFAS solutions on the TSQ Altis Plus. Lower LODs in several cases (e.g., PFOA) were limited by contamination in solvent blanks.
 - PFAS calibration curves from 0.5 – 1000 ng/L, using internal calibration, yielded linear regression calibrations with $r^2 > 0.995$.
 - Nearly all targeted PFAS compounds were able to be measured in spiked surface water samples at 1-2 ng/L (after correcting for 1:1 dilution with MeOH).
- Part 2: A two-method strategy for direct injection to extend testing for both short and longer chain PFAS- dealing with aqueous samples and samples containing a high concentration of organic solvent (>95%)
 - LC system set-up with large sample loop and solvent sandwich injections for short and long chain chromatographic optimization
 - LC-Vanquish Duo configurations available for obtaining high throughput and improved sensitivity

Experimental – Liquid Chromatography


Thermo Scientific Vanquish Flex Binary UHPLC System

- PFAS Delay Column: 3.0 x 50 mm, 1.9 um Hypersil GOLD (Thermo Scientific)
- Analytical Column: 2.1 x 100 mm, 2.2 um Acclaim RSLC C18 (Thermo Scientific)
- Column Temp: 40 C
- Mobile Phase A: H₂O with 2% MeOH + 2 mM Am. Acetate + 0.1% HOAc
- Mobile Phase B: MeOH with 2% H₂O + 2 mM Am. Acetate + 0.1% HOAc
- Gradient: see table
- Injection Volume: 25 uL
- Sample Temp: 20 C

No	Time	Flow [mL/min]	%B	Curve	
1	0.000	Run			
2	0.000	0.400	20.0	5	
3	1.000	0.400	50.0	5	
4	15.000	0.400	100.0	5	
5	17.000	0.400	100.0	5	
6	17.200	0.400	20.0	5	
7	22.500	0.400	20.0	5	
8	New Row				
9	22.500	Stop Run			



Experimental – Liquid Chromatography

PFAS Kit Retrofit & Delay Column setup



Experimental – Liquid Chromatography

Strong Solvent Loop added in autosampler: Large volume direct injections



Peak shape for PFBA is maintained up to 25 uL inj. of 50% MeOH solution

Experimental – Mass Spectrometry

Thermo Scientific TSQ Altis Plus

- Ionization Mode: HESI, Negative ion mode
 - Source Parameters: see figure at right
- MS Acquisition Mode: Timed Selected Reaction Monitoring (SRM)
- SRM Cycle Time: 0.4 s
- Quad Isolation Width: Q1, Q3 = Unit (0.7 Da FWHM)

Ion Source Properties	
Ion Source Type	H-ESI
Spray Voltage	Static
Positive Ion (V)	3500
Negative Ion (V)	1000
Current LC Flow ($\mu\text{L}/\text{min}$)	0
Get Defaults	
Sheath Gas (Arb)	50
Aux Gas (Arb)	10
Sweep Gas (Arb)	1.5
Ion Transfer Tube Temp ($^{\circ}\text{C}$)	175
Vaporizer Temp ($^{\circ}\text{C}$)	250

Experimental – Sample Prep (1)

PFAS Calibration Standard Solutions

- Standard solutions were provided by Wellington Labs. Standards were stored at 4 C until needed.
- PFAS calibration standards were prepared from Wellington PFAC-30PAR stock solution. Final calibration standard solutions were prepared over a concentration range 0.5-1000 ng/L in 50% MeOH.
- The calibration standards were spiked with isotopically-labeled standards (MPFAC-HIF-ES-Wellington Labs) to a final concentration of 50-400 ng/L.
- All calibration solutions were prepared in amber glass autosampler vials with polypropylene caps to prohibit PFAS contamination.
- Final PFAS calibration standards were analyzed by LC-MS/MS shortly after preparation to limit sample adsorption losses.

Experimental – Sample Prep (2)

Surface Water Samples

- Sample volumes of surface water samples were determined by subtracting the mass of provided empty polypropylene centrifuge tubes from mass of each water containing tube. (Based on density = 1.0 g/mL H₂O).
- Prepared 25 mL MeOH solution with 100-800 ng/L isotopically-labeled solution. [Note, this is 2X the concentrations used for calibration standard solutions.]
- This methanolic solution was added to polypropylene tubes at an equal volume to the surface water samples.
- After thoroughly vortexing surface water solutions, ~0.5 mL was transferred to amber glass autosampler vials with polypropylene caps for LC-MS/MS analyses.

PFAS Method


- Example PFAS data via direct injection on TSQ Altis Plus
 - Solvent Blank PFAS contamination
 - Effect of Ion Transfer Tube temperature on PFAS




Full-scan MS & MS/MS (m/z 263, PFPeA isobar)

Background ions detected: ESI(-) at 90% Mobile Phase B

FullMS_90percent_MPB #28 RT: 0.29 AV: 1 NL: 3.36E6
T: - c ESI Q1MS [150.000-800.000]




MSMS_263 #58-72 RT: 0.25-0.31 AV: 15 NL: 2.28E5
T: - c ESI Full ms2 263.000 [10.000-270.000]



Background has ethylene glycol signature; MS2 of m/z 263 shows loss of CO_2

PFAS Contamination – Solvent Blanks



50% MeOH



For some PFAS, a “null” response was not achieved (e.g., PFBA, PFHxA)

Effect of TSQ Ion Transfer Tube Temperature

100 ng/L PFAS: Perfluorinated ether acids, sulfonates and carboxylic acids



All PFAS classes show comparable or improved response at lower ITT temperature


Quantitation Results – PFAS Stds.

- Examples of LC-MS/MS data for direct injection of PFAS standards
 - Chromatogram displays near LODs
 - Calibration curves
 - Final table of estimated LODs



Example Chromatograms near LOD


PFBA (in 50% MeOH, 25 uL inj.)



Est. LOD PFBA is 2 ng/L (50 fg on-column); limited by blank contamination

Calibration Curve: PFBA


2 – 1000 ng/L, Linear, 1/x weighting



PFBA (PAR30) has linear regression $R^2 = 0.9986$ over 3 decade dynamic range

Example Chromatograms near LOD


4:2 FTS (in 50% MeOH, 25 uL inj.)



Est. LOD 4:2FTS is 0.25 ng/L (6.25 fg on-column)

Calibration Curve: 4:2FTS


0.5 – 1000 ng/L, Linear, 1/x weighting



4:2FTS (PAR30) has linear regression $R^2 = 0.9981$ over 3.5 decade dynamic range

Example Chromatograms near LOD


HFPO-DA (in 50% MeOH, 25 uL inj.)



Est. LOD HFPO-DA is 0.25 ng/L (6.25 fg on-column)

Calibration Curve: HFPO-DA

0.5 – 1000 ng/L, Linear, 1/x weighting




HFPO-DA (PAR30) has linear regression $R^2 = 0.9985$ over 3.5 decade dynamic range

Example Chromatograms near LOD


PFOA (in 50% MeOH, 25 uL inj.)

Extraction
Blank




NL: 1.45E3
m/z= 368.50-369.50 F: - c ESI
SRM ms2 413.000
[168.999-169.001,
368.999-369.001] MS ICIS
extraction_std_101

0.5 ng/L



NL: 1.45E3
m/z= 368.50-369.50 F: - c ESI
SRM ms2 413.000
[168.999-169.001,
368.999-369.001] MS ICIS
PAR30_0pt5ppt_102

1 ng/L




NL: 1.45E3
m/z= 368.50-369.50 F: - c ESI
SRM ms2 413.000
[168.999-169.001,
368.999-369.001] MS ICIS
par30_1ppt_103

Est. LOD PFOA is 0.25 ng/L (6.25 fg on-column); limited by blank contamination

Calibration Curve: PFOA

0.5 – 1000 ng/L, Linear, 1/x weighting




PFOA (PAR30) has linear regression $R^2 = 0.9985$ over 3.5 decade dynamic range

Example Chromatograms near LOD


9CI-PF3ONS, (in 50% MeOH, 25 uL inj.)

Extraction
Blank




NL: 6.06E2
m/z= 350.45-351.45 F: - c
ESI SRM ms2 530.900
[350.949-350.951] MS
extraction_std_101

0.5 ng/L



NL: 6.06E2
m/z= 350.45-351.45 F: - c
ESI SRM ms2 530.900
[350.949-350.951] MS
ICIS PAR30_0pt5ppt_102

1 ng/L




NL: 6.06E2
m/z= 350.45-351.45 F: - c
ESI SRM ms2 530.900
[350.949-350.951] MS
ICIS par30_1ppt_103

Time (min)

Est. LOD 9CI-PF3ONS is 0.25 ng/L (6.25 fg on-column)

Calibration Curve: 9CI-PF3ONS


0.5 – 1000 ng/L, Linear, 1/x weighting



9CI-PF3ONS (PAR30) has linear regression $R^2 = 0.9972$ over 3.5 decade dynamic range

Example Chromatograms near LOD


PFTeDA (in 50% MeOH, 25 uL inj.)



Est. LOD PFTeDA is 0.25 ng/L (6.25 fg on-column); limited by I.S. contribution

Calibration Curve: PFTeDA

0.5 – 1000 ng/L, Linear, 1/x weighting



PFTeDA (PAR30) has linear regression $R^2 = 0.9993$ over 3.5 decade dynamic range

Table of PFAS Compounds – LODs

Analyte	Acronym	LOD (ng/L)	LOD (fg OC)
N-ethyl perfluorooctanesulfonamidoacetic acid	NEtFOSAA	1	25
N-methyl perfluorooctanesulfonamidoacetic acid	NMeFOSAA	2	50
Perfluorobutanoic acid	PFBA	2 (est.)	50
Perfluorobutanesulfonic acid	PFBS	0.5	12.5
1H,1H,2H,2H-Perfluorodecane sulfonic acid	8:2FTS	1	25
Perfluorodecane sulfonic acid	PFDS	0.5	12.5
Perfluorodecanoic acid	PFDA	0.25 (est.)	6.25 (est.)
Perfluorododecanoic acid	PFDoA	0.25 (est.)	6.25 (est.)
Perfluoroheptane sulfonic acid	PFHpS	0.5	12.5
Perfluoroheptanoic acid	PFHpA	0.5	12.5
1H,1H,2H,2H-Perfluorohexane sulfonic acid	4:2FTS	0.25 (est.)	6.25 (est.)
Perfluorohexanesulfonic acid	PFHxS	0.5	12.5
Perfluorohexanoic acid	PFHxA	1	25
Perfluorononane sulfonic acid	PFNS	0.5	12.5
Perfluorononanoic acid	PFNA	0.25 (est.)	6.25 (est.)
1H,1H,2H,2H-Perfluorooctane sulfonic acid	6:2FTS	0.5	12.5
Perfluoro-1-butanesulfonamide	FBSA	0.25 (est.)	6.25 (est.)
Perfluoro-1-hexanesulfonamide	FHxSA	0.25 (est.)	6.25 (est.)
Perfluoro-1-octanesulfonamide	FOSA	0.5	12.5
Perfluorooctanesulfonic acid	PFOS	0.5	12.5
Perfluorooctanoic acid	PFOA	0.25 (est.)	6.25 (est.)
Perfluoropentanoic acid	PPPeA	5	125
Perfluoropentanesulfonic acid	PPPeS	0.5	12.5
Perfluorotetradecanoic acid	PFTeDA	0.25 (est.)	6.25 (est.)
Perfluorotridecanoic acid	PFTrDA	0.25 (est.)	6.25 (est.)
Perfluoroundecanoic acid	PFUnA	0.25 (est.)	6.25 (est.)
Hexafluoropropylene oxide dimer acid	HFPO-DA	0.25 (est.)	6.25 (est.)
11-chloroeicosfluoro-3-oxaundecane-1-sulfonic acid	11Cl-PF2OUdS	0.25 (est.)	6.25 (est.)
9-chlorohexadecafluoro-3-oxanone-1-sulfonic acid	9Cl-PF3ONS	0.25 (est.)	6.25 (est.)
4,8-dioxa-3H-perfluorononanoic acid	ADONA	N.A.	N.A.

- SRM data from 25 uL direct injection of PFAS standards in 50% MeOH
- LODs are defined as lowest concentration where the compound is observed or expected to be observed based on response difference between the 0.5 ng/L standard and solvent blank.
- Concentrations for PFAS Sulfonates are not corrected for salt form or for presence of branched isomers
- ADONA was not measured, as the Timed SRM window was incorrect


Results

Measurement of PFAS in Surface Water samples



Example Chromatograms – Spiked Surface Water Samples


PFMOAA (in 50% MeOH, 25 uL inj.)



PFMOAA was observed in 2 of 3 surface water samples

Example Chromatograms – Surface Water Samples


PFO2HxA (in 50% MeOH, 25 uL inj.)



PFO2HxA is observed in 3 of 3 surface water samples

Example Chromatograms – Surface Water Samples


HFPO-DA (in 50% MeOH, 25 uL inj.)



HFPO-DA is observed in all surface water samples, including blank

Example Chromatograms – Surface Water Samples

PFNA (in 50% MeOH, 25 uL inj.)



PFNA is observed in all surface water samples, including blank

Conclusions

- Thermo Scientific TSQ Altis Plus MS system provides excellent quantitative performance for measuring PFAS in surface water samples via 25 µL direct injection LC-MS/MS down to low single and or sub- ng/L concentrations.
- A Single LC-MS/MS method was used to measure wide variety PFAS classes, including perfluorinated ether acids.
- Nearly all (40 of 43) targeted PFAS compounds have LODs at or below 1 ng/L for neat PFAS solutions on the TSQ Altis Plus. Lower LODs in several cases (e.g., PFOA) were limited by contamination in solvent blanks.
- PFAS calibration curves from 0.5 – 1000 ng/L, using internal calibration, yielded linear regression calibrations with $r^2 > 0.995$.
- Surface water samples were analyzed by adding equal volume of methanol containing isotopically-labeled PFAS standards. Nearly all targeted PFAS compounds were able to be measured in spiked surface water samples at 1-2 ng/L levels (after correcting for 1:1 dilution).

Addressing the Main PFAS Analytical Challenges

How do I optimize LC injection volumes to obtain highest MS sensitivity ?



- **Heterogenous group of alkyl compounds:**

- Short-chain (C4-C8)
- Long-chain (C8-C18)

→ solubility, adsorption, LC retention



- **Ubiquitous occurrence:**

- Nature, clothes
- Laboratory equipment
- (U)HPLC system

→ consumable selection, blank controls



- **Guidelines:**

- Sample preparation
- Sample solvents
- Limit of quantitation

→ different flavours of PFAS analysis

An Alternate Liquid Chromatography Set-Up for Large Volume Injections

ThermoFisher
SCIENTIFIC

Vanquish Core HPLC system setup



Custom Injection Program

- Sample loop: 1000 μ L
- Sandwich injection
- In needle mixing

Strong Solvent Loop


- Volume: 46.2 μ L

Capillary

- Sampler-column
- Viper MP35N

Delay Column

- Accucore™ aQ, 3.0 x 50 mm, 2.6 μ m)



LC Method

LC System: **Vanquish Core Binary HPLC System**

Column: Acclaim™ Polar Advantage, 2.1 mm x 100 mm, 2.2 µm

Eluent A: Water, 5 mM ammonium formate, 0.1% formic acid

Eluent B: Methanol, 5 mM ammonium formate, 0.1% formic acid

Flow rate: 400 µL/min

Injection volume: **Custom Injection Program (user defined)**

No	Time	Flow	%B	Curve	
1	-2.000	Equilibration			
2	-2.000	0.400	30.0	5	
3	<i>New Row</i>				
4	0.000	Run			
5	0.000	0.400	30.0	5	
6	0.100	0.400	30.0	5	
7	1.100	0.400	50.0	5	
8	15.100	0.400	100.0	5	
9	17.600	0.400	100.0	5	
10	17.800	0.400	30.0	5	
11	22.100	0.400	30.0	5	
12	<i>New Row</i>				
13	22.100	Stop Run			



Custom Injection Program- High Solvent Content


Vanquish Split SamplerCT

Temperature: 5 °C

Wash solvent: methanol/water (95/5; v/v)


Program: **Solvent sandwich injection- 3 x 25 µL sample**

General Settings			User Defined Program	Method Transfer	Temperature Control
<input checked="" type="radio"/> Replace normal injection					
<input type="radio"/> Normal injection with liquid handling					
No	Command	Parameters			
1	UDP_PreparesLiquidHandling	Volume=250 [µL]			
2	UDP_NeedleWash	Duration=5 [s]			
3	UDP_Draw	Position=SR:1, Volume=System.Injection.CustomVariables.UDP_vol_01, Speed=10 [µL/s], NeedleHeight=2000 [µm]			
4	UDP_NeedleWash	Duration=5 [s]			
5	UDP_Draw				
6	UDP_NeedleWash	Duration=5 [s]			
7	UDP_Draw	Position=SR:1, Volume=System.Injection.CustomVariables.UDP_vol_02, Speed=10 [µL/s], NeedleHeight=2000 [µm]			
8	UDP_NeedleWash	Duration=5 [s]			
9	UDP_Draw				
10	UDP_NeedleWash	Duration=5 [s]			
11	UDP_Draw	Position=SR:1, Volume=System.Injection.CustomVariables.UDP_vol_01, Speed=10 [µL/s], NeedleHeight=2000 [µm]			
12	UDP_NeedleWash	Duration=5 [s]			
13	UDP_Draw				
14	UDP_NeedleWash	Duration=5 [s]			
15	UDP_Draw	Position=SR:1, Volume=System.Injection.CustomVariables.UDP_vol_02, Speed=10 [µL/s], NeedleHeight=2000 [µm]			
16	UDP_NeedleWash	Duration=5 [s]			
17	UDP_InNeedleMix	Volume=10 [µL], DrawSpeed=10 [µL/s], DispenseSpeed=10 [µL/s], Cycles=5			
18	UDP_Wait	10 [s]			
19	UDP_PrepareInject				



Chromatogram- Solvent Sandwich Injection


Acclaim™ Polar Advantage, 2.1 mm x 100 mm, 2.2 μ m



Long-chain PFAS approach- Requires higher percent organic in sample

Calibration Curves- Sandwich Injection

PFTreA at low end of calibration range



Long-chain approach

Injection volume: Total 75 μ L (96% methanol)

Theo Amount ngL	Amnt.Dev. %	Cal Point Status	Mean - AMT	% RSD-AMT
1.20	2.28	OK		
1.20	14.26	OK		
1.20	13.47	OK	1.32	6.09
11.60	-14.56	OK		
11.60	-5.84	OK	10.42	6.87
116.00	-0.03	OK		
116.00	-4.10	OK	113.60	2.93
2.30	12.91	OK		
2.30	1.34	OK		
2.30	9.23	OK	2.48	5.48
23.20	-3.78	OK		
23.20	-5.28	OK	22.15	1.11
232.00	4.48	OK		
232.00	3.60	OK	241.37	0.60
5.80	-4.29	OK		
5.80	-6.90	OK	5.48	1.96
58.00	-9.20	OK		
58.00	-7.58	OK	53.13	1.26


Summary Table- Sandwich Injection

Standards (96% methanol): 1-200 ng/L

Compound name	Retention time[min]	Est. limit of detection [ng/L]	Limit of quantitation [ng/L]	Calibration curve
PFBA Perfluorobutanoic acid	6.51	10	20	0.8384
PPPeA Perfluoropentanoic acid	8.86	10	20	0.9963
FtS 4:2	9.82	5	10	0.9944
PFB _S Perfluorobutanesulfonate	10.09	< 1	1	0.9956
PFHxA Perfluorohexanoic acid	10.73	1	2	0.9942
PPPe _S Perfluoropentansulfonate	11.75	< 1	1	0.9981
PFHpA Perfluoroheptanoic acid	12.18	1	2	0.9962
PFHx _S Perfluorohexanesulfonate	13.04	< 1	1	0.9976
PFOA Perfluorooctanoic acid	13.34	2	5	0.9971
PFHp _S Perfluoroheptanesulfonate	14.11	1	2	0.9977
PFNA Perfluorononanoic acid	14.31	1	2	0.9965
PFOS Perfluorooctanesulfonate	15.01	< 1	1	0.9946
NMeFOSAA	15.10	2	5	0.9952
PFDA Perfluorodecanoic acid	15.13	< 1	1	0.9916
PFNS _S Perfluorononanesulfonate	15.79	< 1	1	0.9964
PFUdA Perfluoroundecanoic acid	15.83	< 1	2	0.9946
PFDoA Perfluorododecanoic acid	16.45	< 1	1	0.9956
PFDS _S Perfluorodecanesulfonate	16.46	1	2	0.9981
PFTriA Perfluorotridecanoic acid	16.98	< 1	1	0.9961
PFTreA Perfluorotetradecanoic acid	17.44	< 1	1	0.9972

Chromatogram- Short Chain Aqueous Injection


Acclaim™ Polar Advantage, 2.1 mm x 100 mm, 2.2 μ m



Short-chain PFAS approach- Highly aqueous sample, longer chains compounds less soluble

Calibration Curves

PFBS at low end of calibration range



Short-chain approach



Injection volume: 225 µL (Mobile phase A sample composition)

Theo Amount ngL	Amnt.Dev. %	Cal Point Status	Mean - AMT	% RSD-AMT
1.00	42.35	OK		
1.00	27.16	OK		
1.00	22.42	OK	1.31	7.97
10.00	-0.18	OK		
10.00	3.87	OK	10.18	2.81
100.00	2.02	OK		
100.00	-0.60	OK	100.71	1.84
2.00	1.43	OK		
2.00	10.41	OK		
2.00	10.06	OK	2.15	4.74
20.00	-6.99	OK		
20.00	-2.00	OK	19.10	3.70
5.00	6.20	OK		
5.00	6.86	OK		
5.00	6.41	OK	5.32	0.32
50.00	-0.68	OK		
50.00	-3.96	OK	48.84	2.38

* Triplicates were performed for all calibration levels. However, sequence interruption led to n=2 for calibration levels 10-100 ng/L

Limitation

Short-Chain Approach



Expected Vial adsorption of mid- and long-chain PFAS using mobile phase A.

Summary Table- Short Chain Approach

PFAS Standards (mobile phaseA): 1-100 ng/L

Compound name	Retention time[min]	Est. limit of detection [ng/L]	Limit of quantitation [ng/L]	Calibration curve
PFBA Perfluorobutanoic acid	6.51	2.0-5.0	5	0.9975
PPPeA Perfluoropentanoic acid	8.86	2.0-5.0	5	0.9976
FtS 4:2	9.82	< 2	2	0.9963
PFBS Perfluorobutanesulfonate	10.09	< 1	1	0.9993
PFHxA Perfluorohexanoic acid	10.73	< 1	1	0.9979
PPPeS Perfluoropentan sulfonate	11.75	< 1	1	0.9991
PFHpA Perfluoroheptanoic acid	12.18	< 1	1	0.9976
PFHxS Perfluorohexane sulfonate	13.04	< 1	1	0.9980
PFOA Perfluorooctanoic acid	13.34	< 2	2	0.9916
PFHpS Perfluoroheptane sulfonate	14.11	< 1	1	0.9910
PFNA Perfluorononanoic acid	14.31	< 2	2	0.9850
PFOS Perfluorooctanesulfonate	15.01	< 2	2	0.9780
NMeFOSAA	15.10	2-5	5	0.9715
PFDA Perfluorodecanoic acid	15.13	< 2	2	0.9828
PFNS Perfluorononane sulfonate	15.79	< 1	1 (out of range: Amnt.Dev. RSD-AMT)	0.9683
PFUdA Perfluoroundecanoic acid	15.83	< 1	1 (out of range: Amnt.Dev. RSD-AMT)	0.9724
PFDoA Perfluorododecanoic acid	16.45	< 1	1 (out of range: Amnt.Dev. RSD-AMT)	0.9782
PFDS Perfluorodecanesulfonate	16.46	< 1	1 (out of range: Amnt.Dev. RSD-AMT)	0.9799
PFTriA Perfluorotridecanoic acid	16.98	< 1	1 (out of range: Amnt.Dev. RSD-AMT)	0.9782
PFTreA Perfluorotetradecanoic acid	17.44	< 1	1 (out of range: Amnt.Dev. RSD-AMT)	0.9778

Comparison of Short- and Long-Chain Approach

Standards Only. Limit of Quantitation.


Compound name	Retention time[min]	Limit of quantitation [ng/L] (long-chain approach)	Limit of quantitation [ng/L] (short-chain approach)
PFBA Perfluorobutanoic acid	6.51	20	5
PPeA Perfluoropentanoic acid	8.86	20	5
FtS 4:2	9.82	10	2
PFBS Perfluorobutane sulfonate	10.09	1	1
PFHxA Perfluorohexanoic acid	10.73	2	1
PPeS Perfluoropentan sulfonate	11.75	1	1
PFHpA Perfluoroheptanoic acid	12.18	2	1
PFHxS Perfluorohexane sulfonate	13.04	1	1
PFOA Perfluorooctanoic acid	13.34	5	2
PFHpS Perfluoroheptane sulfonate	14.11	2	1
PFNA Perfluorononanoic acid	14.31	2	2
PFOS Perfluorooctane sulfonate	15.01	1	2
NMeFOSAA	15.10	5	5
PFDA Perfluorodecanoic acid	15.13	1	2
PFNS Perfluorononane sulfonate	15.79	1	1 (out of range: Amnt.Dev. RSD-AMT)
PFUdA Perfluoroundecanoic acid	15.83	2	1 (out of range: Amnt.Dev. RSD-AMT)
PFDoA Perfluorododecanoic acid	16.45	1	1 (out of range: Amnt.Dev. RSD-AMT)
PFDS Perfluorodecane sulfonate	16.46	2	1 (out of range: Amnt.Dev. RSD-AMT)
PFTriA Perfluorotridecanoic acid	16.98	1	1 (out of range: Amnt.Dev. RSD-AMT)
PFTreA Perfluorotetradecanoic acid	17.44	1	1 (out of range: Amnt.Dev. RSD-AMT)

Mobile phase A
225 µL

96% methanol
75 µL

PFAS Decision Tree – Vanquish UHPLC Systems

Vanquish Duo UHPLC System for Tandem LC-MS




Vanquish Online SPE HPLC and UHPLC-Systems

Minor changes in configuration

1. Subtract column #2
2. Subtract valve 2-p 6-p #2
3. Add trap column

Adapt capillaries



Thank you

