

Water in Liquefied Petroleum Gas

Allison Mason
Shimadzu Scientific Instruments
NEMC 2022

Water in Petrochemical Feedstock

● Water Damage

- Freezing of pipeline and valves
- Blockages by gas hydrates
- Vibration Trouble
- Catalyst Poisoning

Focus is on Liquefied
Petroleum Gases (LPG)



Qualitative Water Analysis Methods

Method	Range(MIN)	Range(MAX)	Solid	Liquid	Gas	Online
GC-BID/TCD/VUV	1ppm	100%	x	x	x	x
Karl Fischer Titration	10ppm	100%	x	x	x	x
Quartz Crystal Microbalance	1ppm	1000ppm			x	x
FTIR Spectroscopy	10ppb	%	x	x	x	x
TDLA Spectroscopy	10ppb	-			x	x
CRDS spectroscopy	0.1ppb	-			x	x
Chilled Mirrors(Dew point)	ppm	%			x	x
Electrolytic	ppm	0.10%			x	
NIR Spectroscopy	0.10%	100%	x	x		x
Colorimetry	0.10%	%	x	x	x	x
Dielectric Constant	1%	10%		x	x	
Loss on Drying	0.01%	-	x	x		
Electric Resistance	0.30%	%		x	x	
Distillation(Azeotropic or not)	0.05%	-		x		
Neutron Scattering	%	%	x			
Freeze Valve	10ppm?	-		x		x
Centrifuging	%	%		x		

Pass
(no frost)



Fail
(frost)



Easy to use, but releases propane as part of the test.

ASTM D2713 (Propane only)

Quantitative Water Analysis Methods

Method	Range(MIN)	Range(MAX)	Solid	Liquid	Gas	Online
GC-BID/TCD/VUV	1ppm	100%	x	x	x	x
Karl Fischer Titration	10ppm	100%	x	x	x	x
Quartz Crystal Microbalance	1ppm	1000ppm			x	x
FTIR Spectroscopy	10ppb	%	x	x	x	x
TDLA Spectroscopy	10ppb	-			x	x
CRDS spectroscopy	0.1ppb	-			x	x
Chilled Mirrors(Dew point)	ppm	%			x	x
Electrolytic	ppm	0.10%			x	
NIR Spectroscopy	0.10%	100%	x	x		x
Colorimetry	0.10%	%	x	x	x	x
Dielectric Constant	1%	10%		x	x	
Loss on Drying	0.01%	-	x	x		
Electric Resistance	0.30%	%		x	x	
Distillation(Azeotropic or not)	0.05%	-		x		
Neutron Scattering	%	%	x			
Freeze Valve	10ppm?	-		x		x
Centrifuging	%	%	x			

ASTM D7995

Common for
Breathing Air

KFT
www.metrohm.com



Shimadzu
Loss on Drying

Karl Fischer Titration

- **Karl Fischer Titration is the leading method for moisture analysis**

- Wide dynamic range

- **Challenges of Karl Fischer Titration**

- Susceptible to side reactions from sulfur, ketones, aldehydes, amides, and siloxanes
- Need to control ambient moisture in solvents
- Not accurate for small sample sizes or samples with low water content
- Need an idea of the range of water content in the sample
- Samples should be soluble in Karl Fischer medium
- Labor intensive and not high throughput
- Requires harsh/halogenated chemicals

Green Water Analysis by Gas Chromatography


- **Gas Chromatography with Barrier Discharge Ionization Detector**

- Detection down to 1 ppm water
- No unwanted side reactions
- No solubility issues
- Less hands-on time
- Minimal chemical waste generated
- Reduced amount of sample required
- Less worker exposure to harmful chemicals
- Reduced cost



Barrier Discharge Ionization Detector (BID)

- **Barrier Discharge Ionization Detector**
- **Virtually universal detection (can detect everything except Ne and He)**
- **Low Energy and Low Maintenance**
- **Only requires Helium for operation → No flammable fuel gases required**
- **100x more sensitive than a TCD**





Three electrodes are placed on the quartz tube. High voltage is applied to electrodes, and plasma is generated

Compounds are ionized by plasma, and signal is collected at collection electrode




Ionic Liquid Column

Example Structure of an Ionic Liquid Phase




PEG/Wax
280 °C

Conventional Columns





Ionic Liquid Columns





Analytical Conditions

Model	GC-2030 Nexis 6 port gas sample valve, Split block and BID
Column	Supelco WaterCol 1910 30 m x 0.25 mm x 0.2 um df
Column Temp	Gradient 35 °C (2.0 min) → 5 °C/min → 90 °C
Injection Mode	Split 7:1
Carrier Gas Controller	Constant Linear Velocity Mode (He)
Linear Velocity	38.6 cm/sec
Discharge Gas	50 mL/min (He)
Injection Volume	1 mL




Calibration With Gases





Level	Conc. (ppm)	Mean Area	SD (n=5)	%RSD
1	5.411	146425	1.61E+03	1.10%
2	27.11	215770	1.89E+03	0.88%
3	98.9	421896	4.47E+03	0.58%

Comparison of Gas to Liquid



Standard Type	Calculated Conc. (ppm)	Conc. Recovery	Average S/N	Average LOQ (ppm)	Average LOD (ppm)
Gas	28.2 (± 0.6)	112.80%	140 (± 32)	2.12 (± 0.6)	0.70 (± 0.2)
LPG	75.9 (± 0.8)	75.90%	775 (± 132)	1.01(± 0.2)	0.33(± 0.06)


Repeatability



Standard Type	Standard Conc. (ppm)	Retention Time (min)	Average Area	%RSD (n=5)
Gas	25 ppm	10.514 (± 0.001)	215769 (± 1891)	0.88%
LPG	100 ppm	10.486(± 0.001)	10.486(± 0.001)	0.67%

Selectivity

- Methanol is commonly present in LPG
- Separation indicates that this setup is transferable to the analysis of water in breathing air



Summary

- Water present in petrochemical feedstocks is responsible for billions of dollars of damage and environmental releases every year
- Active ASTM Work Item for new method
- Environmentally friendly method
- Highly sensitive, repeatable, robust analysis
- Easy quantification of water in petrochemical streams
 - Eliminates the need for environmentally hazardous analysis techniques
- No interference of Sulfurs or alcohols in the analysis
- Simple sample introduction with the vaporizer
 - No sample preparation or hazardous chemicals required

Questions?

- **Allison Mason**

- ammason@shimadzu.com
- 410-910-0912