

Evaluation of New Anion Exchange and Synthetic Carbon Sorbents for the Determination of PFAS in Solid Samples Following EPA Method 1633

Matthew Giardina, Meg Juck

Agilent Technologies, Inc., 2850 Centerville Road, Wilmington, DE 19808, USA

matthew_giardina@agilent.com

NEMC 2023



DE77042345

Introduction

US EPA draft method 1633 was developed to consolidate procedures for the extraction and quantitation of PFAS in aqueous (nonpotable water), solids (soil, biosolids, and sediment) and tissue samples.¹ Principally, the method utilizes polymeric weak anion exchange (WAX) solid phase extraction (SPE) for the selective extraction of target analytes in addition to matrix removal using graphitized carbon black (GCB). The target analytes are extracted along with isotopically labeled standards followed by separation and detection using liquid chromatography/tandem quadrupole (LC/TQ) mass spectrometry. To date, the draft method contains validation results for solids based on a single laboratory study for a total of 40 target PFAS across nine compound classes.

The draft method contains rigorous quality control procedures to ensure optimal data reliability. The requirements are described in Section 9¹ and include: the initial demonstration of precision, accuracy, and method detection limits; the recovery of extracted internal standards and non-extracted internal standards; method blank determination; instrument calibration verification and maintenance; laboratory duplicates; analysis of field replicates when necessary; and analysis of matrix spikes when necessary.


The performance of the extraction and analysis procedures for solid matrices was verified following the draft method quality control protocols using Bond Elut PFAS WAX SPE cartridges, Carbon S (Figure 1) as a replacement for GCB, and the Agilent Infinity II 1290 LC and Agilent 6470B triple quadrupole LC/MS. The results were compared to the US EPA draft method 1633 for the single lab validation study.²




Figure 1. Agilent Bond Elut PFAS WAX and Carbon S

Experimental

Sample Workflow



Sample Preparation Procedure



Results and Discussion

Initial Method Detection Limits

- Seven replicates Ottawa sand spikes

Analyte	Bond Elut PFAS WAX Initial MDL (ng/g)	EPA Draft 1633 Aq. MDL (ng/g)*
PFBA	0.111	0.401
PFPeA	0.022	0.021
PFHxA	0.037	0.020
PFHpA	0.039	0.029
PFOA	0.017	0.037
PFNA	0.021	0.086
PFDA	0.029	0.031
PFUnA	0.023	0.033
PFDoA	0.027	0.059
PFTDA	0.031	0.038
PFTeDA	0.030	0.032
PFBS	0.034	0.014
PFPeS	0.051	0.015
PFHxS	0.030	0.018


Analyte	Bond Elut PFAS WAX Initial MDL (ng/g)	EPA Draft 1633 Aq. MDL (ng/g)*
PFHpS	0.050	0.057
PFOS	0.033	0.067
PFNS	0.043	0.046
PFDS	0.040	0.040
PFDoS	0.064	0.038
PFNA	0.087	0.282
6:2FTS	0.239	0.116
8:2FTS	0.122	0.225
PFOSA	0.061	0.068
NMeFOSA	0.045	0.049
NEtFOSA	0.099	0.038
NMeFOSAA	0.039	0.030
NEtFOSAA	0.029	0.044
NMeFOSE	0.267	0.203

Analyte	Bond Elut PFAS WAX Initial MDL (ng/g)	EPA Draft 1633 Aq. MDL (ng/g)*
NEtFOSE	0.255	0.203
HFPO-DA	0.194	0.136
ADONA	0.084	0.057
PFMPA	0.040	0.033
PFMBA	0.033	0.029
NFDHA	0.067	0.084
9CI-PF3ONS	0.020	0.038
11CI-PF3OUdS	0.058	0.071
PFEESA	0.026	0.018
3:3 FTCA	0.066	0.060
5:3 FTCA	0.101	0.363
7:3 FTCA	0.283	0.308

* EPA 1633 3rd Draft – Table 6, single-lab validation results


Accuracy

- Four midlevel sand spikes, average recovery 99.8%




Precision

- Four midlevel sand spikes, average RSD 3.5%




Topsoil Target Recovery

- Two midlevel soil spikes, average recovery 96.8%



Extracted Internal Standard Recovery

- Four soil spikes, average recovery 96.8%



Conclusions

- Fit for purpose consumables, supplies, hardware and instrumentation are the keys to successful PFAS analysis
- Agilent's Bond Elut PFAS WAX and Carbon S provide outstanding performance for PFAS application

References

- 3rd Draft Method 1633, *United States Environmental Protection Agency*, December 2022.
- Giardina, M., Juck, M., Analysis of PFAS in Solid Samples, *Agilent Technologies application note*, publication number 5994-5667EN, 2023.