

Development of a Semi-Volatile Column Optimized for the Analysis of Hazardous Waste

Chris English, Christopher Rattray, Corby Hilliard, Kristi Sellers,
Jason Herrington, Whitney Dudek-Salisbury, Brian Salisbury.

Restek Corporation, 110 Benner Circle, Bellefonte, PA 16823, United States

National Bureau of Standards
Certificate of Analysis
Standard Reference Material 16
Sulfur in Residual Fuel Oil

National Institute of Standards & Technology
Certificate of Analysis
Standard Reference Material® 1975
Diesel Particulate Extract


This Standard Reference Material (SRM) is intended for use in evaluating analytical methods for the determination of selected polycyclic aromatic hydrocarbons (PAHs) and nitro-substituted PAHs in diesel particulate extracts and

National Institute of Standards & Technology
Certificate of Analysis
Standard Reference Material® 1941b
Organics in Marine Sediment

This Standard Reference Material (SRM) is marine sediment collected at the mouth of the Baltimore (MD) Harbor. SRM 1941b is intended for use in evaluating analytical methods for the determination of selected polycyclic aromatic

National Institute of Standards & Technology
Certificate of Analysis
Standard Reference Material® 1991
Mixed Coal Tar/Petroleum Extract in Methylene Chloride

This Standard Reference Material (SRM) is intended for use in evaluating analytical methods for the determination of selected polycyclic aromatic hydrocarbons (PAHs), particularly in determining relative response factors for ASTM Method D7363-11 "Standard Test Method for Determination of Parent and Alkyl Polycyclic Aromatics in Sediment Pore Water Using Solid-Phase Microextraction and Gas Chromatography/Mass Spectrometry in Selected Ion Monitoring Mode" [1]. All of the constituents for which certified, reference, and information values are provided are naturally present in the extract. A unit of SRM 1991 consists of five ampoules each containing 1.2 mL of a methylene chloride solution.





Development of a Semi-Volatile Column Optimized for the Analysis of Hazardous Waste

Chris English, Christopher Rattray, Corby Hilliard, Kristi Sellers,
Jason Herrington, Whitney Dudek-Salisbury, Brian Salisbury.



Restek Corporation, 110 Benner Circle, Bellefonte, PA 16823, United States

**Review Details of the Lifetime Study Results
Using Contaminated Extracts**



Highly Contaminated Samples: Sulfur & Hydrocarbons




What is Dirt?




What is Dirt?




GC/MS Tune Check Solution: DFTPP, DDT, PCP & Benzidine




Sources of activity in the sample pathway




What causes activity in the column?



Isolated
Silanol



Vicinal Silanol
(H-bond)



Geminal Silanol

Increasing silanol acidity

Shifting retention times

Activity in the column will act as a 2nd retention mechanism

- *Retention (and RI) depend on degree of activity*
- *If not carefully controlled, each column can show a different RI*

Impact:


- Retention time increases Peak start to tail
- Peakheight reduces

Retention time becomes a function of mass injected on the column: Lower levels elute later

Impact:

- Retention time increases Peak start to tail
- Peakheight reduces

The consequence of strong interactions with active sites



Retention times shift as column activity increases

- Retention times can shift on active columns as the mass on column drops

Impact:

- Peaks cannot be “found”
- Misinterpretation
- Incorrect data

Other Options: Extract and Analyze Dirt



National Institute of Standards & Technology

Certificate of Analysis

Standard Reference Material[®] 1622e

Sulfur in Residual Fuel Oil
(Nominal Mass Fraction 2 %)



Lifetime test w/ Sulfurs

Other Possible Standards





National Institute of Standards & Technology

Certificate of Analysis

Standard Reference Material® 1991

Mixed Coal Tar/Petroleum Extract in Methylene Chloride



This Standard Reference Material (SRM) is intended for use in evaluating analytical methods for the determination of selected polycyclic aromatic hydrocarbons (PAHs), particularly in determining relative response factors for ASTM Method D7363-11 "Standard Test Method for Determination of Parent and Alkyl Polycyclic Aromatics in Sediment Pore Water Using Solid-Phase Microextraction and Gas Chromatography/Mass Spectrometry in Selected Ion Monitoring Mode" [1]. All of the constituents for which certified, reference, and information values are provided are naturally present in the extract. A unit of SRM 1991 consists of five ampoules each containing 1.2 mL of a methylene chloride solution.





Are fatty acids overwhelming your QuEChERS dSPE PSA cleanup and causing issues in your GC analysis? Get more cleanup capacity with cartridge SPE cleanup!

March 4th, 2014 by Michelle Misselwitz

Fatty acids are important molecules in the human body because they are used as a source of fuel. There are many food sources of both "healthy" and "unhealthy" fatty acids. Many sources of dietary fatty acids come from fruits, vegetables, seeds, nuts, and animal fats. The QuEChERS methodology was developed to analyze



Testing Contamination in Columns *and* Liners



National Institute of Standards & Technology

Certificate of Analysis

Standard Reference Material® 1975

Diesel Particulate Extract

This Standard Reference Material (SRM) is intended for use in evaluating analytical methods for the determination of selected polycyclic aromatic hydrocarbons (PAHs) and nitro-substituted PAHs in diesel particulate extracts and similar matrices. In addition to certified, reference, and information values for selected PAHs and nitro-substituted




https://www-s.nist.gov/srmors/view_cert.cfm?srm=1975

<https://www-s.nist.gov/srmors/certificates/1975.pdf>


SRM 1975 Lifetime Test

GC Inlet Maintenance


- After the 30 SRM Injections, a segment of column was cut from the head of the column (highlighted on the right), and the liner, septum, and gold seal was replaced
 - We've been saving the segments to ensure that the same length is removed each day
- Ideally, the column performance should be restored by these maintenance actions.




The ageing process during routine analysis of challenging samples is normal



Trimming alone shouldn't cause the elution order switch



10 m of column removed
Same head pressure



Column: Rxi-SVOCms, 30.00 m, 0.25 mm ID, 0.25 μ m (cat.# 16623)
Carrier Gas: Helium, Constant Flow @ 1.20 mL/min
Average Velocity: 39.95 cm/sec
Outlet Pressure (abs): 0.00 psi (Vacuum)
Oven Temp.: 60 °C (hold 0.5 min) to 285 °C @ 25 °C/min

Peaks	t_R (min)	R_s	Peak Width (min)	T_{peak} (°C)
1. 4-Nitrophenol	6.43	5.4	0.021	208.2
2. Dibenzofuran	6.54	5.4	0.022	211.0


Column: Rxi-SVOCms, 20.00 m, 0.25 mm ID, 0.25 μ m
Carrier Gas: Helium, Constant Flow @ 1.80 mL/min
Average Velocity: 59.92 cm/sec
Outlet Pressure (abs): 0.00 psi (Vacuum)
Oven Temp.: 60 °C (hold 0.5 min) to 285 °C @ 25 °C/min

Peaks	t_R (min)	R_s	Peak Width (min)	T_{peak} (°C)
1. 4-Nitrophenol	5.20	1.7	0.018	177.6
2. Dibenzofuran	5.24	1.7	0.019	178.4

4-Nitrophenol is a good example of how column activity can change relative retention times.

Adsorption at active sites


- Peak height reduces (sensitivity)
- Peaks become broader (separation)




New Column Inertness Performance

Injection of 12, 8, 4, 2, 1, 0.5 and 0.1 ng on the column

Pyridine



Pentachlorophenol



Inertness Facilitates Linear Calibrations

Average %RSDs for 6 Rxi-SVOCms columns (16623) calibrated from 1 ng/μL to 120 ng/μL

Compound	Calibration Range (μg/mL)	Average %RSD of Response Factors
N-Nitrosodimethylamine	1 - 120	4.70%
Pyridine	1 - 120	6.10%
(SS) 2-Fluorophenol	1 - 120	1.70%
(SS) Phenol-d6	1 - 120	2.10%
Phenol	1 - 120	3.20%
Aniline	1 - 120	3.10%
Bis(2-chloroethyl)ether	1 - 120	2.40%
2-chlorophenol	1 - 120	2.80%
1,3-dichlorobenzene	1 - 120	2.60%
1,4-Dichlorobenzene	1 - 120	2.10%
Benzyl alcohol	1 - 120	3.30%
1,2-Dichlorobenzene	1 - 120	2.70%
2-Methylphenol	1 - 120	3.30%
Bis(2-chloroisopropyl)ether	1 - 120	2.40%
4-Methylphenol/3-methylphenol	1 - 120	3.30%
N-nitroso-di-n-propylamine	1 - 120	3.80%
Hexachloroethane	1 - 120	3.00%
(SS) Nitrobenzene-D5	1 - 120	1.60%
Nitrobenzene	1 - 120	2.60%
Isophorone	1 - 120	3.40%
2-Nitrophenol	1 - 120	7.00%
2,4-Dimethylphenol	1 - 120	3.70%
Benzoic acid	2.5 - 120	25.00%

Compound	Calibration Range (μg/mL)	Average %RSD of Response Factors
Bis(2-chloroethoxy)methane	1 - 120	3.60%
2,4-Dichlorophenol	1 - 120	4.10%
1,2,4-Trichlorobenzene	1 - 120	2.80%
Naphthalene	1 - 120	3.20%
4-Chloroaniline	1 - 120	3.90%
Hexachlorobutadiene	1 - 120	3.70%
4-Chloro-3-methylphenol	1 - 120	4.40%
2-Methylnaphthalene	1 - 120	3.40%
1-Methylnaphthalene	1 - 120	3.60%
Hexachlorocyclopentadiene	1 - 120	6.00%
2,4,6-Trichlorophenol	1 - 120	5.90%
2,4,5-Trichlorophenol	1 - 120	6.00%
(SS) 2-Fluorobiphenyl	1 - 120	3.10%
2-Chloronaphthalene	1 - 120	2.80%
2-Nitroaniline	1 - 120	7.80%
1,4-Dinitrobenzene	1 - 120	11.10%
Dimethyl phthalate	1 - 120	3.40%
1,3-Dinitrobenzene	1 - 120	10.80%

Compound	Calibration Range (μg/mL)	Average %RSD of Response Factors
2,6-Dinitrotoluene	1 - 120	7.80%
Acenaphthylene	1 - 120	4.10%
1,2-Dinitrobenzene	1 - 120	8.10%
3-Nitroaniline	1 - 120	5.80%
Acenaphthene	1 - 120	3.30%
2,4-Dinitrophenol	2.5 - 120	17.30%
4-Nitrophenol	1 - 120	7.90%
Dibenzofuran	1 - 120	3.50%
2,4-Dinitrotoluene	1 - 120	11.60%

Benzo[ghi]perylene	1 - 120	6.40%
Average %RSD:		6.00%

Fluorene	1 - 120	4.40%
4-Nitroaniline	1 - 120	9.10%
4,6-Dinitro-2-methylphenol	2.5 - 120	15.10%
N-nitrosodiphenylamine	1 - 120	4.60%
Diphenylhydrazine	1 - 120	4.60%
(SS) 2,4,6-Tribromophenol	1 - 120	5.50%
4-Bromophenyl phenyl ether	1 - 120	5.50%
Hexachlorobenzene	1 - 120	4.30%
Pentachlorophenol	1 - 120	10.60%
Phenanthrene	1 - 120	3.70%

Compound	Calibration Range (μg/mL)	Average %RSD of Response Factors
Anthracene	1 - 120	4.80%
Carbazole	1 - 120	5.30%
di-n-Butyl phthalate	1 - 120	7.90%
Fluoranthene	1 - 120	5.10%
Benzidine	1 - 120	9.30%
(SS) Pyrene-D10	1 - 120	1.50%
Pyrene	1 - 120	4.30%
(SS) p-Terphenyl-d14	1 - 120	1.80%
3,3'-Dimethylbenzidine	1 - 120	9.50%
Butyl benzyl phthalate	1 - 120	8.60%
Bis(2-ethylhexyl)adipate	1 - 120	10.50%
3,3'-Dichlorobenzidine	1 - 120	8.50%
Benz[a]anthracene	1 - 120	3.20%
Chrysene	1 - 120	3.70%
Bis(2-ethylhexyl)phthalate	1 - 120	10.40%
Di-n-octyl phthalate	1 - 120	13.20%
Benzo[b]fluoranthene	1 - 120	5.60%
Benzo[k]fluoranthene	1 - 120	4.90%
Benzo[a]pyrene	1 - 120	6.30%
Indeno[123-cd]pyrene	1 - 120	7.20%
Dibenz[a,h]anthracene	1 - 120	7.50%
Benzo[ghi]perylene	1 - 120	6.40%
Average %RSD:		6.00%

Benzoic Acid is the only compound > 20% RSD

New Surface Deactivation: Rugged Enough for the Harshest Matrices


300 injections of a diesel particulate extract (NIST SRM 1975)

- Calibration Verification (CCV) performed every 10 SRM injections
- Maintenance (2 loops removed, new liner, new seal) every 30 SRMinjections

After 300 injections, CCV still meets expected performance targets



- 6.4% of targets analytes fall outside the $\pm 20\%$ recovery window

Standard inlet maintenance every 30 sample injections restores column performance.




New Polymer, Optimized PAH Selectivity

Improved Benzofluoranthene and Indeno[123-cd]pyrene – Dibenz[ah]anthracene Separations




New Polymer, Same Chromatogram

Typical Silarylene 5-type





New Polymer




Benzofluoranthene isomers are still sufficiently resolved after removing ten meters of from a 30 m x 0.25 mm x 0.25 μm

Benzo[B]fluoranthene - Benzo[K]fluoranthene




10m column removed
300 SRM 1975 samples analyzed

Benzo[B]fluoranthene - Benzo[K]fluoranthene




No Changes to temperature program or flow rate

Calibration Stability Compared to the Typical Columns



There is an extensive PAH library in EZGC

Extended PAH List - Rxi-SVOCms



Click and drag to zoom; double-click to reset.

Column: Rxi-SVOCms, 30.00 m, 0.25 mm ID, 0.25 µm (cat.# 16623)
Carrier Gas: Helium, Constant Flow @ 1.40 mL/min
Average Velocity: 42.74 cm/sec
Outlet Pressure (abs): 0.00 psi
Oven Temp.: 40 °C (hold 0.5 min) to 280 °C @ 20 °C/min to 330 °C @ 6 °C/min (hold 6 min)

<https://ez.restek.com/proezgc/en>

Peaks	t _r (min)	R _s	Peak Width (min)	T _m (°C)	Peaks	t _r (min)	R _s	Peak Width (min)	T _m (°C)
1. Naphthalene-D8	6.15	1.2	0.023	153.0	50. 2-Methylfluoranthene	11.31	11.6	0.019	256.3
2. Naphthalene	6.18	1.2	0.023	153.6	51. Fluoranthene	11.54	9.8	0.025	260.8
3. Quinoline-d7	6.37	1.1	0.023	157.4	52. Pyrene-d10	11.79	0.8	0.026	265.8
4. Quinoline	6.40	1.1	0.023	157.9	53. Pyrene	11.81	0.8	0.025	266.2
5. 2-Methylnaphthalene	7.00	4.9	0.023	169.9	54. p-Terphenyl	12.06	9.8	0.024	271.2
6. 1-Methylnaphthalene	7.11	4.9	0.024	172.3	55. 1-Methylpyrene	12.31	0.1	0.025	276.1
7. 2-Fluorobiphenyl	7.44	2.3	0.023	178.8	56. Benzo[c]fluorene	12.31	0.1	0.025	276.2
8. 2-Ethyl-naphthalene	7.50	0.6	0.023	179.9	57. 2-(Acetylamino)fluorene	12.92	5.1	0.027	282.5
9. 1-Ethyl-naphthalene	7.51	0.6	0.023	180.2	58. Benzo[ghi]fluoranthene	13.06	--	0.028	283.3
10. Biphenyl	7.56	0.5	0.023	181.2	59. Benzo[c]phenanthrene	13.06	--	0.029	283.4
11. 2,6-Dimethylnaphthalene	7.58	0.3	0.023	181.5	60. Cyclopenta[cd]pyrene	13.34	0.3	0.030	285.0
12. 2-Chloronaphthalene	7.58	0.1	0.024	181.7	61. Benz[a]anthracene	13.35	0.3	0.030	285.1
13. 2,7-Dimethylnaphthalene	7.59	0.1	0.023	181.8	62. Triphenylene	13.38	0.5	0.030	285.3
14. 1-Chloronaphthalene	7.61	0.8	0.024	182.1	63. Chrysene	13.40	0.5	0.030	285.4
15. 1,3-Dimethylnaphthalene	7.67	1.1	0.023	183.4	64. 1-Methylbenz[a]anthracene	13.67	1.7	0.031	287.0
16. 1,6-Dimethylnaphthalene	7.70	1.1	0.023	184.0	65. 3-Methylchrysene	13.72	1.4	0.032	287.3
17. 1,4-Naphthoquinone	7.79	0.9	0.023	185.8	66. 2-Methylchrysene	13.77	1.4	0.032	287.6
18. 1,4-Dimethylnaphthalene	7.81	--	0.023	186.2	67. 4-Methylchrysene	13.83	--	0.032	288.0
19. 2,3-Dimethylnaphthalene	7.81	--	0.023	186.2	68. 6-Methylchrysene	13.83	--	0.032	288.0
20. 1,2-Dimethylnaphthalene	7.90	4.1	0.023	188.1	69. 5-Methylchrysene	14.06	4.1	0.034	289.4
21. 1,8-Dimethylnaphthalene	8.04	1.4	0.023	190.8	70. Benz[a]anthracene-7,12-dione	14.21	4.1	0.034	290.3
22. Acenaphthylene	8.07	1.4	0.024	191.5	71. 7,12-Dimethylbenzo[a]anthracene	15.05	--	0.040	295.3
23. Acenaphthene-d10	8.25	1.6	0.024	195.0	72. Benzo[b]fluoranthene	15.05	--	0.040	295.3
24. Acenaphthene	8.29	1.6	0.024	195.8	73. Benzo[j]fluoranthene	15.07	0.3	0.040	295.4
25. 1,3,7-Trimethylnaphthalene	8.34	2	0.023	196.7	74. Benzo[k]fluoranthene	15.10	0.9	0.040	295.6
26. 1,4,6-Trimethylnaphthalene	8.47	1	0.023	199.5	75. Benzo[a]fluoranthene	15.24	3.2	0.041	296.4
27. Dibenzofuran	8.50	0.3	0.024	200.0	76. Benzo[e]pyrene	15.54	1.2	0.043	298.2
28. 2,3,6-Trimethylnaphthalene	8.51	0.3	0.023	200.1	77. Benzo[a]pyrene-d12	15.59	0.9	0.043	298.5
29. 2,3,5-Trimethylnaphthalene	8.58	0.1	0.023	201.6	78. Benzo[a]pyrene	15.63	0.9	0.043	298.8
30. 1-Naphthalamine	8.58	0.1	0.024	201.7	79. Perylene-D12	15.74	1.1	0.044	299.4
31. 1,2,6-Trimethylnaphthalene	8.60	0.5	0.023	201.9	80. Perylene	15.78	1.1	0.044	299.7
32. 2-Naphthalamine	8.69	1.6	0.024	203.7	81. p-Quaterphenyl	17.02	8	0.048	307.1
33. 2,4,5-Trimethylnaphthalene	8.72	0.4	0.024	204.5	82. Dibenz[a,j]acridine	17.41	3	0.051	309.4
34. 1,2,5-Trimethylnaphthalene	8.74	0.4	0.023	204.7	83. Dibenzo[a,j]anthracene	17.57	3	0.052	310.4
35. 1,2,3-Trimethylnaphthalene	8.84	0.6	0.023	206.8	84. Indeno[1,2,3-cd]pyrene	17.75	1.1	0.053	311.5
36. 1,4,5-Trimethylnaphthalene	8.86	0.6	0.023	207.1	85. Dibenz[a,h]anthracene	17.81	--	0.053	311.8
37. Fluorene	8.92	0.4	0.024	208.3	86. Dibenz[a,c]anthracene	17.81	--	0.053	311.8
38. 2-Methyl-5-nitroaniline	8.93	0.4	0.023	208.5	87. Benzo[b]chrysene	17.99	1.2	0.053	313.0
39. 1,4,6,7-Tetramethylnaphthalene	9.29	12.4	0.023	215.8	88. Picene	18.06	1.2	0.054	313.4
40. 1,2,5,6-Tetramethylnaphthalene	9.58	12.4	0.024	221.6	89. Benzo[ghi]perylene	18.24	3.4	0.055	314.5
41. Dibenzothiophene	9.93	5.1	0.025	228.6	90. Anthanthrene	18.51	4.7	0.056	316.0
42. Phenanthrene-D10	10.06	1.1	0.025	231.1	91. 6,13-Pentacenedione	19.32	14.4	0.057	320.9
43. Phenanthrene	10.08	1.1	0.025	231.7	92. Dibenzo[a,l]pyrene	20.48	0.4	0.060	327.9
44. Anthracene	10.15	2.6	0.025	233.0	93. Dibenzo[b,k]fluoranthene	20.51	0.4	0.060	328.1
45. Carbazole	10.35	8.1	0.024	237.0	94. Dibenzo[a,e]pyrene	21.16	0.1	0.064	330.0
46. 2-Methylanthracene	10.58	2.1	0.024	241.7	95. Coronene	21.17	0.1	0.065	330.0
47. 1-Methylanthracene	10.63	0.7	0.024	242.7	96. Dibenzo[a,l]pyrene	21.42	2	0.067	330.0
48. 1-Methylphenanthrene	10.65	0.7	0.024	243.1	97. Dibenzo[a,h]pyrene	21.56	2	0.069	330.0
49. 4H-Cyclopenta[def]phenanthrene	10.81	6.4	0.025	246.2	98. p-Quinquephenyl	25.51	57.4	0.109	330.0

Method Development Tools

Restek EZGC Method Translator and Flow Calculator

EZGC™ Method Translator

Carrier Gas

Original	Translation
Helium	Helium

Column

Length	30.00	30.00	m
Inner Diameter	0.25	0.32	mm
Film Thickness	0.25	0.25	µm
Phase Ratio	250	320	

Control Parameters

Outlet Flow	1.20	1.54	mL/min
Average Velocity	39.57	44.77	cm/sec
Holdup Time	1.26	1.12	min
Inlet Pressure (gauge)	9.48	1.99	psi
Outlet Pressure (abs)	0.00	0.00	psi

Oven Program

Ramp (°C/min)	Temp (°C)	Hold (min)	Ramp (°C/min)	Temp (°C)	Hold (min)
40	280	0	26.3	280	0
6	325	0	8.3	325	0

Control Method

Constant Flow

Results

Run Time	25.50	18.40	min
Speed		1.39	x

EZGC™ Flow Calculator

Carrier Gas

Helium

Column

Length	30.00	m
Inner Diameter	0.25	mm
Film Thickness	0.25	µm
Temperature	40	°C

Control Parameters

Outlet Flow	1.40	mL/min
Average Velocity	42.74	cm/sec
Holdup Time	1.17	min
Inlet Pressure (gauge)	11.42	psi
Outlet Pressure (abs)	0.00	psi

Inlet

Temperature	250	°C
Liner Volume	1.00	mL
Flow	1.40	mL/min
Splitless Valve Time	1.1 to 1.5	min

RESTEK

For software updates, web-based application, and technical support, visit us at... www.restek.com/ezgc-mtfc

©2022 Restek Corporation. All rights reserved.

Compounds Conditions My EZGC

Search by Phase Search by Name or CAS >>

Phase: Rxi-SVOCms

Library: Choose a compound library

Semivolatiles

Pesticides & Herbicides

Pro EZGC Chromatogram Modeler

Send Us Your Feedback ?

Compounds Conditions My EZGC

Search by Phase Search by Name or CAS >>

Phase: Rxi-SVOCms

Library: Semivolatiles

Deselect All Filter Compounds: X

Compound Name	CAS
1,2,3-Trimethylnaphthalene	879-12-9
1,2,4,5-Tetrachlorobenzene	95-94-3
1,2,4-Trichlorobenzene	120-82-1
1,2,5,6-Tetramethylnaphthalene	2131-43-3
1,2,5-Trimethylnaphthalene	641-91-8
1,2,6-Trimethylnaphthalene	3031-05-8
1,2-Dichlorobenzene	95-50-1
1,2-Dimethylnaphthalene	573-98-8
1,2-Dinitrobenzene	528-29-0
1,3,5-Trinitrobenzene	99-35-4
1,3,7-Trimethylnaphthalene	2131-38-6
1,3-Dichlorobenzene	541-33-1


213 Selected Targeted to resolve


Need more help? If you need further assistance, feel free to send this compound list and any questions you may have to Restek Technical Service.

Column: Rxi-SVOCms, 30 m, 0.25 mm ID, 0.25 µm (cat.# 16623)
 Carrier Gas: Helium, Constant Flow @ 1.40 mL/min
 Average Velocity: 42.74 cm/sec
 Outlet Pressure (abs): 0.00 psi
 Oven Temp.: 40 °C (hold 1 min) to 120 °C @ 16 °C/min to 330 °C @ 4 °C/min (hold 3 min)

Peaks	t _r (min)	R _s	Peak Width (min)	T _{max} (°C)	Peaks	t _r (min)	R _s	Peak Width (min)	T _{max} (°C)
1, 1,4-Dioxane-d8	2.71	1	0.023	67.4	109, 2,4,5-Trimethylnaphthalene	15.19	2.1	0.069	156.8
2, 1,4-Dioxane	2.73	1	0.023	67.8	110, 1,2,5-Trimethylnaphthalene	15.34	2.1	0.068	157.3
3, N-Nitrosodimethylamine	2.97	2.2	0.023	71.5	111, 1,2,3-Trimethylnaphthalene	15.73	0.9	0.069	158.9
4, Pyridine	3.02	2.2	0.024	72.4	112, 1,4,5-Trimethylnaphthalene	15.80	0.3	0.069	159.2
5, Ethyl methacrylate	3.37	12.6	0.024	77.9	113, Fluorene	15.82	0.3	0.071	159.3
6, 2-Picoline	3.67	3.4	0.025	82.8	114, Diethyl Phthalate	15.86	0.5	0.067	159.4
7, N-Nitrosomethylethylamine	3.76	3.4	0.025	84.2	115, 2-Methyl-5-nitroaniline	16.04	0.5	0.070	160.2
8, Methyl methanesulfonate	4.07	7.2	0.025	89.2	116, 4-Nitroaniline	16.08	0.1	0.070	160.3
9, 2-Fluorophenol	4.25	6.4	0.025	92.0	117, 4-Chlorophenyl phenyl ether	16.09	0.1	0.070	160.3

1 ng 8270 MegaMix on the Rxi-SVOCms 20 m x 0.15 mm x 0.15 μm





Thanks to:

Christopher Rattray

Corby Hilliard

Roy Lautamo

Brian Salisbury

Kristi Sellers

Applications, Development

Applications, Development

R&D Polymer Chemist

Product Development

GC Marketing

1. Phase for semivolatiles developed with unique characteristics
 2. Designed for improved lifetime based on testing using NIST SRM 1975
 3. Targeted for PAH resolutions
- *Each column is tested with semivolatiles at 1 ng level & 340°C*

Pro EZGC Chromatogram Modeler

ChromaBLOGraphy



DOWNLOAD

EZGC Method Translator and Flow
Calculator

For Windows 10/8/7/Vista/XP

Questions?