

# GC/MS/MS Determination of Semivolatiles Using Nitrogen Carrier, Reagent and CID Gas

Douglas Stevens  
Stuart Oehrle, Frank Dorman

*Waters Corporation*



# Hurdles and Highlights for Helium Carrier Gas

- Non-renewable
- Petroleum production by-product
- **Supply instability**
- **Pricing instability**
- **Inertness**
- Chromatographic resolving power
- Compatibility with MS vacuum
- Compatibility with MS ionization
- Broad applications/publications

**Will I be able to deliver results on time?**



# Alternatives to Helium



## ▪ Nitrogen

- Local generation possible
- Inertness
- Slower optimum linear velocity
- Requires column scaling
- Adverse effects on sensitivity for EI
  
- How does N<sub>2</sub> effect GC-APCI?


## ▪ Hydrogen

- Local generation possible
- Undesirable reactivity
- Higher optimum linear velocity
- Requires column scaling
- Adverse effects on sensitivity for EI
  
- Not supported for GC-APCI

# GC-APCI MS: Gas Chromatography Atmospheric Pressure Chemical Ionization Mass Spectrometry




# Multiple Reaction Monitoring: MRM (aka SRM)



- **MS1 filters precursor ions**
- **Gas filled cell with fragmenting voltage applied**
- **Ions undergo collision induced dissociation (CID)**
- **MS2 passes only specific product ions**

# Why MS/MS?


- Higher confidence
  - Decrease/eliminate false positives/negatives
- Ease of integration
  - Easier to automate
  - Less manual reprocessing
- Time Savings
  - Allows faster run times
  - Reduce/eliminate sample clean up steps
- Reduced/eliminated chemical noise
  - Increases accuracy
  - Improves precision
  - Improves dynamic range



## Specificity & Sensitivity

# PAHs: N<sub>2</sub> v He, Without Column Scaling

- N<sub>2</sub> ~10% longer runtime
- He v. N<sub>2</sub> functionally equivalent separation
- Rxi-35Sil MS 30m x .25 x.25
  - No commercially available scaled column
- [Nitrogen Carrier Gas White Paper, 2023](#)



# Translating Methods from He to N<sub>2</sub> Carrier Gas

- Define gas type
- Scale column dimensions
  - Keep phase ratio constant
- Outlet Pressure for APGC (Atm) or EI (Vacuum)


▪ [Link to EZGC Method Translator, \(Restek\)](#)

Carrier Gas	Original	Translation
	Helium	Nitrogen
Column		
Length	30.00	20.00 m
Inner Diameter	0.25	0.15 mm
Film Thickness	0.25	0.15 μm
Phase Ratio	250	250
Control Parameters		
Column Flow	1.40	0.38 mL/min
Average Velocity	42.74	22.40 cm/sec
Holdup Time	1.17	1.49 min
Inlet Pressure	11.42	17.93 psi
Outlet Pressure (abs)	0.00	14.70 psi
	Atm Vacuum	Atm Vacuum
Oven Program		
<input type="radio"/> Isothermal		
<input checked="" type="radio"/> Ramps		
	Ramp Rate (°C/min)	Temp (°C) Hold Time (min)
Ramps (1-4)		
1	8.5	330 1
Control Method		
	Constant Flow	
Results		
	Solve for	Efficiency Speed Translate Custom




# OCPs: N<sub>2</sub> v He, Scaled Column

- **N<sub>2</sub> 40m** x .18 x .18  
Rxi-5Sil MS
- **He 60m** x .25 x .25  
Rxi-5Sil MS
- 50 peaks (39/11)
- +/- 1.1s average
- He at 1.5 mL/min
- N<sub>2</sub> at 0.45 mL/min



# Phthalates: N<sub>2</sub> v He, Scaled Column



Rxi-5Sil MS **30 m** x 0.25 mm x 0.25  $\mu\text{m}$  for **He** and **20 m** x 0.15 mm x 0.15  $\mu\text{m}$  for **N<sub>2</sub>**

Average  $r^2$  0.996 for 16 of 17 analytes. Method includes 11 IS/Surr.

# Multi-Class Semivolatiles (SVOCs)

- Single class SVOC analyses successfully translated to N<sub>2</sub> carrier GC-APCI
- Multi-class SVOC methods present additional challenges
  - Wider range of ionization characteristics
  - Different, specific requirements for long term stability assessment
- With the change to N<sub>2</sub> carrier, evaluation of N<sub>2</sub> collision gas in place of Ar was performed




# GC/MS Conditions


MS	Xevo™ TQ-Absolute	GC	Agilent 8890
Source Type	APGC™, dry source	Column	Rxi®-SVOCms <b>20m x 0.15 x 0.15</b>
Source Temp	150°C	Outlet	Atm
Transfer Line Temp	320°C	Injection	SSL, 310°C, Split 10:1, 4 mm pkd liner
Corona Current	2.0 µA	Carrier Gas	<b>N<sub>2</sub>, 0.25 mL/min</b>
Auxiliary Gas	200 L/hr	Temp Program	40°C 1min, to 120°C at 15°C/min, to 335°C at 30°C/min hold 7 min. 20.5min run time
Cone Gas	270 L/hr		
Make Up Gas	350 mL/min		
CID Gas	<b>N<sub>2</sub>, 0.4 mL/min</b>		

# Translating From Ar CID Gas to N<sub>2</sub>

- Dry source APGC method **developed on TQ 2** and **transferred to TQ 1**
- **TQ 2, N<sub>2</sub> carrier gas and N<sub>2</sub> CID gas**
- **TQ 1, He carrier gas and Ar CID gas**
- Identical MRM Transitions
- **TQ 2 N<sub>2</sub> CEs modified for Ar on TQ 1**
  - Based on guidance in graph



# Estimated v. Optimized, N<sub>2</sub> CID to Ar CID Gas





CE estimate  
from graph =  
8eV

Optimized CE  
= 7eV

Multiple  
analytes  
classes  
showed  
similar  
agreement


# Comparison Between Different Configurations



76 analytes  
in 8270  
Megamix  
were within  
+/- 2X

500pg in-  
vial  
50pg on  
column



# Comparison Between Different Configurations






# Comparison Between Different Configurations

**He carrier  
Ar CID**




**N<sub>2</sub> carrier  
N<sub>2</sub> CID**

# All N<sub>2</sub> System Quantitation, TQ2




**Avr  $r^2 > 0.992$   
for all analytes  
100fg to 100pg  
on-column**


# All N<sub>2</sub> System Quantitation, TQ2



# All N<sub>2</sub> System Quantitation, TQ2




# All N<sub>2</sub> System Quantitation, TQ2



**6 IStd**  
**<10% RSD**

# All N2 System, TQ2 v. He/Ar System TQ1

- Coefficient of determination > 0.992 v. 0.994
- %RSD IStds <10% v < 7%
- 65 of 67 analytes gave same LOQ on both



# Conclusions

- GC-APCI enables translation of SVOC methods to N<sub>2</sub> carrier gas without loss of separation or sensitivity
- N<sub>2</sub> CID gas, in place of Ar, allows a GC/MS/MS system to be operated using a single gas supply
- Dry source (charge exchange) ionization mode demonstrated compatibility with multi-class SVOCs
  
- However...
- N<sub>2</sub> carrier gas not allowed for EPA SVOC methods (8270, 525, 1625)
- Ability to meet long term stability requirements CCV of some methods requires further evaluation

# Acknowledgements

- Quebec Laboratory for Environmental Testing (QLET)
  - PAH, OCP, phthalate samples, standards and methods
- Gordon Fujimoto, Sarah Dowd, Lindsay Hatch – Waters Corporation





# Thank You for Your Attention. Questions?



**Douglas Stevens**  
**Principal Scientist**  
**[doug\\_stevens@waters.com](mailto:doug_stevens@waters.com)**