

# Determination of Microplastics in Water and Wastewater Using IR Spectroscopy



William Lipps

Shimadzu Scientific Instruments, Inc.

NEMC 2023




# Why should we measure microplastics?



Chemical contaminants are trapped on the microplastics, then can be consumed and adsorbed into fish, animals, or humans

# ASTM D19 has developed and is developing several guides, practices, and methods for microplastics in water



Reference  
materials  
D8402



Sample  
Collection  
D8332



Sample  
Preparation  
D8333




Analysis

DIA  
D8489

PY-GCMS  
WK67788

IR  
Microscopy  
WK67565

# Sample collection by D8332



# D8332 Sampling at a WWTP



The key is a large volume through sieves

# Sample preparation by D8333



• H<sub>2</sub>O<sub>2</sub>

• Schweizer 5 minutes

• Ammonium Hydroxide 5 minutes

• Buffer

• Protease, buffer, lipase  
• React 20 hours



Final sample = plastics suspended  
in methanol

# Further prep may be needed after D8333, depending on the method



# Once the final sample is prepared it is ready for measurement



Pyrolysis GCMS  
WK67788



Dynamic Imaging  
D8489



WK67565



FTIR Microscopy





Manual IR  
Spectrometry




# Smaller particles can be identified and counted using software and the IR Microscope




Microplastics on a filter



Figs. 3-(a) PE Distribution





Figs. 3-(b) PP Distribution



Figs. 3-(C) PET Distribution

# But what if you don't have an IR Microscope?



Microplastics Contamination in a High Population Density Area of the Chao Phraya River, Bangkok  
June 2020 *Journal of Engineering and Technological Sciences* 52(4):534-545  
DOI:10.5614/j.eng.technol.sci.2020.52.4.6

# Unlike the ISO TC147 workgroup on microplastics, ASTM WK67565 will include optical microscopy and manual IR

ASTM WK7565 is being written in collaboration with The Southern California Coastal Water Research Project (SCCWRP) to help ensure consistency with CA SOP, with special thanks to:

Charles S. Wong, Ph.D. Southern California Coastal  
Water Research Project Authority 3535 Harbor Blvd. Suite  
110 Costa Mesa CA 92626 [charlesw@sccwrp.org](mailto:charlesw@sccwrp.org)

Scott Coffin, Ph.D. State Water Resources Control Board  
1001 I St. Sacramento, CA 95814  
[scott.coffin@waterboards.ca.gov](mailto:scott.coffin@waterboards.ca.gov)

# Scope:

---


- identification and quantification of microplastic particles in
  - municipal wastewater influent,
  - treated effluent from a wastewater treatment plant,
  - ambient waters,
  - finished drinking water,
  - and bottled water.
- 20  $\mu\text{m}$ -5 mm
- using visual microscopy for particle counts and infrared (IR) spectroscopy for chemical identification of counted particles.
  - Fourier transform IR (FTIR),
  - laser direct infrared (LDIR) imaging,
  - and other techniques capable of measuring IR spectra from particles in this size range.

# Method Summary


- samples collected using Practice D8332
- prepared for analysis using Practice D8333
- Samples are viewed using a stereomicroscope and microplastic particles are identified.
- A representative subsample of particles is selected and prepared for IR spectroscopy.
- Each subsampled particle is measured and optionally photographed
- identified individually using IR spectroscopy and spectral matching to a spectral reference library or standards of known composition.
- This method is applicable for the positive identification of common plastics including polyethylene (HDPE and LDPE), polypropylene (PP), polyvinyl chloride (PVC), polyurethane (PUR), polyethylene terephthalate (PET), and polystyrene (PS). (PVC may not be detectable if using a specific gravity separation process during sample preparation.)

# Interferences

- Undigested organic matter and inorganics
- Extra digestion or removal after D8333 may be needed
- May be limited to particles > 20  $\mu\text{m}$  in diameter
- Surface contamination, such as oxide coatings
- Thermal and UV damage
- Other techniques, such as Pyrolysis GC-MS may be needed
- Visual microscopy limited to particles  $\geq 106 \mu\text{m}$



# Stereo Microscope Counting and IR Analysis of Sieve Fractions <500 $\mu\text{m}$



- place a 20  $\mu\text{m}$  filter onto the filter holder and secure the filtering funnel on top
- pour the prepared sieve fraction through the filtration system and rinse
- carefully remove the filter, avoiding sample loss, transfer into a clean petri dish
- Repeat for all sieve fractions

# Use Visual Microscopy with appropriate magnification and illumination settings

- In each sieve fraction  $> 212 \mu\text{m}$ :
- Identify similar particles
- Count the similar particles
- Characterize by color and morphology
- Select a minimum of 30 each and transfer to the IR for reflectance or ATR






# IR Microscopy Analysis of Samples on slides or filters < 212 $\mu\text{m}$

- For slides, transfer a representative sample aliquot from each prepared sieve fraction <212  $\mu\text{m}$  to a mirrored microscope slide dropwise and evaporate to dryness.
- For filters, transfer each filter containing particulates from a prepared sieve fraction <212  $\mu\text{m}$  to a glass petri dish, cover, and let dry.
- Dry each slide or filter completely before analysis and, if facilitated by heating, do not heat in excess of 50°C to prevent possible degradation of the plastic material.
- Initiate optical image of the slide containing the sample for obtaining the location and size distribution of particles.



# Identify the various plastics by comparing spectrum to a library




- Generate a spectrum of each of the above selected particles.
- Use a commercial or laboratory generated library to compare spectrum with reference spectrum. Positive identification is made with a  $\geq 60\%$  match.
- Less than 60 % may mean interferences are present.
- An automated system may be used

# What are the next steps

---

- Use D8402 to create reference materials
- Prepare a D8333 mix containing various plastics at different counts
- Analyze to obtain precision of counts and accuracy of ID
  
- Prepare samples for an Inter-lab study with a minimum of six labs

# In addition, develop new techniques and methods that can be used to measure and ID microplastic pollution



Questions?

wclipps@shimadzu.com