

# Microplastics Analysis using Accelerated Solvent Extraction (ASE) and Pyrolysis Gas Chromatography / Mass Spectrometry (Pyr-GC/MS)

**Chris Shevlin**

Scientific Affairs Manager

NEMC 2023

August 1, 2023


 The world leader in serving science



# The microplastic problem

## Microplastics detected in meat, milk and blood of farm animals

Particles found in supermarket products and on Dutch farms, but human health impacts unknown



Scientists found microplastics in 75% of meat and milk in their pilot study. Photograph: David Kelly

Microplastic contamination has been reported in the blood of cows and pigs

## Microplastics found in human blood for first time

Exclusive: The discovery shows microplastics in the body and may lodge in organs



Microplastics cause damage to human health. Photograph: David Kelly

## Microplastics found in freshly fallen Antarctic snow for first time


New Zealand researchers identified tiny plastics, which are toxic to plants and animals, in 19 snow samples



Research published in The Cryosphere journal identified microplastics in Antarctica for the first time. Photograph: Alex Aves

## Viruses survive in fresh water by 'hitchhiking' on plastic, study finds

Intestinal viruses such as rotavirus were found to be infectious for up to three days by attaching to microplastics, research shows



Human enteric, or intestinal, viruses survive for days on plastic particles easily swallowed by swimmers, Stirling University researchers found. Photograph: GJ Flick and DD Kuhn/Virginia Tech


# Analysis of microplastics: challenges

## Sample collection and preparation

- Critical to understand distribution, composition, size
- Automation where feasible to reduce costs and accelerate sample throughput
- Include sub 1µm size (nanoplastic)

Technique	Attribute
FTIR	<ul style="list-style-type: none"> <li>• Widely used particle number concentration</li> <li>• Polymer identification (&gt; 10 µm)</li> <li>• Spectra match to library</li> </ul>
Raman micro-spectroscopy	<ul style="list-style-type: none"> <li>• Size and shape to 1 µm</li> </ul>
Py-GC-MS	<ul style="list-style-type: none"> <li>• Provides molecular level information to determine chemical composition/structure <u>regardless of size</u></li> <li>• Polymer determination</li> <li>• Contaminants on MP surface</li> <li>• Quantitative and unknowns</li> </ul>


# Microplastics analysis



# Microplastics analysis

Analysis Goal

Target analysis



- Detection of target polymers
- Quantitation

Non-targeted




- Detection of non-target polymers
- Identification of unknowns



# Accelerated solvent extraction (ASE) Pyr-GC-MS

Sample prep is extremely important!



# Accelerated Solvent Extraction

- Automates sample preparation for solid and semisolid samples using solvents at elevated temperatures and pressure.
- Operates above the boiling point of extraction solvents by using pressurized sealed extraction cells.
- Ensures complete dissolution of all polymer materials



**Thermo Scientific™ EXTREVA™ ASE™  
Accelerated Solvent Extractor**



**Thermo Scientific™ Dionex™ ASE™ 350  
Accelerated Solvent Extractor system**

# Extraction and analysis of rice samples

## Sample Preparation Procedure




### Matrix spiked store-bought rice

- Polyethylene (PE)
- Polystyrene (PS)
- Polypropylene (PP)
- Polyethylene Terephthalate (PET)
- Polycarbonate (PC)
- Polyvinyl Chloride (PVC)

Samples freeze dried  
then milled to fine  
powder for 30 min.

Glass fiber filter



Glass fiber filter

10mL ASE Cell



# ASE extraction conditions

Parameter	Extraction parameters
Cell type	Stainless steel
Extraction solvent	Dichloromethane (DCM)
Extraction/Oven temperature (°C)	180
Static time	5
Cycles	2
Rinse volume (%)	80
Purge time (s)	75
System rinse volume (mL)	9
Heating time (min)	9
Pressure (psi)	1,500



After Extraction

- 80  $\mu$ L extract transferred into a pyrolysis cup
- Evaporated for 30 min in fume hood
- Loaded onto autosampler for Pyr-GC-MS

**The six polymers investigated remained sufficiently dissolved  
for over 2.5 hours post extraction**

# Pyrolysis GC-MS parameters

Parameter	Setting
<b>Micro-furnace Multi-Shot Pyrolyzer™ (Double-Shot analysis) EGA/PY-3030D (Frontier Lab)</b>	
First-shot furnace temperature (thermal desorption)	Ramped from 100 °C → 20 °C /min → 300 °C (1 min)
Second-shot furnace temperature (pyrolysis)	650 °C
Interface temperature	320 °C
Pyrolysis time	0.20 min (12 s)
<b>GC</b>	
Column	Ultra Alloy™ 5 capillary column (30 m, 0.25 mm i.d., 0.25 µm film thickness) (Frontier Lab)
Injector port temperature	300 °C
Column oven temperature program	40 °C (2 min) → (20 °C /min) → 320 °C (14 min)
Injector mode	Split/splitless (split 50:1)
Carrier gas	Helium, 1.0 mL/min, constant linear velocity
<b>MS</b>	
Ion source temperature	250 °C
Ionization energy	Electron ionization (EI); 70 eV
Scan range	40 to 600 <i>m/z</i>

# Polymer targets and recoveries

Plastic type	Pyrolysis product	Indicator ions ( <i>m/z</i> )	Molecular ion ( <i>m/z</i> )	LOQ (µg/g)
PP	2,4-dimethyl-1-heptene	70, 83, <b>126</b>	126	1.25
PS	5-hexene-1,3,5-triyltribenzene (styrene trimer)	<b>91</b> , 117, 194, 312	312	0.94
PET	Dimethyl terephthalate* Vinyl benzoate	194, <b>163</b> <b>105</b> , 77, 148, 51	194 148	2.86
PC	Bisphenol A (BA)	<b>213</b> , 119, 91, 165, 228	242	1.73
PE	n-alkene (C <sub>10</sub> , C <sub>12</sub> , C <sub>14</sub> )	<b>83</b> , 111, 140	140, 196	3.95
PVC	Naphthalene	<b>128</b> , 132, 146, 116, 102	128	3.97
Internal standard				
Polystyrene-d <sub>5</sub>	Styrene monomer	<b>109</b> , 82, 54, 107, 108		

\*Only after TMAH treatment

	PE	PP	PET (VB)	PET*	PS	PVC	PC
Spike 1	95	96	99	91	86	80	128
Spike 2	75	75	97	66	96	73	131
Spike 3	86	78	97	87	90	93	132
<b>Avg</b>	<b>85</b>	<b>83</b>	<b>97</b>	<b>81</b>	<b>90</b>	<b>82</b>	<b>130</b>
<b>St Dev</b>	<b>10</b>	<b>11</b>	<b>1</b>	<b>13</b>	<b>5</b>	<b>10</b>	<b>2</b>

\*After TMAH treatment, VB: quantification using vinyl benzoate

# Real world rice samples

Sample	Concentration of polyethylene (µg/g dw)				Concentration of polypropylene (µg/g dw)				Concentration of polyethylene terephthalate (µg/g dw)				Overall (µg/g dw)				
	Not shaken		Shaken		Not shaken		Shaken		Not shaken		Shaken		Not shaken		Shaken		
	Not washed	Washed	Not washed	Washed	Not washed	Washed	Not washed	Washed	Not washed	Washed	Not washed	Washed	Not washed	Washed	Not washed	Washed	
1	317	143	207	166	3	3	10	<LOD	<LOD	<LOD	<LOD	<LOD	<LOD	320	146	217	166
2	314	132	260	209	3	<LOD	3	<LOD	4	<LOD	8	17	322	132	272	226	
3	64	47	55	51	105	<LOD	<LOD	3	<LOD	<LOD	<LOD	<LOD	<LOD	168	47	55	54
4	94	51	56	51	14	10	<LOD	<LOD	<LOD	<LOD	<LOD	<LOD	108	61	56	51	
5	60	45	56	45	3	<LOD	<LOD	<LOD	<LOD	<LOD	<LOD	<LOD	64	45	56	45	
6	57	48	67	53	3	<LOD	7	<LOD	<LOD	<LOD	<LOD	<LOD	61	48	74	53	
7	55	58	55	54	3	<LOD	<LOD	<LOD	<LOD	<LOD	<LOD	<LOD	59	58	55	54	
8	60	46	57	46	3	<LOD	<LOD	<LOD	<LOD	<LOD	<LOD	<LOD	63	46	57	46	
9	60	54	57	54	3	<LOD	3	<LOD	<LOD	<LOD	<LOD	<LOD	63	54	61	54	
10	63	53	54	50	0	3	<LOD	<LOD	<LOD	<LOD	<LOD	<LOD	63	56	54	50	
11	79	51	47	45	3	3	<LOD	<LOD	<LOD	<LOD	<LOD	<LOD	82	54	47	45	
12	68	52	63	56	3	<LOD	<LOD	<LOD	<LOD	<LOD	<LOD	<LOD	71	52	63	56	
13	51	51	48	60	3	<LOD	<LOD	<LOD	<LOD	<LOD	<LOD	<LOD	54	51	48	60	

dw = dry weight

# Find out more:



## Extraction and analysis of plastics in rice samples using accelerated solvent extraction and pyrolysis-gas chromatography-mass spectrometry

### Authors

Jake O'Brien, Elvis Okofio, Cassie Rauert,  
Kevin Thomas  
University of Queensland,  
Brisbane, Australia

### Goal

To describe a new method for the extraction and analysis of microplastics using accelerated solvent extraction

### Introduction

Increased production and use of plastics have resulted in growth in the amount of plastic debris accumulating in the environment. This plastic debris can potentially fragment into smaller pieces, with particles <5 mm and <0.1 µm defined as microplastics (MPs) and nanoplastics (NPs), respectively.<sup>1</sup> Over the past decades, an increasing number of studies have reported the occurrence of MPs/NPs in the aquatic and terrestrial environments, including oceans, rivers, lakes, air, soil, and dust.<sup>2-6</sup> Most of the previously reported studies have typically relied on visual inspection and spectroscopic imaging approaches, reporting data on the size, shape, color, number, and polymer type of particles. These measurements may not reflect the total mass concentration of polymers in samples because the approaches are typically limited by size. The number of studies now reporting mass-based concentrations of common plastics is growing, as these methods can report the total plastic content of samples, including plastics in both the

### Keywords

Accelerated solvent extraction (ASE),  
Dionex ASE 350 Accelerated Solvent  
Extractor, microplastics, nanoplastics

thermo scientific




- [www.thermofisher.com/microplastics](http://www.thermofisher.com/microplastics)

# Microplastics analysis

Analysis Goal

Target analysis



- Detection of target polymers
- Quantitation

Non-target analysis




- Detection of non-target polymers
- Identification of unknowns



# Polymer standards & samples

Supplied by NILU & University of Queensland



- Polymer standards supplied by NILU
- Enables pyrolysis products to be determined
- HRAM spectrum obtained and added to library
- Two types of samples:
  - Stormwater samples (Filtration)
  - Milk & meat samples (ASE)

# Preparation of milk and meat samples



Extraction parameter	Setting
Cell type	5 mL Stainless Steel
Extraction solvent	Dichloromethane (DCM)
Extraction temperature	180 °C
Static time	5 min.
Extraction cycles	3
Pressure	1,500 psi



- 80 µL extract transferred into a pyrolysis cup
- Evaporated in fume hood
- Loaded onto autosampler for Pyr-GC-HRAM



# Method parameters

## Multi-Shot Pyrolyzer EGA/PY-3030D parameters

Analysis type Double-shot analysis

### Thermal desorption

Initial (°C) 100

Initial (min) 0

Rate (°C /min) 20

Final (°C) 300

Final (min) 1

Total time (min) 11

### Pyrolysis

Initial (°C) 650

Initial (min) 0.2

Interface temperature °C 320

## Trace 1310 GC System parameters

Injector type SSL with an adapter kit for gas injection

Injection mode Split

Temperature (°C) 300

Split ratio 200:1

Carrier gas (mL/min) He, 1

### Oven temperature program

Temperature 1 (°C) 40

Hold time (min) 2

Rate (°C /min) 20


Temperature 2 (°C) 320

Hold time (min) 14

# Method parameters

## Orbitrap Exploris GC 240 MS parameters

Transfer line temperature (°C)	300
Ionization type	EI
Ion source temperature (°C)	280
Electron energy (eV)	70
Emission current (μA)	50
Acquisition mode	Full scan
Mass range ( <i>m/z</i> )	40–600
Resolving power setting	60,000
Lock masses ( <i>m/z</i> )	133.01356; 207.03235; 225.04292; 281.05114; 299.06171; 355.06993





Thermo Scientific™ Orbitrap Exploris™ GC 240


**Orbitrap Exploris GC 240 with Frontier multi shot pyrolizer installed on GC**

# Interferences

## Single quadrupole vs Orbitrap







Single quad selectivity




Orbitrap GC selectivity

# Chromatography and recovery




Spiked fishmeal  
Calculated amount  
= 2.2 mg PMMA  
(spiked at 2.5 mg)



# Pyrolysis products


Pattern diagram of  
typical polymeric material



● : Additives

— : Polymer

● : Inorganics




X: Additives

Y: Polymer

Z: Inorganic




Pyrolysis (Single-Shot)



Mixed Information of polymer and volatiles

**Sometimes difficult to interpret results !!**

# Total ion chromatogram vs blank after TD




Total ion current chromatogram (m/z 40–600) obtained for a milk sample (black chromatogram) compared with a solvent standard of a mix of polymers (red chromatogram) after the thermal desorption (TD) step.

# Polymers found in real samples


Polymer	Pyrolysis products
Polystyrene (PS)	Styrene; styrene dimer; styrene trimer; allylbenzene; $\alpha$ -methylstyrene; toluene
Polypropylene (PP)	2,4-dimethyl-1-heptane; 3-5-dimethyl-1-hexane
Polyvinyl chloride (PVC)	Benzene, naphthalene, fluorene
Polymethyl methacrylate (PMMA)	Methyl methacrylate
Polycarbonate (PC)	Bisphenol A
Polyethylene terephthalate (PET)	Vinyl benzoate

# High resolution spectra

Matching pyrolysis products with NIST 2020 & HRAM contaminants library



- NIST nominal mass library




- Internal Library




# Milk samples

## Standard




## Milk sample




# Environmental samples

## Standard



## Storm water sample



# Conclusions

## Targeted Analysis:

- Accelerated solvent extraction can extract a variety of MPs/NPs including PS, PP, PE, PP, PC, and PVC from rice samples.
- Reduces processing and labor time needed to prepare MP/NP particles from samples
- Effective extraction of MPs/NPs from a wide range of environmental samples
- Rice samples that were found to contain between 17 and 317  $\mu\text{g/g}$  dw of MPs/NPs and this approach has been applied elsewhere

## Untargeted Analysis:

- Positive confirmation of the presence and identity of microplastics in different sample types.
- High selectivity and sensitivity were achieved by using the unique characteristics of the Orbitrap mass spectrometer, in combination with a targeted screening
- Combination of automated sample prep and automated analysis using the pyrolyzer and targeted data processing enables an automated analysis of environmental samples.



# Acknowledgements

## University of Queensland

- Dr. Kevin Thomas, Dr. Jake O'Brien, Dr. Elvis Okoffo, and Dr. Cassie Rauert
- An extra thank you to Dr Cassie Rauert for supplying milk and meat samples



## NILU

- Dr. Vladimir Nikiforov and Dr. Dorte Herzke



**Thank you!**  
**Are there any questions?**

 The world leader in serving science

