

Analysis of Semivolatile compounds by EPA 8270 with Shimadzu GCMS Triple Quadrupole Mass Spectrometer

Yoshiro Hiramatsu

GCMS Product Liaison

Shimadzu Scientific Instrument

A decorative graphic element consisting of a red curved line that starts at the bottom left and curves upwards and to the right, ending at the top right. The line is surrounded by a soft, grey-to-white gradient shadow.


Today's presentation

- **What are Semivolatile compounds?**
- **GC Triple quadrupole Mass Spectrometer**
- **SVOC application**
- **Non-Target analysis in Triple quadrupole**
- **Conclusion**
- **Q&A**


What are Semivolatile compounds?

SVOCs (**SemiVolatile Organic Compounds**)


Phenols




Nitro-Aromatics



Phthalates



PAHs



Others (Organic acids, Ketones, Amines, etc.)

How to analyze SVOCs?

GC/MS technology is most used to acquire chromatograms, identify and quantitate for SVOCs



Search EPA.gov

Environmental Topics ▾ Laws & Regulations ▾ Report a Violation ▾ About EPA ▾

Environmental Sampling and Analytical Methods (ESAM) Program [CONTACT US](#)

- ESAM Home
- Sample Collection Information Documents (SCIDs)
- Sample Collection Procedures and Strategies
- Selected Analytical Methods (SAM)
- Analytical Methods and Protocols

EPA Method 8270E (SW-846): Semivolatile Organic Compounds by Gas Chromatography/Mass Spectrometry (GC-MS)

This document is included in [Selected Analytical Methods for Environmental Remediation and Recovery \(SAM\)](#).



EPA methods: Update SW-846

GC-MS/MS detector was allowed to use for the analysis of SVOCs by EPA 8270E, updated in 2018


6.1.3.3 An MS/MS detector may be used if the detector has the necessary pumps, collision cell, collision gases, and high-vacuum system capable of performing transitions in product ion scan mode or the selected reaction monitoring mode (SRM) for the target analytes of interest. Recommendations for specific precursor and product ions in SRM are available for some target analytes from the manufacturers of the equipment. When analysis is performed using product ions for quantitation, it is not an appropriate verification of the system to perform DFTPP analysis and meet the criteria outlined in Sec. 11.3.1. The system, however, must be capable of documenting the performance of both MSs against manufacturer specifications for mass resolution, mass assignment, and sensitivity using the internal calibrant (e.g., Perfluorotributylamine). The performance of the system should be checked at least weekly, or at a frequency appropriate to meet the needs of the project. At a minimum, the performance of the system must be checked just prior to the initial calibration (ICAL).

Today's presentation

- **What are Semivolatile compounds?**
- **GC Triple quadrupole Mass Spectrometer**
- **SVOC application**
- **Non-Target analysis in Triple quadrupole**
- **Conclusion**
- **Q&A**


Triple quadrupole GCMS (GC-MS/MS)

GC-MS/MS can achieve great selectivity using Selected Reaction Monitoring, also referred to as Multiple Reaction Monitoring (MRM)






What are the benefits of GC-MS/MS?

GC-MS/MS can achieve great selectivity using MRM



What are the benefits of GC-MS/MS?

GC-MS/MS can perform high sensitivity Scan analysis with MRM analysis



Today's presentation

- **What are Semivolatile compounds?**
- **GC Triple quadrupole Mass Spectrometer**
- **SVOC application**
- **Non-Target analysis in Triple quadrupole**
- **Conclusion**
- **Q&A**

SVOC Application following 8270E

GC and MS conditions

| GC conditions | |
|-------------------|---|
| Injection Temp. | 275 ° C |
| Injection Mode | Split (1:10) |
| Column Flow | 1.2 mL/min |
| Flow Control Mode | Linear Velocity |
| Column | Rxi-5MS (30m x 0.25 mm, 0.25 µm) |
| Oven Temp. | 40 ° C (1 min) ->20 ° C/min to 280 ° C -> 5 ° C to 320 ° C (1 min) Total run: 22 min |
| Carrier Gas | Helium |
| MS conditions | |
| Solvent Cut Time | 2 min |
| Ion Source Temp. | 200 ° C |
| Interface Temp. | 320 ° C |
| Detector Voltage | 1.3 kV (Absolute) |

Compound list and MSMS parameter


| Serial# | Type | Acq. Mode | ISTD Group | Level1 Conc (1s) | Method No. | Compound Name (E) | Ret. Index 1 | Ret. Index 2 | Ret. Index 3 | Ret. Time | Cl |
|---------|--------|-----------|------------|------------------|------------|---|--------------|--------------|--------------|-----------|----|
| 1 | Target | MRM | 1 | | 1 | 1,4-Dioxane | 718 | | | 2.353 | |
| 2 | Target | MRM | 1 | | 1 | N-Nitrosodimethylamine | 745 | | | 2.603 | |
| 3 | Target | MRM | 1 | | 1 | Pyridine | 750 | | | 2.649 | |
| 90 | Target | MRM | 1 | | 1 | 2-Fluorophenol | 876 | | | 3.817 | |
| 4 | Target | MRM | 1 | | 1 | Benzaldehyde | 970 | | | 4.667 | |
| 5 | Target | MRM | 1 | | 1 | Phenol | 979 | | | 4.747 | |
| 91 | Target | MRM | 1 | | 1 | Phenol-d6 | 979 | | | 4.747 | |
| 6 | Target | MRM | 1 | | 1 | Aniline | 983 | | | 4.783 | |
| 7 | Target | MRM | 1 | | 1 | Bis(2-chloroethyl)ether | 990 | | | 4.846 | |
| 8 | Target | MRM | 1 | | 1 | 2-Chlorophenol | 997 | | | 4.909 | |
| 9 | Target | MRM | 1 | | 1 | 1,3-Dichlorobenzene | 1017 | | | 5.078 | |
| 10 | Target | MRM | 1 | | 1 | 1,4-Dichlorobenzene | 1026 | | | 5.154 | |
| 84 | ISTD | MRM | 1 | | 1 | 1,4-Dichlorobenzene-d4 | 1026 | | | 5.154 | |
| 11 | Target | MRM | 1 | | 1 | Benzyl Alcohol | 1041 | | | 5.280 | |
| 12 | Target | MRM | 1 | | 1 | 1,2-Dichlorobenzene | 1046 | | | 5.322 | |
| 14 | Target | MRM | 1 | | 1 | 2-Methylphenol(o-Cresol) | 1055 | | | 5.397 | |
| 13 | Target | MRM | 1 | | 1 | Bis(2-chloro-1-methylethyl) ether | 1059 | | | 5.431 | |
| 15 | Target | MRM | 1 | | 1 | N-nitroso-di-n-propylamine | 1076 | | | 5.574 | |
| 16 | Target | MRM | 1 | | 1 | 3-Methylphenol(m-Cresol)/4-Methylphenol(p-Cresol) | 1076 | | | 5.574 | |
| 17 | Target | MRM | 1 | | 1 | Hexachloroethane | 1091 | | | 5.700 | |
| 18 | Target | MRM | 2 | | 1 | Nitrobenzene | 1099 | | | 5.767 | |
| 92 | Target | MRM | 2 | | 1 | Nitrobenzene-d5 | 1099 | | | 5.767 | |
| 19 | Target | MRM | 2 | | 1 | Isophorone | 1134 | | | 6.042 | |
| 20 | Target | MRM | 2 | | 1 | 2-Nitrophenol | 1147 | | | 6.144 | |
| 22 | Target | MRM | 2 | | 1 | 2,4-Dimethylphenol | 1152 | | | 6.183 | |
| 21 | Target | MRM | 2 | | 1 | Benzoic Acid | 1160 | | | 6.246 | |
| 23 | Target | MRM | 2 | | 1 | Bis(2-chloromethoxy)methane | 1172 | | | 6.340 | |
| 24 | Target | MRM | 2 | | 1 | 2,4-Dichlorophenol | 1182 | | | 6.418 | |

Commercially available libraries or database can accelerate method optimization when new compounds need to be added

Quantitation Result

All 83 compounds passed guidance criteria outlined in EPA 8270E

>60% compounds presented %RSD of RF \leq 20% and excellent linearity
with wide calibration ranges (0.01-50ppm or 0.01-20ppm)




Acenaphthene


Quantitation Result

All 83 compounds passed guidance criteria outlined in EPA 8270E

Calibration curve criteria was met for challenging compounds, such as benzoic acid and nitro-group containing compounds, from 0.2 to 50 ppm when using a quadratic fitting




Benzoic Acid




Quantitation Result

At Lowest calibration point good peak shape and good repeatability were obtained


1,3-Dichlorobenzene
(%RSD=3.39%)



Diphenylamine
(%RSD=4.20%)



Pyrene
(%RSD=3.96%)




Chromatogram of 10 ppb standard (n=6)

Quantitation Result


Good repeatability was also achieved for challenging compounds at lowest calibration point

2-Nitrophenol
(%RSD=19.67%)




20 ppb Chromatogram (n=6)

2-Nitroaniline
(%RSD=12.39%)



20 ppb Chromatogram (n=6)

Benzoic Acid
(%RSD=16.9%)





200 ppb Chromatogram (n=6)

Qualification Result


Scan/MRM simultaneous analysis support unknown compound identification

Black:
Total Scan ion chromatogram
Other color:
MRM chromatogram



Library Search

Propyzamide




Application Result

- **All 83 compounds pass calibration curve criteria outlined 8270E with wide range**
- **All 83 compounds show good chromatography and repeatability using MSMS**
- **Unknown peaks can be identified using library search with MRM analysis**

Today's presentation


- What are Semivolatile compounds?
- GC Triple quadrupole Mass Spectrometer
- SVOC application
- **Non-Target analysis in Triple quadrupole**
- Conclusion
- Q&A

Non-Target analysis in Triple quadrupole



AIQS-DB (Automated Identification and Quantification System)

Theory



Sample prep

- Sample 1 L
 - ← Phosphate buffer (1 mol/L, pH 7.0, 1 mL)
 - ← Sodium chloride (50 g)
- Liquid-liquid extraction
 - Shaking with 50 mL dichloromethane for 10 min (twice)
- Dehydration (anhydrous sodium sulfate)
- Concentration (5 mL, rotary evaporator)
 - ← 1 mL hexane
- Concentration (1 mL, nitrogen gas)
 - ← Internal standard solution (10 µg/mL, 100 µL)
- Measurement sample

| Category 1 Substances | Number of Substances | Category 2 Substances | Number of Substances |
|---------------------------------|----------------------|-----------------------|----------------------|
| CH-containing substances | 194 | Polycyclic aromatic | 79 |
| | | PCBs | 62 |
| | | Other | 53 |
| Oxygen-containing compounds | 150 | Phenols | 50 |
| | | Other | 100 |
| Nitrogen-containing compounds | 113 | Aromatic amines | 43 |
| | | Nitrogen compounds | 42 |
| | | Other | 28 |
| Sulfur-containing compounds | 12 | | 12 |
| Phosphorus-containing compounds | 8 | | 8 |
| | | PPCPs | 14 |
| Pesticides | 451 | Insecticides | 184 |
| | | Herbicides | 118 |
| | | Fungicides | 116 |
| | | Other | 33 |
| Total | 942 | | |
| Internal Standards | 8 | | |

Evaluation of AIQS-DB

Using 4 GCMS at 4 location (Column: DB-5ms 30mx0.25mm i.d. 0.25um)

| Reproducibility | %RSD (Number of Substances) |
|-----------------|-----------------------------|
| Under 10% | 41%(47) |
| 10-20% | 37%(42) |
| Over 20%(~70%) | 22%(25) |

Many of the substances having no or one highly polar functional group.

Polycyclic aromatic hydrocarbons with high boiling point, substances with relatively high water solubility containing a high polarity function group


At least 2 high-polarity functional groups such as 2,6-diaminotoluene, m-aminophenol, etc.



When we try to improve those compounds, we should consider about polar column. But as screening purpose, this method is acceptable.


Demonstrate the use of AIQS-DB

TIC of River water sample (Q3 Scan)




Detected pollutants example


n-C₉H₂₀(Nonane)




Naphthalene




Bis(2-ethylhexyl)phthalate




Event#1: Q3 Scan Ret.Time: [6.698 > 6.702] - [6.665 <> 6.735] Scan#: [680 > 682] - [667 <> 695] Ret.Index: 899!




ID#10 Standard Spectrum




Event#1: Q3 Scan Ret.Time: [12.610 > 12.620] - [12.530 <> 12.680] Scan#: [3045 > 3049] - [3013 <> 3073] Ret.Index: 1189




ID#95 Standard Spectrum



Event#1: Q3 Scan Ret.Time: [30.645 > 30.655] - [30.555 <> 30.810] Scan#: [22058 > 22070] - [21950 <> 22256] Ret.Index: 2530




ID#811 Standard Spectrum




Detected pollutants example





Butoxyethanol



Quinoline



Benzothiazole



Demonstrate Results

Semi-Quantitative Results for detected environmental pollutants using AIQS-DB

| Categ | Detected Compounds | Conc. (ng/L in water) | Categ | Detected Compounds | Conc. (ng/L in water) | Categ | Detected Compounds | Conc. (ng/L in water) |
|-------|-------------------------|--------------------------|-------|-------------------------------|--------------------------|-------|---------------------------------------|--------------------------|
| 1 | n-C11H24 | 62.7 | 1 | Biphenyl | 160.7 | 2 | Coprostanol | 14121.5 |
| 1 | n-C12H26 | 120.1 | 1 | Fluorene | 44.6 | 2 | Cyclohexanol | 899.7 |
| 1 | n-C13H28 | 241.3 | 1 | Naphthalene | 270.4 | 2 | Ethanol, 2-phenoxy- | 1989.0 |
| 1 | n-C14H30 | 190.5 | 1 | Phenanthrene | 126.4 | 2 | Phenylethyl alcohol | 101.2 |
| 1 | n-C15H32 | 188.6 | 2 | Diphenyl ether | 213.6 | 2 | Stigmasterol | 1697.2 |
| 1 | n-C16H34 | 215.3 | 2 | Isophorone | 8601.1 | 2 | 1,3-Dichloro-2-propanol | 304.5 |
| 1 | n-C17H36 | 419.8 | 2 | 2-Methylphenol | 506.7 | 3 | 2-Naphthylamine | 93.4 |
| 1 | n-C18H38 | 420.6 | 2 | 2-Naphthol | 126.0 | 3 | Acetamide, N-phenyl- | 571.0 |
| 1 | n-C19H40 | 435.1 | 2 | 2-Phenylphenol | 52.4 | 3 | 2-Chloroaniline | 533.5 |
| 1 | n-C20H42 | 541.7 | 2 | 4-Methyl-2,6-di-t-butylphenol | 167.6 | 3 | 3,4-Dichloroaniline | 959.2 |
| 1 | n-C21H44 | 677.3 | 2 | 4-tert-Octylphenol | 155.4 | 3 | Quinoline | 175.4 |
| 1 | n-C22H46 | 794.4 | 2 | Bisphenol A | 3169.6 | 4 | 2(3H)-Benzothiazolone | 2544.3 |
| 1 | n-C23H48 | 1118.2 | 2 | Nonylphenol | 9845.0 | 4 | 2-(Methylthio)-benzothiazol | 277.7 |
| 1 | n-C24H50 | 1505.7 | 2 | 2,4,5-Trichlorophenol | 38.1 | 4 | 2-Acetyl-5-methylthiophene | 65.4 |
| 1 | n-C25H52 | 2383.4 | 2 | 2,4,6-Tribromophenol | 3292.6 | 4 | Benzothiazole | 177.0 |
| 1 | n-C26H54 | 1299.5 | 2 | 2,4,6-Trichlorophenol | 44.7 | 5 | Tris(1,3-dichloro-2-propyl) phosphate | 677.0 |
| 1 | n-C27H56 | 1392.3 | 2 | Triclosan | 211.3 | 6 | Caffeine | 1736.5 |
| 1 | n-C28H58 | 966.2 | 2 | Bis(2-ethylhexyl)phthalate | 11036.3 | 6 | Diethyltoluamide | 345.0 |
| 1 | n-C29H60 | 1445.2 | 2 | Diethyl phthalate | 2306.6 | 6 | Ibuprofen | 1596.3 |
| 1 | n-C30H62 | 1625.2 | 2 | Diisobutyl phthalate | 914.1 | 6 | L-Menthol | 4281.0 |
| 1 | n-C31H64 | 2393.8 | 2 | Dimethyl phthalate | 261.8 | 6 | Thymol | 223.1 |
| 1 | n-C32H66 | 979.3 | 2 | Di-n-butyl phthalate | 807.1 | 6 | Nicotine | 846.2 |
| 1 | n-C33H68 | 863.2 | 2 | 2-Butoxyethanol | 4013.9 | 7 | Fenobucarb | 551.5 |
| 1 | n-C9H20 | 152.0 | 2 | 2-Ethyl-1-hexanol | 1628.3 | 7 | Permethrin 1 | 661.7 |
| 1 | 4-Cymene | 351.2 | 2 | alpha-Terpineol | 900.9 | 7 | Permethrin 2 | 114.2 |
| 1 | 1,3-Dimethylnaphthalene | 618.1 | 2 | Benzyl alcohol | 203.7 | 7 | Piperonyl butoxide | 85.7 |
| 1 | 2,6-Dimethylnaphthalene | 403.9 | 2 | beta-Sitosterol | 3799.5 | 7 | 2-Phenylphenol (OPP) | 63.9 |
| 1 | 2-Methylnaphthalene | 100.9 | 2 | Cholesterol | 14408.1 | 7 | Biphenyl | 34.6 |

Non-target and target analysis with GC-MS/MS




Using 1 method we can acquire both data of non-target result and accurate quant result following 8270E

Create Calibration curve for Target analysis using MRM

Running sample with Scan/MRM method

Report semi-quant result using Scan Data

Report quant result using MRM data



Conclusion




GC-MSMS is Ready to use solution for SVOC analysis following 8270E guidance

GC-MSMS can do both quantitation for analyte and qualification for unknown compounds at the same time

GC-MSMS can do non-target semi-quant and accurate target analysis at the same time

Q&A



For any questions, contact:
Yoshiro Hiramatsu
yohiramatsu@shimadzu.com

For more information, visit:
www.OneLabOneEarth.com

Connect with us:

Twitter - [@shimadzussi](https://twitter.com/shimadzussi)

Instagram - [@shimadzussi](https://www.instagram.com/shimadzussi)

LinkedIn - [/company/shimadzu-scientific-instruments/](https://www.linkedin.com/company/shimadzu-scientific-instruments/)

YouTube - [@ShimadzuScientificInstruments](https://www.youtube.com/@ShimadzuScientificInstruments)