


Applying new ion chromatography technology to enable faster determinations of inorganic anions and disinfection byproducts in drinking water

Jingli Hu, Hua Yang, Jeff Rohrer, Chris Shevlin, and **Carl Fisher**

August 4, 2023

 The world leader in serving science



Agenda

1 Drinking water analysis

2 Inorganic anion determinations

3 Shortening runtime

4 IC innovations

5 Resources



Drinking water analysis

- National Primary Drinking Water Regulations (NPDWR)
 - Microorganisms
 - Disinfectants
 - Disinfection byproducts
 - Inorganic chemicals
 - US EPA 300
 - Organic chemicals
 - Radionuclides



Drinking water analysis

- National Primary Drinking Water Regulations (NPDWR)
 - Microorganisms
 - Disinfectants
 - Disinfection byproducts
 - Inorganic chemicals
 - US EPA 300
 - Organic chemicals
 - Radionuclides



EPA Method 300.0 and 300.1


Outlines the method for determination of inorganic anions by ion chromatography (IC).

- Specifies use of suppressed conductivity for determination of:

- Applies to:
 - Drinking water
 - Ground and Surface water
 - Wastewater (domestic and industrial)
 - Raw water (unfinished drinking water)

Part A- Common Anions	Part B- Inorganic Disinfection Byproducts
Bromide	Bromide
Chloride	Bromate
Fluoride	Chlorite
Nitrate	Chlorate
Nitrite	
ortho-Phosphate-P	
Sulfate	


IC system



Inorganic ion determination


Thermo Scientific™ Dionex™ IonPac™ AS4A column

1.7 mM Na₂CO₃ / 1.8 mM NaHCO₃



Inorganic ion determination

Thermo Scientific™ Dionex™ IonPac™ AS4A column
1.7 mM Na₂CO₃ / 1.8 mM NaHCO₃



Carbonate-based eluents

Benefits

- Strong elution potential
- Inexpensive

Limitations


- Conductivity is 10–20x higher than pure water which impacts sensitivity
- Response is nonlinear at lower analyte concentrations
- Gradient separations are difficult
 - Ramping, baseline disturbances, slow buffering

Solution: Hydroxide-based eluent

Inorganic ion determination

Thermo Scientific™ Dionex™ IonPac™ AS4A column

1.7 mM Na₂CO₃ / 1.8 mM NaHCO₃



Carbonate-based eluents

Benefits

- Strong elution potential
- Inexpensive


Limitations

- Conductivity is 10–20x higher than pure water which impacts sensitivity
- Response is nonlinear at lower analyte concentrations
- Gradient separations are difficult
 - Ramping, baseline disturbances, slow buffering

Solution: Hydroxide-based eluent


Inorganic ion determination

Thermo Scientific™ Dionex™ IonPac™ AS18 column
22–40 mM KOH



Inorganic ion determination

Thermo Scientific™ Dionex™ IonPac™ AS18 column
22–40 mM KOH



Hydroxide-based eluents

Benefits


- Suppressed to water, yielding low background, greater sensitivity
- Produces linear gradients

Limitations of manual eluent preparation

- Absorbance of CO_2 from the air causes instability of
 - Retention time, peak area, baseline
- For gradients, a proportioning pump is required

Inorganic ion determination

Thermo Scientific™ Dionex™ IonPac™ AS18 column
22–40 mM KOH



Hydroxide-based eluents

Benefits

- Suppressed to water, yielding low background, greater sensitivity
- Produces linear gradients

Limitations of manual eluent preparation

- Absorbance of CO_2 from the air causes instability of
 - Retention time, peak area, baseline
- For gradients, a proportioning pump is required

Solution: Electrolytic eluent generation (Reagent-free IC (RFIC))

- Precise and accurate eluent concentrations
- An isocratic pump can be used for gradients

Automated eluent generation (EG)

ONE current for
ONE concentration
(of eluent)



Deionized water used as the carrier


- Eliminates need to handle acids or bases
- Decreases pump maintenance

Electrolysis automatically generates high purity acid, base, or carbonate eluents in-line

Delivers concentrations specified in Chromeleon CDS through control of electrical **current**


Gradients created without need for a more costly, proportioning pump

Gradient reproducibility: Thermo Scientific Dionex IonPac AS11 4 mm (n=30)




	Fluoride	Chloride	Nitrite	Sulfate	Bromide	Nitrate	Phosphate
Area (µS·min)	0.110	0.324	0.236	0.137	0.242	0.182	0.217
% RSD Area	0.059	0.075	0.075	0.130	0.061	0.101	0.049
Ret. Time (min)	2.400	5.297	5.755	7.223	10.047	7.394	13.139
%RSD Ret.	0.060	0.020	0.029	0.021	0.014	0.023	0.022

Gradient reproducibility: Thermo Scientific Dionex IonPac AS11 4 mm (n=30)



	Fluoride	Chloride	Nitrite	Sulfate	Bromide	Nitrate	Phosphate
Area (µS·min)	0.110	0.324	0.236	0.137	0.242	0.182	0.217
% RSD Area	0.059	0.075	0.075	0.130	0.061	0.101	0.049
Ret. Time (min)	2.400	5.297	5.755	7.223	10.047	7.394	13.139
%RSD Ret.	0.060	0.020	0.029	0.021	0.014	0.023	0.022

Peak focusing using a KOH eluent generator



Separator: Thermo Scientific™ Dionex™ IonPac™ AG11, AS11, 4-mm

Eluent: A: 15.5 mmol/L KOH
B: 0.5 to 25 mmol/L KOH in 8 min

Flow rate: 2 mL/min

Inj. volume: 25 µL


Detection: Suppressed conductivity

Suppressor: ASRS, AutoSuppression

Recycle Mode

Peaks:	1. Fluoride	2 [mg/L]
	2. Chloride	3
	3. Nitrate	10
	4. Sulfate	15
	5. Orthophosphate	15

Peak focusing using a KOH eluent generator



Separator: Thermo Scientific™ Dionex™ IonPac™ AG11, AS11, 4-mm

Eluent: A: 15.5 mmol/L KOH
B: 0.5 to 25 mmol/L KOH in 8 min

Flow rate: 2 mL/min

Inj. volume: 25 µL

Detection: Suppressed conductivity


Suppressor: ASRS, AutoSuppression

Recycle Mode


Peaks:	1. Fluoride	2 [mg/L]
	2. Chloride	3
	3. Nitrate	10
	4. Sulfate	15
	5. Orthophosphate	15

Ease of use and consistency of RFIC, but how do we shorten run time?

IC column parameters: Speed, capacity, and resolution




IC column parameters: Speed, capacity, and resolution



Column efficiency: Increase with smaller particles (4 μm)

Smaller Particle Columns:



Produce more efficient peaks


Impact chromatographic speed and resolution

Improve peak integration – more accurate and reliable results

Increase sample throughput without compromising data quality

Improve quality of analytical results

4 μm columns: Combining speed and resolution



150mm

- Faster run times
- Higher flow rates

250mm

- Better resolution
- Standard flow rates

Achieving the optimal combination of speed and resolution

4 μm columns: Combining speed and resolution

150mm

- Faster run times
- Higher flow rates



250mm

- Better resolution
- Standard flow rates

Achieving the optimal combination of speed and resolution

Need High Pressure IC (HPIC) system to take full advantage of benefits

Column efficiency: 4 μm vs. 6.5 μm particles



Column A: Thermo Scientific™ Dionex™ IonPac™ AS22-Fast (2 x 150 mm)

Column B: Thermo Scientific™ Dionex™ IonPac™ AS22-Fast-4 μm (2 x 150 mm)

Eluent: 4.5 mM Na₂CO₃ / 1.4 mM NaHCO₃

Flow Rate: 0.3 mL/min



Inj. Volume: 2.5 μL

Temperature: 30°C

Detection: Suppressed Conductivity, Thermo Scientific™ Dionex™ AERS 500 (2 mm) Suppressor, AutoSuppression, Recycle mode

Peaks:	mg/L (ppm)
1. Fluoride	5.0
2. Chloride	10.0
3. Nitrite	15.0
4. Bromide	25.0
5. Nitrate	25.0
6. Phosphate	40.0
7. Sulfate	30.0

Column efficiency: 4 μm vs. 6.5 μm particles




Column A: Thermo Scientific™ Dionex™ IonPac™ AS22-Fast (2 x 150 mm)
 Column B: Thermo Scientific™ Dionex™ IonPac™ AS22-Fast-4 μm (2 x 150 mm)
 Eluent: 4.5 mM Na₂CO₃ / 1.4 mM NaHCO₃
 Flow Rate: 0.3 mL/min
 Inj. Volume: 2.5 μL
 Temperature: 30°C
 Detection: Suppressed Conductivity, Thermo Scientific™ Dionex™ AERS 500 (2 mm) Suppressor, AutoSuppression, Recycle mode

Peaks:	mg/L (ppm)
1. Fluoride	5.0
2. Chloride	10.0
3. Nitrite	15.0
4. Bromide	25.0
5. Nitrate	25.0
6. Phosphate	40.0
7. Sulfate	30.0

Smaller particles deliver more efficient peaks and higher resolution.


Faster runs using smaller particle columns



Thermo Scientific™ Dionex™
IonPac™ AS18-Fast-4µm column
(4 x 150 mm)
15–44 mM KOH

Peaks:	A	B	
1. Fluoride	1	0.63	mg/L
2. Chloride	12	9.7	
3. Nitrite	8	0.023	
4. Carbonate	NA	NA	
5. Bromide	9	NA	
6. Sulfate	10	16.3	
7. Nitrate	9	1.3	
8. Phosphate	6.5	0.085	

Faster runs using smaller particle columns




Thermo Scientific™ Dionex™
IonPac™ AS18-Fast-4 μm column
(4 x 150 mm)
15–44 mM KOH

Peaks:	A	B	
1. Fluoride	1	0.63	mg/L
2. Chloride	12	9.7	
3. Nitrite	8	0.023	
4. Carbonate	NA	NA	
5. Bromide	9	NA	
6. Sulfate	10	16.3	
7. Nitrate	9	1.3	
8. Phosphate	6.5	0.085	

Run time reduced by 25%

Fast determinations using carbonate eluent



Column:

Thermo Scientific™ Dionex™ IonPac™
AS22-Fast-4µm, 4 mm


Eluent:

4.5 mM Na₂CO₃ / 1.4 mM NaHCO₃

Flow rate:

2 mL/min

Fast determinations using carbonate eluent



Column:

Thermo Scientific™ Dionex™ IonPac™
AS22-Fast-4µm, 4 mm

Eluent:


4.5 mM Na₂CO₃ / 1.4 mM NaHCO₃

Flow rate:

2 mL/min

Baseline resolution in under five minutes.

Fast anion determinations: drinking water




Sample: Municipal drinking water

Peaks:


	Min	mg/L
1. Fluoride	1.2	0.6
2. Chloride	1.6	9.5
3. Nitrite	2.0	0.2
4. Unknown	—	—
5. Nitrate	2.6	1.2
6. Phosphate	3.5	0.2
7. Sulfate	4.0	16.5

Thermo Scientific Dionex IonPac anion exchange column development


Carbonate-based



Hydroxide selective




Faster haloacetic acid determinations




Peaks (Standard):	mg/L
1. Monochloroacetate	1.0
2. Monobromoacetate	1.0
3. Bromate	1.0
4. Chloride	316
5. Carbonate	150
6. Sulfate	250
7. Dalapon	1.0
8. Dichloroacetate	1.0
9. Bromochloroacetate	1.0
10. Nitrite	0.25
11. Dibromoacetate	1.0
12. Nitrate	20.0
13. Trichloroacetate	1.0
14. Bromodichloroacetate	1.0
15. Chlorodibromoacetate	1.0
16. Tribromoacetate	1.0

Faster haloacetic acid determinations




Peaks (Standard):	mg/L
1. Monochloroacetate	1.0
2. Monobromoacetate	1.0
3. Bromate	1.0
4. Chloride	316
5. Carbonate	150
6. Sulfate	250
7. Dalapon	1.0
8. Dichloroacetate	1.0
9. Bromochloroacetate	1.0
10. Nitrite	0.25
11. Dibromoacetate	1.0
12. Nitrate	20.0
13. Trichloroacetate	1.0
14. Bromodichloroacetate	1.0
15. Chlorodibromoacetate	1.0
16. Tribromoacetate	1.0

Faster haloacetic acid determinations



Peaks (Standard):	mg/L
1. Monochloroacetate	1.0
2. Monobromoacetate	1.0
3. Bromate	1.0
4. Chloride	316
5. Carbonate	150
6. Sulfate	250
7. Dalapon	1.0
8. Dichloroacetate	1.0
9. Bromochloroacetate	1.0
10. Nitrite	0.25
11. Dibromoacetate	1.0
12. Nitrate	20.0
13. Trichloroacetate	1.0
14. Bromodichloroacetate	1.0
15. Chlorodibromoacetate	1.0
16. Tribromoacetate	1.0

Faster haloacetic acid determinations



Peaks (Standard):	mg/L
1. Monochloroacetate	1.0
2. Monobromoacetate	1.0
3. Bromate	1.0
4. Chloride	316
5. Carbonate	150
6. Sulfate	250
7. Dalapon	1.0
8. Dichloroacetate	1.0
9. Bromochloroacetate	1.0
10. Nitrite	0.25
11. Dibromoacetate	1.0
12. Nitrate	20.0
13. Trichloroacetate	1.0
14. Bromodichloroacetate	1.0
15. Chlorodibromoacetate	1.0
16. Tribromoacetate	1.0

Saving time by altering column specificity

Thermo Scientific Dionex ion chromatography family

Reagent-Free Ion Chromatography (RFIC)

High Pressure Ion Chromatography (HPIC)



Thermo Scientific™
Dionex™ Easion™
- 2 mm/4 mm IC
- Carbonate or MSA eluents
- **Chemical suppression**



Thermo Scientific™
Dionex™ Aquion™
- 2 mm/4 mm IC
- Carbonate or MSA eluents
- **Electrolytic suppression**



Thermo Scientific™
Dionex™ Integrion™
- 2 mm/4 mm HPIC
- Field upgradable options
- **RFIC** or manual eluents
- Tablet control
- Thermo Scientific™ Viper™ fittings



Thermo Scientific™
Dionex™ ICS-6000
- 2 mm/4 mm HPIC
- RFIC, dual EGC and gradient pump
- **Flexible & modular** with tablet
- Versatile detection capabilities
- Performance monitoring/ tracking



- Reagent-free IC
 - Reproducibility
 - Ease of use
 - Reduced exposure to chemicals



- Column Chemistry
 - Capacity
 - Tolerance to matrix loading
 - Resolution
 - Reduced run times
 - Specificity



- Alternative detectors
 - Broaden the range of applications

Extended application range: optional detectors



Conductivity

- Anions and cations



Electrochemical

- Electroactive compounds
(e.g. saccharides)



Spectrophotometric

- UV and visible absorbing compounds
- Post-column and pre-column derivatization techniques
(e.g. Perchlorate)



Mass spectrometry

- Ionic and polar compounds
(e.g. Haloacetic acids (HAAs))

Inductively coupled plasma spectroscopy

- Multiple elemental species
(e.g. Bromate, selenium)



Future directions in IC

- Columns
 - Selectivity, capacity, resolution
- Combining other systems with IC
 - Combustion IC; IC-MS
- Component layout
 - Unobstructed accessibly; tubing management
- Upgradability and customization
 - Base that can easily grow as needs expand
- Software
 - Usability; automation, diagnostics
- Online Support
 - Dedicated to post-sales resources



AppsLab: the latest applications

The screenshot shows the ThermoFisher AppsLab Library website. At the top, there is a navigation bar with the ThermoFisher logo and 'AppsLab Library' text. Below this is a search bar with the placeholder text 'enter compound, matrix or instrument type'. A central banner features three icons: a magnifying glass for 'find a method', a download icon for 'download 1-click workflow', and a chromatogram for 'run, process, report'. The main content area is divided into sections: 'AppsLab Library: Find your Methods, eWorkflows and more' with a descriptive paragraph; 'Latest Methods' featuring three method cards with titles like 'AN72680: Determination of Zinc Oxide in Sunscreen Using Ion Chromatography with Visible Absorbance Detection' and 'NIBRT: Universal methods for biotherapeutics' with a molecular structure image.

- Fully searchable online, analytical method repository
- Latest applications for **IC**, LC, GC, GC-MS, LC-MS, ICP-MS, ICP-OES and DIA instruments
- Download one-click eWorkflows for use with Thermo Scientific Chromeleon Chromatography Data System (CDS) software
- Ever-expanding database of field-tested workflows

www.appslab.com


Resources for IC information



IC Hub

- What's new
- Best practices
- Tips and tricks


thermofisher.com/iceducation



Analyte Guru

- Scientific communities
- Community support
- Knowledgebase
- Blogs
- Events

analyteguru.com



Dionex IC Products

- Systems
- Detectors
- Autosamplers
- Columns
- Consumables

thermofisher.com/IC

Thank you!

Any questions?

