


Polyfluoroalkyl substances (PFAS) in the environment

**An investigation of contamination and recovery of PFAS
in analytical methods that require filter membranes**

Lindsay D. Lozeau, Ph.D., Senior Scientist
MilliporeSigma | 400 Summit Dr. | Burlington MA, 01803

Amy Laws, Ph.D., Senior Portfolio Lead
Maricar Dube, Ph.D., Technical Marketing
MilliporeSigma | 400 Summit Dr. | Burlington MA, 01803




MilliporeSigma is the U.S. and
Canada Life Science business of
Merck KGaA, Darmstadt, Germany.



PFAS (Per- and Polyfluoroalkyl Substances) Perspective on Forever Chemicals

- PFAS is a group of **>4,000** very stable synthetic chemicals [1]
- **Impact on human health** includes cancer, liver damage, weakened immune system, high cholesterol [2]
- Report of **PFAS** found **in the blood of 97% Americans** [1]
- **Stockholm Convention regulations** for **PFOS, PFOA** (subset of PFAS compounds). Effective in **2020** [3]
 - **Exemptions:** medical devices, active negotiations for other applications
- EPA's drinking water health advisory level is **70 parts per trillion** for PFAS
- **EPA** recommends testing wastewater for PFAS (**2020**)
- **EPA PFAS Strategic Roadmap (2021)** outlines actions for 2021-2024
- **US EPA proposed levels of PFOA and PFOS to 4.0 ppt (2022)**



[1] https://www.niehs.nih.gov/health/materials/perfluoroalkyl_and_polyfluoroalkyl_substances_508.pdf

[2] EU/REACH webinar: <https://echa.europa.eu/-/restriction-of-per-and-polyfluoroalkyl-substances-pfas-under-reach>

[3] Stockholm convention regulations: <http://chm.pops.int/Countries/StatusofRatifications/Amendmentstoannexes/tabid/3486/Default.aspx>



Analytical Methods for PFAS Developed at a Rapid Pace

Regulatory Landscape

Method(s)	Update/Revision	Matrix/Matrices	Sample Prep	Analytical Method
EPA 537.1	Jun. 2020	Drinking water	SPE	LC-MS/MS
EPA 533	Dec. 2019	Drinking water	SPE	LC-MS/MS
SW-846 Method 8327*	Jul. 2021	Non-potable groundwater, surface water, wastewater	SPE, filtration	LC-MS/MS
ASTM D7968-17a	Sep. 2017	Environmental solids	Solvent extraction, filtration	LC-MS/MS
ASTM D7979-19	Sep. 2021	Water matrix (no drinking water)	Solvent extraction, filtration	LC-MS/MS
ISO 25101	Mar. 2019	Drinking water, non-potable water	SPE, solvent elution	LC-MS/MS
ISO 21675	Oct. 2019	Drinking, natural and wastewater	SPE, filtration as needed	LC-MS/MS
FDA C-010.02	Dec. 2021	Foods	QuEChERS, SPE, filtration	LC-MS/MS
OTM-45	Jan. 2021	Air Emissions (stationary sources)	Sampling train: filtration , impingers	LC-MS/MS
EPA Draft 1633	Feb. 2022	Aqueous, soil, biosolids, sediment, tissue	SPE, filtration	LC-MS/MS
EPA Draft 1621*	Apr. 2022	Aqueous matrices	TSS, GAC column cleanup	CIC

- Almost exclusively LC-MS/MS based methods
- SPE & filtration are common sample preparation
- Increased focus on high-particulate matrices

**screening method only*

Abbreviations: SPE = solid phase extraction; TSS = total suspended solids determination; GAC = granular activated carbon; CIC = combustion ion chromatography
Selected methods; does not include all drinking water and international methods



Analytical Methods for PFAS Developed at a Rapid Pace

Regulatory Landscape – **Filtration only**

Method(s)	Sample Prep	Matrix/Matrices	Filter(s) Required by Method
SW-846 Method 8327*	SPE, filtration	Non-potable groundwater, surface water, wastewater	0.2µm polypropylene syringe filter (hydrophilic)
ASTM D7968-17a	Solvent extraction, filtration	Environmental solids	0.2µm polypropylene syringe filter (hydrophilic)
ASTM D7979-19	Solvent extraction, filtration	Water matrix (no drinking water)	0.2µm polypropylene syringe filter (hydrophilic)
ISO 21675	SPE, filtration as needed	Drinking, natural and wastewater	1µm GFF syringe filter (rapid) 1-10µm nylon or GFF (if >2g/L suspended matter)
FDA C-010.02	QuEChERS, SPE, filtration	Foods	0.2µm nylon syringe filter
OTM-45	Sampling train: filtration , impingers	Air Emissions (stationary sources)	GFF or QFF (no pore size listed) membrane filter
EPA Draft 1633	SPE, filtration	Aqueous, soil, biosolids, sediment, tissue	0.2µm nylon syringe filter

- Polypropylene, nylon and glass fiber (GFF) are the most common materials
- Filters should not have detectable levels of PFAS compounds above reporting limit (RL) of method

**screening method only*

Abbreviations: SPE = solid phase extraction; TSS = total suspended solids determination; GAC = granular activated carbon; CIC = combustion ion chromatography; philic = hydrophilic; GFF = glass fiber filter; QFF = quartz fiber filter
Selected methods; does not include all drinking water and international methods



Developing Analytical Methods for PFAS

Major concerns about consumables in sample prep

1: Contamination

- Filter may be manufactured in a place with PFAS
 - *Trace amounts observed in nylon, for example [1,2]*
- Filter is fluorinated / contains fluorinated compounds
 - *Generally, want to avoid PVDF, PTFE, but they are common*
- Solvents and/or containers may be contaminated

2: Adsorption

- Filters have different levels of non-specific binding and polarity
 - *In nylon, flushing filters can help desorb PFAS analytes [3]*
- Tubes, filters, containers all demonstrate time-dependent sorption of various PFAS compounds, leading to recovery losses [3]
 - *Glass fiber & quartz fiber sorption debated [4]*



[1] So, M.K., et al. *Archives of Environmental Contamination and toxicology*, **2006**, 50, 240-8.

[2] Yamashita, N., et al. *Marine Pollution Bulletin*, **2005**, 51(8-12), 658-68.

[3] Lath, S., et al. *Chemosphere*, **2019**, 222, 671-8.

[4] Zhou, J. et al. *Environmental Science: Processes Impacts*, **2021**, 23, 580-7.



So, if filters can be “sticky”, then why filter?

Importance of Sample Preparation

- Without sample preparation, data quality can be less than optimal

Clean Sample





Raw Sample




So, if filters can be "sticky", then why filter? Importance of Sample Preparation

- Without sample preparation, data quality can be less than optimal





- Undissolved particles **mobile phase** can clog and reduce the life of the chromatography column




So, if filters can be "sticky", then why filter?

Importance of Sample Preparation


- Without sample preparation, data quality can be less than optimal



- Undissolved particles **mobile phase** can clog and reduce the life of the chromatography column




- Undissolved particles in a **sample** can do the same




Which Membrane, and Why?


Membrane Morphology & Applications




Polyvinylidene Fluoride (PVDF): *Low binding and fast flow for protein sample prep*




Mixed Cellulose Esters (MCE): *Biologically inert, versatile, smooth and uniform*




Polyethersulfone (PES): *Quick flow and high capacity, asymmetric for high-particulate water samples*




PES Cross-section




Polycarbonate (PC): *preferred for microscopy and cell-based applications*



Polytetrafluoroethylene (PTFE): *Low extractables and high chemical compatibility*



Nylon: *Broad compatibility and commonly used for HPLC*




Nylon with a Glass Fiber Prefilter (HPF): *One-step cleanup of large and small particulates without clogging*




Which Membrane, and Why?


Membrane Morphology & Applications




Polyvinylidene Fluoride (PVDF): *Low binding and fast flow for protein sample prep*




Mixed Cellulose Esters (MCE): *Biologically inert, versatile, smooth and uniform*




Polycarbonate (PC): *preferred for microscopy and cell-based applications*




Polytetrafluoroethylene (PTFE): *Low extractables and high chemical compatibility*




Polyethersulfone (PES): *Quick flow and high capacity, asymmetric for high-particulate water samples*



PES Cross-section



Nylon: *Broad compatibility and commonly used for HPLC*




Nylon with a Glass Fiber Prefilter (HPF): *One-step cleanup of large and small particulates without clogging*



Are Membranes Contaminated with PFAS out-of-box? Testing Syringe Filters using EPA 537.1

Overview of Modified EPA 537.1



LC-MS/MS Conditions

Column:	C18, 100 x 2.1mm ID, 2.7 µm superficially porous particles			
Mobile phase:	[A] DI Water, 0.1% (v/v) acetic acid; [B] Methanol (MeOH), 0.1% (v/v) acetic acid			
Gradient:	Time (min)	A %	B %	Flow (mL/min)
	0-0.0	65%	35%	0.4
	0-7.0	0%	100%	0.4
	7.0-10.0	0%	100%	0.7
	10.0-11.0	0%	100%	0.7
	11.0-15.0	65%	35%	0.4
Flow rate:	See gradient table			
Detection:	MS/MS, ESI(-), details of MS/MS conditions can be requested from the author			
Column temp:	50.0 °C			
Injection volume:	3-5 µL autosampler injection			
Sample :	SPE eluate concentrated to 1 mL methanol: water, 96:4% (v/v)			

Testing Syringe Filters using EPA 537.1

Millex® syringe filters do not have detectable levels of PFAS contaminants

28 PFAS analyzed

Compound(s)	RL (ppb)	MDL (ppb)	Millex® Syringe Filter			
			0.22µm PES ^a	0.45µm PES ^a	0.20µm nylon ^a	0.20µm nylon-HPF ^b
Perfluoroalkyl carboxylic acids – PFBA	0.0040	0.0020	ND – not detected			
Perfluoroalkyl carboxylic acids – PFPeA, PFHxA, PFHpA, PFOA, PFNA, PFDA, PFUnDA, PFDoDA, PFTTrDA, PFTeDA	0.0020	0.0010	ND – not detected			
Perfluoroalkyl sulfonic acids – PFBS, PFPeS, PFHxS, PFHpS, PFOS, PFNS, PFDS	0.0020	0.0010	ND – not detected			
Fluorotelomer sulfonic acids – 4:2 FTS, 6:2 FTS, 8:2 FTS	0.0080	0.0020	ND – not detected			
Perfluorooctane sulfonamides – PFOSA	0.0040	0.0020	ND – not detected			
Perfluorooctane sulfonamidoacetic acids – N-MeFOSAA, NEt-FOSAA	0.0040	0.0020	ND – not detected			
Per and Polyfluoroether sulfonic acids – 9Cl-PF3ONS, 11Cl-PF3OUdS	0.0080	0.0020	ND – not detected			
Per and Polyfluoroether carboxylic acids – GenX	0.0040	0.0020	ND – not detected			
Per and Polyfluoroether carboxylic acids – ADONA	0.0080	0.0020	ND – not detected			

A] Values represent average of 3 lots tested, n=3 devices per lot; B] Values represent average of 2 lots tested, n=3 devices per lot.


ABBREVIATIONS: RL = reporting limit, ppb; MDL = minimum detection limit, ppb; PES = polyethersulfone



Testing Syringe Filters using EPA 537.1

There are variations in recovery of C-13 labeled standards

28 PFAS analyzed



- Filtered solvent: **water**
- Nylon demonstrates lower recoveries than PES
- Presence of glass fiber did not make a difference



- **Similar trends reported** previously for PES vs. nylon in aqueous solutions [1]

[1] Lath, S., et al. *Chemosphere*, 2019, 222, 671-8.



Analytical Methods for PFAS Developed at a Rapid Pace

Regulatory Landscape – EPA Draft 1633

Method(s)	Sample Prep	Matrix/Matrices	Filter(s) Required by Method
SW-846 Method 8327*	SPE, filtration	Non-potable groundwater, surface water, wastewater	0.2µm polypropylene syringe filter (hydrophilic)
ASTM D7968-17a	Solvent extraction, filtration	Environmental solids	0.2µm polypropylene syringe filter (hydrophilic)
ASTM D7979-19	Solvent extraction, filtration	Water matrix (no drinking water)	0.2µm polypropylene syringe filter (hydrophilic)
ISO 21675	SPE, filtration as needed	Drinking, natural and wastewater	1µm GFF syringe filter (rapid) 1-10µm nylon or GFF (if >2g/L suspended matter)
FDA C-010.02	QuEChERS, SPE, filtration	Foods	0.2µm nylon syringe filter
OTM-45	Sampling train: filtration , impingers	Air Emissions (stationary sources)	GFF or QFF (no pore size listed) membrane filter
EPA Draft 1633	SPE, filtration	Aqueous, soil, biosolids, sediment, tissue	0.2µm nylon syringe filter

- Polypropylene, nylon and glass fiber (GFF) are the most common materials
- Filters should not have detectable levels of PFAS compounds above reporting limit (RL) of method


**screening method only*

Abbreviations: SPE = solid phase extraction; TSS = total suspended solids determination; GAC = granular activated carbon; CIC = combustion ion chromatography; philic = hydrophilic; GFF = glass fiber filter; QFF = quartz fiber filter
Selected methods; does not include all drinking water and international methods




Are Membranes Contaminated with PFAS out-of-box? Testing Syringe Filters using EPA 537.1

Overview of Modified EPA Draft 1633



- Isolate filtration portion of EPA Draft 1633 for Aqueous samples



- 5 mL methanol sample + extracted & non-extracted internal standards mixes (EIS/NIS)
- Filtration with Millex® syringe filter
- LC-MS/MS, C18 column
- Analysis by internal standards

LC-MS/MS Conditions

Column:	C18, 2.7µm, 100x2.1mm ID, superficially porous column																																										
Mobile phase:	[A] Acetonitrile [B] 2 mM ammonium acetate in 95:5 water: acetonitrile																																										
Gradient:	<table border="1"> <thead> <tr> <th>Time (min)</th> <th>A %</th> <th>B%</th> <th>Flow (mL/min)</th> </tr> </thead> <tbody> <tr><td>0.20</td><td>10.0%</td><td>90.0%</td><td>0.350</td></tr> <tr><td>4.00</td><td>30.0%</td><td>70.0%</td><td>0.350</td></tr> <tr><td>7.00</td><td>55.0%</td><td>45.0%</td><td>0.350</td></tr> <tr><td>9.00</td><td>75.0%</td><td>25.0%</td><td>0.350</td></tr> <tr><td>10.00</td><td>95.0%</td><td>5.0%</td><td>0.400</td></tr> <tr><td>10.30</td><td>95.0%</td><td>5.0%</td><td>0.400</td></tr> <tr><td>10.40</td><td>2.0%</td><td>98.0%</td><td>0.400</td></tr> <tr><td>11.80</td><td>2.0%</td><td>98.0%</td><td>0.400</td></tr> <tr><td>13.00</td><td>2.0%</td><td>98.0%</td><td>0.400</td></tr> </tbody> </table>	Time (min)	A %	B%	Flow (mL/min)	0.20	10.0%	90.0%	0.350	4.00	30.0%	70.0%	0.350	7.00	55.0%	45.0%	0.350	9.00	75.0%	25.0%	0.350	10.00	95.0%	5.0%	0.400	10.30	95.0%	5.0%	0.400	10.40	2.0%	98.0%	0.400	11.80	2.0%	98.0%	0.400	13.00	2.0%	98.0%	0.400		
Time (min)	A %	B%	Flow (mL/min)																																								
0.20	10.0%	90.0%	0.350																																								
4.00	30.0%	70.0%	0.350																																								
7.00	55.0%	45.0%	0.350																																								
9.00	75.0%	25.0%	0.350																																								
10.00	95.0%	5.0%	0.400																																								
10.30	95.0%	5.0%	0.400																																								
10.40	2.0%	98.0%	0.400																																								
11.80	2.0%	98.0%	0.400																																								
13.00	2.0%	98.0%	0.400																																								
Flow rate:	See gradient table																																										
Detection:	MS/MS ESI(-), details of MS/MS conditions can be found at sigmaaldrich.com/pfassamplefiltration																																										
Column temp:	50.0°C																																										
Injection volume:	6-10µL autosampler injection																																										
Sample:	5mL methanol sample + EIS/NIS mixes																																										

Testing Syringe Filters using EPA Draft 1633

Millex® syringe filters do not have levels of PFAS above RL

40 PFAS analyzed

Compound(s)	RL* (ppb)	Syringe Filter Device				
		Millipore® 0.22µm PES ^P	Millipore® 0.20µm Nylon ^P	MFR-5 0.2µm Nylon ^P	Millipore® 0.20µm nylon-HPF	
					Lot1	Lot2 ^Q
Perfluoroalkyl carboxylic acids						
PFBA	0.8 ^A	ND	ND	ND	ND	ND
PFPeA	0.4 ^B	ND	ND	ND	ND	ND
PFHxA, PFHpA, PFOA, PFNA, PFUnDA, PFDODA, PFTrDA, PFTeDA	0.2 ^C	ND	ND	ND	ND ^R	ND
Perfluoroalkyl sulfonic acids						
PFBS, PFPeS, PFHxS, PFHpS, PFOS, PFNS, PFDS, PFDoS	0.2 ^D	ND	ND	ND	ND	ND
Fluorotelomer sulfonic acids						
4:2 FTS, 6:2 FTS, 8:2 FTS	0.8 ^E	ND	ND	ND	ND	ND
Fluorotelomer carboxylic acids						
3:3 FTCA	1 ^F	ND	ND	ND	ND	ND
5:3 FTCA, 7:3 FTCA	5 ^G	ND	ND	ND	ND	ND
Perfluorooctane sulfonamides, Perfluorooctane sulfonamidoacetic acids & Perfluorooctane sulfonoamido ethanols						
FOSA, N-MeFOSA, N-EtFOSA,	0.2 ^H	ND	ND	ND	ND	ND
N-MeFOSAA, N-EtFOSAA	0.2 ^I	ND	ND	ND	ND	ND
N-MeFOSE, N-EtFOSE	2 ^J	ND	ND	ND	ND	ND
Per and Polyfluoroether carboxylic acids & Per and Polyfluoroether sulfonic acids						
GenX, ADONA	0.8 ^K	ND	ND	ND	ND	ND
PFMPA, PFMBA, NFDHA	0.4 ^L	ND	ND	ND	ND	ND
9CI-PF3ONS, 11CI-PF3OUdS	0.8 ^M	ND	ND	ND	ND	ND
PFEESA	0.4 ^N	ND	ND	ND	ND	ND

*Minimum detection limits, MDL (ppb), ranged depending on analyte of interest:

- A] 0.4
- B] 0.1
- C] 0.05-0.081
- D] 0.1-0.11
- E] 0.4-0.41
- F] 0.5
- G] 1
- H] 0.1
- I] 0.1-0.13
- J] 1
- K] 0.2
- L] 0.1-0.12
- M] 0.2
- N] 0.1

P] Average of 2 lots, n=3 devices per lot; Q] Fourth device tested after hits reported in first replicate device and result was ND; R] Out of n=3 devices, 1 demonstrated hits above MDL but below RL

ABBREVIATIONS: RL = reporting limit, ppb; MDL = minimum detection limit, ppb; PES = polyethersulfone; HPF = High particulate filter



Testing Syringe Filters using EPA Draft 1633

Millex® syringe filters do not have levels of PFAS above RL

40 PFAS analyzed

Compound(s)	RL* (ppb)	Syringe Filter Device				
		Millipore® 0.22µm PES ^P	Millipore® 0.20µm Nylon ^P	MFR-5 0.2µm Nylon ^P	Millipore® 0.20µm nylon-HPF	
					Lot1	Lot2
Perfluoroalkyl carboxylic acids						
PFBA	0.8 ^A	ND	ND	ND	ND	ND
PFPeA	0.4 ^B	ND	ND	ND	ND	ND
PFHxA, PFHpA, PFOA, PFNA, PFUnDA, PFDODA, PFTTrDA, PFTeDA	0.2 ^C	ND	ND	ND	ND ^R	ND
Perfluoroalkyl sulfonic acids						
PFBS, PFPeS, PFHxS, PFHpS, PFOS, PFNS, PFDS, PFDoS	0.2 ^D	ND	ND	ND	ND	ND
Fluorotelomer sulfonic acids						
4:2 FTS, 6:2 FTS, 8:2 FTS	0.8 ^E	ND	ND	ND	ND	ND
Fluorotelomer carboxylic acids						
3:3 FTCA	1 ^F	ND	ND	ND	ND	ND
5:3 FTCA, 7:3 FTCA	5 ^G	ND	ND	ND	ND	ND
Perfluorooctane sulfonamides, Perfluorooctane sulfonamidoacetic acids & Perfluorooctane sulfonamido ethanols						
FOSA, N-MeFOSA, N-EtFOSA, N-MeFOSAA, N-EtFOSAA	0.2 ^H	ND	ND	ND	ND	ND
N-MeFOSE, N-EtFOSE	2 ^J	ND	ND	ND	ND	ND
Per and Polyfluoroether carboxylic acids & Per and Polyfluoroether sulfonic acids						
GenX, ADONA	0.8 ^K	ND	ND	ND	ND	ND
PFMPA, PFMBA, NFDHA	0.4 ^L	ND	ND	ND	ND	ND
9CI-PF3ONS, 11CI-PF3OUdS	0.8 ^M	ND	ND	ND	ND	ND
PFEESA	0.4 ^N	ND	ND	ND	ND	ND

Compound	Device1	Device2	Device3
PFHxA	ND	ND	0.0835
PFHpA	ND	ND	0.0688
PFOA	ND	ND	0.0992
PFNA	ND	ND	0.0827
PFDA	ND	ND	0.115
PFUnDA	ND	ND	0.107
PFDODA	ND	ND	0.113
PFTTrDA	ND	ND	0.116
PFTeDA	ND	ND	0.149

*Minimum detection limits, MDL (ppb), ranged depending on analyte of interest:

- A] 0.4
- B] 0.1
- C] 0.05-0.081
- D] 0.1-0.11
- E] 0.4-0.41
- F] 0.5
- G] 1
- H] 0.1
- I] 0.1-0.13
- J] 1
- K] 0.2
- L] 0.1-0.12
- M] 0.2
- N] 0.1

P] Average of 2 lots, n=3 devices per lot; Q] Fourth device tested after hits reported in first replicate device and result was ND; R] Out of n=3 devices, 1 demonstrated hits above MDL but below RL


ABBREVIATIONS: RL = reporting limit, ppb; MDL = minimum detection limit, ppb; PES = polyethersulfone; HPF = High particulate filter



Testing Syringe Filters using EPA Draft 1633

Methanol solvent improves PFAS analyte recoveries

40 PFAS
analyzed



- Filtered solvent: **methanol**
- All filters show similar recoveries
- Presence of glass fiber did not make a difference
- Much improved recovery for Nylon and similar for PES (vs. **water**)



- **Consistent with literature:** a wash with methanol caused PFOS, PFOA and other PFAS to desorb significantly from nylon filter media [1]



[1] So, M.K., et al. Archives of Environmental Contamination and toxicology, 2006, 50, 240-8.



Testing Syringe Filters using EPA Draft 1633 in the ppt range

Millex® syringe filters do not have levels of PFAS above RL

40 PFAS analyzed



- **Full workflow of EPA Draft 1633 tested, using aqueous matrix**
 - Spike: Low- and Mid-level
 - Solvent filtered: **Methanol**
- **Syringe Filters Tested:**
 - Millipore® 0.20µm Nylon (33mm diameter)
 - MFR-5 0.20µm Nylon (25mm diameter)



- No PFAS compounds detected in any device above RL = **no contamination**
- The loss of analyte was similar per cm² of filter area for each syringe filter type = **larger filtration areas do not lead to more analyte loss**

Data collected in collaboration with:




(Mansfield, MA)




Which Membrane, and Why?


Membrane Morphology & Applications




Polyethersulfone (PES)



Nylon



Hydrophilic polypropylene



Hydrophobic polypropylene

Syringe Filter Devices: Millex®



A good option for mobile phase & sample prep because:

- Great chemical compatibility (comparable to hydrophilic PTFE) and low extractables
- Suitable retention of particles in both air and liquid
- High stability to heating
- Hydrophilic & hydrophobic options (most are hydrophobic)
- Suggested in some PFAS methods

Millipore® cut disc membrane filters




Which Membrane, and Why?

Membrane Morphology & Applications



Low Extractables in PFAS relevant solvents



PTFE = polytetrafluoroethylene; PP = polypropylene (hydrophilic)

A good option for mobile phase & sample prep because:

- **Great chemical compatibility (comparable to hydrophilic PTFE) and low extractables**
- Suitable retention of particles in both air and liquid
- High stability to heating
- Hydrophilic & hydrophobic options (most are hydrophobic)
- Suggested in some PFAS methods




Which Membrane, and Why?

Membrane Morphology & Applications




Low Extractables in PFAS relevant solvents



PTFE = polytetrafluoroethylene; PP = polypropylene (hydrophilic)


A good option for mobile phase & sample prep because:

- **Great chemical compatibility (comparable to hydrophilic PTFE) and low extractables**
- Suitable retention of particles in both air and liquid
- High stability to heating
- Hydrophilic & hydrophobic options (most are hydrophobic)
- Suggested in some PFAS methods



Use Swinnex® to test polypropylene according to EPA 537.1 (water)

Use Swinnex® to test polypropylene according to EPA Draft 1633 (methanol)



Swinnex®



Testing Polypropylene Membranes using EPA 537.1

Polypropylene membrane filters (and Swinnex®) do not have detectable levels of PFAS contaminants

28 PFAS analyzed



Compound(s)	RL (ppb)	MDL (ppb)	Millipore® Polypropylene Membranes in Swinnex®			
			Hydrophobic ^a		Hydrophilic ^a	
			0.20µm	0.45µm	0.2µm	0.45µm
Perfluoroalkyl carboxylic acids – PFBA	0.0040	0.0020	ND	ND	ND	ND
Perfluoroalkyl carboxylic acids – PFPeA, PFHxA, PFHpA, PFOA, PFNA, PFDA, PFUnDA	0.0020	0.0010	ND	ND	ND	ND
Perfluoroalkyl carboxylic acids – PFDoDA, PFTTrDA, PFTeDA	0.0020	0.0010	ND ^b	ND ^b	ND	ND
Perfluoroalkyl sulfonic acids – PFBS, PFPeS, PFHxS, PFHpS, PFOS, PFNS, PFDS	0.0020	0.0010	ND	ND	ND	ND
Fluorotelomer sulfonic acids – 4:2 FTS, 6:2 FTS, 8:2 FTS	0.0080	0.0020	ND	ND	ND	ND
Perfluorooctane sulfonamides – PFOSA	0.0040	0.0020	ND	ND	ND	ND
Perfluorooctane sulfonamidoacetic acids – N-MeFOSAA	0.0040	0.0020	ND	ND	ND	ND
Perfluorooctane sulfonamidoacetic acids – NET-FOSAA	0.0040	0.0020	ND ^b	ND ^b	ND	ND
Per and Polyfluoroether sulfonic acids – 9Cl-PF3ONS, 11Cl-PF3OUdS	0.0080	0.0020	ND	ND	ND	ND
Per and Polyfluoroether carboxylic acids – GenX	0.0040	0.0020	ND	ND	ND	ND
Per and Polyfluoroether carboxylic acids – ADONA	0.0080	0.0020	ND	ND	ND	ND

A] Values represent average of 3 devices tested over 1 lot; B] Associated internal standard outside of control limits for all three devices tested: 13C2-PFDoDA, 13C2-PFTeDA, d5-EtFOSAA.
ABBREVIATIONS: RL = reporting limit, ppb; MDL = minimum detection limit, ppb; PP = polypropylene; philic = hydrophilic; phobic = hydrophobic



Testing Hydrophilic Polypropylene Membranes using EPA Draft 1633

Hydrophilic polypropylene membrane filters (and Swinnex®) do not have detectable levels of PFAS contaminants

40 PFAS analyzed



Compound(s)	RL (ppb)	MDL Range* (ppb)	Millipore® Hydrophilic Polypropylene Membranes in Swinnex®, 0.2µm		
			Lot1	Lot2	Lot3
Perfluoroalkyl carboxylic acids					
PFBA	0.8	0.4	ND	ND	ND
PFPeA	0.4	0.1	ND	ND	ND
PFHxA, PFHpA, PFOA, PFNA, PFUnDA, PFDoDA, PFTrDA, PFTeDA	0.2	0.05-0.081	ND	ND	ND
Perfluoroalkyl sulfonic acids					
PFBS, PFPeS, PFHxS, PFHpS, PFOS, PFNS, PFDS, PFDoS	0.2	0.1-0.11	ND	ND	ND
Fluorotelomer sulfonic acids					
4:2 FTS, 6:2 FTS, 8:2 FTS	0.8	0.4-0.41	ND	ND	ND
Fluorotelomer carboxylic acids					
3:3 FTCA	1	0.5	ND	ND	ND
5:3 FTCA, 7:3 FTCA	5	1	ND	ND	ND
Perfluorooctane sulfonamides, Perfluorooctane sulfonamidoacetic acids & Perfluorooctane sulfonoamido ethanols					
FOSA, N-MeFOSA, N-EtFOSA	0.2	0.1	ND	ND	ND
N-MeFOSAA, N-EtFOSAA	0.2	0.1-0.13	ND	ND	ND
N-MeFOSE, N-EtFOSE	2	1	ND	ND	ND
Per and Polyfluoroether carboxylic acids & Per and Polyfluoroether sulfonic acids					
GenX, ADONA	0.8	0.2	ND	ND	ND
PFMPA, PFMBA, NFDHA	0.4	0.1-0.12	ND	ND	ND
9Cl-PF3ONS, 11Cl-PF3OUdS	0.8	0.2	ND	ND	ND
PFEESA	0.4	0.1	ND	ND	ND

- 0.2µm hydrophilic polypropylene tested only

*Minimum detection limits, MDL (ppb), ranged depending on analyte of interest. Refer to [sigmaaldrich.com/pfassampfiltration](https://www.sigmaaldrich.com/pfassampfiltration) for details.

a] Extracted internal standards (EIS): 13C4-PFBA, 13C5-PFPeA, 13C5-PFHxA, 13C4-PFHpA, 13C8-PFOA; b] Values represent the average of n=3 devices for each lot. No PFAS compound detected in any device above RL or MDL.

ABBREVIATIONS: RL = reporting limit, ppb; MDL = minimum detection limit, ppb; PP = polypropylene; philic = hydrophilic; phobic = hydrophobic




Testing Polypropylene Membranes with EPA 537.1 and EPA Draft 1633


Recoveries vary widely from philic to phobic, water to methanol



Average Recovery of C13 labeled standards in Water → lower for hydrophobic polypropylene



Average Recovery of C13 labeled standards in water and methanol for hydrophilic polypropylene → similar overall



Analytical Methods for PFAS Developed at a Rapid Pace

Regulatory Landscape – Air Testing & OTM-45

Method(s)	Sample Prep	Matrix/Matrices	Filter(s) Required by Method
SW-846 Method 8327*	SPE, filtration	Non-potable groundwater, surface water, wastewater	0.2µm polypropylene syringe filter (hydrophilic)
ASTM D7968-17a	Solvent extraction, filtration	Environmental solids	0.2µm polypropylene syringe filter (hydrophilic)
ASTM D7979-19	Solvent extraction, filtration	Water matrix (no drinking water)	0.2µm polypropylene syringe filter (hydrophilic)
ISO 21675	SPE, filtration as needed	Drinking, natural and wastewater	1µm GFF syringe filter (rapid) 1-10µm nylon or GFF (if >2g/L suspended matter)
FDA C-010.02	QuEChERS, SPE, filtration	Foods	0.2µm nylon syringe filter
OTM-45	Sampling train: filtration , impingers	Air Emissions (stationary sources)	GFF or QFF (no pore size listed) membrane filter
EPA Draft 1633	SPE, filtration	Aqueous, soil, biosolids, sediment, tissue	0.2µm nylon syringe filter

- Polypropylene, nylon and glass fiber (GFF) are the most common materials
- Filters should not have detectable levels of PFAS compounds above reporting limit (RL) of method


**screening method only*

Abbreviations: SPE = solid phase extraction; TSS = total suspended solids determination; GAC = granular activated carbon; CIC = combustion ion chromatography; philic = hydrophilic; GFF = glass fiber filter; QFF = quartz fiber filter; OTM = other test method Selected methods; does not include all drinking water and international methods




The First PFAS Method in Air, Focused on Stack Emissions

Other Test Method (OTM)-45



The First PFAS Method in Air, Focused on Stack Emissions

Other Test Method (OTM)-45





Filter for collecting particulate phase:

- 82, 82.6, 90 mm cut disc membrane filter
- Glass fiber without binder
- At least 99.95% efficiency on 0.3- μm dioctyl phthalate smoke particles



AP 40 Glass fiber (0.7 μm)

- TSS
- TCLP
- Hot gases






Are Membranes Contaminated with PFAS out-of-box?

Testing Glass Fiber Membrane Filter Suitability for OTM-45

Cleaning (prior to sample collection)

- Soak filters in 5% (v/v) NH₄OH in methanol (30 min)
- Soak filters in methylene chloride (30 min)

Sampling on Sample Train
(skipped in this study to examine filters alone)



Front-half Extraction

- Addition of isotope dilution internal standards to filters
- Extraction of filters in methanol, 18 hours
- Hot-block concentration and blow-down step
 - Add non-extracted internal standards
 - Final volume: 10mL (80% v/v methanol in water)
- Analyze by internal standards

Concentration & Evaporation

LC-MS/MS

Data collected in collaboration with:





Testing Glass Fiber Membrane Filter Suitability for OTM-45

Glass fiber membrane filters demonstrate similar levels of PFAS after cleaning and front-half extraction

Compound(s)	RL (ng/sample)	MDL** Range (ng/sample)	Millipore® Glass fiber without Binder		MFR-1, Glass Fiber without Binder
			Lot 1, ^a pooled	Lot 2, ^a pooled	Lot 1, ^a pooled
Perfluoroalkyl carboxylic acids					
PFHpA	1.00	0.620	2.04	1.82	1.91 ←
PFBA	2.00	1.30	ND	ND	ND
PFHxA	1.00	0.210	0.258 ^J	0.270 ^{J,I}	ND ←
PFPeA, PFOA, PFNA, PFDA, PFUnDA, PFDoDA, PFTTrDA, PFTeDA, PFHxDA, PFOA	1.00	0.0850-0.650	ND	ND	ND
Perfluoroalkyl sulfonic acids					
PFBS, PFPeS, PFHxS, PFHpS, PFNS, PFDS, PFDoS	1.00	0.0950-0.890	ND	ND	ND
PFOS	1.00	0.450	0.592 ^J	0.565 ^{J,I}	0.566 ^{J,I} ←
Perfluorooctane sulfonamides, Perfluorooctane Sulfonamidoacetic acids, Perfluorooctane sulfonamido ethanols					
PFOSA, N-EtFOSA, N-MeFOSA, N-MeFOSAA, N-EtFOSAA, N-EtFOSE	1.00	0.0880-0.160	ND	ND	ND
N-MeFOSE	5.00	4.90	ND	ND	ND
Fluorotelomer sulfonic acids					
4:2 FTS, 8:2 FTS, 10:2 FTS	1.00	0.0910-0.320	ND	ND	ND
6:2 FTS	5.00	4.70	ND	ND	ND
Per and Polyfluoroether carboxylic acids, Per and Polyfluoroether sulfonic acids, Fluorotelomer carboxylic acids, Next-Generation PFAS analytes					
GenX	5.00	4.70	ND	ND	ND
ADONA, 9CI-PF3ONS, 11CI-PF3OUdS, 6:2 FTUCA, 7:3 FTCA*, 10:2 FTCA, 8:2 FTCA, PFEESA, 8:2 FTUCA, PFMPA, PFMBA, 5:3 FTCA*, 6:2 FTCA, 3:3 FTCA, PFECHS, NFDHA	1.00	0.0980-480	ND	ND	ND


a] n=3 membranes per lot tested and results pooled together. J] less than RL but greater than MDL and concentration is estimated. I] Value is the estimated maximum possible concentration. * LCS and/or LSCD is outside acceptance limits. **MDL ranged depending on the analyte of interest.



Testing Glass Fiber Membrane Filter Suitability for OTM-45

Average recovery of standards fall within acceptable QC range

49 PFAS analyzed




Most recoveries stay from 80-100%



Testing Glass Fiber Membrane Filter Suitability for OTM-45

Average recovery of standards fall within acceptable QC range

49 PFAS analyzed



Most recoveries stay from 80-100%



Investigation of Contamination and Recovery of PFAS in Filter Membranes

Summary & Conclusions

- Filtration is important for high data quality, however, **consumables must not interfere** with sensitive PFAS analysis
- We tested different syringe filters and membrane filters in several PFAS methods
 - Nylon, nylon with a glass fiber prefilter, polyethersulfone (PES), hydrophilic polypropylene, hydrophobic polypropylene, glass fiber without binder
 - **Analyzed extractables** to determine baseline contamination
 - **Analyzed** internal standard and analyte **recoveries** to determine adsorption
- The syringe filters and membrane filters tested were found to be **suitable for sample filtration/preparation** in:
 - Water (via EPA 537.1)
 - Methanol (via EPA Draft 1633)
 - Cleanup/Front-half extraction for stack sampling (via OTM-45)



Some key takeaways:

- PES and philic polypropylene show similar recoveries in water & methanol, while nylon does not
- Rinsing nylon with methanol desorbs certain PFAS
- Larger filtration areas do not lead to loss of analyte (33mm vs. 25mm)
- Polypropylene is a suitable alternative to PTFE for sample prep in many cases



Lindsay D. Lozeau, Ph.D.


EMD Millipore Corporation | 600 Summit Drive
Burlington, MA 01803 | USA
Email: lindsay.lozeau@milliporesigma.com

Amy Laws, Ph.D.

EMD Millipore Corporation | 600
Summit Drive
Burlington, MA 01803 | USA
Email:
amy.laws@milliporesigma.com

Maricar Dube, Ph.D.

EMD Millipore Corporation | 600
Summit Drive
Burlington, MA 01803 | USA
Email:
Maricar.dube@milliporesigma.com



The vibrant M, Millex®, Millipore®, Swinnex® are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2021 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.