

Prepared by Bradley VanMiddlesworth<sup>1</sup>, Emily Parry<sup>2</sup>  
 1 - ePrep, 2820 Shadelands Drive, Walnut Creek, CA, 94598, United States  
 2 - Agilent Technologies, 2850 Centerville Rd, Wilmington, DE, 19808, United States

## Introduction

Per- and Polyfluoroalkyl Substances (PFAS) are of considerable interest to regulatory agencies across the globe. In January 2024, the US EPA released method 1633 to test 40 PFAS compounds in a variety of matrices, including soils. An automated screening method for 51 PFAS compounds in soils has been developed and tested with Ottawa sand, topsoil, and humic peat soil. Sixty samples can be extracted in parallel with a specialized, grinding magnetic stirrer and automatically controlled magnetic stirring plate. Samples are extracted with methanolic ammonium hydroxide, cleaned using a carbon SPE cartridge, then concentrated prior to LC-MS/MS analysis. Extracted internal standard recoveries were within EPA 1633 guidelines. Reproducibility of the automated sample preparation workstation were < 10% RSD in laboratory control samples. Carryover was measured with a solvent blank after a mid-level spike and is generally < 1%.

## Extraction workflow

- 2 g soil, polypropylene (PP) sealed vials, PP stirring magnets (added manually)
- Add Extraction Standard at 10 ppb
- Add PFAS spikes
- Add 15 mL Extraction solvent  
3% ammonium hydroxide 1:1 MeOH:Acetone
- Magnetic Stir 6 mins at 600 RPM
- Pause for 10 min settling time

### Soil Extraction


- Condition carbon SPE cartridge with 5 mL MeOH followed by 5 mL of extraction solvent
- Load 6mL of sample extract onto carbon cartridges, collect in 15 mL centrifuge tubes
- Elute with 3mL of extraction solvent, collect in 15 mL centrifuge tubes

### SPE


- Concentrate sample to 1 mL
- Add Non-extracted standard
- LC-MS/MS Analysis

### LC-MS/MS Analysis


Automated




**ePrep ONE Workstation**




**Agilent 1290 Infinity II LC with PFC Free Kit and 6475 TQ**



**ePrep Software**



**MassHunter 12.1**



Solvent Blank   Ottawa Sand   Humic Peat



**Blank Evaluation**  
Ottawa Sand Blank

PFAS	Conc. (ng/g)	Area	Height	Retention Time (min)
13C4-PFBSA	10	100	100	1.2
13C5-PFPeA	10	100	100	1.5
13C2-4-2-FTSA	10	100	100	1.8
13C3-PFPeA	10	100	100	2.1
13C3-PFBS	10	100	100	2.4
13C3-PFPeDA	10	100	100	2.7
13C4-PFPeA	10	100	100	3.0
13C2-4-2-FTSA	10	100	100	3.3
13C3-PFPeDA	10	100	100	3.6
13C3-PFHUS	10	100	100	3.9
13C3-PFPeA	10	100	100	4.2
13C2-2-2-FTSA	10	100	100	4.5
13C3-PFPeDA	10	100	100	4.8
2H9-EFOSAA	10	100	100	5.1
2H5-N-EFOSAA	10	100	100	5.4
13C3-PFOS	10	100	100	5.7
13C7-PFPeDA	10	100	100	6.0
13C2-PFPeDA	10	100	100	6.3
13C3-PFPeDA	10	100	100	6.6
13C2-PFPeDA	10	100	100	6.9
2H9-EFOSAA	10	100	100	7.2
2H5-N-EFOSAA	10	100	100	7.5

**Carryover Evaluation**  
Solvent Blank extracted after mid-level spike

PFAS	Conc. (ng/g)	Area	Height	Retention Time (min)
13C4-PFBSA	10	100	100	1.2
13C5-PFPeA	10	100	100	1.5
13C2-4-2-FTSA	10	100	100	1.8
13C3-PFPeA	10	100	100	2.1
13C3-PFBS	10	100	100	2.4
13C3-PFPeDA	10	100	100	2.7
13C4-PFPeA	10	100	100	3.0
13C2-4-2-FTSA	10	100	100	3.3
13C3-PFPeDA	10	100	100	3.6
13C3-PFHUS	10	100	100	3.9
13C3-PFPeA	10	100	100	4.2
13C2-2-2-FTSA	10	100	100	4.5
13C3-PFPeDA	10	100	100	4.8
2H9-EFOSAA	10	100	100	5.1
2H5-N-EFOSAA	10	100	100	5.4
13C3-PFOS	10	100	100	5.7
13C7-PFPeDA	10	100	100	6.0
13C2-PFPeDA	10	100	100	6.3
13C3-PFPeDA	10	100	100	6.6
13C2-PFPeDA	10	100	100	6.9
2H9-EFOSAA	10	100	100	7.2
2H5-N-EFOSAA	10	100	100	7.5

## Results



### Reproducibility

Extraction Standard	RSD in Humic Peat (n=6)	RSD in Ottawa Sand (n=7)
13C4-PFBSA	5%	0.2%
13C5-PFPeA	6%	0.1%
13C2-4-2-FTSA	5%	0.1%
13C3-PFPeA	7%	0.0%
13C3-PFBS	6%	0.2%
13C3-PFPeDA	6%	0.1%
13C4-PFPeA	8%	0.1%
13C2-4-2-FTSA	8%	0.1%
13C3-PFPeDA	8%	0.1%
13C3-PFHUS	5%	0.0%
13C3-PFPeA	7%	0.0%
13C2-2-2-FTSA	4%	0.2%
13C3-PFPeDA	6%	0.0%
2H9-EFOSAA	6%	0.2%
2H5-N-EFOSAA	8%	0.2%
13C3-PFOS	6%	0.0%
13C7-PFPeDA	6%	0.0%
13C2-PFPeDA	19%	0.0%
13C3-PFPeDA	7%	0.0%
13C2-PFPeDA	21%	0.0%
2H9-EFOSAA	42%	0.0%
2H5-N-EFOSAA	50%	0.4%
2H9-EFOSAA	31%	1.9%
2H5-N-EFOSAA	64%	0.4%

Likely caused by variability in the Humic Peat

## Conclusion

PFAS in soils screening analysis can be performed in an automated method to ensure reliable data for sound decision making.

Up to 60 soil samples can be processed on an ePrep workstation and an Agilent LC-MS/MS with low carryover and no system interference.

Minimal sampling handling is required. Laboratory control sample spiking, internal standard addition, extraction solvent addition, magnetic stirring, and SPE cleanup is performed by the robot in 8-10 minutes per sample.