


## **Novel developments in Inductively Coupled Plasma Mass Spectrometry: How can the analysis of complex samples be made simple?**

**Andy Fornadel, PhD**

Product Marketing Manager

Thermo Fisher Scientific

 The world leader in serving science



# Typical challenges faced by laboratories analyzing trace elements



Highly diverse matrix samples



Interruptions due to maintenance



Operational complexity





Personnel that operate several different instrument types




Reducing environmental impact




# Application of ICP-MS



Need for method optimization



Increased need of cleaning and system maintenance




Compromised method robustness

**ICP-MS is best suited in terms of performance but with historical challenges**

# EPA methods for environmental analysis


- EPA method 6020B is one of the most applied guidelines in environmental laboratories
- Often used as a starting point for method creation in other industries
- In contrast to EPA method 200.8\*, the use of Collision/Reaction Cells (CRC) is permitted

EPA 200.8 (1994)




21 elements

- ✓ Method for regulatory compliance
- ✓ QC acceptance criteria for lab accreditation



18 common elements

EPA 6020B (2014)




23 elements

- ✓ For guidance purpose, performance-based method

\* For the analysis of drinking water only



# Water Analysis using ICP-MS



# What interferences are there in ICP-MS?

## 1: Non-spectral interferences

- Signal **suppression** and **enhancement**
  - **Suppression** caused by high levels of easily ionised elements (Na, K, Mg) in the sample
  - **Enhancement** caused by the presence of carbon in the samples (As and Se in particular affected)





# What interferences are there in ICP-MS?

## 1: Non-spectral interferences

- Signal **suppression** and **enhancement**
  - **Suppression** caused by high levels of easily ionised elements (Na, K, Mg) in the sample
  - **Enhancement** caused by the presence of carbon in the samples (As and Se in particular affected)
- Can correct for signal **suppression** using internal standards and matrix matching standards to samples
- Can overcome **enhancement** effects also by sample to standard matrix matching
  - By adding an organic solvent such as isopropanol to all blanks, standards and samples



# What interferences are there in ICP-MS?

## 2: Spectral interferences

- Sub-divided into 'polyatomic' and 'isobaric'
- **Polyatomic interferences** formed from combination of plasma gases and sample matrix constituents
  - $^{40}\text{Ar}^{16}\text{O}^+$  on  $^{56}\text{Fe}^+$ ,  $^{40}\text{Ar}^{35}\text{Cl}^+$  on  $^{75}\text{As}^+$ ,  $^{40}\text{Ar}_2^+$  on  $^{80}\text{Se}^+$

The image shows a stylized, tilted periodic table of elements. The elements are arranged in a grid, with their symbols, names, and atomic numbers visible. The table is tilted at an angle, giving it a three-dimensional appearance. The elements shown include Al, Si, P, Ga, Ge, As, Se, Br, Zn, Ga, Ge, As, Se, Br, Cd, In, Sn, Sb, Te, I, Au, Hg, Tl, Pb, Bi, Po, At, Uuh, and Uus. The table is set against a dark background, and the elements are highlighted in a light color.



# What interferences are there in ICP-MS?

## 2: Spectral interferences

- Sub-divided into 'polyatomic' and 'isobaric'
- **Polyatomic interferences** formed from combination of plasma gases and sample matrix constituents
  - $^{40}\text{Ar}^{16}\text{O}^+$  on  $^{56}\text{Fe}^+$ ,  $^{40}\text{Ar}^{35}\text{Cl}^+$  on  $^{75}\text{As}^+$ ,  $^{40}\text{Ar}_2^+$  on  $^{80}\text{Se}^+$
- **Isobaric interferences** caused by overlap of isotopes of different elements that have the same mass
  - $^{48}\text{Ca}^+$  on  $^{48}\text{Ti}^+$ ,  $^{64}\text{Ni}^+$  on  $^{64}\text{Zn}^+$

The image shows a stylized, tilted periodic table of elements. The elements are arranged in rows and columns, with their atomic numbers and symbols visible. The elements shown include Al, Si, P, S, Ga, Ge, As, Se, Br, Zn, In, Sn, Sb, Te, I, Cd, Hg, Tl, Pb, Bi, Po, At, Au, and Uuh. The table is tilted at an angle, giving it a three-dimensional appearance.

# What interferences are there in ICP-MS?

## 2: Spectral interferences

- Sub-divided into 'polyatomic' and 'isobaric'
- **Polyatomic interferences** formed from combination of plasma gases and sample matrix constituents
  - $^{40}\text{Ar}^{16}\text{O}^+$  on  $^{56}\text{Fe}^+$ ,  $^{40}\text{Ar}^{35}\text{Cl}^+$  on  $^{75}\text{As}^+$ ,  $^{40}\text{Ar}_2^+$  on  $^{80}\text{Se}^+$
- **Isobaric interferences** caused by overlap of isotopes of different elements that have the same mass
  - $^{48}\text{Ca}^+$  on  $^{48}\text{Ti}^+$ ,  $^{64}\text{Ni}^+$  on  $^{64}\text{Zn}^+$
- **Doubly charged ion** interferences, formed from elements having a 2<sup>nd</sup> ionization potential lower than the ionization potential of argon (15.8 eV)
  - Type of isobaric interference
  - Appear at half the parent isotope mass
  - $^{150}\text{Nd}^{2+}$  on  $^{75}\text{As}^+$ ,  $^{156}\text{Gd}^{2+}$  on  $^{78}\text{Se}^+$


The image shows a portion of a periodic table, tilted at an angle. The elements visible include:

- Al (Aluminium, 13)
- Si (Silicon, 14)
- P (Phosphorus, 15)
- S (Sulfur, 16)
- Zn (Zinc, 30)
- Ga (Gallium, 31)
- Ge (Germanium, 32)
- As (Arsenic, 33)
- Se (Selenium, 34)
- Br (Bromine, 35)
- Cd (Cadmium, 48)
- In (Indium, 49)
- Sn (Tin, 50)
- Sb (Antimony, 51)
- Te (Tellurium, 52)
- I (Iodine, 53)
- Au (Gold, 79)
- Hg (Mercury, 80)
- Tl (Thallium, 81)
- Pb (Lead, 82)
- Bi (Bismuth, 83)
- Po (Polonium, 84)
- At (Astatine, 85)
- Uuh (Ununhexium, 116)
- Uus (Ununseptium, 117)


# Removing interferences with single quadrupole ICP-MS

5% HNO<sub>3</sub> 5% HCl, 1% IPA 1% H<sub>2</sub>SO<sub>4</sub> 200ppm Sodium, 500ppm Phosphorous, 200ppm Calcium

- **Chlorine** causes severe interferences on vanadium ( $^{35}\text{Cl}^{16}\text{O}^+$ ), chromium ( $^{35}\text{Cl}^{16}\text{O}^1\text{H}^+$ ) and arsenic ( $^{40}\text{Ar}^{35}\text{Cl}^+$ )
- **Carbon** interferes with chromium ( $^{40}\text{Ar}^{12,13}\text{C}^+$ )
- **Sulfur** causes interferences with titanium ( $^{32}\text{S}^{16}\text{O}^+$ ) or vanadium ( $^{32}\text{S}^{18}\text{O}^1\text{H}^+$ )
- **Sodium** may bias results for copper ( $^{40}\text{Ar}^{23}\text{Na}^+$ )



48Ti+




52, 53Cr+



51V+



63Cu+




75As+




# He KED collision cell

- **Cell with ion lenses purged with He reduces polyatomic interferences.**
- Some designs also serve as a “pre-filter” to the quadrupoles
- **Single analysis mode He KED** achieved excellent interference removal and detection of low concentration (below  $1 \mu\text{g}\cdot\text{L}^{-1}$ ) analytes with high ratio of signal/background.



Polyatomic interference of selenium



Symbol	Mass	Abundance	Interferences
74Se	73.9225	0.90	74Ge(36.500%); 16O
76Se	75.9192	9.00	76Ge(7.800%); 36Ar ..
77Se	76.9199	7.60	40Ar + 37Cl(24.133%)
78Se	77.9173	23.60	78Kr(0.350%); 14N +
80Se	79.9165	49.70	80Kr(2.250%); 40Ar +
82Se	81.9167	9.20	82Kr(11.600%); 1H +




# Removing interferences with single quadrupole ICP-MS

5% HNO<sub>3</sub> 5% HCl, 1% IPA 1% H<sub>2</sub>SO<sub>4</sub> 200ppm Sodium, 500ppm Phosphorous, 200ppm Calcium


Same matrix spiked with Li, Be, B, Na, Mg, Al, K, Ca, Sc, Ti, V, Cr, Mn, Fe, Ni, Co, Ni, Cu, Zn, Ga, Ge, As, Se at 10 µg·L<sup>-1</sup> each



# Removing interferences with single quadrupole ICP-MS



**STD mode:** Polyatomic interference leads to poor IDL and elevated BEC




**KED mode:** Polyatomic interference removed - IDL below 5 ppt



# Argon Gas Dilution



- **Auxiliary flow of argon gas added to aerosolized sample**
- Performs **online dilution** in aerosol, reduces time/error associated with manual liquid dilution
- **Reduces matrix load into plasma**
  - Reduces signal suppression from abundant, easily-ionized species (Na, Ca, K, etc.)
  - Reduces matrix/salt build-up in instrument, carryover



# Argon Gas Dilution


- Compatible with the standard sample introduction system components
- PFA-ST Microflow nebulizer is more resistant against blockage when analyzing high salt loads
- Argon humidification beneficial for salt-rich matrices

Sample Introduction System Component	Dilution Level		
	<i>Low</i>	<i>Mid</i>	<i>High</i>
Glass Concentric Nebulizer	✓	✓	✗
Baffled Cyclonic Spray chamber	✓	✓	✓
2.5 mm i.d. Quartz injector	✓	✓	✓
Torch (Quartz, Thermo Scientific™ PLUS Torch)	✓	✓	✓
Skimmer Cone insert	High Matrix		
ESI Pergo Humidifier*			✓
PFA-ST Microflow Nebulizer*			✓

# Argon Gas Dilution


- Tuning of the dilution level is achieved by variation of nebulizer gas flow and additional gas flow added
- As a consequence of the dilution, a lower oxide level is achieved, indicating a more robust plasma
- Tuning all dilution levels is fully automated

Sample Matrices	% TDS Content [%]	Dilution level
<b>Drinking Water and Surface Water</b>	<b>&lt; 0.5</b>	<b>Low</b>
<i>Wastewaters</i>	<i>&lt; 1.0</i>	<i>Mid</i>
<i>Soil digests, geological &amp; mining samples</i>	<i>&lt; 1.0</i>	
<i>Brackish waters, fracking flowback solutions</i>	<i>&lt; 1.5</i>	
<b>Brackish waters, sea water, brine solutions</b>	<b>&lt; 3.0</b>	<b>High</b>
Highly concentrated brine solutions	> 4.0	







# Experimental set up



Parameter	Value
Nebulizer	Micromist Nebulizer (400 $\mu\text{L}\cdot\text{min}^{-1}$ )
Interface cones	Ni – tipped sample and Skimmer
Skimmer cone Insert	High Matrix
Spray chamber	Cyclonic quartz
Injector	Quartz, 2.5 mm ID
Torch	Quartz Torch
RF Power (W)	1550
Number of Replicates	3
Spray Chamber Temp ( $^{\circ}\text{C}$ )	2.7
KED settings (gas flow rate in $\text{mL}\cdot\text{min}^{-1}$ )	4.8 (with a 3V kinetic energy barrier)

# Analysis samples

Item	Place	Category	Note
Tap water 1	Bremen West	Tap water	-
Tap water 2	Bremen South	Tap water	-
Tap water 3	Bremen North	Tap water	-
Surface Water 1	Bremen South	Lake	Sampling location is close to a major highway
Surface Water 2	Bremen North	Lake	Sampling location is close to an area with heavy traffic
Well water	Bremen North	Well water	Ground water sample, no additional treatment
Waste water 1	Bremen	Industrial waste water	Elevated Na, Ca, Fe
Waste water 2	Bremen	Industrial waste water	Elevated Na, Ca, Fe
Brackish water 1	-	Brackish water	Simulated
Brackish water 2	-	Brackish water	Simulated
SLRS-5	Ottawa	River	CRM
NASS-6	Nova Scotia	Seawater	CRM



# Accuracy - river water CRM SLRS-5

Concentration results in  $\mu\text{g}\cdot\text{L}^{-1}$

	CRM values	Measured	Recovery (%)	Result
Al	49.5	50.6	102%	✓
Sb	0.3	0.35	117%	✓
As	0.413	0.478	116%	✓
Ba	14.0	15.2	109%	✓
Cr	0.208	0.216	104%	✓
Co	0.05	0.052	104%	✓
Cu	17.4	18.7	107%	✓
Fe	91.2	91.3	100%	✓
Pb	0.081	0.077	95%	✓
Mn	4.33	4.64	107%	✓

	CRM values	Measured	Recovery (%)	Result
Mo	0.27	0.27	100%	✓
Ni	0.476	0.525	110%	✓
Sr	53.6	55.9	104%	✓
U	0.093	0.092	99%	✓
V	0.317	0.304	96%	✓
Zn	0.845	0.960	114%	✓
Na	5,380	4,890	91%	✓
Mg	2,540	2,450	96%	✓
K	839	823	98%	✓
Ca	10,500	9,900	94%	✓




# Accuracy - Seawater

## NAAS-6: Seawater CRM

Analyte	Certified conc ( $\mu\text{g.L}^{-1}$ )	Obtained conc ( $\mu\text{g.L}^{-1}$ )	% Accuracy
As	1.43	1.64	114.7
Cd	0.0311	0.026	87.5
Co	0.015	0.027	118.6
Mn	0.53	0.598	112.9
Mo	9.89	10.553	106.7
Ni	0.31	0.359	115.9
V	1.46	1.64	112.3

## Spike and accuracy check in seawater

Analytes	Spiked conc ( $\mu\text{g.L}^{-1}$ )
As, Cd, Co, Mn, Ni, Pb & Tl	4
Cr, Cu, Sb, Zn	8
Se, V	20
Ba	40
Al	80
Fe	800





# Spike and recovery – waste and brackish water

Analyte	Observed concentration (mg·L <sup>-1</sup> )		Spiked concentration (mg·L <sup>-1</sup> )	Observed concentration (mg·L <sup>-1</sup> )		% RPD	Average % recovery
	Sample 1 (unspiked)	Sample 2 (unspiked)		Sample 1 (spiked)	Sample 2 (spiked)		
Ag	<0.0001	<0.0001	0.02	0.021	0.021	2.4	103
Al	<0.025	<0.025	10	10.1	9.4	7.2	97
As	<0.001	<0.001	0.2	0.179	0.181	1.1	90
Ba	<0.001	<0.001	0.2	0.185	0.185	0.0	92
Be	<0.001	<0.001	0.2	0.221	0.196	12.0	104
Ca	455.2	451.6	N/A	N/A	N/A	0.9	N/A
Cd	<0.001	<0.001	0.2	0.174	0.171	1.7	86
Co	<0.001	<0.001	0.2	0.177	0.177	0.0	88
Cr	<0.001	<0.001	0.2	0.179	0.176	1.7	88
Cu	<0.001	<0.001	0.2	0.183	0.177	3.3	90
Fe	38.8	39.6	10	47.3	47.1	0.4	80
Hg	<0.0001	<0.0001	0.5	0.473	0.513	8.1	99
K	86.1	85.2	10	95.2	94.4	0.8	92
Mg	90.2	93.1	N/A	N/A	N/A	3.1	N/A
Mn	<0.001	<0.001	0.2	0.185	0.182	1.6	91
Mo	<0.001	<0.001	0.2	0.198	0.2	1.0	99
Na	175	169	N/A	N/A	N/A	3.7	N/A
Ni	<0.001	<0.001	0.2	0.184	0.183	0.5	91
Pb	<0.001	<0.001	0.2	0.176	0.177	0.6	88
Sb	<0.001	<0.001	0.2	0.192	0.188	2.1	95
Se	<0.001	<0.001	0.2	0.178	0.184	3.3	90
Tl	<0.001	<0.001	0.2	0.171	0.17	0.6	85
V	<0.001	<0.001	0.2	0.176	0.178	1.1	88
Zn	<0.001	<0.001	0.2	0.171	0.167	2.4	84

# Robustness

- Daily sample load was approximately 300 samples per day
- Stability of CCV analytes and internal standards as a check for deviation



# Conclusions

- **ICP-MS is a highly sensitive, multi-element technique for the analysis of environmental samples, such as drinking water, surface water, and wastewater**
- **Historically challenging with diverse and complex matrices**
  - Some challenges addressed through careful sample preparation, offline dilution, matrix matching, use of internal standards
- **Modern developments in ICP-MS technology have lessened common challenges**
  - KED collision cell (with low mass cutoff)
  - Argon gas dilution
  - Robust sample introduction and interface components
- **Enables routine, high-throughput analysis of mixed-matrix samples while yielding excellent recoveries and stability**
- **ICP-MS is growing in routine use to benefit from high sensitivity while minimizing operational challenges and spectral interferences.**



# Thank you

The line has been unmuted for questions.

