



# TN using persulfate digestion and DMP



Standard under development

Dr. Christian Prokisch | CProkisch@mn-net.com | 01.08.2024

# Agenda





# Introduction



# Introduction

## Significance of detecting Nitrogen in Water

- Essential for environmental and human health
- Indicator of water quality and ecosystem balance
- Critical for regulatory compliance and pollution control
- Impacts drinking water safety and aquatic life





# Introduction

## Current regulatory background

- EPA sets standards for nitrogen in drinking water
- Clean Water Act regulates nitrogen discharge
- Nutrient Criteria Program guides state-level standards
- Ongoing efforts to adopt numeric nutrient criteria





# Introduction

## Goal of this presentation

- Revisit an established method for total nitrogen detection in water
- Propose standardization for potential nationwide regulation
- Present data supporting method's reliability and applicability
- Discuss implications for future water quality regulations






# Introduction

## Sources and forms of Nitrogen

- Sources: Agriculture, wastewater, atmospheric deposition, urban runoff
- Forms: Organic nitrogen, ammonium ( $\text{NH}_4^+$ ), nitrate ( $\text{NO}_3^-$ ), nitrite ( $\text{NO}_2^-$ )
- Natural vs. anthropogenic sources
- Nitrogen cycle in aquatic environments





# Introduction

## Effects and Global Trends

- Effects of excess Nitrogen
  - Eutrophication and algal blooms
  - Hypoxia and fish kills
- Global trends in Nitrogen pollution
  - Drinking water contamination (e.g., blue baby syndrome)
  - Global increase in nitrogen pollution over past century
  - Regional variations and hotspots





# Introduction

## Existing methods

- ISO 23697-1:2023 and ISO 23697-2:2023:  
Small-scale sealed tube methods
- ISO 29441:2010: UV digestion method
- Standard Methods including but not limited to
  - Flow injection analysis
  - Persulfate digestion followed by any Nitrate determination
  - Conductimetric determination of inorganic nitrogen
  - High temperature combustion





# Introduction

## Limitations of the current methods

- Lack of validation data acceptable for EPA regulation
- Insufficient statistical robustness for US regulatory purposes
- Absence of comprehensive interlaboratory studies for some methods



Gap between existing methods and EPA regulatory requirements



Persulfate digestion followed by detection using DMP



# Persulfate digestion followed by detection using DMP

## Brief overview

- Established method for total nitrogen detection
- Persulfate digestion converts all nitrogen to nitrate
- Detection using 2,6-Dimethylphenol colorimetry
- Small-scale sealed tube method Potential applications





# Persulfate digestion followed by detection using DMP

## Digestion step using Persulfate

- Alkaline oxidation at 120 °C
- Converts organic and inorganic nitrogen to nitrate
- Uses potassium persulfate and NaOH
- Requires heating block and digestion tubes





# Persulfate digestion followed by detection using DMP

## Detection step using Dimethylphenol

- Nitrate detection via colorimetric reaction
- 2,6-Dimethylphenol as the color reagent
- Photometric evaluation at specific wavelengths
- Requires acid mixture and photometer





# Persulfate digestion followed by detection using DMP

## Preparation

- Collection of water samples  
(groundwater, surface water, wastewater)
- Filtration if necessary to remove particulates
- Preparation of calibration standards using nitrate solution
- Creation of digestion check standard using glutamic acid solution
- Preparation of reagent blanks





# Persulfate digestion followed by detection using DMP

## Equipment and reagents

- Equipment: Heating block, digestion tubes, detection tubes, photometer
- Key reagents: Potassium persulfate, NaOH, 2,6-Dimethylphenol
- Acid mixture: Phosphoric acid and sulfuric acid (1:1)
- Ammonia-free and nitrate-free water
- Stock and intermediate solutions for standards
- Collection of water samples (groundwater, surface water, wastewater)





# Persulfate digestion followed by detection using DMP

## Analytical Procedure

- Sample preparation and standard creation
- Digestion: Sample + persulfate reagent, heated at 120°C for 30 min
- Cooling of digested sample
- Detection: Addition of acid mixture and 2,6-Dimethylphenol
- Photometric measurement at 340–365 nm
- Quality control measures throughout





# Persulfate digestion followed by detection using DMP

## Advantages

- High accuracy and precision in nitrogen measurement
- Consistent results across different water types
- Simplified procedure with commercially available test kits
- Reduced need for complex sample preparation





# Persulfate digestion followed by detection using DMP

## Advantages

- Lower operational costs compared to traditional methods
- Reduced reagent consumption and waste generation
- Environmentally friendly with fewer hazardous chemicals
- Supports sustainable water quality monitoring practices

Low operational costs





Data



## ISO interlab trial – design

Sample ID	Matrix	NO <sub>3</sub> -N	NH <sub>4</sub> -N	TN <sub>b</sub>
A	Ground water	1	0,1	n.a.
B	Swimming pool water	7 – 8	0,5	n.a.
C	Waste water 1	40	3	n.a.
D	Waste water 2	8	40	n.a.
E	Waste water 3	n.a.	n.a.	10
F	Waste water 4	n.a.	n.a.	60
G	Waste water 5	n.a.	n.a.	80
H	Surface water	10 – 12	0.8	n.a.
I	QC standard	1	8	13



## ISO interlab trial – Symbols used

Symbols used	
$l$	number of laboratories after outlier rejection
$n$	number of individual test results after outlier rejection
$o$	percentage of outliers
$X$	assigned value
$\bar{\bar{x}}$	overall mean of results (without outliers)
$\eta$	recovery rate
$s_R$	reproducibility standard deviation
$C_{V,R}$	coefficient of variation of reproducibility
$s_r$	repeatability standard deviation
$C_{V,r}$	coefficient of variation of repeatability



## ISO interlab trial – results

Sample	Matrix	l	n	o	X	$\bar{x}$	$\eta$	$s_R$	$C_{V,R}$	$s_r$	$C_{V,r}$
				%	$\mu\text{g/L}$	$\mu\text{g/L}$	%	$\mu\text{g/L}$	%	$\mu\text{g/L}$	%
E	Waste water 3	41	164	10,9	-	16,30	-	0,864	5,3	0,281	1,7
F	Waste water 4	45	180	4,3	-	61,63	-	3,526	5,7	1,530	2,5
G	Waste water 5	45	179	4,8	-	70,42	-	9,996	14,2	2,341	3,3
I	QC standard	44	174	8,4	13,033	13,15	100,9	0,698	5,3	0,260	2,0



## Standard Methods 4500-N F



# Standard Methods 4500-N F

## Current status

- Received charge in January 2024
- First draft created in March 2024
- Comments received in April 2024
- Second draft sent out in July 2024
- Comment period open until July 19, 2024
- Review and incorporation of comments pending





# Standard Methods 4500-N F

## Method validation plans

- Aiming for Tier 3-like validation (nationwide use)
- Nine-laboratory interlaboratory study planned
- Minimum of nine different matrix types to be tested
- Method detection limit (MDL) determination by each lab
- Recovery studies and real sample analysis across application range





## Standard Methods 4500-N F – Current status and next steps

### Facilitating Participation in the Interlaboratory Trial

- Custom-made vials compatible with market-leading instruments
- Fits HACH photometers and heating blocks
- Eliminates need for reagent preparation by participants
- Reduces potential sources of error
- Lowers barrier to participation for laboratories Aiming for Tier 3-like validation (nationwide use)





# Standard Methods 4500-N F – Current status and next steps

## Call to Action – Join Our Interlaboratory Trial

- Seeking laboratories for method validation
- Opportunity to shape future water quality standards
- Easy participation with provided custom vials
- Gain early experience with potential future regulatory method
- Contact information for expressing interest:
  - [cprokisch@mn-net.com](mailto:cprokisch@mn-net.com)





## Summary

# Summary

- Established method for total nitrogen detection revisited
- Persulfate digestion with 2,6-Dimethylphenol colorimetry
- Advantages: accuracy, ease of use, cost-effectiveness
- Potential basis for nationwide regulation
- Ongoing standardization and validation process
- Call for participation in interlaboratory trial



Volunteers needed for total nitrogen method validation!

Dr. Christian Prokisch | [cprokisch@mn-net.com](mailto:cprokisch@mn-net.com) | 01.08.2024

## Image credits

---

© Jürgen Nickel (1, 31), Kekyalyaynen (4), psdesign1 (6), Alexey Stiop (8), [Orlando Florin Rosu](#) (10), Paulista (15), Alexander Raths (17), Uwe Bumann, (18,19), fotomek (25), rdnzl (26), Alex White (27), kerkezz (28), Swapan(30) – Adobe Stock

EPA Logo (5)

[https://de.wikipedia.org/wiki/Datei:Environmental Protection Agency logo.png](https://de.wikipedia.org/wiki/Datei:Environmental_Protection_Agency_logo.png)

ISO Logo (9)

[https://commons.wikimedia.org/wiki/File:ISO\\_Logo\\_%28Red\\_square%29.svg](https://commons.wikimedia.org/wiki/File:ISO_Logo_%28Red_square%29.svg)