


# Redefined Analysis Workflow following Environment method

- How to use latest GC-MS technology to enhance your laboratory operation? -

Yoshiro Hiramatsu - GCMS Product Liaison -  
Shimadzu Scientific Instrument

A decorative graphic element consisting of a red curved line that starts at the bottom left and curves upwards and to the right, ending at the top right. The line is surrounded by a soft, grey gradient shadow.

# What will be big deal to optimize laboratory workflow?





How to use latest GC-MS technology to enhance your laboratory operation?


# Enviromental application common workflow

## Instrument Checkout and report

Tune mix injection



ID#: 1 R. Time: 11.210(Scan#:1843)  
 MassPeaks:451  
 Spectrum Mode:Averaged 11.210-11.220(1843-1845)  
 BG Mode:Calculated from Peak Segment 1 - Event 1




| m/z | Criteria              | Rel. Int. | Abs. Int. | Status |
|-----|-----------------------|-----------|-----------|--------|
| 68  | < 2% of mass 69       | 1.22      | 1072      | Pass   |
| 69  | Present               | 19.57     | 88105     | Pass   |
| 70  | < 2% of mass 69       | 0.56      | 489       | Pass   |
| 197 | < 2% of mass 198      | 0.25      | 1125      | Pass   |
| 198 | Base peak or present  | 100.00    | 450208    | Pass   |
| 199 | 5 - 9% of mass 198    | 6.55      | 29500     | Pass   |
| 365 | > 1% of the base peak | 1.54      | 6941      | Pass   |
| 441 | < 150% of mass 443    | 82.00     | 52358     | Pass   |
| 442 | Base peak or present  | 77.59     | 349315    | Pass   |
| 443 | 15 - 24% of mass 442  | 18.28     | 63853     | Pass   |

# Environmental application common workflow

Run Initial Calibration and ICV(CCV)

TIC (100pg on column) :

8270 Megamix, Benzidine Mix, 1,4-Dioxane, Benzoic Acid, Atrazine, Benzaldehyde, Biphenyl, Caprolactam and ISTD/Surrogate



# Enviromental application common workflow



Check ICAL and ICV result

MC result :

1,4-dioxane

%RSD: 2.53% (5pg on column n=5)



RSE: 8.26% (0.01ppm to 20ppm)



2-Nitroaniline

%RSD: 3.23% (5pg on column n=5)



RSE: 7.96% (0.01ppm to 20ppm)



Di-n-butyl phthalate

%RSD: 1.87% (5pg on column n=5)



RSE: 10.43% (0.01ppm to 10ppm)



Benzo(a)pyrene


%RSD: 2.58% (5pg on column n=5)

RSE: 8.78% (0.01ppm to 10ppm)



# Enviromental application common workflow

## Inject Sample and report



Insight Report  
 Results File: C:\Users\Klomanadu\Desktop\KPA\_Data\124\LabSolutions\Data\Cal\_curve\CCV\Report\prec - 0FAS24\_2.gm  
 Printed at 2023/05/17 21:41:01

CCV2

| Type      | Name                      | Found RT | m/z    | Area   | Conc. Unit | RT Diff | Accuracy (%) | Flag Result | Mode ID | Sealed | Recovery (Enviro) | Comment 1         |
|-----------|---------------------------|----------|--------|--------|------------|---------|--------------|-------------|---------|--------|-------------------|-------------------|
| STD       | Fluorobenzene_Std         | 8.218    | 96.00  | 144480 | 50.0 ug/L  | 0.032   | Auto         | 30.0000     | Auto    | 0.0000 | Auto              |                   |
| Target    | Bromofom                  | 12.535   | 173.00 | 97556  | 40.0 ug/L  | -0.009  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
| Surrogate | 4-Bromofluorobenzene_Sur  | 12.895   | 95.00  | 276013 | 50.0 ug/L  | -0.009  | Auto         | 50.0000     | Auto    | 0.0000 | 101.49            | Flag Range 70-130 |
|           | Dichlorodifluoromethane   | 1.864    | 85.00  | 43278  | 55.0 ug/L  | -0.004  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Chloromethane             | 1.552    | 50.00  | 102145 | 45.0 ug/L  | -0.007  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Vinyl chloride            | 1.640    | 62.00  | 98183  | 49.1 ug/L  | -0.004  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Bromomethane              | 2.094    | 94.00  | 40589  | 45.0 ug/L  | -0.006  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Chloroethane              | 2.391    | 64.00  | 100939 | 48.0 ug/L  | -0.008  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Trichlorofluoromethane    | 2.762    | 101.00 | 132981 | 45.4 ug/L  | -0.011  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | 1,1-Dichloroethene        | 3.750    | 96.00  | 117865 | 49.1 ug/L  | -0.012  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Iodomethane               | 4.031    | 142.00 | 79692  | 72.1 ug/L  | 0.021   | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Carbon disulfide          | 4.137    | 76.00  | 344993 | 44.6 ug/L  | -0.013  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Allyl chloride            | 4.507    | 76.00  | 89601  | 49.1 ug/L  | -0.012  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Methylene chloride        | 4.737    | 84.00  | 146739 | 41.2 ug/L  | -0.012  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Acrylonitrile             | 5.176    | 53.00  | 159381 | 46.1 ug/L  | -0.012  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | trans-1,2-Dichloroethene  | 5.232    | 96.00  | 149968 | 48.9 ug/L  | -0.011  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Methyl tert-butyl ether   | 5.240    | 73.00  | 528206 | 43.0 ug/L  | -0.012  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | 1,1-Dichloroethane        | 5.937    | 83.00  | 156730 | 45.0 ug/L  | -0.011  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | cis-1,2-Dichloroethene    | 6.756    | 98.00  | 110233 | 46.0 ug/L  | 0.033   | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | 2,2-Dichloropropane       | 6.756    | 77.00  | 162752 | 35.8 ug/L  | -0.011  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Propionitrile             | 6.831    | 54.00  | 54114  | 50.1 ug/L  | -0.012  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Methyl acrylate           | 6.878    | 55.00  | 359860 | 46.2 ug/L  | -0.011  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Bromochloromethane        | 7.066    | 128.00 | 63238  | 45.6 ug/L  | -0.011  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Methyl acrylonitrile      | 7.091    | 67.00  | 140205 | 47.3 ug/L  | -0.012  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Tetrahydrofuran           | 7.079    | 71.00  | 41666  | 49.1 ug/L  | -0.011  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Chloroform                | 7.219    | 83.00  | 258868 | 45.9 ug/L  | -0.011  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | 1,1,1-Trichloroethane     | 7.406    | 97.00  | 139955 | 47.3 ug/L  | 0.031   | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Carbon tetrachloride      | 7.597    | 117.00 | 149015 | 47.3 ug/L  | 0.021   | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | 1-Chlorobutane            | 7.567    | 56.00  | 496602 | 49.0 ug/L  | 0.034   | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Benzene                   | 7.847    | 78.00  | 746521 | 44.9 ug/L  | -0.012  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | 1,2-Dichloroethane        | 7.919    | 62.00  | 255970 | 49.1 ug/L  | 0.034   | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Trichloroethene           | 8.636    | 95.00  | 134002 | 43.0 ug/L  | -0.011  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | 1,2-Dichloropropane       | 8.952    | 63.00  | 230289 | 44.3 ug/L  | 0.031   | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Dibromomethane            | 9.029    | 93.00  | 89920  | 49.1 ug/L  | 0.031   | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Methyl methacrylate       | 8.999    | 89.00  | 200384 | 49.1 ug/L  | -0.011  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Bromodichloromethane      | 9.232    | 83.00  | 205943 | 46.7 ug/L  | 0.029   | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | 2-Nitropropane            | 9.470    | 41.00  | 124536 | 49.1 ug/L  | 0.032   | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Chloroacetonitrile        | 9.715    | 48.00  | 8071   | 49.1 ug/L  | 0.031   | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | cis-1,3-Dichloropropene   | 9.715    | 75.00  | 287309 | 44.8 ug/L  | 0.026   | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Toluene                   | 10.094   | 92.00  | 449097 | 40.8 ug/L  | 0.031   | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | trans-1,3-Dichloropropene | 10.336   | 75.00  | 254137 | 44.5 ug/L  | -0.010  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |
|           | Ethyl methacrylate        | 10.386   | 69.00  | 327140 | 49.1 ug/L  | -0.010  | Auto         | 0.0000      | Auto    | 0.0000 | Auto              |                   |

# Analytical Throughput

## Conceptualizing Analytical Throughput

Instruments  
per laboratory



Injections  
per Time

# Instruments per laboratory area


---

How to increase number of instruments in limited laboratory space?



# Instruments per laboratory area

How to increase number of instruments in limited laboratory space?  
→ Choose Small Foot print GC-MS system




# Instruments per laboratory area


You can double your productivity by using two instruments instead of one.



SVOC



SVOC




# Instruments per laboratory area


You can double your productivity by using two instruments instead of one.



SVOC



Other application



# Sample per time

---

How can we increase the number of injections given in a specific time slot?


# Sample per time

How can we increase the number of injections given in a specific time slot?  
→ **Choose narrow bore column to run sample**




30m x 0.25mm (I.D)

R.T: 19.8 min



20m x 0.18mm (I.D)  
20m x 0.15mm (I.D)


R.T: 10.7 min



# Sample per time

When use Narrow bore column, run time is 50% and injection to injection is 65% compare with conventional column.

TIC (500pg on column) :  
8270 Megamix, Benzidine Mix, 1,4-Dioxane, Benzoic Acid, Atrazine, Benzaldehyde, Biphenyl, Caprolactam



# Sample per time

What should we keep in mind when using a narrow bore column?  
 → **Column Capacity** and Scan rate

30m x 0.25mm (I.D)

20m x 0.15mm (I.D)




20m x 0.15mm (I.D)

Split 10:1

Split 10:1

Split **30:1**

Benzo(a)pyrene (10ppm on vial)



# Sample per time

What should we keep in mind when using a narrow bore column?  
 → Column Capacity and **Scan rate(Scan Speed)**

30m x 0.25mm (I.D)

20m x 0.15mm (I.D)




20m x 0.15mm (I.D)

1,2-Dichlorobenzene

Scan/SIM: 5,000u/sec  
 10 points

Scan/SIM: 5,000u/sec  
 6 points

Scan/SIM: **10,000u/sec**  
 10 points






# Sample per time


How does the latest GC-MS technology support the use of narrow bore columns?  
→ **High sensitivity and Fast Scan/SIM analysis**



DuraEase ion source (section)




Quadrupole rods with Prerod




Ultra-fast scan performance without sacrificing accuracy.

# Sample per time

Using narrow bore column method, we can meet 8270E criteria for all compounds




N-Nitrosodimethylamine (1ppm with Split 30)




Calibration (0.05ppm to 10ppm)

RSE: 6.581% (0.05ppm to 10ppm)




# Sample per time

Using narrow bore column method, we can meet 8270E criteria for all compounds




Benzo(b)fluoranthene (1ppm with Split 30)



Calibration (0.05ppm to 10ppm)


RSE: 5.897% (0.05ppm to 10ppm)



# Sample per time

Using a faster scan speed, the library search scores for all compounds are similar across all methods.

| Method                                  | Scan only n=3 |        | Scan/SIM n=3<br>0.25 sec or less<br>(Scan: 0.1 sec,<br>SIM 0.15 sec or less) |
|-----------------------------------------|---------------|--------|------------------------------------------------------------------------------|
|                                         | 0.3sec        | 0.1sec |                                                                              |
| Loop Time (Sampling rate)               | 0.3sec        | 0.1sec |                                                                              |
| Scan range(m/z)                         | 40-450        | 40-600 | 40-450                                                                       |
| Scan Speed(u/sec)                       | 1,428         | 10,000 | 4,285                                                                        |
| Average Similarity Score (85 compounds) | 92.62         | 91.79  | 91.49                                                                        |



● Scan Speed 1,428u/sec   ● Scan Speed 10,000u/sec   ● Scan Speed 4,285u/sec  
 ※Some compounds are co-eluting and show lower similarity score

# Post Analysis

Conceptualizing post processing in laboratory



# Pre-set report template

## Pre-set report template makes easier ICAL/CCV/sample report

ICAL report



Quant Report (CCV, Sample, Duplicate)

**Insight Report** Printed at 2023/05/17 16:15:04

Result File: C:\Users\Shimadzu\Desktop\EPA Data\5248\LabSolutions\Data\ICAL\_curve\ICAL\_Report\ipoc - EPA524\_2.gm

| Sample Name | Data File          | Date Acquired       | Date Modified |
|-------------|--------------------|---------------------|---------------|
| 0.5 ppb     | 2018/12/28 2:53:08 | 2023/04/24 13:48:11 |               |
| 1.0 ppb     | 2018/12/28 3:13:30 | 2023/04/24 13:48:14 |               |
| 2.0 ppb     | 2018/12/28 4:18:16 | 2023/04/24 13:48:19 |               |
| 5.0 ppb     | 2018/12/28 4:46:41 | 2023/04/24 13:48:21 |               |
| 10.0 ppb    | 2018/12/28 5:15:05 | 2023/04/24 13:48:24 |               |
| 20.0 ppb    | 2018/12/28 5:43:36 | 2023/04/24 13:48:26 |               |
| 50.0 ppb    | 2018/12/28 6:12:03 | 2023/04/24 13:48:29 |               |
| 100.0 ppb   | 2018/12/28 6:40:20 | 2023/04/24 13:48:33 |               |
| 200.0 ppb   | 2018/12/28 7:08:41 | 2023/04/24 13:48:36 |               |


| Name                     | 0.25 ppb | 0.5 ppb | 1.0 ppb |
|--------------------------|----------|---------|---------|
| Dichlorodifluoromethane  | Pass     | Pass    | Pass    |
| Chloromethane            | Pass     | Pass    | Pass    |
| Vinyl chloride           | Pass     | Pass    | Pass    |
| Bromomethane             | Pass     | Pass    | Pass    |
| Chloroethane             | Pass     | Pass    | Pass    |
| Trichlorofluoromethane   | Pass     | Pass    | Pass    |
| 1,1-Dichloroethane       | Pass     | Pass    | Pass    |
| Iodomethane              | Pass     | Pass    | Pass    |
| Carbon disulfide         | Pass     | Pass    | Pass    |
| Aral chloride            | Fail     | Fail    | Fail    |
| Methylene chloride       | Pass     | Pass    | Pass    |
| Acrylonitrile            | Pass     | Pass    | Pass    |
| trans-1,2-Dichloroethane | Pass     | Pass    | Pass    |
| Methyl tert-butyl ether  | Pass     | Pass    | Pass    |
| Acrylonitrile            | Pass     | Pass    | Pass    |
| cis-1,2-Dichloroethane   | Pass     | Pass    | Pass    |
| 2,2-Dichloropropane      | Pass     | Pass    | Pass    |
| Propionitrile            | Pass     | Pass    | Pass    |
| Methyl acrylate          | Pass     | Pass    | Pass    |
| Bromochloromethane       | Pass     | Pass    | Pass    |
| Methyl acrylonitrile     | Pass     | Pass    | Pass    |
| Tetrahydrofuran          | Fail     | Fail    | Fail    |
| Chloroform               | Pass     | Pass    | Pass    |
| 1,1,1-Trichloroethane    | Pass     | Pass    | Pass    |
| Carbon tetrachloride     | Pass     | Pass    | Pass    |
| 1-Chlorobutane           | Fail     | Fail    | Fail    |
| Benzene                  | Pass     | Pass    | Pass    |
| 1,2-Dichloroethane       | Pass     | Pass    | Pass    |
| Fluorobenzene (STD)      | Pass     | Pass    | Pass    |
| Trichloroethene          | Pass     | Pass    | Pass    |
| 1,2-Dichloropropane      | Pass     | Pass    | Pass    |
| Bromomethane             | Fail     | Fail    | Fail    |
| Methyl methacrylate      | Fail     | Fail    | Fail    |
| Bromochloromethane       | Pass     | Pass    | Pass    |
| 2-Nitropropane           | Fail     | Fail    | Fail    |
| Chloroacetonitrile       | Fail     | Fail    | Fail    |
| cis-1,3-Dichloropropene  | Pass     | Pass    | Pass    |
| Toluene                  | Pass     | Pass    | Pass    |



**Insight Report** Printed at 2023/05/17 2:14:01



Result File: C:\Users\Shimadzu\Desktop\EPA Data\5248\LabSolutions\Data\ICAL\_curve\CCV\_Report\ipoc - EPA524\_2.gm

| Type      | Name                      | Found RT | m/z    | Area   | Conc. Unit | RT Diff | Accuracy(%) | Flag Result | Mode ID   | Spiked | Recovery (Rvino)  | Comment 1 |
|-----------|---------------------------|----------|--------|--------|------------|---------|-------------|-------------|-----------|--------|-------------------|-----------|
| STD       | Fluorobenzene (STD)       | 8.213    | 96.00  | 744480 | 50.0 ug/L  | 0.012   | 100.00      | Auto        | 0.000000  | Auto   | 80.000000         |           |
| Target    | Bromofom                  | 12.825   | 173.00 | 97936  | 40.0 ug/L  | -0.009  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
| Surrogate | 4-Bromofluorobenzene_Sur  | 12.895   | 95.00  | 276013 | 50.7 ug/L  | -0.009  | 100.00      | Auto        | 50.000000 | 101.49 | Flag Range 70-130 |           |
|           | Dichlorodifluoromethane   | 1.364    | 85.00  | 43278  | 55.5 ug/L  | -0.004  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Chloroethane              | 1.552    | 50.00  | 150145 | 65.8 ug/L  | -0.007  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Vinyl chloride            | 1.640    | 62.00  | 98183  | 49.1 ug/L  | -0.004  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Bromomethane              | 2.094    | 94.00  | 46989  | 45.7 ug/L  | -0.006  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Chloroethane              | 2.261    | 54.00  | 100939 | 48.6 ug/L  | -0.008  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Trichlorofluoromethane    | 2.762    | 101.00 | 132081 | 45.4 ug/L  | -0.011  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | 1,1-Dichloroethene        | 3.750    | 96.00  | 117865 | 49.1 ug/L  | -0.012  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Iodomethane               | 4.021    | 141.00 | 76652  | 72.5 ug/L  | 0.021   | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Carbon disulfide          | 4.137    | 76.00  | 364893 | 44.6 ug/L  | -0.013  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Allyl chloride            | 4.587    | 76.00  | 88601  | 44.4 ug/L  | -0.012  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Methylene chloride        | 4.737    | 84.00  | 146738 | 41.2 ug/L  | -0.012  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Acrylonitrile             | 5.176    | 53.00  | 193981 | 48.1 ug/L  | -0.012  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | trans-1,2-Dichloroethane  | 5.232    | 96.00  | 143968 | 48.8 ug/L  | -0.011  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Methyl tert-butyl ether   | 5.240    | 73.00  | 528206 | 43.0 ug/L  | -0.012  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | 1,1-Dichloroethane        | 5.937    | 63.00  | 356770 | 45.7 ug/L  | -0.011  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | cis-1,2-Dichloroethane    | 6.756    | 98.00  | 110233 | 46.8 ug/L  | 0.033   | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | 2,2-Dichloropropane       | 6.756    | 77.00  | 162752 | 35.8 ug/L  | -0.011  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Propionitrile             | 6.831    | 54.00  | 64114  | 50.1 ug/L  | -0.012  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Methyl acetate            | 6.878    | 55.00  | 309860 | 46.2 ug/L  | -0.011  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Bromochloromethane        | 7.066    | 128.00 | 63238  | 45.6 ug/L  | -0.011  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Methyl acrylonitrile      | 7.091    | 67.00  | 140705 | 47.3 ug/L  | -0.012  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Tetrahydrofuran           | 7.079    | 71.00  | 41466  | 47.8 ug/L  | -0.011  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Chloroform                | 7.219    | 83.00  | 258868 | 45.9 ug/L  | -0.011  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | 1,1,1-Trichloroethane     | 7.266    | 97.00  | 139955 | 47.2 ug/L  | 0.011   | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Carbon tetrachloride      | 7.597    | 177.00 | 149075 | 47.1 ug/L  | 0.021   | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | 1-Chlorobutane            | 7.567    | 56.00  | 496602 | 44.8 ug/L  | 0.034   | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Benzene                   | 7.847    | 78.00  | 746531 | 44.8 ug/L  | -0.012  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | 1,2-Dichloroethane        | 7.919    | 62.00  | 255970 | 48.1 ug/L  | 0.034   | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Trichloroethene           | 8.636    | 95.00  | 134002 | 43.9 ug/L  | -0.011  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | cis-1,2-Dichloropropene   | 8.932    | 63.00  | 232389 | 44.2 ug/L  | 0.031   | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Dibromomethane            | 9.029    | 93.00  | 89920  | 40.7 ug/L  | 0.031   | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Methyl methacrylate       | 9.999    | 69.00  | 203094 | 40.9 ug/L  | -0.011  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Bromodichloromethane      | 9.232    | 83.00  | 209483 | 46.7 ug/L  | 0.029   | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | 2-Nitropropane            | 9.470    | 43.00  | 124536 | 49.1 ug/L  | 0.032   | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Chloroacetonitrile        | 9.715    | 48.00  | 8071   | 49.1 ug/L  | 0.031   | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | cis-1,3-Dichloropropene   | 9.715    | 75.00  | 287399 | 44.8 ug/L  | 0.036   | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Toluene                   | 10.084   | 92.00  | 449097 | 40.8 ug/L  | 0.031   | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | trans-1,3-Dichloropropene | 10.336   | 75.00  | 254127 | 44.6 ug/L  | -0.010  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |
|           | Methyl chloroacetate      | 10.386   | 69.00  | 327480 | 40.7 ug/L  | -0.010  | 100.00      | Auto        | 0.000000  | Auto   | 0.000000          |           |




# Report export

Directly exporting table from the processing software helps integrate the LIMS system.





| Name           | Found RT | Conc.  | Accuracy | Area    | Recovery | Spiked | RF %RSD  | ISTD %Dif | ISTD %Dif | Status  | Actual Ret Unit | m/z | ISTD Area | Area Ratio | Type   | %RSD (CoI | ISTD Grou | Tailing F. | Asymmetri S/ |
|----------------|----------|--------|----------|---------|----------|--------|----------|-----------|-----------|---------|-----------------|-----|-----------|------------|--------|-----------|-----------|------------|--------------|
| 1 Chlorome     | 1.738    | 8.805  | 44.02    | 853217  | 44.02    | 20     | 152.988  | 100.58    | 100       | Pending | 0 ppb           | 50  | 3374229   | 0.253      | Target | 1         | 0         |            |              |
| 2 Vinyl Chlc   | 1.823    | 9.386  | 46.93    | 679438  | 46.93    | 20     | 150.7404 | 100.58    | 100       | Pending | 0 ppb           | 62  | 3374229   | 0.201      | Target | 1         | 0         |            |              |
| 3 Bromome      | 2.074    | 11.238 | 56.19    | 242271  | 56.19    | 20     | 146.087  | 100.58    | 100       | Pending | 0 ppb           | 94  | 3374229   | 0.072      | Target | 1         | 0         |            |              |
| 4 Chloroeth    | 2.153    | 9.55   | 47.75    | 373766  | 47.75    | 20     | 152.1311 | 100.58    | 100       | Pending | 0 ppb           | 64  | 3374229   | 0.111      | Target | 1         | 0         |            |              |
| 5 Trichlorof   | 2.331    | 9.402  | 47.01    | 435694  | 47.01    | 20     | 148.4978 | 100.58    | 100       | Pending | 0 ppb           | 101 | 3374229   | 0.129      | Target | 1         | 0         |            |              |
| 6 Acrolein     | 2.614    | 8.718  | 43.59    | 158415  | 43.59    | 20     | 157.5276 | 100.58    | 100       | Pending | 0 ppb           | 56  | 3374229   | 0.047      | Target | 1         | 0         |            |              |
| 7 1,1-Dichlo   | 2.695    | 5.992  | 47.96    | 867141  | 47.96    | 20     | 155.1256 | 100.58    | 100       | Pending | 0 ppb           | 61  | 3374229   | 0.257      | Target | 1         | 0         |            |              |
| 8 Acetone      | 2.699    | 11.146 | 55.73    | 326257  | 55.73    | 20     | 160.2486 | 100.58    | 100       | Pending | 0 ppb           | 43  | 3374229   | 0.097      | Target | 1         | 0         |            |              |
| 9 Methylen     | 3.023    | 9.217  | 46.09    | 877011  | 46.09    | 20     | 152.1939 | 100.58    | 100       | Pending | 0 ppb           | 49  | 3374229   | 0.26       | Target | 1         | 0         |            |              |
| 10 Acrylonitr  | 3.166    | 8.771  | 43.86    | 344527  | 43.86    | 20     | 158.691  | 100.58    | 100       | Pending | 0 ppb           | 53  | 3374229   | 0.102      | Target | 1         | 0         |            |              |
| 11 trans-1,2-l | 3.208    | 9.33   | 46.65    | 929175  | 46.65    | 20     | 154.6388 | 100.58    | 100       | Pending | 0 ppb           | 61  | 3374229   | 0.275      | Target | 1         | 0         |            |              |
| 12 1,1-Dichlo  | 3.503    | 9.252  | 46.26    | 1073240 | 46.26    | 20     | 153.1664 | 100.58    | 100       | Pending | 0 ppb           | 63  | 3374229   | 0.318      | Target | 1         | 0         |            |              |
| 13 Ethyl Acet  | 3.867    | 7.321  | 36.61    | 1085266 | 36.61    | 20     | 154.1255 | 100.58    | 100       | Pending | 0 ppb           | 43  | 3374229   | 0.322      | Target | 1         | 0         |            |              |
| 14 Chlorofori  | 4.109    | 8.26   | 41.3     | 641997  | 41.3     | 20     | 151.2245 | 100.58    | 100       | Pending | 0 ppb           | 83  | 3374229   | 0.19       | Target | 1         | 0         |            |              |
| 15 1,1,1-Tric  | 4.254    | 9.32   | 46.6     | 463597  | 46.6     | 20     | 150.7924 | 100.58    | 100       | Pending | 0 ppb           | 97  | 3374229   | 0.137      | Target | 1         | 0         |            |              |
| 16 Carbon Te   | 4.367    | 9.757  | 48.78    | 374614  | 48.78    | 20     | 151.5158 | 100.58    | 100       | Pending | 0 ppb           | 117 | 3374229   | 0.111      | Target | 1         | 0         |            |              |
| 17 Isopropyl   | 4.475    | 2.487  | 12.43    | 1401670 | 12.43    | 20     | 149.2    | 100.58    | 100       | Pending | 0 ppb           | 43  | 3374229   | 0.415      | Target | 1         | 0         |            |              |
| 18 1,2-Dichlo  | 4.479    | 29.567 | 98.56    | 957214  | 98.56    | 30     | 0.518219 | 100.58    | 100       | Pending | 0 ppb           | 65  | 3374229   | 0.284      | Target | 1         | 0         |            |              |
| 19 Benzene     | 4.51     | 5.461  | 27.31    | 2095380 | 27.31    | 20     | 196.674  | 100.58    | 100       | Pending | 0 ppb           | 78  | 3374229   | 0.621      | Target | 1         | 0         |            |              |
| 20 1,2-Dichlo  | 4.535    | 9.268  | 46.34    | 688200  | 46.34    | 20     | 152.6257 | 100.58    | 100       | Pending | 0 ppb           | 62  | 3374229   | 0.204      | Target | 1         | 0         |            |              |
| 21 Fluorober   | 4.712    | 30     | 100      | 3374229 |          | 0      |          | 100.58    | 100       |         | 0 ppb           | 96  | 3374229   |            | ISTD   | 1         | 0         |            |              |
| 22 Trichloroe  | 4.981    | 9.172  | 45.86    | 509510  | 45.86    | 20     | 150.1063 | 101.45    | 100       | Pending | 0 ppb           | 130 | 1744723   | 0.292      | Target | 2         | 0         |            |              |
| 23 1,2-Dichlo  | 5.183    | 1.011  | 5.05     | 658430  | 5.05     | 20     | 0        | 101.45    | 100       | Pending | 0 ppb           | 63  | 1744723   | 0.377      | Target | 2         | 0         |            |              |
| 24 Bromodici   | 5.371    | 9.666  | 48.33    | 460888  | 48.33    | 20     | 152.5379 | 101.45    | 100       | Pending | 0 ppb           | 83  | 1744723   | 0.264      | Target | 2         | 0         |            |              |
| 25 2-Chloroe   | 5.567    | 8.18   | 40.9     | 397757  | 40.9     | 20     | 138.6593 | 101.45    | 100       | Pending | 0 ppb           | 43  | 1744723   | 0.228      | Target | 2         | 0         |            |              |
| 26 cis-1,3-Dic | 5.715    | 9.447  | 47.23    | 668830  | 47.23    | 20     | 148.4615 | 101.45    | 100       | Pending | 0 ppb           | 75  | 1744723   | 0.383      | Target | 2         | 0         |            |              |
| 27 Toluene-C   | 5.944    | 29.521 | 98.4     | 3374795 | 98.4     | 30     | 1.444597 | 101.45    | 100       | Pending | 0 ppb           | 98  | 1744723   | 1.934      | Target | 2         | 0         |            |              |
| 28 Toluene     | 6.002    | 9.337  | 46.68    | 2187715 | 46.68    | 20     | 157.3203 | 101.45    | 100       | Pending | 0 ppb           | 91  | 1744723   | 1.254      | Target | 2         | 0         |            |              |
| 29 trans-1,3-l | 6.174    | 9.566  | 47.83    | 543315  | 47.83    | 20     | 149.0686 | 101.45    | 100       | Pending | 0 ppb           | 75  | 1744723   | 0.311      | Target | 2         | 0         |            |              |
| 30 2-Bromo-    | 6.303    | 30     | 100      | 1744723 |          | 0      |          | 101.45    | 100       |         | 0 ppb           | 41  | 1744723   |            | ISTD   | 2         | 0         |            |              |
| 31 1,1,2-Tric  | 6.348    | 6.111  | 30.55    | 413880  | 30.55    | 20     | 201.7038 | 94.55     | 100       | Pending | 0 ppb           | 97  | 2204193   | 0.188      | Target | 3         | 0         |            |              |
| 32 Tetrachlor  | 6.443    | 9.565  | 47.82    | 602621  | 47.82    | 20     | 153.0763 | 94.55     | 100       | Pending | 0 ppb           | 166 | 2204193   | 0.273      | Target | 3         | 0         |            |              |
| 33 Dibromoc    | 6.689    | 10.144 | 50.72    | 352658  | 50.72    | 20     | 156.1117 | 94.55     | 100       | Pending | 0 ppb           | 129 | 2204193   | 0.16       | Target | 3         | 0         |            |              |



# Speedup data processing

Manual peak identification/integration sometimes need to finish processing.





# AI-Assisted peak integration





AI assisted peak integration accelerate data processing

**Conventional Shimadzu Algorithm**

The simultaneous analysis of hundreds of regulated compounds. Adjusting the peak integration parameters for each compound is a tremendous task.

**Peakintelligence for GCMS**

Accurate peak integration is achieved without configuring the parameter settings.




# Focusing data

By flagging the data, the focus data becomes clearer, which improves processing efficiency.

| #                                   | Fla... | Flag ID   | Type      | Name                        | Conc.        | ISTD %Diff (Selected) Area | ISTD %Diff (CC) Area | Recovery (Enviro) | RPD    | Flag Result |
|-------------------------------------|--------|-----------|-----------|-----------------------------|--------------|----------------------------|----------------------|-------------------|--------|-------------|
| <input checked="" type="checkbox"/> |        |           |           |                             |              |                            |                      |                   |        |             |
| <input checked="" type="checkbox"/> |        |           | Target    | 1,1,1,2-Tetrachloroethane   | 44.822076... | 93.77                      | 96.75                | 89.64             | 3.45   | Pass        |
| <input checked="" type="checkbox"/> |        |           | Target    | Xylene Total                | 127.18399... | 93.77                      | 96.75                | 84.79             | 2.32   | Pass        |
| <input checked="" type="checkbox"/> |        | RPD       | Target    | m/p-Xylene                  | 120.05264... | 93.77                      | 96.75                | 120.05            | 36.68  | Fail        |
| <input checked="" type="checkbox"/> |        |           | Target    | o-Xylene                    | 42.925714... | 93.77                      | 96.75                | 85.85             | 3.75   | Pass        |
| <input checked="" type="checkbox"/> |        | rEnv, RPD | Target    | Styrene                     | 14.893732... | 93.77                      | 96.75                | 29.79             | 95.73  | Fail        |
| <input checked="" type="checkbox"/> |        |           | Target    | Bromoform                   | 40.945274... | 93.77                      | 96.75                | 81.56             | 8.53   | Pass        |
| <input checked="" type="checkbox"/> |        |           | Target    | Isopropylbenzene            | 41.973292... | 93.77                      | 96.75                | 83.95             | 2.88   | Pass        |
| <input checked="" type="checkbox"/> |        | rEnv      | Surrogate | 4-Bromofluorobenzene_Su...  | 50.743044... | 93.77                      | 96.75                | 169.14            | 0.16   | Fail        |
| <input checked="" type="checkbox"/> |        |           | Target    | 1,1,2,2-Tetrachloroethane   | 44.296847... | 93.77                      | 96.75                | 88.59             | 7.62   | Pass        |
| <input checked="" type="checkbox"/> |        |           | Target    | Bromobenzene                | 43.310274... | 93.77                      | 96.75                | 86.62             | 4.91   | Pass        |
| <input checked="" type="checkbox"/> |        |           | Target    | trans-1,4-Dichloro-2-butene | 40.251729... | 93.77                      | 96.75                | 80.50             | 13.03  | Pass        |
| <input checked="" type="checkbox"/> |        |           | Target    | 1,2,3-Trichloropropane      | 44.596189... | 93.77                      | 96.75                | 89.19             | 9.61   | Pass        |
| <input checked="" type="checkbox"/> |        |           | Target    | n-Propylbenzene             | 41.823714... | 93.77                      | 96.75                | 83.65             | 3.77   | Pass        |
| <input checked="" type="checkbox"/> |        |           | Target    | 2-Chlorotoluene             | 42.453903... | 93.77                      | 96.75                | 84.91             | 4.66   | Pass        |
| <input checked="" type="checkbox"/> |        |           | Target    | 4-Chlorotoluene             | 44.608940... | 93.77                      | 96.75                | 89.22             | 4.14   | Pass        |
| <input checked="" type="checkbox"/> |        |           | Target    | 1,2,4-Trimethylbenzene      | 42.331516... | 93.77                      | 96.75                | 84.66             | 5.20   | Pass        |
| <input checked="" type="checkbox"/> |        | RPD       | Target    | tert-Butylbenzene           | 41.404369... | 93.77                      | 96.75                | 82.81             | 120.69 | Fail        |
| <input checked="" type="checkbox"/> |        |           | Target    | 1,3,5-Trimethylbenzene      | 42.672844... | 93.77                      | 96.75                | 85.35             | 4.97   | Pass        |
| <input checked="" type="checkbox"/> |        |           | Target    | Pentachloroethane           | 46.455152... | 93.77                      | 96.75                | 92.91             | 4.85   | Pass        |
| <input checked="" type="checkbox"/> |        |           | Target    | sec-Butylbenzene            | 41.726976... | 93.77                      | 96.75                | 83.45             | 5.58   | Pass        |

# Instrument downtime

Conceptualizing downtime in laboratory



# Keep running instrument

How can we keep the instrument running without interruption?  
 → **Use BFB/DFTPP tuning algorithm with proven robustness**

Tuning Information

Target Condition

Adjust Resolution  
 FWHM of Peak Profile: 0.60

Adjust Sensitivity  
 Target Mass: 264



Calibrate Mass

Adjust Mass Pattern

| m/z                                     | Inten. Ratio(%) | m/z                                     | Inten. Ratio(%) |
|-----------------------------------------|-----------------|-----------------------------------------|-----------------|
| <input checked="" type="checkbox"/> 50  | 1.00            | <input checked="" type="checkbox"/> 69  | 100.00          |
| <input checked="" type="checkbox"/> 131 | 50.00           | <input checked="" type="checkbox"/> 219 | 80.00           |
| <input checked="" type="checkbox"/> 414 | 6.00            | <input checked="" type="checkbox"/> 502 | 9.00            |

Initialize

OK Cancel




Even using 524.2 BFB criteria which know well as strict criteria, BFB tuning algorithm works over 800 injections (Standards and samples)

# Reduce maintenance frequency


How to reduce maintenance frequency?

→ **Use robust MS instrument without compromising effectiveness**

## DuraEase ion source (section)



In the continuous analysis of VOCs in water, **the effect of elution of the column liquid phase was eliminated to the extreme, and very stable measurement was realized.**





# Reduce maintenance time

Intuitive/less step operation reduce Instrument downtime



# Reduce maintenance time

Intuitive/less step operation reduce Instrument downtime



# Reduce maintenance time

Intuitive/less step operation reduce Instrument downtime







# Reduce Human Error

A simplified design makes maintenance easier for operators, and accessible instructions from anywhere help reduce human error.


|                | Latest GC-MS   | Conventional GC-MS |
|----------------|----------------|--------------------|
| Parts (PC)     | 7              | 10                 |
| Tools (PC)     | 3              | 5                  |
| Assemble (Min) | 1              | 10                 |
| Cleaning (Min) | 0 (Disposable) | 15                 |



Easy access to maintenance help



# Take home message




Integrating a compact design with rapid Scan/SIM analysis improves injection efficiency in terms of area and time.

AI-assisted peak integration and visualization software streamline data processing, reducing manual tasks and enhancing processing efficiency

Highly durable instruments minimize maintenance needs, while intuitive maintenance procedures prevent human error-related issues

# Q&A

|               |                                           |
|---------------|-------------------------------------------|
| HPLC          | Carbamate<br>Diquat<br>Glyphosate         |
| IC            | Anions<br>Chromium VI<br>Ammonia          |
| TOC           | Nitrogen<br>Organic Carbon<br>Phosphorous |
| FTIR & Raman  | <b>Microplastics</b>                      |
| UV-Vis        | Chlorine<br>UV254<br>Others               |
| Particle Size | <b>Solids</b>                             |



|                                                                                    |       |
|------------------------------------------------------------------------------------|-------|
| Cyanotoxins<br>PFAS<br>Pesticides<br>Emerging Contaminants<br>Unknowns (QTOF)      | LCMS  |
| HAA's<br>Herbicides<br>Pesticides                                                  | GC    |
| Phenols<br>Other Organics<br>PCBs                                                  | GCMS  |
| 1,4-dioxane<br>Dioxins & Furans<br>Nitrosamines<br>Other Organics<br>Microplastics | GCMS  |
| Semivolatiles<br>Taste & Odor<br>Volatiles<br>THMs                                 | ICPMS |
| <b>Metals</b>                                                                      | ICPMS |
| <b>Metals</b>                                                                      | ICP   |
| <b>Metals</b>                                                                      | AA    |

For any questions, contact:  
**Yoshiro Hiramatsu**  
[yohiramatsu@shimadzu.com](mailto:yohiramatsu@shimadzu.com)

For more information, visit:  
[www.OneLabOneEarth.com](http://www.OneLabOneEarth.com)

Connect with us:

X - [@shimadzussi](https://twitter.com/shimadzussi)

Instagram - [@shimadzussi](https://www.instagram.com/shimadzussi)

LinkedIn - [/company/shimadzu-scientific-instruments/](https://www.linkedin.com/company/shimadzu-scientific-instruments/)

YouTube - [@ShimadzuScientificInstruments](https://www.youtube.com/@ShimadzuScientificInstruments)