


**Simplification of laboratory workflows for the analysis of common contaminants as required per EPA methods 8270E, 625, Appendix IX, 8081B, 8141B by the use of gas chromatography mass spectrometry**

**Andy Fornadel**

Regional Product Marketing Manager  
Thermo Fisher Scientific

 The world leader in serving science




# Overview


- **Challenges in the laboratory**
- **EPA 8081:**
  - Advantages of GC-MS/MS over GC ECD
  - How GC-MS/MS improves selectivity
  - Data examples including Method detection limits
- **EPA 8270E**
  - Compliant tuning
  - Critical pair separation
- **Consolidating EPA methods:**
  - Single run covering multiple methods
  - Helium saving technology
  - Single software across platforms
- **Conclusions**



# Challenges faced by environmental testing laboratories



# Modernizing workflows onto a common platform



- **Single platform for modernized workflow**

- ✓ **Same GC-MS/MS Platform**

TSQ9610 NV-AEI

TRACE 1610 GC w/ SSL

TriPlus RSH SMART or AS1610 Autosampler

- ✓ **Same Column**

TG-5SilMS 30 m × 0.25 mm ID x 0.25 μm

w/ 5 m × 0.25 mm ID GuardGold

- ✓ **Simplified software, operation, data processing and reporting**

Chromeleon 7.3.2

- ✓ **Flexibility to choose how and where to run your samples**





# Single platform and configuration



- **Safer Lab Hardware**
  - No more ECDs
  - Forget about licensing fees, training, and safety concerns
- **Modern Methods**
  - Ditch the dual column methods
  - Reduce consumables cost and complexity
  - Faster calibration, analysis, and data processing
  - Optimized each SVOC method for seamless method transfer
- **Confident Results**
  - Fast, confident identification at ppb sensitivities

## Study of **EPA 8081** on GC-ECD versus GC-MS/MS



- Dual Column
- Dual Detector



- Single Column
- MS/MS Detector

# Improving Productivity through Data Quality

## Interference on GC-ECD

- Leads to systematic errors
  - Biases
- Can raise LODs
- Reduce signal/noise




## Selectivity to resolve interferences on TSQ 9610

- MS mass filtering to target specific precursor ions
- Fragmentation
- MS mass filtering of product ions
- All yields higher selectivity than ECD, reducing LODs, and increasing S/N

# How does GC-MS/MS generate better data?

## Single Quadrupole Detection


- Sensitivity in SIM (Selected Ion Monitoring)



**SIM on  $m/z$  182**


## Triple Quadrupole Detection

- Sensitivity and Selectivity in SRM (Selected Reaction Monitoring)




**SRM on  $m/z$  182 > 153**

Single Quadrupole GC-MS



Triple Quadrupole GC-MS/MS






# Improving Productivity through Data Quality

## Conventional Detector Confirmation

- Retention Time on two columns of different phase



*Methoxychlor*




## Mass Spectrometer Detector Confirmation

- Retention Time
- Unique Quantitation Ion
- Unique Confirmation ion(s)
- Ion ratio (Quantitation to Confirmation ion(s))
- Library Match in Full Scan



*Methoxychlor*



# EPA 8081B Chromatography




# EPA 8081B Method detection limits




**10 replicates for Aldrin at 0.5 µg/L**

# EPA 8081B Lowest calibration point





# EPA 8270E: Built in tuning compliance



# EPA 8270E: Built in tuning compliance

### MS Tune Check Report (DFTPP)

Sequence Details					
Sequence Name:	20240221 - A Tune and ISTD Replicates	Created On:	21/Feb/24 17:53:39		
Directory:	App Note ENV	Created By:	Thermo		
Data Vault:	2024	Updated On:	26/Mar/24 11:07:17		
No. of Injections:	35	Updated By:	sarah.crumlett		
Calibration Source:	20240221 - A Tune and ISTD Replicates				

Injection Details					
Injection Name:	Tune Check	Injection Volume:			
Analysis Type:	MS Tune Check	Dilution Factor:			
Injection Type:	Unknown	Vial Number:			
Injection Date/Time:	21/Feb/24 18:50	Processing Method:			
		Instrument Method:			

Apex Decafluorotriphenylphosphine (surr) Scan: #1849 AV: 6.58 - 6.58 min (3) NL: 6.33E+008

No.	Name	Eval. Result	Operator	Ref. Value 1	Ref. V
1	m/z 68 - Less than 2% of m/z 69	0.3	<	2	
2	m/z 69 - Present	51184357.7	>	0	
3	m/z 70 - Less than 2% of m/z 69	0.6	<	2	
4	m/z 197 - Less than 2% of m/z 198	0.0	<	2	
5	m/z 198 - Base Peak or Present	632495297.7	>	0	
6	m/z 199 - 5 to 9% of m/z 198	5.7	between	5	
7	m/z 365 - Greater than 1% of m/z 198	2.1	>	1	
8	m/z 441 - Less than 150% of m/z 443	121.5	<	150	
9	m/z 442 - Base Peak or Present	611513468.7	>	0	
10	m/z 443 - 15 to 24% of m/z 442	18.1	between	15	

### Breakdown Report

Sequence Details					
Sequence Name:	20240221 - A Tune and ISTD Replicates	Created On:	21/Feb/24 17:53:39		
Directory:	App Note ENV	Created By:	Thermo		
Data Vault:	2024	Updated On:	24/Feb/24 09:14:25		
No. of Injections:	35	Updated By:	sarah.crumlett		
Calibration Source:	20240221 - A Tune and ISTD Replicates				

Injection Details					
Injection Name:	Tune Check	Injection Volume:	1.00		
Analysis Type:	n.a.	Dilution Factor:	1.0000		
Injection Type:	Unknown	Vial Number:	2		
Injection Date/Time:	21/Feb/24 18:50	Processing Method:	EPA 8270E - DFTPP & Tailing V30		
		Instrument Method:	DFTPP		

#### DDT Breakdown

Compound	Type	Response	%Breakdown	Max %Breakdown	Result
4,4'-DDT	Native	81546587	3.345	20	Pass
4,4'-DDE	Breakdown	130821			
4,4'-DDD	Breakdown	2691157			

#### Endrin Breakdown

Compound	Type	Response	%Breakdown	Max %Breakdown	Result
Endrin	Native	72779601	5.933	20	Pass
Endrin aldehyde	Breakdown	509389			
Endrin ketone	Breakdown	4080643			

# EPA 8270E: Built in tuning compliance

### MS Tune Check Report (DFTPP)

Sequence Details					
Sequence Name:	20240221 - A Tune and ISTD Replicates	Created On:	21/Feb/24 17:53:39		
Directory:	App Note ENV	Created By:	Thermo		
Data Vault:	2024	Updated On:	26/Mar/24 11:07:17		
No. of Injections:	35	Updated By:	sarah.crumlett		
Calibration Source:	20240221 - A Tune and ISTD Replicates				

Injection Details					
Injection Name:	Tune Check	Injection Volume:			
Analysis Type:	MS Tune Check	Dilution Factor:			
Injection Type:	Unknown	Vial Number:			
Injection Date/Time:	21/Feb/24 18:50	Processing Method:			

Apex Decafluorotriphenylphosphine (surr) Scan: #1849 AV: 6.58 - 6.58 min (3) NL: 6.33E+008

No.	Name	Eval. Result	Operator	Ref. Value 1	Ref. V
1	m/z 68 - Less than 2% of m/z 69	0.3	<	2	
2	m/z 69 - Present	51184357.7	>	0	
3	m/z 70 - Less than 2% of m/z 69	0.6	<	2	
4	m/z 197 - Less than 2% of m/z 198	0.0	<	2	
5	m/z 198 - Base Peak or Present	632495297.7	>	0	
6	m/z 199 - 5 to 9% of m/z 198	5.7	between	5	
7	m/z 365 - Greater than 1% of m/z 198	2.1	>	1	
8	m/z 441 - Less than 150% of m/z 443	121.5	<	150	
9	m/z 442 - Base Peak or Present	611513468.7	>	0	
10	m/z 443 - 15 to 24% of m/z 442	18.1	between	15	

### Breakdown Report

Sequence Details					
Sequence Name:	20240221 - A Tune and ISTD Replicates	Created On:	21/Feb/24 17:53:39		
Directory:	App Note ENV	Created By:	Thermo		
Data Vault:	2024	Updated On:	24/Feb/24 09:14:25		
No. of Injections:	35	Updated By:	sarah.crumlett		
Calibration Source:	20240221 - A Tune and ISTD Replicates				

Injection Details					
Injection Name:	Tune Check	Injection Volume:			
Analysis Type:	n.a.	Dilution Factor:			
Injection Type:	Unknown	Vial Number:			
Injection Date/Time:	21/Feb/24 18:50	Processing Method:			

### Peak Tailing Report

Sequence Details					
Sequence Name:	20240221 - A Tune and ISTD Replicates	Created On:	21/Feb/24 17:53:39		
Directory:	App Note ENV	Created By:	Thermo		
Data Vault:	2024	Updated On:	24/Feb/24 09:14:25		
No. of Injections:	35	Updated By:	sarah.crumlett		
Calibration Source:	20240221 - A Tune and ISTD Replicates				


  

Injection Details					
Injection Name:	Tune Check	Injection Volume:	1.00		
Analysis Type:	n.a.	Dilution Factor:	1.0000		
Injection Type:	Unknown	Vial Number:	2		
Injection Date/Time:	21/Feb/24 18:50	Processing Method:	EPA 8270E - DFTPP & Tailing V30		
		Instrument Method:	DFTPP		

Compound	Ret Time min	Asymmetry (AIA)	Limit	Outside Limit?
TIC	TIC	TIC	TIC	
Pentachlorophenol	6.201	0.97	2.00	
Benzidine	7.466	1.00	2.00	

# EPA 8270E separation of structural isomers







# EPA 8270E Critical Pair Resolution

Isomers must be at least 50% resolved on the mid-point of the ICAL


- Benzo(a)anthracene and Chrysene



- Benzo(b)fluoranthene and Benzo(k)fluoranthene




- Phenanthrene and Anthracene




# Consolidating methods Improving Productivity

## Better results, less time

- Consolidate methods into single injection
  - In some cases, common SVOC workflows can be consolidated into a single injection
  - 336 SVOCs analyzed in a single run
  - PCBs, pesticides, and PAHs/Phenols



Timed SRM



# Combining EPA methods Improving Productivity


Better results, less time

Peak Name	Cal.Type	R <sup>2</sup> RSD (%)	R <sup>2</sup> / %RSD Failed?	Manually Integrated?
MS Quantitation	MS Quantitation	MS Quantitation Peak	MS Quantitation	
Hexachlorocyclopentadiene	AvgCalFact	5.85 (%RSD)		
TCMX (S)	AvgCalFact	7.55 (%RSD)		
alpha-BHC	AvgCalFact	7.21 (%RSD)		
Hexachlorobenzene	AvgCalFact	6.74 (%RSD)		
beta-BHC	AvgCalFact	6.68 (%RSD)		
gamma-BHC	AvgCalFact	6.73 (%RSD)		
delta-BHC	AvgCalFact	5.82 (%RSD)		
Heptachlor	AvgCalFact	7.07 (%RSD)		
Aldrin	AvgCalFact	8.75 (%RSD)		
Heptachlor Epoxide	AvgCalFact	7.71 (%RSD)		
gamma-Chlordane	AvgCalFact	8.72 (%RSD)		
alpha-Chlordane	AvgCalFact	6.31 (%RSD)		
Endosulfan I	AvgCalFact	4.58 (%RSD)		
4,4-DDE	AvgCalFact	5.15 (%RSD)		
Dieldrin	AvgCalFact	8.38 (%RSD)		
Endrin	AvgCalFact	5.22 (%RSD)		
4,4-DDD	AvgCalFact	5.55 (%RSD)		
Endosulfan II	AvgCalFact	10.35 (%RSD)		
Endrin Aldehyde	AvgCalFact	12.57 (%RSD)		
4,4-DDT	AvgCalFact	8.57 (%RSD)		
Endosulfan Sulfate	AvgCalFact	8.22 (%RSD)		
Methoxychlor	AvgCalFact	9.99 (%RSD)		
Endrin Ketone	AvgCalFact	12.39 (%RSD)		
DCBP (S)	AvgCalFact	18.38 (%RSD)		
		Number of failed tests:	0	0



Compound	Ret. Time min	Quantitation Peak (min)	Response	Cal Amount	Amount
MS Quantitation	MS Quantitation	MS Quantitation	MS Quantitation	MS Quantitation	MS Quantitation
Hexachlorocyclopentadiene	2.828	143.000	519	1.00	0.97
alpha-BHC	4.138	183.000	513	1.00	0.93
Hexachlorobenzene	4.189	248.800	1450	1.00	1.02
beta-BHC	4.299	145.000	741	1.00	1.07
gamma-BHC	4.394	145.000	773	1.00	0.94
delta-BHC	4.576	147.000	346	1.00	0.99
Heptachlor	5.003	237.000	817	1.00	0.93
Aldrin	5.316	191.000	239	1.00	0.94
Heptachlor Epoxide	5.501	283.000	195	1.00	0.88
gamma-Chlordane	5.778	237.000	270	1.00	0.90
alpha-Chlordane	5.875	266.000	272	1.00	0.97
Endosulfan I	5.988	206.000	104	1.00	0.99
4,4-DDE	5.977	248.000	596	1.00	0.98
Dieldrin	6.077	193.000	212	1.00	1.13
Endrin	6.228	173.000	63	1.00	1.07
4,4-DDD	6.266	199.000	514	1.00	1.04
Endosulfan II	6.288	123.000	58	1.00	1.01
Endrin Aldehyde	6.385	138.000	68	1.00	1.24
4,4-DDT	6.534	165.000	1837	1.00	1.14
Endosulfan Sulfate	6.545	237.000	478	1.00	1.02
Methoxychlor	6.903	169.000	2867	3.00	3.25
Endrin Ketone	6.922	281.000	50	1.00	1.13

Calibration results from 1-60 ppb using TSQ 9610 (left) and 1 ppb check (right)

- Unambiguous results for known and unknowns
- GC-MS/MS provides equivalent or better data than GC-ECD without the challenges
- Positive identification in a single run
- Better data reduces time spent on data review




Calibration curve for heptachlor



Comparison of 60 ppb analytes using ECD (left) and MS (right)

# Sustainable operation




PCBs analyzed using TSQ 9610 with HeSaver/H2Safer at 0.2 ng/mL, extending lifetime of He cylinder by ~4.5x

## Reduce carrier gas usage

- A unique solution for reducing helium usage, conversion to hydrogen
- Conserve carrier gas during analysis and standby
- Reduce helium consumption by up to 70% without revalidating methods
- Safely convert to hydrogen



	Helium consumption with standard SSL	Helium consumption with HeSaver SSL	Gas savings
Daily He usage	44.0 L	12.5 L	-72%
He cylinder life	~7.5 months	~2 years	



Deltamethrin full scan spectra using He (top) and H<sub>2</sub> (bottom)



# Simplified data review: Single software

Same Data Management Workflows for IC, LC, GC, and related MS workflows



**Interactive Results**

Peak Name	B	C	D	E	F	G	H	I	J	K	L	M
434 flutianil	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ
435 Etrimphos	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ
436 Fluthiacet-methyl	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ
437 Benzovindiflupyr	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ
438 Fipronil sulfide	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	n.a.	<LQ
439 <b>Tebuconazole</b>	<b>8.9</b>	<b>40.5</b>	<b>9.0</b>	<b>40.7</b>	<b>41.3</b>	<b>42.6</b>	<b>25.1</b>	<b>26.5</b>	<b>41.4</b>	<b>43.2</b>	<b>41.4</b>	<b>41.9</b>
440 penthiopyrad	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.
441 Benalaxyl	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ
442 Fenoprop	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.
443 Fonofos	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ
444 Propiconazole	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ
445 Prochloraz Metabolite BTS 4	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ
446 mandestrobin	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ
447 Isofenphos	<LQ	<LQ	<LQ	n.a.	<LQ	<LQ	n.a.	<LQ	n.a.	<LQ	n.a.	<LQ
448 Anilofos	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.
449 Coumaphos	n.a.	n.a.	n.a.	n.a.	n.a.	<LQ	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.
450 Pirimiphos-Methyl	17.5	82.8	17.8	88.5	85.6	90.0	53.8	5				
451 Haloxyfop	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.
452 Cyanofenphos	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.	n.a.
453 Carpropamid	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ	<LQ

**Chromatogram Data:**  
 Tebuconazole  
 RT=14.55  
 AA=26832  
 [C]=8.9 ppb  
 IR=11%\_Exp(9%\_14%)

**Processing Steps:**

- Step 1 - Retention review 2.0
- Step 2 - Check ISTD 2.0
- Step 3 - Calibration curve\_linearity\_R<sup>2</sup>\_%diff 2.0
- Step 4 - Ion Ratio Adjustment 2.0
- Step 5 - QC\_%diff 2.0
- Step 6 - Quan & Ion Ratio check 2.0

**Logic:**  
 If RT OK  
 If ion ratio OK  
 If > LQ  
 Then give Quan value

# Moving environmental methods to GC-MS/MS



## NeverVent technology

- Increasing instrument uptime




## Off-axis ion guide pre-filter

- Eliminates the neutral noise

## XLXR detector

- Extended dynamic range and lifetime

## Excellent sensitivity



8 x 1 fg on-column OFN injections with %RSD of 4.1%. IDL is 0.12 fg

## Modern GC design

- Unique modular injector and detector design
- Easy to use touch screen
- Real time instrument monitoring and video guides



## Software productivity tools

- Compliant ready software
- SmartStatus and SmartTune

## Sustainable operation


- Partnership with MyGreenLab
- Helium saving capabilities



## Proven robustness

- Over multiple applications

900 continuous injection  
no maintenance



- **The modern laboratory demands the next generation of analytical capability**
- **Modernizing common SVOC workflows can be realized using the TSQ 9610 GC-MS/MS**
  - EPA 8270E, EPA 625, Appendix IX, EPA 8081B, EPA 8141B,
  - Covering SVOCs, organochlorine pesticides, organophosphorus compounds, PCBs
- **Simplifies and accelerates routine operation for an increasingly stringent market**
  - Generate better data to reduce common challenges
  - Improve sample throughput
    - Automated sample preparation, faster calibrations and chromatography, single-run confirmation, unique features of TSQ 9610
  - Streamline operations
    - Simpler training, fewer specialized operators, harmonized consumables and laboratory supplies, unified chromatography software (IC, GC, GCMS, LC, LCMS)
  - Potential to consolidate of methods for multiple assays in a single injection
  - Peace-of-mind for future operations as methods become more stringent

# Thank you

