


Emerging Methods for the Analysis of Volatile PFAS in Water: HS-SPME-GC/MS

Andy Sandy, Yoshiyuki Okamura, Evelyn Wang,
Alan Owens, Ruth Marfil-Vega



The PFAS universe




* Strictly speaking, these substances are not fluorotelomers, as they are not derived from the telomerization process. Despite this, they are termed here "n:1 fluorotelomer-based" substances for readability. Future work may consider to identify more proper terminology for this group of PFASs.
 ** Note that for many compounds such as HFP and TFE, there are different synthesis routes with different starting materials, and here shows only one of them.
 *** Note that there are three synthesis routes shown here for manufacturing of PFCAs, from PACFs, PFAs and n:2 FTIs. Note that different synthesis routes may generate PFCAs with different perfluorocarbon chain lengths.


Sources: (1) Siegmund G, Schwertfeger W, Feiring A, Smart B, Behr F, Vogel H, McKusick B. *Fluorine Compounds, Organic*, 3rd ed., Wiley-VCH Verlag GmbH & Co. KGaA: Weinheim, Germany, 2000; Vol. 33. (2) Banks RE, Smart BE, Tatlow JC. *Organofluorine Chemistry: Principles and Commercial Applications*. New York: Plenum, 1994. (3) Buck RC, Franklin J, Berger U, Conder JM, Cousins IT, De Voogt P, Jensen AA, Kannan K, Mabury SA, van Leeuwen SPJ. *Perfluoroalkyl and polyfluoroalkyl substances in the environment: terminology, classification, and origins*. *Integr Environ Assess Manag* 2011, 7 (4), 513–541. (4) Wang Z, Cousins IT, Scheringer M, Buck RC, Hungerbühler K. *Global emission inventories for C4–C14 perfluoroalkyl carboxylic acid (PFCA) homologues from 1951 to 2030. Part I: production and emissions from quantifiable sources*. *Environ Int* 2014, 70, 62–75. (5) Moffett RH, Howell JL, Hoerter JM, Shtarov AB, Jannerfeldt G, Johnston SB, Keenan J, Warriner C, Closser DM. *Perfluoroalkylpolyethers in Synthetics, Mineral Oils, and Bio-based Lubricants: Chemistry and Technology* (third edition), Edited by Rudnick LR. 2020. CRC Press. ISBN: 978-1-138-06821-6. (6) Grot W. *Fluorinated Ionomers*. William Andrew 2011. ISBN: 978-1-437-74457-6.

Figure 10. An overview of some common synthesis routes of different individual or groups of PFASs based on publicly accessible source

The analytical instrumentation universe for PFAS



Historically, focus on ionic PFAS



Growing interest in neutral and volatile PFAS

- Diverse sample types
- New methods being published by EPA and other standardization organizations
- LCMS can be used... but GCMS is most suitable



Sample introduction for GCMS

Sample Introduction Technique	Sensitivity (GC-MS)	Extraction Mode
SHS	ppb-ppm level	Static equilibrium gas extraction
DHS	ppt-ppb level	Dynamic non-equilibrium gas extraction
SPME	ppt-ppb level	Sorptive extraction
TD	ppt-ppb level	Sorptive extraction
Direct-TD	ppt-ppm level	Direct thermal extraction
Liq	ppb level	-
Py	µg level	Destructive thermal decomposition
DI	ng level	-

SHS	Static Headspace	Direct TD	Direct Thermal Desorption
DHS	Dynamic Headspace (Purge & Trap)	Liq	Liquid Injection
SPME	Solid Phase Microextraction	Py	Pyrolysis
TD	Thermal Desorption	DI	Direct Injection

Sample introduction for GCMS

Environmental Compliance

Sample Introduction Technique	Sensitivity (GC-MS)	Extraction Mode
SHS	ppb-ppm level	Static equilibrium gas extraction
DHS	ppt-ppb level	Dynamic non-equilibrium gas extraction
SPME T&O	ppt-ppb level	Sorptive extraction
TD	ppt-ppb level	Sorptive extraction
Direct-TD	ppt-ppm level	Direct thermal extraction
Liq	ppb level	-
Py	µg level	Destructive thermal decomposition
DI	ng level	-

SHS	Static Headspace	Direct TD	Direct Thermal Desorption
DHS	Dynamic Headspace (Purge & Trap)	Liq	Liquid Injection
SPME	Solid Phase Microextraction	Py	Pyrolysis
TD	Thermal Desorption	DI	Direct Injection

Sample introduction for GCMS


Environmental Compliance

Emerging


Sample Introduction Technique	Sensitivity (GC-MS)	Extraction Mode
SHS	ppb-ppm level	Static equilibrium gas extraction
DHS	ppt-ppb level	Dynamic non-equilibrium gas extraction
SPME T&O	ppt-ppb level	Sorptive extraction
TD	ppt-ppb level	Sorptive extraction
Direct-TD	ppt-ppm level	Direct thermal extraction
Liq	ppb level	-
Py	µg level	Destructive thermal decomposition
DI	ng level	-

SHS	Static Headspace	Direct TD	Direct Thermal Desorption
DHS	Dynamic Headspace (Purge & Trap)	Liq	Liquid Injection
SPME	Solid Phase Microextraction	Py	Pyrolysis
TD	Thermal Desorption	DI	Direct Injection

Comparison of techniques



Comparison of techniques for PFAS analysis



Most commonly used

- Liquid injection
- Thermal Desorption


What about Headspace and SPME?

SHS

Static Headspace


Passive extraction of analytes into gas phase

Pre-heating and Equilibrium



>>>

Sampling and Injection




SPME

Solid Phase Microextraction


Inserting extraction sorbent into a sample vial

Sorption and Extraction



>>>

Injection



Successful analysis of PFAS by HS-GCMS

Henry's Law constants of 15 per- and polyfluoroalkyl substances determined by static headspace analysis




Fig. 1. Schematic of the method of preparing and analyzing headspace samples by GC-MS/MS for the determination of the k_H of PFAS.

Successful analysis of PFAS by HS-GCMS

PFAS	k_H	Slope	Intercept	R^2	pH
4:2 FTOH	0.31 ± 0.07	$2.77 \pm 0.39 \times 10^{-5}$	$8.91 \pm 1.70 \times 10^{-5}$	0.920	6.59 ± 0.29
6:2 FTOH	1.26 ± 0.40	$5.18 \pm 0.29 \times 10^{-6}$	$4.12 \pm 1.20 \times 10^{-6}$	0.986	6.55 ± 0.08
8:2 FTOH*	1.98 ± 0.69	$1.51 \pm 0.06 \times 10^{-5}$	$7.66 \pm 2.60 \times 10^{-6}$	0.993	6.74 ± 0.17
10:2 FTOH*	2.82 ± 1.12	$6.92 \pm 0.22 \times 10^{-4}$	$2.45 \pm 0.96 \times 10^{-4}$	0.995	6.69 ± 0.26
4:2 FTS	0.09 ± 0.02	$4.11 \pm 0.88 \times 10^{-4}$	$4.33 \pm 0.38 \times 10^{-3}$	0.833	5.66 ± 0.02
6:2 FTS	0.16 ± 0.01	$2.37 \pm 0.17 \times 10^{-4}$	$1.48 \pm 0.08 \times 10^{-3}$	0.976	5.70 ± 0.01
8:2 FTS	0.18 ± 0.02	$8.55 \pm 0.65 \times 10^{-5}$	$4.75 \pm 0.28 \times 10^{-4}$	0.975	5.84 ± 0.03
PFHxI	1.01 ± 0.46	$4.70 \pm 0.47 \times 10^{-7}$	$4.66 \pm 2.01 \times 10^{-7}$	0.956	6.90 ± 0.01
6:2 FTUI	0.48 ± 0.10	$4.41 \pm 0.39 \times 10^{-6}$	$9.27 \pm 1.70 \times 10^{-6}$	0.966	6.39 ± 0.01
6:2 FTI	0.30 ± 0.07	$1.32 \pm 0.19 \times 10^{-6}$	$4.39 \pm 0.85 \times 10^{-6}$	0.910	6.71 ± 0.08
N-EtFOSA-M	0.43 ± 0.12	$1.18 \pm 0.15 \times 10^{-3}$	$2.77 \pm 0.68 \times 10^{-3}$	0.927	6.33 ± 0.16
N-MeFOSA-M	0.92 ± 0.48	$2.03 \pm 0.26 \times 10^{-3}$	$2.21 \pm 1.11 \times 10^{-3}$	0.933	6.24 ± 0.07
6:2 FTO	3.86 ± 1.49	$4.40 \pm 0.18 \times 10^{-7}$	$1.14 \pm 0.43 \times 10^{-7}$	0.987	6.23 ± 0.02
8:2 FTCA	0.69 ± 0.07	$8.18 \pm 0.26 \times 10^{-5}$	$1.19 \pm 0.11 \times 10^{-4}$	0.995	5.66 ± 0.28
8:2 FTAC	0.32 ± 0.11	$3.09 \pm 0.06 \times 10^{-3}$	$9.75 \pm 2.61 \times 10^{-4}$	0.858	6.58 ± 0.05






Fig. 3. k_H of 4:2 FTOH and N-EtFOSA-M in varying water matrixes.

- Contribution to closing the mass balance of PFAS in the environment:
 - Simplified analytical method
 - Determination of physico-chemical properties

More about SPME




- Direct immersion
- Head Space



HS-SPME-GCMS

Occurrence of volatile PFAS in **liquid** and solids samples of diverse origin (**environmental**, food, consumer products) with same set-up



Instrumentation and targets



Shimadzu's GCMS QP-2020NX



Shimadzu's GCMS TQ-8040

Chemical Class	Compound	Acronym	CAS Number
Perfluoroalkyl iodides (PFIs)	Perfluorohexyl iodide	PFHxI	355-43-1
	Perfluorooctyl iodide	PFOI	507-63-1
(n:2) Fluorotelomer iodides (FTIs)	4:2 Fluorotelomer iodide	4:2 FTI	2043-55-2
	6:2 Fluorotelomer iodide	6:2 FTI	2043-57-4
	8:2 Fluorotelomer iodide	8:2 FTI	2043-53-0
(n:2) Fluorotelomer acrylates (FTACs)	6:2 Fluorotelomer acrylate	6:2 FTAC	17527-29-6
	8:2 Fluorotelomer acrylate	8:2 FTAC	27905-45-9
	<i>1H,1H,2H,2H-Perfluoro-n-octyl acrylate-d3</i>	<i>6:2 FTAC d3</i>	<i>7527-29-6</i>
(n:2) Fluorotelomer methacrylates (FTMACs)	6:2 Fluorotelomer methacrylate	6:2 FTMAC	2144-53-8
	8:2 Fluorotelomer methacrylate	8:2 FTMAC	1996-88-9
(n:2) Fluorotelomer alcohols (FTOHs)	8:2 Fluorotelomer alcohol	8:2 FTOH	678-39-7
	10:2 Fluorotelomer alcohol	10:2 FTOH	865-86-1
	<i>2-perfluorooctyl-[1,1-2H2-1,2-13C2]-ethanol</i>	<i>8:2 FTOH 13C2</i>	<i>872398-73-7</i>
	<i>2-perfluorodecyl-[1,1-2H2-1,2-13C2]-ethanol</i>	<i>10:2 FTOH 13C2</i>	<i>865-86-1</i>
Perfluoroalkane sulfonamides (FASAs)	N-Methyl perfluorooctane sulfonamide	MeFOSA	31506-32-8
	N-Ethyl perfluorooctane sulfonamide	EtFOSA	4151-50-2
	<i>n-ethyl-d5-perfluoro-1-octanesulfonamide</i>	<i>EtFOSA d5</i>	<i>936109-40-9</i>


Optimized method conditions

Gas Chromatography	Nexis GC-2030
Injection mode	Splitless
Carrier gas	Helium
Injection port temperature (°C)	240
High pressure injection	Auto, 250 kPa, 1 min
Column	SH-I-624Sil MS Capillary, 30 m x 0.25 mmID x 1.40 µm
Flow control mode (cm/sec)	Linear velocity, 44.4
Total flow (mL/min)	50
Oven temperature	40°C (7 min.), 5°C/min. to 188°C (0 min.), 40°C/min. to 300°C, (5 min.)
SPME analysis	AOC-6000 Plus
SPME Fiber	50/30 µm DVB/CAR/PDMS
Incubation time (min)	5
Extraction time (min)	30
Desorption time (min)	7
Agitation speed (rpm)	300
Extraction temperature (°C)	50
Sample volume (mL)	10
Desorption temperature (°C)	240
Sampling salinity	2 % NaCl (w/v)


Mass Spectrometer	Common Parameters
Interface temperature (°C)	280
Ion source temperature (°C)	200
Detector voltage (kV)	Relative to Tune 0.4
Threshold	0
Mass Spectrometer	QP-2020NX
Acquisition mode	Qualitative analysis: Full scan: m/z 50 to 600 Quantitative analysis: SIM, Event time 0.3 sec.
Tuning mode	High Sensitivity
Mass Spectrometer	TQ-8040
Acquisition mode	Acquisition mode: MRM, Loop time: 0.3 sec.
Tuning mode	Normal mode

Experimental plan

- Chromatographic and MS methods, by liquid injection
 - Single Quadrupole
 - Triple Quadrupole
- SPME parameters - optimization
- Calibration curve and linear range
- Carry-over
- PFAS in the background
- Analysis of environmental samples
 - Low system background
 - Precision and accuracy
 - Matrix effects




Results - Chromatography



TIC chromatogram of the 13 targeted PFAS compounds at 5 mg/L

SIM and MRM




Calibration curve

Compound	QP-2020NX			TQ-8040		
	Calibration Range (ng/L)	R2	RF (Response Factor) %RSD	Calibration Range (ng/L)	R ²	RF (Response Factor) %RSD
PFHxI	2.5-2000	0.993	10.89	2.5-2000	0.999	13.68
PFOI	2.5-2000	0.997	10.26	2.5-1000	0.998	18.94
4:2 FTI	2.5-800	0.993	8.28	2.5-2000	0.997	9.30
6:2 FTI	25-800	0.994	13.53	1-2000	0.998	17.18
8:2 FTOH	25-2000	0.997	5.37	2.5-2000	>0.999	6.31
6:2 FTAC	25-2000	0.998	19.87	2.5-2000	0.998	4.03
8:2 FTI	2.5-800	0.996	13.59	2.5-2000	0.999	9.05
10:2 FTOH	2.5-2000	0.999	10.38	2.5-2000	>0.999	6.45
6:2 FTMAC	2.5-800	0.995	12.43	2.5-2000	0.998	10.41
8:2 FTAC	5-250	0.995	14.81	2.5-2000	0.999	11.32
8:2 FTMAC	2.5-250	0.998	19.51	2.5-2000	0.999	9.98
MeFOSA	5-2000	>0.999	17.79	2.5-2000	0.999	6.85
EtFOSA	10-2000	0.999	11.40	1-2000	>0.999	7.17

The linear range of each PFAS target includes at least seven calibration levels
 Results showed a good linear fit for all compounds with $R2 \geq 0.993$ and $RF \%RSD < 20$.

%Accuracy and precision with TQ-8040



What about PFAS in the background?



Consumables and reagents

- SPME vials
- SPME fibers
- Salt
- Solvents
- Commercial standards


No fluorinated components
in sample flow path

Background and carryover


- None of the target PFAS in the laboratory blank samples showed quantifiable results
- The area of peaks of in the blanks were less than 1/5 of the lowest calibration standard
- The carryover effect was evaluated by analyzing a blank immediately after the highest calibration standard: <0.2%






More results - Samples






More results - Samples



Take home messages



Take home messages




Take home messages

- Demonstrated the suitability of HS-SPME-GCMS for the analysis of selected classes of volatile PFAS in environmental liquid samples, without the need of extensive sample preparation
 - Sensitivity:
 - SQ: 2.5 - 25 ppt
 - TQ: 1 – 2.5 ppt
 - Minimal sources of PFAS in background
- More work on-going...



- HPLC** Carbamate
Diquat
Glyphosate
- IC** Anions
Chromium VI
Ammonia
- TOC** Nitrogen
Organic Carbon
Phosphorous
- FTIR & Raman** **Microplastics**
- UV-Vis** Chlorine
UV254
Others
- Particle Size** **Solids**



- LCMS** Cyanotoxins
PFAS
Pesticides
Emerging Contaminants
Unknowns (QTOF)
- GC** HAAs
Herbicides
Pesticides
Phenols
Other Organics
PCBs
- GCMS** 1,4-dioxane
Dioxins & Furans
Nitrosamines
Other Organics
Microplastics
Semivolatiles
Taste & Odor
Volatiles
THMs
- ICPMS** **Metals**
- ICP** **Metals**
- AA** **Metals**

For any questions, contact:
Ruth Marfil-Vega
rmmarfilvega@shimadzu.com

For more information, visit:
www.OneLabOneEarth.com

Connect with us:

X - [@shimadzussi](https://twitter.com/shimadzussi)

Instagram - [@shimadzussi](https://www.instagram.com/shimadzussi)

LinkedIn - [/company/shimadzu-scientific-instruments/](https://www.linkedin.com/company/shimadzu-scientific-instruments/)

YouTube - [@ShimadzuScientificInstruments](https://www.youtube.com/@ShimadzuScientificInstruments)