

Orbitrap Workflows for Non-targeted Analysis Using Dispersive Liquid-liquid Microextraction (DiLLME) Sample Preparation


Ed George

Vertical Marketing Manager

NEMC 2024

August 8, 2024

 The world leader in serving science




History of PFAS and where we are today

PFAS →

“fluorinated substances that contain at least one fully fluorinated methyl or methylene carbon atom (without any H/Cl/Br/I atom attached to it), i.e. with a few noted exceptions, any chemical with at least a perfluorinated methyl group (–CF₃) or a perfluorinated methylene group (–CF₂–)”

Organization for Economic Co-operation and Development (OECD)



Current

WATERSHED | WHYY NEWS CLIMATE DESK

EPA finalizes new rules for companies producing toxic PFAS chemicals

Companies must now notify the EPA if they want to resume the production of certain PFAS chemicals.

What the Hazardous Substance Designation of PFAS Chemicals Means for Local Governments
National League of Cities

CLIMATE
Biden administration sets first-ever limits on 'forever chemicals' in drinking water
AP News

(P)FASten your seatbelts
Harvard Law Today
By Scott Young

PFAS transport through the environment



Aqueous firefighting foam

Airports
Military bases



Household products

Cleaning products
Personal care products
Paints, waxes, polishes



Manufacturing processes

Commercial products
Electronics
Petrochemicals




Fabrics

Stain resistant carpets and fabrics
Waterproof or resistant clothing



Food contact materials

Non-stick cookware
Microwave popcorn bags
Fast food wrappers




>10,000 possible PFAS compounds present in the environment


Workflow strategies



Source of PFAS in the United States




PFAS contamination of drinking water in the United States




Workflows for Targeted and Non-targeted Analysis

Thousands of possible PFAS compounds present in the environment



Workflows for Targeted and Non-targeted Analysis

Thousands of possible PFAS compounds present in the environment




Objective

Demonstrate the comprehensive workflow for PFAS in drinking water samples with a fully automated novel extraction technique followed by targeted quantitation and non-targeted analysis.

Accelerated Solvent Extractor


Thermo Scientific™ Vanquish UHPLC + TSQ Plus Series Triple Quadrupole Mass Spectrometer

Key method conditions for NTA LC-Orbitrap Workflow




Automate DLLME with the Thermo Scientific™ TriPlus RSH SMART

- Samples were obtained from multiple drinking water sources, including tap water and drinking water from the area of Villebon, France.
- Prior to DLLME samples were spiked with 56 PFAS compounds at two different levels: 1 and 75 ng/L



Thermo Scientific™ Vanquish Flex UHPLC


- Delay Column: Thermo Scientific™ Hypersil Gold C18
- Analytical column: Thermo Scientific™ Acclaim™ 120 C18 column
- Flow rate: 0.4 mL/min



Thermo Scientific™ Orbitrap Exploris™ 240

- Heated electrospray ionization (HESI)
- Negative mode
- Combination of Full Scan (MS1) and data-dependent MS2 (ddMS2)
- Mass resolution:
 - Full Scan = 240,000
 - ddMS2 Scan = 30,000
- Stepped HCD fragmentation at 10, 25, 55% (normalized collision energy)

Dispersive Liquid-Liquid Microextraction (DLLME) for PFAS



Water containing dispersion and extraction solvent




Injection from extract layer at the bottom of the vial

Targeted quantification results


Figure of merit for 56 target PFAS compounds

- The **mean accuracy and precision** of the workflow was evaluated for both spike levels (1 and 75 ng/L) in both bottled and tap water from 6 injections over 3 days.
- Majority of target PFAS compounds exhibited **spike recoveries** within 70–130% and **precision** (coefficient of variation, CV) within 20%.

Spike Recoveries




Spike Precision



Targeted quantification results

Figure of merit for 56 target PFAS compounds

- **Limit of quantification (LOQ)** values were obtained based on the concentration level for which both reproducibility (<30%) and accuracy (70-130%) criteria were met on 6 injections over 3 days in neat solution.
- Majority of target PFAS compounds (~80%) exhibited LOQs of ≤ 0.5 ng/L, with the remaining ~20% exhibiting LOQs within the range of 1–5 ng/L.



Non-targeted analysis (NTA) of PFAS: The Wild West

Lots of great ideas and ambitions.....but no rules to follow




**Goal of PFAS
NTA workflows**

Detect unknown PFAS compounds and annotate them as confidently as possible using multiple layers of mass spectrometry evidence

Increasing confidence in identifying PFAS compounds


Harnessing the Power of High Resolution Accurate Mass and MSⁿ Analysis

Increasing Confidence in Compound Identification



Workflow for non-targeted PFAS analysis

Compound Discoverer™ version 3.3 SP3



Initial results – All annotated compounds

Starting out with ~4,000 detected and annotated compounds, but how many PFAS?


Name	Formula	F	RT [min]	m/z	Area (Max)	Annot. Source	mzVault Library Matches	Mass List Matches	mzCloud Best Match	mzVault Best Match	Class Coverage	Mass Defect	MS2	Reference Ion
Embelin	C17 H26 O4	0	11.202	293.17564	6206955746	■■■■■	■■■■■	■■■■■			0.00 0.00	0.183 622 0.202 0.202	■	[M-H]-1
2,2,3,3,4,4,5,5,6,6,7,7-Dodecafluoroheptanoic acid	C7 H2 F12 O2	12	8.571	344.97796	2723233336	■■■■■	■■■■■	■■■■■		98.1	0.37 6.25	-0.015 -43 0.007 0.007	■	[M-H]-1
Sebacic acid	C10 H18 O4	0	8.231	201.11318	2557665245	■■■■■	■■■■■	■■■■■		74.2	0.00 0.00	0.12 596 0.133 0.133	■	[M-H]-1
Sebacic acid	C10 H18 O4	0	8.610	201.11317	2532548468	■■■■■	■■■■■	■■■■■		76.7	0.00 0.00	0.12 596 0.133 0.133	■	[M-H]-1
Perfluoro-1-nonanesulfonate	C9 H F19 O3 S	19	12.750	548.92651	2528470664	■■■■■	■■■■■	■■■■■	93.4		0.99 31.25	-0.066 -121 -0.031 -0.031	■	[M-H]-1
1-(2-NITROPHENOXY)OCTANE	C14 H21 N O3	0	13.848	250.14459	2443776168	■■■■■	■■■■■	■■■■■			0.00 0.00	0.152 605 0.168 0.168	■	[M-H]-1
((3,5-Bis(trifluoromethyl)phenyl)ethyl)amino)methylen	C14 H9 F6 N3	6	14.600	332.06237	2441071481	■■■■■	■■■■■	■■■■■			0.50 6.25	0.07 209 0.091 0.091	■	[M-H]-1
Perfluoro-1-hexanesulfonamide	C6 H2 F13 N O2 S	13	11.828	397.95152	2395022278	■■■■■	■■■■■	■■■■■	99.5	95.4	0.75 31.25	-0.041 -103 -0.016 -0.016	■	[M-H]-1
2,2,3-Trifluoro-3-(1,1,2,2,3,3-hexafluoro-3-trifluoromethyl)heptanesulfonate	C7 H2 F12 O4	12	10.374	376.96779	2389386699	■■■■■	■■■■■	■■■■■	89.7	96.0	0.25 6.25	-0.025 -66 -0.001 -0.001	■	[M-H]-1
Perfluoro-1-heptanesulfonate	C7 H F15 O3 S	15	11.301	448.93240	2349891893	■■■■■	■■■■■	■■■■■	97.8		0.87 31.25	-0.06 -134 -0.032 -0.032	■	[M-H]-1
Perfluoro-1-decanesulfonic acid (PFDS)	C10 H F21 O3 S	21	13.337	598.92313	2345931354	■■■■■	■■■■■	■■■■■	94.7		1.37 25.00	-0.07 -116 -0.031 -0.031	■	[M-H]-1
Perfluorohexanoic acid	C6 H F11 O2	11	9.290	312.97248	2338026334	■■■■■	■■■■■	■■■■■		95.3	1.24 18.75	-0.02 -64 0 0	■	[M-H]-1
Perfluoroheptanoic acid	C7 H F13 O2	13	10.284	362.96852	2204485204	■■■■■	■■■■■	■■■■■		91.9	0.87 18.75	-0.024 -66 -0.001 -0.001	■	[M-H]-1
Dibutyl phthalate	C16 H22 O4	0	12.774	277.14443	2204379622	■■■■■	■■■■■	■■■■■	87.6		53.0 0.00 0.00	0.152 545 0.169 0.169	■	[M-H]-1
Ricinoleic acid	C18 H34 O3	0	14.379	297.24321	2165582651	■■■■■	■■■■■	■■■■■		86.1	0.00 0.00	0.251 840 0.27 0.27	■	[M-H]-1
Perfluoro-1-octanesulfonic acid (PFOS)	C8 H F17 O3 S	17	12.080	498.92939	2092317947	■■■■■	■■■■■	■■■■■	91.1	56.2	1.86 56.25	-0.063 -127 -0.031 -0.031	■	[M-H]-1
9-Chlorohexadecafluoro-3-oxanonane-1-sulfonate	C8 H Cl F16 O4 S	16	12.455	530.89499	2086059949	■■■■■	■■■■■	■■■■■	99.2	94.5	0.50 18.75	-0.098 -184 -0.064 -0.064	■	[M-H]-1
Perfluoro-1-butanedisulfonamide	C4 H2 F9 N O2 S	9	9.716	297.95864	2062831194	■■■■■	■■■■■	■■■■■	99.1		25.6 0.87 31.25	-0.034 -114 -0.015 -0.015	■	[M-H]-1
Diffuoro[1,1,2,2-tetrafluoro-2-(trifluoromethoxy)ethoxy	C5 H F9 O4	9	9.094	200.97920	1954326070	■■■■■	■■■■■	■■■■■		99.8	0.37 6.25	-0.014 -67 -0.001 -0.001	■	[M-H]-1
Perfluorooctanoic acid	C8 H F15 O2	15	11.301	412.96555	1932879256	■■■■■	■■■■■	■■■■■		96.4	0.87 31.25	-0.027 -66 -0.001 -0.001	■	[M-H]-1
11-Chlorohexadecafluoro-3-oxanonane-1-sulfonate	C10 H Cl F20 O4 S	20	13.607	630.88856	1908371270	■■■■■	■■■■■	■■■■■	94.5	88.4	0.62 18.75	-0.104 -165 -0.064 -0.064	■	[M-H]-1
Perfluoropentanoic acid	C5 H F9 O2	9	7.902	262.97576	1894512695	■■■■■	■■■■■	■■■■■		92.7	1.24 31.25	-0.017 -64 0 0	■	[M-H]-1
Perfluoro-3-ethoxypropanoic acid	C5 H F9 O3	9	8.224	278.97074	1782839110	■■■■■	■■■■■	■■■■■		28.6	0.37 12.50	-0.022 -79 -0.004 -0.004	■	[M-H]-1
Perfluorodecanoic acid (PFDA)	C10 H F19 O2	19	12.788	512.95949	1780720217	■■■■■	■■■■■	■■■■■	82.7	87.5	1.49 31.25	-0.033 -65 0 0	■	[M-H]-1
Perfluorononanoic acid (PFNA)	C9 H F17 O2	17	12.107	462.96246	1703398034	■■■■■	■■■■■	■■■■■	75.6	84.6	1.12 31.25	-0.03 -65 -0.001 -0.001	■	[M-H]-1
Cyclohexanesulfonic acid, nonafluorobis(trifluoromethyl)	C6 H F15 O3 S	15	11.154	460.93250	1692089139	■■■■■	■■■■■	■■■■■		52.3	0.50 12.50	-0.06 -130 -0.031 -0.031	■	[M-H]-1
10-Hydroxy-2-decanoic acid	C10 H18 O3	0	9.197	185.11810	1582476257	■■■■■	■■■■■	■■■■■		32.9	0.00 0.00	0.125 674 0.137 0.137	■	[M-H]-1
Perfluorododecanoic acid (PFDoA)	C12 H F23 O2	23	13.903	612.95295	1555517238	■■■■■	■■■■■	■■■■■	83.6	91.6	1.61 31.25	-0.04 -65 -0.001 -0.001	■	[M-H]-1
Laurylsulfuric acid	C12 H26 O4 S	0	13.311	265.14770	1518837203	■■■■■	■■■■■	■■■■■		97.2	0.25 6.25	0.155 582 0.172 0.172	■	[M-H]-1
Perfluoro-1-dodecanesulfonate	C12 H F25 O3 S	25	14.302	698.91628	1443433368	■■■■■	■■■■■	■■■■■	98.4		1.49 31.25	-0.076 -109 -0.032 -0.032	■	[M-H]-1
Perfluorooctylsulfonamide (PFOSA)	C8 H2 F17 N O2 S	17	13.301	497.94546	1336086426	■■■■■	■■■■■	■■■■■	98.7		1.37 62.50	-0.047 -95 -0.015 -0.015	■	[M-H]-1
7-Hydroxy-3-(2-hydroxyethyl)-4-methyl-2H-chromen-2-yl	C12 H12 O4	0	9.030	201.05574	1228450551	■■■■■	■■■■■	■■■■■		66.5	0.12 0.00	0.063 312 0.076 0.076	■	[M-H]-1
Laurylsulfuric acid	C12 H26 O4 S	0	12.835	265.14771	1148524771	■■■■■	■■■■■	■■■■■		95.8	0.75 18.75	0.155 582 0.172 0.172	■	[M-H]-1
Saccharin	C7 H5 N O3 S	0	4.672	181.99154	1132965009	■■■■■	■■■■■	■■■■■	99.9		0.12 0.00	-0.001 -6 0.011 0.011	■	[M-H]-1
Saccharin	C7 H5 N O3 S	0	4.541	181.99154	1101742601	■■■■■	■■■■■	■■■■■	99.9		0.12 0.00	-0.001 -6 0.011 0.011	■	[M-H]-1
Ricknic acid A	C15 H22 O3	0	13.702	249.14934	1093438074	■■■■■	■■■■■	■■■■■		94.5	0.00 0.00	0.157 626 0.173 0.173	■	[M-H]-1
2-[(4-Aminobenzoyl)amino]benzoic acid	C14 H12 N2 O3	0	9.324	211.08749	1089696762	■■■■■	■■■■■	■■■■■		67.8	0.00 0.00	0.095 447 0.108 0.108	■	[M-H]-1
2-(2-Trifluoromethyl-benzimidazol-1-yl)-ethanol	C10 H9 F3 N2 O	3	11.539	229.05917	1067219280	■■■■■	■■■■■	■■■■■		40.2	0.00 0.00	0.074 310 0.09 0.09	■	[M-H]-1
6-Isopropyl-3-methyl-2,4-dinitrophenol	C10 H12 N2 O5	0	12.115	239.06715	1083937330	■■■■■	■■■■■	■■■■■		84.7	1.12 25.00	-0.036 -64 0 0	■	[M-H]-1
Perfluoroundecanoic acid (PFUDA)	C11 H F21 O2	21	13.382	562.95644	1045300000	■■■■■	■■■■■	■■■■■	5		0.75 31.25	-0.057 -143 -0.032 -0.032	■	[M-H]-1
Perfluoro-1-hexanesulfonic acid (PFHxS)	C6 H F13 O3 S	13	10.342	398.93566	1034990000	■■■■■	■■■■■	■■■■■	2		0.75 31.25	-0.057 -143 -0.032 -0.032	■	[M-H]-1
Perfluorotridecanoic acid (PFTrDA)	C13 H F25 O2	25	14.358	662.94974	1014910000	■■■■■	■■■■■	■■■■■	3		1.37 25.00	-0.043 -65 -0.001 -0.001	■	[M-H]-1
Perfluoro-3-methoxypropanoic acid	C4 H F7 O3	7	7.178	228.97389	999322296	■■■■■	■■■■■	■■■■■		25.3	0.50 12.50	-0.019 -82 -0.004 -0.004	■	[M-H]-1
3-[β-(2-Hydroxypropan-2-yl)-3-methoxycarbonyl-3,4a,4b,8a-tetrahydro-2H-chromen-2-yl]-2-(2-Trifluoromethyl-benzimidazol-1-yl)-ethanol	C26 H38 O8	0	16.084	381.23733	986917737	■■■■■	■■■■■	■■■■■		21.8	0.00 0.00	0.245 875 0.268 0.268	■	[M-H]-1
2-(2-Trifluoromethyl-benzimidazol-1-yl)-ethanol	C10 H9 F3 N2 O	3	11.606	229.05917	964975463	■■■■■	■■■■■	■■■■■		40.2	0.12 6.25	0.086 289 0.081 0.081	■	[M-H]-1
4-Bromo-4'-dodecyloxybiphenyl	C24 H33 Br O	0	12.537	415.16372	852243417	■■■■■	■■■■■	■■■■■		0.00	0.00 0.00	0.171 411 0.188 0.188	■	[M-H]-1
3,5-Di-tert-butyl-2-hydroxybenzoic acid	C15 H22 O3	0	12.045	249.14935	824548158	■■■■■	■■■■■	■■■■■		97.4	0.00 0.00	0.167 638 0.179 0.179	■	[M-H]-1

We need to narrow this list down to only PFAS, and then determine the confidence levels of each annotation

Thousands more

Data reduction approach

Moving towards a finalized list of compounds annotated with Level 1 confidence




Annotation Criteria	Annotation Confidence Level				
	5	4	3	2	1
Measured mass ± 2 ppm of at least one PFAS Mass List	✓	✓	✓	✓	✓
Standard mass defect is between -0.11–0.12		✓	✓	✓	✓
Isotopic pattern match		✓	✓	✓	✓
≥ 1 MS ² fragment with match to FluoroMatch database and/or >50% similarity match to in-silico PFAS spectral libraries			✓		
>50% similarity match to mzCloud™ or 2023 NIST MS/MS spectral libraries				✓	✓
Retention time match to Reference Standard					✓

Annotation of PFAS in DLLME extracts of spiked water samples

Evaluating the detection and annotation efficiency of 56 spiked PFAS compounds


Annotation Criteria	Annotation Confidence Level			
	5	4	3	2
Measured mass ± 2 ppm of at least one PFAS Mass List	✓	✓	✓	✓
Standard mass defect is between -0.11–0.12		✓	✓	✓
Isotopic pattern match		✓	✓	✓
≥ 1 MS ² fragment with match to FluoroMatch database and/or >50% similarity match to in-silico PFAS spectral libraries			✓	
>50% similarity match to mzCloud™ or 2023 NIST MS/MS spectral libraries				✓
Retention time match to Reference Standard				




Overlaid Extracted Ion Chromatograms of annotated PFAS

Efficient detection and annotation from unspiked to 75 ppt


Spiked at
75 ppt



Spiked at
1 ppt




Unspiked



Confidently annotating unknown PFAS compounds



Bringing together multiple resources to annotate at Level 2 confidence



Confidently annotating unknown PFAS compounds


Bringing together multiple resources to annotate at Level 3 confidence

Tags	Name	Formula	F	RT [min]	m/z	Area (Max)	Annot. Source	mzVault Library Matches	Mass List Matches	mzCloud Best Match	mzVault Best Match	Class Cov
43	N-Methylperfluoro-1-octanesulfonamidoacetic acid (N	C11 H6 F17 N O4 S	17	13.437	569.96704	214202729				93.7	79.1	1.37 31.25





Diving deeper in the interpretation of MS2 spectra

Utilizing the SledgeHammer tool in Mass Frontier to better understand MS² fragmentation



Leverages >54,000 fragmentation schemes that cover mechanisms published in major journals



The screenshot shows the 'SledgeHammer 3' interface. On the left, there is a table of predicted fragments with their m/z values and counts. On the right, a reaction mechanism is shown for the fragment at m/z 469.99581, illustrating the fragmentation of the polyglycerol molecule.

m/z	Count
435.973932	14
436.981757	
440.956685	
441.964510	
442.972335	
449.989582	
450.997407	59
453.964510	4
454.972335	4
455.960174	1
455.980160	15
457.975824	1
466.972335	3
467.980160	2
468.987985	6
469.995811	25
472.962914	1
473.970739	3
480.987985	33
481.978564	9
499.986389	32
500.994214	61
504.949155	1
506.964805	1
517.976967	6
518.984792	8
519.992617	38
521.951895	1
522.959720	7
524.955384	1

Reaction mechanism diagram showing the fragmentation of polyglycerol (M 889.97594) into a fragment at m/z 888.96867, which then undergoes dissociation (diss) to form a fragment at m/z 888.96867, and finally a fragment at m/z 469.99581.

Generates list of predicted fragments from polyglycerol


For each predicted fragment a full reaction mechanism is provided

Evaluating overall PFAS composition across samples


Built-in plotting tools to visually investigate difference in PFAS composition

Differential analysis plot


Spiked (1 ppt) vs. Unspiked



Spiked (75 ppt) vs. Unspiked




Orthogonal MS1 plot*




Why PFAS by GCMS ?

- Added value to LCMS
- Complementary results, alternative for LC difficult compounds
- Cross confirmation of results
- Target quantitation, screening, and non-targeted (unexpected) analyses all available
- High throughput
- System always ready- EI, PCI modes and easy change between them
- Obtain Solvent Quality information in mobile phase or sample prep



Tools added for GC-Orbitrap NTA analysis- CD version 3.4

- Dedicated GC-Orbitrap PFAS workflow template
- Updated PFAS compound database
- Adding more chemicals to GC Orbitrap Contaminants library
- CSC and collaborator labs completing full methodology



Filename	Description	File Size	Uploaded	Updated	Control
Arta Lab 6549 Flavonoid Structure Database	Contributed by Arta Laboratory, National Institute of Genetics, Japan.	4,405 KB	6/11/2024 8:52 AM	5/25/2018 1:32 AM	
Chemical List PRASSTRUCT-2022-04-20	Total 10737 compounds from PFAS Structure lists. https://comptox.epa.gov/dashboard/chemical-lists/PRASSTRUCT	11,480 KB	6/11/2024 8:52 AM	1/16/2023 11:36 AM	
EPF HRAM Compound Database		288 KB	6/11/2024 8:52 AM	9/26/2017 3:00 PM	
Example Mass List	Example mass list for the metabolomics tutorial	156 KB	6/11/2024 8:52 AM	8/26/2021 8:11 AM	
Extractables and Leachables HRAM Compound Database	Total of 2862 compounds pertaining to EML analyses, expanded on a prior version with contributions from Gyroq/Vet, Intertek	1,704 KB	6/11/2024 8:52 AM	2/7/2023 8:18 AM	
FCCDB_2022	Food Contact Chemicals database. Kansas Groh, Brigit Geweke, & Jane Muncie (2022). FCCdb: Food Contact Chemicals database. [Data set]. Zenodo. https://doi.org/10.5281/zenodo.3240108 . Groh, K	1,896 KB	6/11/2024 8:52 AM	11/3/2023 11:17 AM	
GC Orbitrap Contaminants Library Compound Database	Total 888 compounds including pesticides, environmental and EML compounds	616 KB	6/11/2024 8:52 AM	11/3/2023 10:32 AM	
GC Orbitrap Flavor and Fragrances Compound Database	Total 411 flavor and fragrances in 11 categories	336 KB	6/11/2024 8:52 AM	1/30/2024 8:23 AM	
GC Orbitrap Metabolomics Library Compound Database	Total 1015 NIST (msl) and TRMAS deconvoluted metabolites with some synonyms	912 KB	6/11/2024 8:52 AM	11/3/2023 10:31 AM	
GC Orbitrap Off Color Compound Database	Total 49 terpenes and off color compounds	172 KB	6/11/2024 8:52 AM	1/29/2024 8:30 AM	
GC Orbitrap PFAS Compound Database	Contributed by Jeremy Kodiveli, from Innovative Clinics. The PFAS in the list are all acquired on GC Orbitrap. Their mass spectra are hosted in this website https://innovativeclinics.com/software/	216 KB	6/14/2024 7:53 PM	5/28/2024 2:16 PM	
HMDB All metabolites (v5.217719 cpdb)	Wilhart DS, Guo AC, Clew E, et al. HMDB 5.0: The Human Metabolome Database for 2022. Nucleic Acids Res. 2022. Jan 7;50(1):D402-31. https://pubmed.ncbi.nlm.nih.gov/34886597/	339,638 KB	6/11/2024 8:51 AM	11/3/2023 8:49 AM	
HMDB C5F metabolites (v5.402 cpdb)	Wilhart DS, Guo AC, Clew E, et al. HMDB 5.0: The Human Metabolome Database for 2022. Nucleic Acids Res. 2022. Jan 7;50(1):D402-31. https://pubmed.ncbi.nlm.nih.gov/34886597/	548 KB	6/11/2024 8:51 AM	12/8/2023 1:24 PM	
HMDB Feces metabolites (v5.6790 cpdb)	Wilhart DS, Guo AC, Clew E, et al. HMDB 5.0: The Human Metabolome Database for 2022. Nucleic Acids Res. 2022. Jan 7;50(1):D402-31. https://pubmed.ncbi.nlm.nih.gov/34886597/	20,404 KB	6/11/2024 8:51 AM	12/8/2023 1:47 PM	
HMDB Saliva metabolites (v5.1197 cpdb)	Wilhart DS, Guo AC, Clew E, et al. HMDB 5.0: The Human Metabolome Database for 2022. Nucleic Acids Res. 2022. Jan 7;50(1):D402-31. https://pubmed.ncbi.nlm.nih.gov/34886597/	1,772 KB	6/11/2024 8:51 AM	12/8/2023 1:55 PM	
HMDB Serum metabolites (v5.4000 cpdb)	Wilhart DS, Guo AC, Clew E, et al. HMDB 5.0: The Human Metabolome Database for 2022. Nucleic Acids Res. 2022. Jan 7;50(1):D402-31. https://pubmed.ncbi.nlm.nih.gov/34886597/	58,716 KB	6/11/2024 8:51 AM	11/3/2023 9:49 AM	
HMDB Sweet metabolites (v5.89 cpdb)	Wilhart DS, Guo AC, Clew E, et al. HMDB 5.0: The Human Metabolome Database for 2022. Nucleic Acids Res. 2022. Jan 7;50(1):D402-31. https://pubmed.ncbi.nlm.nih.gov/34886597/	216 KB	6/11/2024 8:51 AM	12/8/2023 1:57 PM	
HMDB Urine metabolites (v5.4641 cpdb)	Wilhart DS, Guo AC, Clew E, et al. HMDB 5.0: The Human Metabolome Database for 2022. Nucleic Acids Res. 2022. Jan 7;50(1):D402-31. https://pubmed.ncbi.nlm.nih.gov/34886597/	6,200 KB	6/11/2024 8:51 AM	12/8/2023 1:59 PM	
LipidMaps Structure Database 2023-01-11	Lipid MAPS Structure Database. https://www.ebi.ac.uk/chem/summary/17098933	79,076 KB	6/11/2024 8:51 AM	3/6/2023 6:44 PM	
Natural Products Atlas 2023_06	van Santen et al. The Natural Products Atlas 2.0: a database of microbially-derived natural products. Nucleic Acids Research. 2022. 50, D1, 11. D1317-D1323. 10.1093/nar/gkac041	58,496 KB	7/8/2024 9:26 PM	12/8/2023 3:21 PM	
PFAS_NEG	Novel PFAS adapted from Barzen-Hanson et al. Discovery of 40 Classes of Per- and Polyfluoroalkyl Substances in Historical Aqueous Film-Forming Foams (AFFF) and AFFF-impacted Groundwater. En	164 KB	6/11/2024 8:51 AM	1/16/2023 11:36 AM	
PFAS_NIST	Suspect List of Possible Per- and Polyfluoroalkyl Substances (PFAS). https://data.nist.gov/od/id/mdb-2387	7,028 KB	6/11/2024 8:51 AM	1/23/2023 8:00 AM	
PFAS_suspectDB_Duke	PFAS compound database. Getzinger, G. J., Higgins, C. P., Ferguson, R. L. Structure Database and In Silico Spectral Library for Comprehensive Suspect Screening of Per- and Polyfluoroalkyl Substances (P	51,800 KB	6/11/2024 8:51 AM	11/1/2023 10:17 PM	

PFAS

Core GC PCI - Before Reprocessing.cdProcessingWF

- PFAS GC PCI w Molecular Networks and unknown ID using Local and Online Databases.cdProcessingWF
- PFAS GC PCI w unknown ID using Local and Online Databases.cdProcessingWF
- PFAS GC PCI w unknown ID using Local Databases.cdProcessingWF

Some recent publications PFAS using GC-HRMS

Journal of Chromatography A 1693 (2023) 463884

Contents lists available at ScienceDirect

 **Journal of Chromatography A**

journal homepage: www.elsevier.com/locate/chroma

The use of gas chromatography – high resolution mass spectrometry for suspect screening and non-targeted analysis of per- and polyfluoroalkyl substances

Jonathan S Casey^a, Stephen R Jackson^b, Jeff Ryan^b, Seth R Newton^{b,*}

^aORISE, Office of Research & Development, U.S. Environmental Protection Agency, Research Triangle Park, NC 27711, United States
^bOffice of Research and Development, U.S. Environmental Protection Agency, Research Triangle Park, NC 27711, United States



ENVIRONMENTAL
Science & Technology

pubs.acs.org/est Article

Characterization of Per- and Polyfluorinated Alkyl Substances Present in Commercial Anti-fog Products and Their *In Vitro* Adipogenic Activity

Nicholas J. Herkert, Christopher D. Kassotis, Sharon Zhang, Yuling Han, Vivek Francis Pulikkal, Mei Sun, P. Lee Ferguson, and Heather M. Stapleton*

 Cite This: *Environ. Sci. Technol.* 2022, 56, 1162–1173  Read Online



Conclusions

- Automated dispersive liquid-liquid microextraction (DLLME) is a promising technique to extract and concentrate PFAS from drinking water samples. The resulting extract from a single sample preparation can then be measured by both quantitative targeted and non-targeted analysis in separate runs, as well as further exploration by GC-Orbitrap analysis.
- Quantitative targeted analysis showed good accuracies at low (1 ng/L) and high (75 ng/L) spiking level, as well as reproducibility (<30%, n=7) over several days. The use of internal standards for selected PFAS classes can additionally improve obtained results in terms of accuracy.
- The non-targeted analysis workflow in Compound Discoverer 3.3 SP3 provides a comprehensive package enabling confident annotation of unknown PFAS compounds through the use of multiple spectral libraries, fragmentation databases and as well as understanding differences in PFAS composition across samples.
- CD version 3.4 will greatly expand the workflow capabilities for NTA using GC-Orbitrap

Acknowledgements

- **Collaborators**
 - Dr. Lee Ferguson, Duke University
- **Thermo Fisher Scientific**
 - Cynthia Grim
 - Valérie Thibert
 - Aristide Ganci
 - Bénédicte Gauriat
 - Jean-François Garnier