

Comparison of the Hydrolnert and EI Source using US EPA Method 8260 with the Tekmar Atomx XYZ P&T and the Agilent 7890B GC and 5977B MS System with Hydrogen and Helium Carrier Gas

Introduction

The United States Environmental Protection Agency (US EPA) last revised its 8260D method, revision 4, in February 2017, which differs from its 8260A-C methods with an updated Bromofluorobenzene (BFB) tune check acceptance criteria to help with interactions between the carrier gas and water vapor. This application is a demonstration of hydrogen carrier gas for the quantitation of volatile organic compounds (VOCs) with the method requirements of US EPA Method 8260D in conjunction with US EPA Method 5030 for the aqueous samples and US EPA Method 5035 for the solid waste samples.

The Agilent Hydrolnert source is a GC/MS ion source that improves chromatographic performance when using hydrogen as carrier gas. Its inertness minimizes hydrogenation and dechlorination reactions in the MS source, avoiding loss of sensitivity and spectral anomalies while offering a high-boiler peak shape. The Teledyne LABS Tekmar Atomx XYZ combined Purge and Trap (P&T) autosampler is compatible with hydrogen carrier gas and allows the use of nitrogen as purge gas, avoiding the need for a helium supply.

The Tekmar Atomx XYZ was paired with an Agilent 7890B GC and 5977B MS and the Hydrolnert ion source (Figure 1) to demonstrate the method requirements of 8260D: . This work also provides a comparison of the stated method requirements on the Hydrolnert source and the Electron Ionization (EI) source with hydrogen and helium carrier gas.

Experimental

Working 5, 50, and 250 parts per million (ppm) or milligram per liter (mg/L) calibration standards were prepared in methanol from the following commercially available standards: 8260B MegaMix, 8260B Acetate, California Oxygenates, VOA (Ketones), 502.2 Calibration Mix, 2-Chloroethyl Vinyl Ether, and Hexachloroethane. In total, the standards contained 97 compounds.

Nine-point average response factor (%RSD) calibration curves for both water and soil were prepared from 0.5 ppb to 200 parts per billion (ppb) for all compounds. The %RSD was calculated for each compound using four internal standards. Internal and surrogate standards were prepared in methanol from commercially available standards at a concentration of 25 ppm, after which 5 microliters (µL) was then mixed with each 5 milliliter (mL) sample for a resulting concentration of 25 ppb.

Seven 0.5 ppb water standards and seven 1 ppb soil standards were prepared to calculate the MDL calculations. Seven 20 ppb water and soil standards were also prepared for the accuracy and precision calculations of the mid-point calibration check.

All calibration, MDL, and mid-point calibration check standards were analyzed with the Tekmar Atomx XYZ conditions in Table I (water) and Table II (soil). GC-MS conditions are shown in Table III.

Experimental

Table I: Tekmar Atomx XYZ Water Method Conditions				
Standby	Variable	Desorb	Variable	
Valve Oven Temp	140°C	Methanol Needle Rinse	Off	
Transfer Line Temp	140°C	Water Needle Rinse Vol	7.00 mL	
Sample Mount Temp	90°C	Sweep Needle Time	0.25 min	
Water Heater Temp	90°C	Desorb Preheat Temp	245°C	
Sample Cup Temp	20°C	Desorb Temp	250°C	
Soil Valve Temp	50°C	Desorb Time	2.00 min	
Standby Flow	10 mL/min	Drain Flow	300 mL/min	
Purge Ready Temp	40°C	GC Start Signal	Begin Desorb	
Purge	Variable	Bake	Variable	
Sample Equib Time	0.00 min	Methanol Glass Rinse	Off	
Pre-sweep Time	0.25 min	Water Bake Rinses	1	
Prime Sample Fill Vol	3.00 mL	Water Bake Rinse Vol	7.00 mL	
Sample Volume	5.00 mL	Bake Rinse Sweep Time	0.25 min	
Sweep Sample Time	0.25 min	Bake Rinse Sweep Flow	100 mL/min	
Sweep Sample Flow	100 mL/min	Bake Rinse Drain Time	0.40 min	
Sparge Vessel Heater	Off	Bake Time	2.00 min	
Purge Time	11.00 min	Trap Bake Temp	270°C	
Purge Flow	40 mL/min	MCS Bake Temp	180°C	
Purge Temp	20°C	Bake Flow	200 mL/min	
MCS Purge Temp	20°C			
Dry Purge Temp	20°C	Trap	9	
Dry Purge Time	0.50 min	Chiller Tray	Off	
Dry Purge Flow	100 mL/min	Purge Gas	Nitrogen	

Table II: Tekmar Atomx XYZ Soil Method Conditions				
Standby	Variable	Purge	Variable	
Valve Oven Temp	140°C	Purge Temp	20°C	
Transfer Line Temp	140°C	MCS Purge Temp	20°C	
Sample Mount Temp	90°C	Dry Purge Time	2.00 min	
Water Heater Temp	90°C	Dry Purge Flow	100 mL/min	
Sample Cup Temp	40°C	Dry Purge Temp	20°C	
Soil Valve Temp	100°C	Desorb	Variable	
Standby Flow	10 mL/min	Methanol Needle Rinse	Off	
Purge Ready Temp	40°C	Water Needle Rinse Vol	7.00 mL	
Purge	Variable	Sweep Needle Time	0.25 min	
Pre-purge Time	0.00 min	Desorb Preheat Temp	245°C	
Pre-purge Flow	0 mL/min	Desorb Time	2.00 min	
Pre-heat Mix Speed	Slow	Desorb Temp	250°C	
Sample Pre-heat Time	0.00 min	Drain Flow	300 mL/min	
Pre-sweep Time	0.25 min	GC Start Signal	Begin Desorb	
Water Volume	10.00 mL	Bake	Variable	
Sweep Water Time	0.25 min	Bake Time	2.00 min	
Sweep Water Flow	100 mL/min	Bake Flow	200 mL/min	
Sparge Vessel Heater	Off	Bake Temp	270°C	
Purge Mix Speed	Medium	MCS Bake Temp	180°C	
Purge Time	11.00 min	Trap	9	
Purge Flow	40 mL/min	Purge Gas	Nitrogen	

Table III: Agilent 7890B GC and 5977A MSD System Conditions			
Agilent 7890B GC Conditions			
	EI Source	Hydrolnert Source	
Column	DB-624 Ultra Inert	DB-624 Ultra Inert	
	20m x 0.18 mm, 1 µm Film	20m x 0.18 mm, 1 µm Film	
	Carrier gas (He or H2) 1.0 mL/min	Hydrogen carrier gas 1.0 mL/min	
Oven Profile	35 °C hold 2 min, 15°C/min to 100°C, 30°C/min to 230°C, hold 1 min, Run Time 11.67 min		
Inlet	220°C, 80:1 Split, 19.752 psi		200°C, 75:1 Split, 8.5032 psi
Agilent 5977A MSD Conditions			
	EI Source	Hydrolnert Source	
Temp	Transfer Line 250°C; Source 250°C; Quad 150°C	Transfer Line 250°C; Source 230°C; Quad 150°C	
Scan	Range 35 m/z to 270 m/z, Solvent Delay 0.75 min, Normal Scanning	Range 35 m/z to 270 m/z, Solvent Delay 0.50 min, Normal Scanning	
Gain	Gain Factor 5.00, BFB Auto tune	Gain Factor 10.00, Auto tune	

Results and Discussion


Table IV: US EPA Method 8260D Calibration Data in Water and Soil, 0.5-100 ppb													
Compound	Hydrolnert Source		EI Source, Helium		EI Source, Hydrogen		Compound	Hydrolnert Source		EI Source, Helium		EI Source, Hydrogen	
	Water	Soil	Water	Soil	Water	Soil		Water	Soil	Water	Soil	Water	Soil
	RRF (≤% RSD)	RRF (≤% RSD)	RRF (≤% RSD)	RRF (≤% RSD)	RRF (≤% RSD)	RRF (≤% RSD)		RRF (≤% RSD)	RRF (≤% RSD)	RRF (≤% RSD)	RRF (≤% RSD)	RRF (≤% RSD)	RRF (≤% RSD)
Pentafluorobenzene (IS 1)							Propyl Acetate	4.9	7.1	6.5	10.4	10.5	5.8
Dibromofluoromethane (SS 1)	8.5	7.1	2.2	1.8	10.3	3.3	2-Nitropropane	5.9	19.2	11.5	9.7	15.0	0.999 ²
Dichlorodifluoromethane	7.5	12.3	16.7	16.1	8.3	9.2	2-Chloroethyl Vinyl Ether	11.2	12.1	5.6	16.7	8.1	8.5
Chloromethane	8.9	8.3	8.0	10.5	9.2	5.1	cis-1,3-Dichloropropene	7.1	9.6	2.9	9.3	16.9	6.2
Vinyl Chloride	6.2	8.9	10.1	9.4	5.4	8.1	4-Methyl-2-Pentanone	13.3	12.3	4.6	7.7	11.3	14.2
Bromomethane	14.0	10.4	0.998 ²	0.995 ²	0.999 ²	13.4	Toluene	7.9	10.5	4.0	5.3	8.1	11.8
Chloroethane	11.5	9.2	9.4	7.9	7.0	9.5	trans-1,3-Dichloropropene	6.8	6.2	3.7	10.0	15.6	6.8
Trichlorofluoromethane	5.1	8.9	8.5	7.3	4.5	7.6	Ethyl Methacrylate	8.1	9.4	4.1	13.5	12.4	13.7
Diethyl Ether	8.9	3.7	6.3	5.9	3.6	10.3	1,1,2-Trichloroethane	8.3	5.5	4.5	6.5	9.2	6.4
1,1-Dichloroethene	6.8	8.8	7.1	8.1	5.3	11.2	Tetrachloroethylene	9.0	9.5	18.0	12.3	7.1	9.8
1,1,2-Trichlorotrifluoroethane	6.4	10.7	7.2	10.2	4.5	6.7	1,3-Dichloropropane	5.0	9.4	2.7	6.6	5.3	8.0
Acetone	0.996 ¹	0.999 ²	0.999 ¹	0.998 ¹	0.999 ¹	0.998 ¹	2-Hexanone	16.4	7.7	5.9	16.5	10.3	12.6
Iodomethane	17.2	11.7	0.999	20.0	0.999	0.999	Dibromochloromethane	8.5	8.3	4.8	7.6	14.2	6.3
Carbon Disulfide	11.2	0.999 ¹	15.6	12.5	12.3	19.9	Butyl Acetate	8.0	3.7	2.9	8.2	6.5	5.4
Allyl Chloride	11.6	9.8	7.8	10.9	12.0	14.8	1,2-Dibromoethane	6.4	12.2	4.2	16.1	13.7	15.5
Methyl Acetate	5.6	10.1	9.8	14.0	7.3	9.3	Chlorobenzene-d5 (IS 3)						
Methylene Chloride	16.1	0.999 ¹	15.7	0.999 ¹	7.5	0.999 ¹	Toluene-d8 (SS 3)	1.7	1.4	1.5	1.6	1.4	2.0
tert-Butyl alcohol	9.7	7.8	10.8	15.0	5.0	13.4	Chlorobenzene	5.4	7.0	6.2	6.1	7.5	5.3
trans-1,2-Dichloroethene	9.2	11.1	9.5	9.0	4.3	7.0	1,1,1,2-Tetrachloroethane	5.8	9.7	8.0	7.6	0.999 ¹	9.3
Methyl tert-butyl ether	8.4	10.9	7.0	6.9	5.1	12.1	Ethylbenzene	7.6	11.3	5.3	6.3	11.4	6.9
Acetonitrile	10.5	10.9	8.6	10.2	5.1	7.7	m, p-Xylene	9.4	14.4	7.7	6.9	15.4	11.7
1,1-Dichloroethane	8.5	9.5	7.3	7.1	4.4	10.7	o-Xylene	10.2	12.8	8.0	6.3	14.4	9.7
Acrylonitrile	6.1	6.4	3.7	8.9	5.2	12.5	Styrene	6.1	13.3	8.7	9.0	9.4	13.7
Chloroprene	6.1	9.5	3.7	8.9	5.3	5.3	Bromoforn	6.3	8.6	5.3	8.3	0.999 ¹	5.8
Vinyl Acetate	5.5	7.5	12.6	7.5	11.4	6.1	Amyl Acetate	8.2	10.4	11.2	13.7	6.7	14.8
Diisopropyl ether	5.5	10.7	8.5	7.4	4.4	9.3	Isopropylbenzene	7.0	10.8	7.3	7.3	11.9	7.3
Ethyl tert-butyl ether	8.7	12.3	7.1	10.2	5.9	12.8	cis-1,4-Dichloro-2-Butene	7.1	10.3	10.9	14.2	14.1	9.8
2,2-Dichloropropane	8.2	8.5	5.5	7.7	4.6	8.2	Bromobenzene	7.8	11.6	8.4	7.4	5.5	7.8
cis-1,2-Dichloroethene	5.9	13.0	5.1	5.7	3.5	9.5	1,1,2,2-Tetrachloroethane	11.7	11.7	8.5	12.4	17.6	4.2
2-Butanone	11.1	8.0	10.1	7.5	9.3	10.7	1,2,3-Trichloropropane	6.3	9.5	7.2	10.1	7.6	4.1
Propionitrile	10.2	13.8	11.9	5.1	15.4	14.2	Trans-1,4-dichloro-2-butene	7.7	10.4	8.1	12.9	4.2	9.2
Ethyl Acetate	7.8	8.9	7.6	6.0	10.2	9.3	n-Propylbenzene	10.3	14.1	9.6	7.0	14.1	9.2
Methyl Acrylate	8.7	6.9	8.3	8.8	12.3	11.8	2-Chlorotoluene	7.8	11.0	9.0	6.6	12.2	8.9
Bromochloromethane	6.3	12.0	6.7	4.6	7.0	9.7	1,3,5-Trimethylbenzene	6.5	11.5	7.9	9.2	14.8	11.3
Methacrylonitrile	5.7	9.0	10.2	7.2	6.9	9.4	4-Chlorotoluene	7.1	10.2	12.4	7.9	14.8	10.7
Chloroform	4.2	7.3	7.4	5.2	3.8	7.2	tert-Butylbenzene	8.1	8.8	7.9	10.3	13.0	10.4
Tetrahydrofuran	10.5	11.7	10.5	9.6	11.3	16.9	1,2,4-Trimethylbenzene	7.3	9.5	7.6	8.6	14.3	11.8
1,1,1-Trichloroethane	5.5	5.7	4.4	4.9	5.7	5.8	sec-Butylbenzene	8.9	10.9	9.2	8.1	14.2	10.3
Trichloroethylene	6.5	12.7	6.9	7.0	2.7	7.2	1,3-Dichlorobenzene	8.8	14.3	16.5	13.2	11.3	11.9
Isobutyl alcohol	6.4	12.2	7.4	12.0	9.0	7.7	p-Isopropyltoluene	8.8	8.6	8.6	7.0	16.3	11.3
Isopropyl Acetate	6.4	9.4	7.4	11.8	10.7	6.8	1,4-Dichlorobenzene-d4 (IS 4)						
tert-Amyl methyl ether	8.7	8.6	7.1	12.0	7.9	13.1	Bromofluorobenzene (SS 4)	2.5	2.8	3.8	4.6	5.3	4.8
1,2-Dichloroethane	6.2	10.3	5.8	4.1	8.6	9.1	1,4-Dichlorobenzene	6.1	13.5	14.8	15.2	8.7	8.0
Benzene	4.3	3.8	2.6	3.5	7.9	5.8	1,2-Dichlorobenzene	4.6	8.8	8.7	12.3	6.5	6.1
1,1-Dichloropropene	8.5	10.2	2.7	6.1	4.9	8.8	n-Butylbenzene	6.9	8.3	8.3	11.4	11.3	6.1
Carbon Tetrachloride	4.1	6.2	4.9	7.5	9.4	4.6	Hexachloroethane	8.4	13.2	6.1	11.9	9.2	9.2
1,4-Difluorobenzene (IS 2)							1,2-Dibromo-3-Chloropropane	6.9	7.3	4.0	10.9	12.2	10.0
1,2-Dichloroethane-d4 (SS 2)	3.7	5.2	1.7	1.6	3.7	9.6	Nitrobenzene	8.8	0.999 ²	16.7	0.999 ²	14.6	0.999 ²
Dibromomethane	6.5	7.4	7.7	9.0	8.6	13.7	1,2,4-Trichlorobenzene	4.7	12.8	14.1	16.7	4.5	5.9
1,2-Dichloropropane	9.8	12.4	4.8	5.4	3.9	9.6	Hexachlorobutadiene	7.4	8.4	14.9	11.8	0.999 ¹	7.2
Methyl Methacrylate	5.9	7.4	6.0	8.8	13.6	9.3	Naphthalene	7.4	10.8	4.7	15.9	9.5	6.6
Bromodichloromethane	3.6	12.7	2.7	5.0	10.9	10.4	1,2,3-Trichlorobenzene	4.8	10.7	8.2	14.9	5.1	5.5

¹Compound used a linear regression calibration curve fit, ²Compound used a quadratic regression calibration curve fit

Tom Hartlein¹, Amy Nutter¹, Tyler Trent¹, Simone Novaes-Card²

¹Teledyne LABS, Mason, Ohio, USA; ²Agilent Technologies, Inc.

EMS 2025



Agilent



TELEDYNE

Teledyne LABS




Figure 2. Comparison of average response factor calibration curves in water

Water Method Results