




Environmentally Friendly Implementation of EPA 8270E by Solid-Phase Extraction and Hydrogen Carrier Gas GC/MS Analysis

Arielle Coccozza | UCT, Inc. | 2731 Bartram Rd, Bristol, PA 19007 | (800)-385-3153 | info@unitedchem.com




INTRODUCTION

EPA Method 8270E is a standardized protocol for analyzing semivolatile organic analytes (SVOCs) in various matrices, including water, using gas chromatography/mass spectrometry (GC/MS). Using UCT's proprietary 8270 cartridge blend for solid phase extraction (SPE), this method effectively isolates acidic, basic, and neutral analytes, offering analytical coverage of a wide range of analytes in an aqueous matrix. The method's use of activated carbon cartridges improves the retention of highly polar analytes like 1,4-dioxane. In the wake of the recent restrictions on methylene chloride usage, reducing solvent usage by using SPE is paramount, while adopting hydrogen as a carrier gas provides a sustainable and economical alternative to helium. This alternative procedure achieves recoveries using SPE that are similar to traditional liquid-liquid extraction procedures, with 98% of analytes meeting EPA's recommended 50-150% recovery criteria.

References:

- [1] <https://www.epa.gov/sites/default/files/2015-12/documents/8000d.pdf>
- [2] https://www.epa.gov/sites/default/files/2020-10/documents/method_8270e_update_vi_06-2018_0.pdf
- [3] <https://www.epa.gov/sites/default/files/2015-12/documents/3535a.pdf>

SPE PROCEDURE



CHROMATOGRAPHY




Figure 1. Total Ion Chromatogram (TIC) of a method blank with surrogates at 10 and 20 µg/L and internal standard at 40 µg/L.




Figure 2. Total Ion Chromatogram (TIC) of the 1L fortified sample at 20 µg/L (with internal standards at 40 µg/L). (133 Compounds)

DISCUSSION & CONCLUSION


Initial demonstration of proficiency (IDP) was performed by fortifying 4 replicates of reagent water samples (133 analytes) at different concentration levels (1 µg/L, 5 µg/L, and 20 µg/L). The IDP data shown represents the 20 µg/L study (Figure 3). At this concentration, recovery values for 131 analytes ranged from 56% to 121%, with RSD values ≤20%. The analytes 1,4-phenylenediamine (34%) and hexachloropropene (41%) did not meet acceptable recovery criteria. These analytes are either known to pose challenges in extraction methods or lack historical data for aqueous extraction techniques.

Additional low-level concentration extractions were performed over three days to generate seven replicate data points per analyte. Using statistical analysis, these data were used to calculate the limit of detection (LOD) for each analyte. The LOD data shown represents the 5 µg/L study (Figure 4). At this concentration, LODs ranged from 0.25 µg/L to 6.28 µg/L, with a median LOD of 0.78 µg/L.

According to EPA Method 8270E, recovery values between 50–150% with ≤20% RSD are considered acceptable to validate the extraction method; laboratories are encouraged to establish in-house recovery limits based on sufficient ongoing data, as outlined in Method 8000D Section 9.6. This study demonstrates that the SPE method using EC8270 and EU521 cartridges, in conjunction with analysis by GC/MS with hydrogen carrier gas, provides accurate and precise recovery of a wide range of acidic, basic, and neutral semi volatile analytes. The results met these criteria across multiple concentration levels for 131 out of 133 compounds.

This approach offers a reliable alternative to traditional aqueous extraction methods, with added benefits of increased sample throughput, reduced solvent usage, and lower analyst exposure to hazardous chemicals, enhancing both laboratory efficiency and environmental safety. Furthermore, using hydrogen as the carrier gas in GC/MS analysis reduces operational costs and offers a more environmentally sustainable alternative to traditional helium-based methods.

RESULTS



GC/MS PARAMETERS

| | |
|-------------------|--|
| Instrument | Agilent 6890/5975C |
| Column | Rxi-5Sil MS, 30m x 250 µm x 0.25 µm |
| Inlet Temperature | 280 °C |
| Injection Volume | 1 µL; Split 10:1 |
| Liner | Restek Topaz Splitless, 4mm |
| Oven Program | 40°C (1 min) to 280 °C at 20 °C/min, 300 °C at 5 °C/min |
| Carrier Gas | Hydrogen at 1.2 mL/min constant flow rate |
| MS Parameters | Source at 230 °C; Quad at 150 °C; Transfer line at 280 °C; Electron energy; at 70 eV |
| Detector | MSD Quadrupole, full scan mode (EI), 50-500 amu |
| Tune File | DFTPP.U |

RESULTS

