

Detection of VOCs by Agilent 8697 Headspace with 7010D GC/MS/MS Using Hydrogen Carrier Gas

ASMS 2025
TP-###



Introduction

Volatile organic compounds (VOCs) are present in the environment and can alter the chemistry of the atmosphere. Current methods for analyzing these compounds use gas chromatography with a single quadrupole detector and the sample is introduced using a purge and trap¹. With a triple quadrupole mass spectrometer (GC/TO) as a detector, we can acquire in dynamic MRM (dMRM) mode which enables lower detection levels by minimizing interference from co-eluting compounds. Along with this TQ detector, we employed headspace as our sample introduction method and hydrogen carrier gas.

Table 1. List of volatile organic compounds (VOCs) analyzed.

Compound Name			
1,1,1-Tetrachloroethane	2-Propene nitrate (acrylonitrile)	Isonone	
1,1,1-Trichloroethane	4-Chlorobutene	Isopropylbenzene (Cumene)	
1,1,2-Dichloroethane	Bromobenzene	Methylcyclopentadiene	
1,1,2-Trichloroethane	Benzene	Methylstyrene	
1,1-dichloroethane	Bromobenzene	Methyl methacrylate	
1,1-dichloroethane	Bromobenzene	Norbornene	
1,1-dichloropropane	Bromodichloromethane	MIBK	
1,2,3-Trichlorobenzene	Bromoforn	MTBE	
1,2,3-Trichloropropane	Bromonitroethane	Naphthalene	
1,2,4-Trichlorobenzene	Carbon tetrachloride	Norbornene	
1,2,4-Trimethylbenzene	Chlorobenzene	n-Propylbenzene	
1,2-Dibromobenzene (EDB)	Chloroethane	o-Xylene	
1,3-Dichlorobenzene	Cis-1,2-Dichloroethane	p-Toluenesulfonic acid	
1,3-Dichlorobenzene	Chlorobenzene	p-Isopropylbenzene (Xyrene)	
1,3-Dichloropropane	cis-1,2-Dichloroethane	sec-Butylbenzene	
1,3-trimethylbenzene (cs) or 1,3-Dichloropropane	Chlorobenzene	Symne	
1,3,4-Trichlorobenzene	Dibromochloromethane	Tetrahydrofuran	
1,3,5-Trichlorobenzene	Dibromochloromethane	Tetralin	
1,4-dichlorobenzene	Dibromomethane	Toluene	
1,4-dichlorobenzene	Dichloromethane	trans-1,3-Dichlorobenzene	
2,2-dichloropropane	Diethyl ether	trans-1,4-Dichlorobenzene	
2-Butanone (MEK)	Ethyl methacrylate	Trichloroethene	
2-Chlorobutanol	Ethylbenzene	Trichloromethyl isopropylmethane	



Figure 1. Agilent 8890 GC with 8697 headspace and 7010D TQ was used for this analysis.

Experimental

Standard mixtures containing VOCs were prepared at 0.05–200 ppb in water to a final volume of 10 mL with 5 g sodium sulfate present to decrease the partition coefficient and thus improve sensitivity. Internal standard was spiked into each vial at a final concentration of 10 ppb. Matrix spikes were prepared using an artificial wastewater mimic containing surfactant, nutrients, and other organics. The matrix spikes were injected from the headspace into a split/splitless inlet (SSI) in pulsed split mode with a 1 mm direct liner. A 20 m \times 180 μ m \times 1 μ m DB-624 column was used with a constant flow rate of 0.6 mL/min of hydrogen. Data was acquired in dMRM mode at 6 cycles per second with a gain of 5. The method required 203 MRM transitions across 110 MRM groups. 72 analytes elute in less than seven minutes, making this a fast method while maintaining good resolution.

Table 2. Agilent 8890 GC and 7010D TQ Parameters

GC and MS Conditions:

Inlet	MMI, pulsed split mode 200°C, split 30:1 10 psi pulse until 0.3 min
Column	DB-624 UI Column (20 m x 180 µm, 1 µm)
Carrier gas	0.6 mL/min constant flow (hydrogen)
Oven (8.45 min run)	35°C hold 0.25 min Ramp 25°C/min → 240°C, no hold
MSD Transfer Line	250°C
Gain	5
Total MRMs	293

Table 3. Agilent 8697 HS Parameters

HS Conditions:	
Injection vol.	1 mL
Temperatures	Oven 75°C Loop 75°C Transfer line 115°C
Times	12 min vial equilibration 0.5min injection duration
Shaking	High, 250 shakes/min

Results and Discussion

Most compounds are detected at 0.05 ppb.

Under our method conditions, 64 compounds are detected at 0.05 ppb with a signal to noise ratio greater than 3. Of the 72 compounds analyzed, and 51 compounds calibrate over the full 0.05-200 ppb range, with another 9 calibrating from 0.1-200 ppb, 4 calibrating 0.2-200 ppb, and the final 8 calibrating 0.5-200 ppb. Average RF curve fits were prioritized, and 45 compounds had an Avg RF RSD less than 20%. Linear and quadratic fits were used when applicable where R^2 was greater than 0.99 and accuracies were +/- 30%.

Table 4. Full calibration data for VOCs analyzed in this study, showing retention time (minutes), curve fit type, curve fit value (Avg RF RSD, linear R², or quadratic R²), calibrated range (ppb), and what level the analyte was detected at (ppb) with S/N > 3.

Compound Name	RT (min)	C ₂ Type	FR	Range	Detectable	Compound Name	RT (min)	C ₂ Type	FR	Range	Detectable
Dichloro-Fluoromethane	0.83	Quadrupole	0.98	0.5-2.00	0.05	Ethyl methylacetate	3.60	Quadrupole	0.99	0.1-2.00	0.05
Chloromethane	0.94	Avy	16.0	0.5-2.00	0.05	1,2-Dichloromethane	3.66	Avy	14.9	0.5-2.00	0.05
Linear	0.97	Linear	0.97	0.5-2.00	0.05	Tetrahydrofuran	3.84	Avy	18.4	0.5-2.00	0.05
Bromomethane	1.20	Avy	35.9	0.5-2.00	0.05	1,1-Dichloroethane	3.74	Avy	14.9	0.5-2.00	0.05
Chloroethane	1.25	Avy	15.5	0.5-2.00	0.05	2-Hexanone	3.78	Quadrupole	0.99	0.1-2.00	0.05
1,2-Dichloroethane	1.38	Quadrupole	0.98	0.5-2.00	0.05	Dibromomethane	3.84	Avy	18.4	0.5-2.00	0.05
Ethyl ether	1.48	Avy	17.4	0.5-2.00	0.05	1,2-Dichloroethane (E08)	3.93	Quadrupole	0.98	0.5-2.00	0.05
1,1-Dichloroethane	1.58	Linear	0.91	0.5-2.00	0.05	Chlorobenzene	4.18	Avy	13.9	0.5-2.00	0.05
Acetone	1.60	Linear	0.98	0.5-2.00	0.05	1,1,1,2-Tetrachloroethane	4.21	Avy	17.5	0.5-2.00	0.05
Avy	1.64	Avy	34.4	0.5-2.00	0.05	Ethyl ethylacetate	4.28	Avy	17.9	0.5-2.00	0.05
Methylene chloride	1.77	Quadrupole	0.99	0.5-2.00	0.05	m-Xylene	4.29	Avy	17.3	0.5-2.00	0.05
2-Propanethiol	1.86	Linear	0.98	0.5-2.00	0.1	n-Xylene	4.51	Avy	13.5	0.5-2.00	0.05
1,2-Dichlorobenzene	1.87	Avy	51.3	0.5-2.00	0.05	Styrene	4.53	Avy	11.9	0.5-2.00	0.05
MTBE	1.88	Avy	39.3	0.5-2.00	0.05	Isononyl	4.62	Avy	16.3	0.5-2.00	0.05
1,1-Dichlorobenzene	2.04	Avy	31.1	0.5-2.00	0.05	Isopropylbenzene (cumene)	4.70	Avy	18.2	0.5-2.00	0.05
Avy	2.07	Avy	36.0	0.5-2.00	0.05	1,1,1-Trichloroethane	4.74	Avy	18.2	0.5-2.00	0.05
Methylchloride	2.27	Linear	0.93	0.1-2.00	0.1	1,2-Trichloropropane	4.86	Linear	0.98	0.5-2.00	0.05
1,2-Dichloropropane	2.27	Linear	0.98	0.5-2.00	0.05	Bromobenzene	4.88	Linear	11.4	0.5-2.00	0.05
2-Butanol (MBA)	2.37	Quadrupole	0.98	0.5-2.00	0.05	n-Propylbenzene	4.82	Avy	18.2	0.5-2.00	0.05
Chloroethane	2.38	Avy	39.7	0.5-2.00	0.05	2-Chloropropane	4.99	Avy	16.8	0.5-2.00	0.05
Chloroform	2.38	Linear	0.97	0.5-2.00	0.05	3,3,5-Trimethylbenzene	5.03	Avy	18.2	0.5-2.00	0.05
1,1,1-Trichloroethane	2.47	Linear	0.97	0.5-2.00	0.05	4-Chlorotoluene	5.05	Avy	13.9	0.5-2.00	0.05
Linear	2.50	Linear	0.99	0.5-2.00	0.05	Tetrahydrofuran	5.17	Avy	18.1	0.5-2.00	0.05
Calcium tetrachloride	2.54	Linear	0.97	0.2-2.00	0.05	Pentachloroethane	5.23	Avy	9.3	0.5-2.00	0.05
Benzene	2.62	Avy	17.1	0.5-2.00	0.05	2,4-Trimethylbenzene	5.25	Avy	16.3	0.5-2.00	0.05
1,2-Dichlorobenzene	2.62	Avy	35.3	0.5-2.00	0.05	sec-Butylbenzene	5.28	Avy	19.9	0.5-2.00	0.05
Trichloroethane	2.89	Avy	36.1	0.1-2.00	0.05	1,3-Dichlorobenzene	5.42	Avy	14.7	0.5-2.00	0.05
1,2-Dichloropropane	2.99	Avy	33.9	0.1-2.00	0.1	n-Isopropylacetate (cumene)	5.43	Avy	15.9	0.5-2.00	0.05
Methyl methylacetate	2.99	Linear	0.98	0.5-2.00	0.05	1,4-Dichlorobenzene	5.48	Avy	15.9	0.5-2.00	0.05
Avy	3.00	Avy	39.9	0.5-2.00	0.05	n-Butylbenzene	5.67	Quadrupole	0.99	0.5-2.00	0.05
1,4-Dioxane	3.05	Linear	0.99	0.5-2.00	0.1	1,2-Dichlorobenzene	5.69	Avy	14.0	0.5-2.00	0.05
Bromodichloromethane	3.10	Avy	39.5	0.5-2.00	0.05	DBP	6.15	Avy	17.6	0.5-2.00	0.05
Dichloromethane	3.30	Quadrupole	0.99	0.1-2.00	0.05	1,2,4-Trichlorobenzene	6.1	Avy	18.8	0.5-2.00	0.05
MBK	3.38	Quadrupole	1.0	0.1-2.00	0.1	Hexachlorobutadiene	6.74	Quadrupole	0.99	0.5-2.00	0.05
Toluene	3.47	Quadrupole	0.99	0.5-2.00	0.05	Naphthalene	6.79	Avy	12.5	0.5-2.00	0.05
1,2-Dichlorobenzene	3.47	Avy	34.1	0.5-2.00	0.05						




Figure 2. Chromatograms and calibration curves for several VOC compounds.

Results and Discussion

All VOC gases calibrated more than 2 orders of magnitude.

Early eluting VOC gases co-elute with water, resulting in distorted peak shapes. By running in dMRM mode, some of this distortion can be reduced, but not fully eliminated. We observed that at low ppb concentrations, vinyl chloride and bromomethane exhibit better peak shape with hydrogen carrier gas than with helium, but chloromethane peak shape is slightly better with helium carrier. Even though some of the gases exhibit distorted peak shapes at low levels, all six gases still calibrated over a range spanning at least 2.5 orders of magnitude.




Figure 3. Chromatograms and calibration curves of the six VOC gases.

Wastewater matrix spikes maintained sensitivity and reproducibility.

Eight replicate matrix spikes were prepared in artificial wastewater for each of the following levels: 0.05 ppb, 0.5 ppb, and 50 ppb. For compounds that calibrated down to 0.05 ppb, the average response RSD value at 0.05 ppb was 19.8 and the average concentration RSD was 15.9. The average concentration RSD for the 50 ppb replicates was 12.2 and the average response RSD was 14.9. The addition of matrix did not appear to negatively impact peak shape, sensitivity, or reproducibility.




Figure 4. Total ion chromatogram of a 50 ppb wastewater matrix spike.

Conclusions

All 72 compounds successfully calibrated over multiple orders of magnitude using a combination of new techniques for VOC analysis.

- Hydrogen carrier gas does not drastically diminish sensitivity.
- The 7010D TQ in dMRM mode enables sensitivity to detect low ppb levels.
- Static headspace is a reproducible sample introduction technique for VOCs.
- Further method development may be necessary to improve peak shape of the gases.

References

¹Quimby, B., et al. "Volatile organic compounds analysis in drinking water with headspace GC/MSD using hydrogen carrier gas and hydroinert source." *Agilent Technologies application note*, publication number 5994-4963EN, 2022