



Low Flow, High Gains—Using Microflow LC Techniques for Environmental Applications

Karl Oetjen, Diana Tran, Craig Butt, Kendra Adams; SCIEX, USA

- INTRODUCTION
- Liquid chromatography (LC) has been applied to a wide range of environmental samples; combining this with tandem mass spectrometry (MS/MS) allows for highly sensitive and accurate measurements
 - Electrospray ionization (ESI) is susceptible to matrix effects, including ion suppression
 - Microflow LC has been shown to achieve sensitivity gains and operates at significantly lower flow rates (up to 100x lower) compared to traditional analytical high-performance LC systems. These systems operate at flow rates in the range 1–200 µL/min and the droplets created have a diameter of only a few microns; generating more ions and minimizing ion suppression effects.
 - We compare microflow LC and analytical flow LC for the analysis of 69 frequently analyzed pesticides was performed. Both methods use the same SCIEX QTRAP 6500+ LC-MS/MS System. Sensitivity gains of up to 240x were observed for selected pesticides. The sensitivity gains from microflow LC compared to analytical flow LC enabled simpler sample preparation procedures without sacrificing limits of quantitation (LOQs) while also reducing solvent consumption for better sustainability.

ANALYSIS DETAILS

Microflow versus Analytical Flow: Stationary phase kept consistent

Luna Omega 3 µm Polar C18 100 Å 100 mm length

Analytical flow 800 µL/min on a 4.6 mm ID Column

Microflow: 15 µL/min on a 0.5 mm ID column

Mass Spectrometric Analysis:

Parameter	High flow	Microflow
Curtain Gas (CUR)	45 psi	20 psi
Ionspray Voltage (IS)	4500 V	4000 V
Heater Temperature (TEM)	300 ° C	200 ° C
Gas 1	60 psi	15 psi
Gas 2	50 psi	60 psi

Compound parameters like collision energy (CE) were kept consistent between both analytical flow and microflow LC analysis.

Case Study :Analysis of pesticides in 13 varieties of beer. Samples were sonicated for 10 minutes to remove carbonation, followed by direct LC-MS/MS injection.




Figure 1: Yearly solvent cost of micro flow versus high flow assuming 8 min run times and \$240 USD for 4L bottle of solvent

