

Membrane Filter Choice for Preparing PFAS Testing Samples

Lindsay D. Lozeau¹, Gunjan Shah¹, Ranjani Muralidharan², Amy Laws¹

¹EMD Millipore Corporation, Burlington, MA 01803, USA

²Sigma Aldrich Inc., St. Louis, MO 63178, USA



Scan for more
PFAS data!

Introduction

Highly sensitive, LC-MS/MS-based analytical methods for measuring perfluoroalkylated substances (PFAS) in complex matrices, such as wastewater, are becoming more prevalent in today's regulatory landscape. To improve data quality and preserve both instrument and column lifetime, analysts should carefully consider their sample preparation steps.

Filtration with a disc-type membrane filter or a syringe filter (**Fig. 1A-B**) is a simple, cost-effective way to clean up samples and mobile phases¹, and has been increasingly included in PFAS analytical methods for capture of particulates in liquid and air matrices (selected methods in **Table 1**).



Fig. 1: **A**, Syringe filters; **B**, disc-type membrane filters

With consumables like filters, major concerns about their use are: **1**, contamination of samples and **2**, loss of analytes due to unanticipated binding²; however, there are **many considerations** for choosing the right filter for a particular analytical method.

Table 1: Selected PFAS analytical methods that require filters

Method(s)	Matrix/Matrices	Sample Preparation	Analytical Method
ASTM D7968-17a	Environmental solids	Solvent extraction, filtration	LC-MS/MS
ASTM D7979-19	Water matrix (no drinking water)	Solvent extraction, filtration	LC-MS/MS
FDA C-010.03	Foods	QuEChERS, filtration	LC-MS/MS
OTM-45 (OTM-50)	Stationary sources (Air emissions)	Sampling train: filtration (particulates); Impingers (gaseous)	LC-MS/MS (GC-MS/MS)
ASTM D8535-23	Soil, Biosolids	Solvent extraction, filtration	LC-MS/MS
EPA 1633	Aqueous, soil, biosolids, sediment, tissue	SPE, filtration	LC-MS/MS
ASTM 8421-22	Aqueous matrices	Cosolvation, filtration	LC-MS/MS

The **goal of this series of studies** was to draw parallels between the performance of membrane filters in PFAS analytical methods and their physicochemical characteristics, in order to recommend the best filters to use depending on the method parameters.


The life science business of Merck KGaA, Darmstadt, Germany operates as MilliporeSigma in the U.S. and Canada

EMDMillipore.com


© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. MilliporeSigma, Millipore, and the vibrant M are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources. Lit. No. PBXXXXXXX XXX

Consideration #1: Filter Material



There are many types of filter materials




Filter material influences PFAS binding



Consideration #3: Diameter & Format



Syringe filters and cut discs perform similarly



Effect of Filter Diameter on Recovery of PFAS: 25mm vs. 33mm


For nearly every compound, a Millipore® nylon (33mm diameter) filter demonstrated better average recovery per cm² versus a device of 25mm. Thus, larger filtration areas do not necessarily lead to more analyte loss.

Effect of Filter Format on Recovery of PFAS: Disc-type vs. syringe filter

Recovery performance of 0.2µm PES is similar for most PFAS compounds in both disc-type and syringe filter, regardless of solvent used (left: water; right: methanol)

Consideration #2: Pore Size

Smaller pore sizes may lead to higher recovery



However, not all syringe filters with the same pore size retain particles in the same way

MFR	Material	% Retention
MilliporeSigma	Nylon	100.0 ± 0.10
	PES	69.4 ± 28.1
	PP	57.1 ± 3.5
MFR-5	Nylon	98.0 ± 1.1
	PP	25.3 ± 0.90
MFR-3	PP	7.9 ± 1.1
MFR-2	RC	15.8 ± 2.2

Table 2: Percent retention of 0.24 µm diameter polystyrene beads (0.05% solution) by various 0.2 µm pore size syringe filters

Recovery of PFCAs by MCE of various pore sizes

There may be a slightly negative relationship in recovery with increasing pore size, from 0.2 to 0.8µm, in water.


This trend was seen up through 5.0µm pores and in both methanol and water (not shown) and may be caused by changes in diffusional gradients of small molecules within filter pores.

Not all syringe filters perform the same way

Why?
• Pore size definition and type of testing
• Force involved in pressure-based filtration devices
• Manufacturing approaches and lot-to-lot variation
• Chemical interactions with particles or aggregation at the surface

Consideration #4: Solvent

ASTM 8421 demonstrates that PFAS Recovery increases with methanol content for many filter types



Adding even a small amount of solvent increases recovery dramatically, even for non-polar filter types. PES recovery of long chain PFCAs was the highest in a mixture of polar and non-polar solvents.

Excellent recovery of even the longest chain PFCAs in 50% (v/v) methanol (via ASTM 8421) seen for all filter types. Similar results observed for philic polypropylene.

References

- [1] Lozeau, L.D. & Dube, M., *Analytix Reporter*, **2024**, 16, 28-31
- [2] Lath, S., et al. *Chemosphere*, **2019**, 222, 671-8.

Acknowledgements

SGS Orlando, Alpha Analytical (Pace), Pace Analytical, David Bell, Dave Brewster, Vivek Joshi, Taylor Reynolds, and Maricar Dube

Millipore®

Preparation, Separation, Filtration & Monitoring Products