

Analysis of PFAS in Drinking Water: Validation Studies of EPA Method 537.1 Using StayClean™ QSight LC-MS/MS



Sheng-Suan (Victor) Cai
Senior Field Application Scientist

August 5th, 2025



HUMAN HEALTH • ENVIRONMENTAL HEALTH


© 2015 PerkinElmer

StayClean™ QSight LC-MS/MS

- Dual Binary LC Pump
 - Two Sets of Mobile Phases
 - A1 and B1 for ESI
 - A2 and B2 for APCI
- Dual ESI/APCI Source
 - ESI for Polars
 - APCI for Non-polars



Robustness: PFAS in Food Matrices (19,100 Continuous Injections)



Food Matrices:
salmon/avocado/tomato mix,
spice powder, dog/cat feed

Preparation using QuEChERS¹
method developed for PFAS
analysis in food, but the final
SPE step was omitted to
maximize sample matrix
components in the final matrix
extracts

Instrument robustness
evaluated by intermittently
monitoring solvent QC samples
between large blocks of 100
consecutive matrix injections


Salmon, Avocado, Tomato Mix, Spice Powder, Dog/Cat Feed

Source Images: PFAS in Food Matrices

Before




After



After 19,100 Continuous Injections of Salmon, Avocado, Tomato Mix, Spice Powder, Dog/Cat Feed

Robustness: 25,500 Continuous Injections of Fetal Bovine Serum

Peak Area Ratios (All QC Injections)



● Sulfamethoxazole ● Alprazolam ● Diazepam

- Protein precipitation of FBS with methanol followed by a 50% dilution with water prior to injecting 5 μ L for analysis
- Monitoring 3 analytes and their deuterated standards (1 ng/mL)
- Instrument robustness evaluated by intermittently monitoring SST samples between large blocks of FBS matrix injections

Source Images



View from Source Window



Sample Cone




Plenum Chamber

After 25,500 Continuous Injections of Fetal Bovine Serum

Summary of PFAS EPA Drinking Water Methods

- EPA Method 537.1
 - Selected linear PFAS
 - 18 PFAS with chain lengths C₄-C₁₄
 - Isotopic internal standards and reverse phase SPE
- EPA Method 533
 - Short-chain & polar PFAS
 - 25 PFAS with chain length C₄-C₁₂
 - Addition of more polar fluorotelomers and ether carboxylic acids
 - Isotopic dilution and ion exchange SPE
- EPA Method 1633:
 - Broad range of PFAS in multiple matrices.
 - 40 PFAS in aqueous, solid, biosolids and tissue samples.
 - Isotopic dilution and multiple sample preparation techniques.

EPA Method 537.1 Summary






Major Challenges of EPA Method 537.1

- PFAS ubiquitous, leading to elevated background level
 - Found in clothing, carpet, food packaging, etc.
 - LC-MS and SPE systems constructed from PTFE
 - Glass containers absorb PFAS. Use polyethylene/polypropylene containers
 - Wear nitrile gloves to reduce contamination





Remediation Steps to Reduce PFAS Contamination

Source of Contamination	Remediation
Mobile phases	<ul style="list-style-type: none">• Purchase LC-MS grade solvents• Use a delay column
PTFE parts & tubing in HPLC pump	Use a delay column
PTFE Tubing in HPLC autosampler	Replace with PEEK tubing
Vials and PTFE lined caps	Use only polyethylene vials and caps
PTFE tubing in SPE apparatus	Replace with polyethylene tubing

Reducing Background from Pump and Mobile Phases




Reducing Background from Pump and Mobile Phases (Cont'd)



Use a delay column to retain background PFAS interferences from instrument related sources


Replace PTFE Autosampler Tubing with PEEK

- Tubing kit available
- Can be easily installed by chemist or install engineer




Manual Large Volume SPE Apparatus

- Eliminate any PTFE or fluoropolymer components that will contribute to PFAS background.
- Replace PTFE transfer tubing with either Linear Low-Density Polyethylene (LLDPE) or PEEK tubing.
- Use stopcocks constructed of polyethylene (PE).
- Collect samples in PE centrifuge tubes.



Sample Preparation Procedure (Manual)



Automated SPE Device

- Pre-programmed for EPA methods 537.1, 533 and 1633, etc.
- SPE-03 Automated SPE with MOD-004 for PFAS and applications that require automated rinsing of up to 250mL bottles



Modified “Fast” HPLC Method

Original EPA 537.1 LC Method

Time (min)	% 20 mM ammonium acetate	% Methanol
Initial	60.0	40.0
1.0	60.0	40.0
25.0	10.0	90.0
32.0	10.0	90.0
32.1	60.0	40.0
37.0	60.0	40.0

Column: Waters Atlantis® C₁₈, 2.1 x 150 mm, 5.0 µm

Flow Rate: 0.3 mL/min

Injection Volume: 10 µL

Runtime: 37 min

Modified PerkinElmer EPA 537.1 LC Method

Time (min)	% 10 mM ammonium acetate	% Methanol
Initial	95.0	5.0
0.7	95.0	5.0
1.0	55.0	45.0
7.0	2.0	98.0
8.0	2.0	98.0
8.1	95.0	5.0
10.0	95.0	5.0

Column: Brownlee SPP® C₁₈, 4.6 x 75 mm, 2.7 µm

Flow Rate: 0.8 mL/min


Injection Volume: 10 µL

Runtime: 10 min


- Runtime reduced from 37 min to 10 min (73% reduction).
- Same retention time order
- Excellent resolution of branched isomers

Original Method vs. Modified Method






Original EPA method
Total runtime: 37 min.



Mod. PerkinElmer method
Total runtime: 10 min.



Baseline Resolution Between Branched & Linear Chain Isomers



QSight 220 MS Parameters

- Two MRM experiments per analyte (only quantifier shown here)
- Quantifier/Qualifier ion ratios used to confirm identity
- Dwell times optimized with MRM Time Manager
- Minimum of 10 data points across all peaks

ESI Source Parameters

Parameter	Value
Polarity	Negative
Drying Gas	110.0
HSID Temperature (°C)	280.0
Nebulizer Gas 1	400.0
Electrospray V1 Negative	-3500.00
Source 1 Temperature (°C)	350.0

Quantitative MRM Experiments





Analyte	Q1 Mass	Q2 Mass	CE	Voltages	
				EV	CCL2
PFBS	299.5	79.8	59	-40	250
13C2-PFHxA (SS)	315.0	270.0	13	-10	48
PFHxA	313.0	269.1	13	-10	52
13C3-HFPO-DA (SS)	286.9	168.9	12	-5	44
HFPO-DA	285.0	168.9	14	-5	40
PFHpA	363.0	319.0	14	-10	56
PFHxS	399.0	80.0	91	-45	120
ADONA	377.0	251.1	17	-10	64
PFOA	413.0	368.9	14	-10	68
13C2-PFOA (IS)	421.0	376.0	15	-5	76
PFNA	463.0	419.0	16	-10	76
13C4-PFOS (IS)	507.0	79.9	103	-45	124
PFOS	499.1	79.9	100	-50	300
9Cl-PF3ONS	530.9	350.9	35	-30	112
13C2-PFDA (SS)	515.0	469.9	16	-13	84
PFDA	513.0	468.9	16	-10	84
d3-NMeFOSAA (IS)	573.0	419.0	27	-25	104
NMeFOSAA	570.0	419.0	27	-20	108
PFUnA	562.9	518.9	17	-10	96
d5-NEtFOSAA (SS)	589.0	419.0	28	-20	112
NEtFOSAA	584.0	418.9	27	-20	96
11Cl-PF3OUdS	630.9	450.9	36	-40	176
PFDaA	612.9	568.9	17	-10	104
PFTTrDA	662.9	618.9	18	-11	104
PFTA	712.9	668.8	17	-10	116

Calibration Linearity

All calibration curves forced through zero as required in method

Compound	Instrument Calibration Range (ng/L) ^a	Method Calibration Range (ng/L) ^b	R ²
PFBS	16.4 - 26287	0.07 - 105.1	0.99944
PFHxA	5.5 - 29703	0.02 - 118.8	0.99870
¹³ C ₂ -PFHxA	4.6 - 24752	0.02 - 99.0	0.99890
¹³ C ₃ -HFPO-DA	67.5 - 24752	0.27 - 99.0	0.99923
HFPO-DA	18.5 - 29703	0.07 - 118.8	0.99848
PFHpA	5.5 - 29703	0.02 - 118.8	0.99840
PFHxS	5.2 - 28218	0.02 - 112.9	0.99976
ADONA	5.2 - 28218	0.02 - 112.9	0.99904
PFOA	5.5 - 29703	0.02 - 118.8	0.99983
PFOS	5.3 - 28515	0.02 - 114.1	0.99740
PFNA	18.5 - 29703	0.07 - 118.8	0.99926
9CI-PF3ONS	5.1 - 27772	0.02 - 111.1	0.99979
PFDA	81.0 - 29703	0.32 - 118.8	0.99900
¹³ C ₂ -PFDA	4.6 - 24752	0.02 - 99.0	0.99880
NMeFOSAA	5.5 - 29703	0.02 - 118.8	0.99983
PFUnA	18.5 - 29703	0.07 - 118.8	0.99680
NEtFOSAA	5.5 - 29703	0.02 - 118.8	0.99680
d5-NEtFOSAA	18.3 - 99010	0.07 - 396.0	0.99620
11CI-PF3OUdS	5.2 - 28069	0.02 - 112.3	0.99972
PFDoA	18.5 - 29703	0.07 - 118.8	0.99630
PFTTrDA	5.5 - 29703	0.02 - 118.8	0.99590
PFTA	5.5 - 29703	0.02 - 118.8	0.99670



- a) Instrument calibration range is the actual concentration range of calibration standards used to determine calibration curves.
- b) Method calibration range is determined by multiplying the instrument calibration range by 1/250 to account for the SPE sample preparation/concentration.



Instrument Sensitivity (LOQ & LOD)

- LOQ and LOD values estimated from the S/N of calibration standards and adjusted for 250-fold sample concentration
- Most LODs below < 0.08 ng/L demonstrating the MS system exceeds the sensitivity required for this method

Analyte	LOD (ng/L)	LOQ (ng/L)
PFBS	0.012	0.040
PFHxA	0.015	0.049
HFPO-DA	0.029	0.098
PFHpA	0.012	0.041
PFHxS	0.003	0.009
ADONA	0.002	0.007
PFOA	0.016	0.053
PFNA	0.015	0.048
PFOS	0.004	0.012
PFNA	0.014	0.004
9Cl-PF3ONS	0.004	0.013
PFDA	0.042	0.141
NMeFOSAA	0.002	0.007
PFUnA	0.016	0.054
NEtFOSAA	0.002	0.007
11Cl-PF3OUdS	0.006	0.019
PFDoA	0.009	0.030
PFTTrDA	0.018	0.060
PFTA	0.325	1.084



Method Detection Limits and LCMRLs

➤ Detection limits (DL) in this study well below detection limits reported in EPA 537.1

➤ $DL = s \times t(n-1, 1-\alpha=0.99)$

➤ LCMRL = Lowest Concentration Minimum Reporting Limit

➤ LCMRL calculated using EPA calculator

<https://www.epa.gov/dwanalyticalmethods/lowest-concentration-minimum-reporting-level-lcmrl-calculator>

Analyte	Experimental DL (ng/L) ^a	EPA 537.1 DL (ng/L) ^b	Experimental LCMRL (ng/L) ^c	EPA 537.1 LCMRL (ng/L) ^d	Experimental MRL (ng/L) ^e
PFBS	1.1	6.3	0.720	1.800	1.42
PFHxA	1.6	1.7	0.930	1.000	1.00
HFPO-DA	0.91	4.3	0.570	1.900	1.60
PFHpA	0.41	0.63	0.100	0.710	1.60
PFHxS	0.85	2.4	0.600	1.400	0.29
ADONA	0.10	0.55	0.10	0.880	0.28
PFOA	0.34	0.82	0.340	0.530	0.30
PFOS	1.5	2.7	1.000	1.100	1.00
PFNA	0.98	0.83	0.500	0.700	1.60
9Cl-PF3ONS	0.95	1.8	0.680	1.400	1.50
PFDA	0.82	3.3	0.4	1.600	0.30
NMeFOSAA	1.9	4.3	0.220	2.400	1.00
PFUnA	0.58	5.2	0.300	1.600	1.60
NEtFOSAA	1	4.8	0.730	2.800	1.60
11Cl-PF3OUdS	0.45	1.5	0.390	1.500	0.28
PFDoA	0.3	1.3	0.190	1.200	0.30
PFTTrDA	1.2	0.53	0.820	0.720	4.00
PFTA	2.1	1.2	1.500	1.100	4.00

a) Experimental DL was determined from ten LFB replicates fortified at ~1.6 ng/L measured over three days and calculated according to section 9.2.8 in EPA Method 537.1

b) Reference DL values from EPA Method 537.1 rev 2.0 (Table 5) determined from seven LFB replicates fortified at 4.0 ng/L measured over three days and calculated according to section 9.2.8

c) Experimental LCMRLs were determined from ten replicates each at five fortification levels ranging from ~0.2 – 80 ng/L using the EPA LCMRL Calculator.^[11]

d) Reference LCMRL values from EPA Method 537.1 rev 2.0 (Table 5).

e) Experimental MRLs were determined from ten LFBs fortified at concentrations ranging from ~0.2 to 4.0 ng/L according to section 9.2.6 of EPA Method 537.1 rev 2.0 using the Half Range prediction interval method with confirmed upper and lower Prediction Interval Results (PIR) ≤150% and ≥50%, respectively.

% Recovery at Low Concentration

- LFRBs fortified with surrogates and analytes at a conc. of ~1.6 ng/L & ~4.0 ng/L
- Experimental recovery results for ten replicate LFRBs (n=10)
- Experimental recoveries at ~1.6 ng/L range from 86-154%
- Experimental recoveries at ~4.0 ng/L range from 86-108%

Analyte	Experimental Recovery			Experimental Recovery		
	Fortified Conc. (ng/L)	Mean % Recovery	% RSD	Fortified Conc. (ng/L)	Mean % Recovery	% RSD
PFBS	1.4	98	1.6	3.5	106	3.6
PFHxA	1.6	115	9.6	4.0	109	4.1
HFPO-DA	1.6	86	6.7	4.0	86	4.2
PFHpA	1.6	116	6.3	4.0	112	1.9
PFHxS	1.5	100	4.0	3.8	105	1.3
ADONA	1.5	112	4.3	3.8	110	2.4
PFOA	1.6	103	2.8	4.0	107	1.4
PFNA	1.6	96	10.3	4.0	114	1.9
PFOS	1.5	117	6.8	3.8	109	1.7
9Cl-PF3ONS	1.5	98	3.0	3.7	99	3.0
PFDA	1.6	107	5.6	4.0	101	5.8
NMeFOSAA	1.6	108	5.0	4.0	103	2.4
PFUnA	1.6	117	5.9	4.0	108	2.3
NEtFOSAA	1.6	111	5.9	4.0	105	3.2
11Cl-PF3OUdS	1.5	97	3.3	3.8	97	4.3
PFDoA	1.6	108	5.8	4.0	104	2.7
PFTTrDA	1.6	127	16.0	4.0	108	3.5
PFTA	1.6	154	24.1	4.0	106	2.8
¹³ C ₂ -PFHxA	40.0	93	4.6	40.0	104	3.4
¹³ C ₃ -HFPO-DA	40.0	110	6.7	40.0	121	3.3
¹³ C ₂ -PFDA	40.0	102	5.2	40.0	105	1.9
d ₅ -NEtFOSAA	160	99	8.3	160	106	4.2

%Recovery at Medium & High Concentrations

Analyte	Experimental Recovery (n=10)			EPA 537.1 Recovery (n=7)		
	Fortified Conc. (ng/L)	Mean % Recovery	% RSD	Fortified Conc. (ng/L)	Mean % Recovery	% RSD
PFBS	14.1	105	4.4	16.0	90.8	6.8
PFHxA	16.0	114	2.4	16.0	101	8.0
HFPO-DA	16.0	94	1.0	16.0	97.8	1.8
PFHpA	16.0	118	1.9	16.0	105	3.3
PFHxS	15.2	107	1.4	16.0	109	6.7
ADONA	15.1	116	4.0	16.0	108	1.3
PFOA	16.0	109	4.4	16.0	106	1.8
PFNA	16.0	108	1.7	16.0	110	2.6
PFOS	15.4	116	4.6	16.0	111	4.7
9CI-PF3ONS	15.0	105	3.8	16.0	108	8.8
PFDA	16.0	115	4.3	16.0	111	2.4
NMeFOSAA	16.0	111	4.9	16.0	104	5.2
PFUnA	16.0	115	2.7	16.0	107	2.8
NEtFOSAA	16.0	111	3.5	16.0	97.7	6.8
11CI-PF3OUdS	15.1	102	3.6	16.0	109	3.4
PFDaA	16.0	112	3.2	16.0	101	7.2
PFTTrDA	16.0	112	4.7	16.0	108	2.6
PFTA	16.0	109	4.4	16.0	110	0.9
¹³ C ₂ -PFHxA	40.0	100	5.9	40.0	88.5	6.4
¹³ C ₃ -HFPO-DA	40.0	95	5.6	40.0	94.5	3.2
¹³ C ₂ -PFDA	40.0	104	7.2	40.0	99.1	3.4
d ₅ -NEtFOSAA	160	103	4.2	160	90.0	2.6


Experimental recoveries at 16 ng/L range from 94-116%

Analyte	Experimental Recovery (n=10)			EPA 537.1 Recovery (n=7)		
	Fortified Conc. (ng/L)	Mean % Recovery	% RSD	Fortified Conc. (ng/L)	Mean % Recovery	% RSD
PFBS	70.8	101	4.2	80.0	85.1	6.7
PFHxA	80.0	109	6.5	80.0	97	4.6
HFPO-DA	80.0	90	6.8	80.0	96.8	5.1
PFHpA	80.0	111	7.3	80.0	104	2.7
PFHxS	76.0	103	3.9	80.0	107	4.4
ADONA	75.6	110	6.8	80.0	106	3.6
PFOA	80.0	103	4.6	80.0	104	3.1
PFNA	80.0	102	4.7	80.0	104	3.6
PFOS	76.8	109	7.4	80.0	107	4.8
9CI-PF3ONS	74.8	99	4.1	80.0	101	3.8
PFDA	80.0	109	7.2	80.0	107	3.6
NMeFOSAA	80.0	101	3.8	80.0	102	5.4
PFUnA	80.0	108	7.0	80.0	101	1.3
NEtFOSAA	80.0	104	4.0	80.0	101	2.5
11CI-PF3OUdS	75.6	97	3.2	80.0	103	6.1
PFDaA	80.0	105	7.6	80.0	107	3.7
PFTTrDA	80.0	103	6.5	80.0	99.1	3.6
PFTA	80.0	101	6.6	80.0	97.2	3.6
¹³ C ₂ -PFHxA	40.0	103	6.2	40.0	97.0	4.9
¹³ C ₃ -HFPO-DA	40.0	94	6.0	40.0	101	9.9
¹³ C ₂ -PFDA	40.0	104	6.4	40.0	106	2.7
d ₅ -NEtFOSAA	160	100	3.8	160	99.5	4.8

Experimental recoveries at 80 ng/L range from 90-111%

PFAS Spiked Recovery Summary

- Spike level at 0.3, 4, 16 and 80 ppt
- Acceptance criteria: 70–130%



Field Sample Analysis

- Field samples of tap water were collected from three different municipalities in the Southeast US, and were designated M1, M2 and M3.
- Four field samples and one FRB were collected at each location
- Ongoing QC requirements
 - Laboratory reagent blank (LRB)
 - Continuing calibration check (CCC)
 - Laboratory fortified blank (LFB)
 - Surrogate recovery
 - Laboratory fortified sample matrix (LFSM) and duplicate (LFSMD)
 - Field sample duplicate (FD)

Field Sample Analyte Test Results

Analyte	Average FD Conc (ng/L)			Average LFSM % Recovery ^a			LFSM RPD ^b		
	M1	M2	M3	M1	M2	M3	M1	M2	M3
PFBS	2.0	14.9	<MRL	120	100	119	5.6	16.0	1.4
PFHxA	1.8	2.0	<MRL	101	95	120	2.1	2.8	6.0
HFPO-DA	<MRL	<MRL	<MRL	116	90	108	4.1	18.0	1.1
PFHpA	<MRL	<MRL	<MRL	103	88	99	2.3	1.2	0.4
PFHxS	0.32	0.56	<MRL	89	75	81	5.3	0.4	0.0
ADONA	<MRL	<MRL	<MRL	114	107	111	2.5	6.8	1.3
PFOA	1.1	1.9	0.39	88	78	88	3.6	8.9	7.2
PFOS	<MRL	2.0	<MRL	129	111	126	0.1	7.0	2.9
PFNA	<MRL	<MRL	<MRL	90	82	92	9.1	12.8	0.1
9CI-PF3ONS	<MRL	<MRL	<MRL	118	97	115	6.2	0.2	2.2
PFDA	0.35	0.37	<MRL	82	128	121	2.1	3.3	1.0
NMeFOSAA	<MRL	<MRL	<MRL	96	85	94	1.7	6.5	0.7
PFUnA	<MRL	<MRL	<MRL	75	120	139	0.2	1.4	5.0
NEtFOSAA	<MRL	<MRL	<MRL	98	84	97	6.3	6.6	0.3
11CI-PF3OUdS	<MRL	<MRL	<MRL	57	86	100	9.0	2.3	4.3
PFDaA	<MRL	<MRL	<MRL	124	118	129	0.2	2.1	0.3
PFTTrDA	<MRL	<MRL	<MRL	120	106	113	2.4	0.7	9.2
PFTA	<MRL	<MRL	<MRL	94	83	92	6.9	1.6	19.3
				% Recovery acceptance criteria: 70-130% or 50-150% for low level spike			% RPD acceptance criteria: ≤30%		

Field Sample Surrogate %Recovery

Surrogates	Average FD %Recovery			Average LFSM %Recovery		
	M1	M2	M3	M1	M2	M3
$^{13}\text{C}_2$ -PFHxA	104	112	127	100	106	117
$^{13}\text{C}_3$ -HFPO-DA	106	93	103	104	94	97
$^{13}\text{C}_2$ -PFDA	76	81	79	81	76	73
d5-NEtFOSAA	110	106	106	111	106	102

% Recovery acceptance criteria: 70-130%

Conclusions

- EPA Method 537.1 requires **special equipment and sample handling** considerations to minimize the PFAS background to measure analytes at ppt (ng/L) levels.
- **Modifications of the QSight LC-MS/MS system and the SPE system** are required to perform this method.
- An **improved LC method** reduces total runtime by 73% (10 min. vs. 37 min.).
- **Method detection limits (DL) well below** those reported in the EPA method.
- **Most Instrument LOQs are <0.060 ppt**, confirming that the QSight 220 LC-MS/MS exceeds the sensitivity required to perform Method 537.1.
- **Analyte recoveries** in LFRBs fortified at 16 and 80 ng/L were all within the range of **94-118%** as required in the method with **most analyte recoveries between 94-111%**.
- **Field samples were analyzed** with good analyte and surrogate spiked recoveries.
- This method has been implemented and validated at multiple customer sites.
- Instrument robustness with **StayClean™ technology** would allow injections of thousands of samples with no instrument cleaning and maintenance.



Thank you very much for attention! Any question?


Victor Cai

Senior Field Application Scientist

sheng-suan.cai@perkinelmer.com

Mobile 951-258-2470

Innovative QSight Design



First vertical MS system

- Dimensions: 50x50x110 cm
- Smallest footprint
- Saves customer lab space

Modularized design

- Improved serviceability
- Fast recovery, low down time
- Quick replacement parts

Smart

- Full diagnostics
- Remote controlling



Plug and play design

- Require minimum tooling
- Reduced cabling by 80%
- Self connecting modules

How StayClean™ Works?



- Heated Surface Induced Desolvation (HSID)
- Self Clean 24/7
 - Better desolvation
 - Less Maintenance
 - Increased uptime

Self-cleaning Heated Surface Induced Desolvation (HSID)



Laminar Flow Ion Guide (Patent Technology)

- No voltage needed to guide ions
- Ions go with the flow, no mass discrimination
- No voltage re-optimization needed (Never re-tune ion guide)
- ~98% ion transmission rate



Summary of PFAS EPA Drinking Water Methods

- EPA Method 537.1 (March 2020): Selected Linear PFAS Compounds
 - Targeted analysis of **18 PFAS** compounds with **chain lengths C₄-C₁₄** including carboxylic acids, sulfonic acids, sulfonamidoacetic acids and GenX compounds (HFPO-DA)
 - LC/MS/MS method using isotopic **internal standards** and reverse phase SPE sample preparation
- EPA Method 533 (December 2019): Short-Chain & Polar PFAS Compounds
 - Targeted analysis of **25 PFAS** compounds **C₄-C₁₂** expanded with the addition of more **polar** fluorotelomers and ether carboxylic acids
 - LC/MS/MS method using **isotopic dilution** and ion exchange SPE sample preparation
- EPA Method 1633 (Draft 2021): Broad range of PFAS in multiple matrices.
 - Targeted analysis of **40 PFAS** compounds in aqueous, **solid, biosolids and tissue** samples.
 - LC/MS/MS using **isotopic dilution** and multiple sample preparation techniques.




EPA Method 537.1 Summary

- A 250-mL water sample is fortified with surrogates
- Polystyrenedivinylbenzene (SDVB) phase is required for SPE
- No deviations allowed for SPE extraction
- SPE extracts dried under N₂ and reconstituted with 1-mL of 96:4 (vol/vol) MeOH/Water (**250-fold concentration**)
- Internal standards are added
- A 10-μL injection is made into an LC-MS/MS instrument with a C18 column
- Analytes are separated and identified by comparing the acquired mass spectra and retention times to reference spectra and retention times for calibration standards acquired under identical LC-MS/MS conditions
- The concentration of each analyte is determined by using the internal standard technique

Original EPA 537.1 Chromatogram – 37 min Runtime

Table 3. Method Analytes, Retention Times (RT) and Suggested IS References


Analyte	Peak # (Fig. 1)	RT (min)	IS# Ref
PFBS	1	7.62	2
PFHxA	2	10.42	1
HFPO-DA	4	11.38	1
PFHpA	6	13.40	1
PFHxS	7	13.58	2
ADONA	8	13.73	1
PFOA	9	15.85	1
PFOS	11	17.91	2
PFNA	13	17.92	1
9Cl-PF3ONS	14	18.91	2
PFDA	15	19.69	1
NMeFOSAA	17	20.50	3
PFUnA	19	21.21	1
NEtFOSAA	20	21.26	3
11Cl-PF3OUdS	22	21.84	2
PFDoA	23	22.52	1
PFTtDA	24	23.66	1
PFTA	25	24.64	1
¹³ C ₂ -PFHxA	3	10.42	1
¹³ C ₃ -HFPO-DA	5	11.40	1
¹³ C ₂ -PFDA	16	19.69	1
d ₅ -NEtFOSAA	21	21.24	3
¹³ C ₂ -PFOA- IS#1	10	15.85	-
¹³ C ₄ -PFOS- IS#2	12	17.91	-
d ₃ -NMeFOSAA-IS#3	18	20.49	-



Modified EPA 537.1 Chromatogram – 10 min Runtime

Analyte	Peak #	RT (min)	IS# Ref
PFBS	1	3.54	2
PFHxA	2	4.15	1
HFPO-DA	4	4.34	1
PFHpA	6	4.78	1
PFHxS	7	4.77	2
ADONA	8	4.84	1
PFOA	9	5.30	1
PFOS	11	5.73	2
PFNA	13	5.74	1
9CI-PF3ONS	14	5.93	2
PFDA	15	6.13	1
NMeFOSAA	17	6.31	3
PFUnA	19	6.45	1
NEtFOSAA	20	6.47	3
11CI-PF3OUdS	22	6.56	2
PFDoA	23	6.72	1
PFTTrDA	24	6.96	1
PDTA	25	7.16	1
¹³ C ₂ -PFHxA SS#1	3	4.15	1
¹³ C ₃ -HFPO-DA SS#2	5	4.34	1
¹³ C ₂ -PFDA SS#3	16	6.12	1
d ₅ -NEtFOSAA SS#4	21	6.46	3
¹³ C ₂ -PFOA IS#1	10	5.29	-
¹³ C ₄ -PFOS IS#2	12	5.74	-
d ₃ -NMeFOSAA IS#3	18	6.30	-

2020-10-22: LFRB-High6
Method TIC (47 pairs) Total Number of Experiments: 14



Manual SPE System Tubing Test

- Do not use PTFE tubing for sample transfer tubing.
- We used 1/8" OD x 1/16" ID linear low density polyethylene tubing (Freelin-Wade).
- Tested 1/8" OD x 1/16" ID PEEK tubing for transfer tubing.
- LLDPE tubing was chosen over PEEK due to flexibility.

