

Analysis of Polycyclic Aromatic Hydrocarbons in Polyurethane Foam (PUF) and Diethylene Glycol/Water Samples using Automated Pressurized Liquid Extraction, Cleanup and Solid Phase Extraction

R. Addink, T. Hall
Fluid Management Systems
Billerica MA



Introduction

- PAHs analyzed worldwide: air, sediment, soil, water.
- Industrial sources: coal gasification, municipal solid waste incineration, petroleum refining.
- PAHs result of incomplete combustion.



- Can be found in incinerator emissions, dust, fly ash, bottom ash
- Subject to monitoring and regulations in most countries
- Stack emissions can be collected with sampling using poly-urethane foam (PUF) and solvent impingers (diethylene glycol-water)



Pressurized Liquid Extraction

- An automated Extraction technique can be used for PAHs collected on PUF
- The Technique Incorporates:
 - Solvent
 - Pressure
 - Heat
 - Time



Why is PLE so effective?

- Performed near the solvent's supercritical region
- Under Programmable Pressure
- Creates a high degree of analyte solubility releasing PAHs from the solid matrix

Extraction

- A solid or semi-solid sample is placed in the Pressurized Extraction Cell 5mL to 200mL
- The Extraction cell is capped and placed into the extraction device which can be pressurized to up 2500psi

Extraction

- The Extraction cell is filled with the extraction solvent put under pressure and depressurized
 - PAHs (suitable for lots of other analytes too)
- The Extract is flushed with Nitrogen into a collection vessel



The PLE®

Pressurized Liquid Extraction

PLE - Pressurized Liquid Extraction


- High Speed
- Modular and expandable from 1 to 8
- Process 1 to 8 samples in 10 to 15 min
- Extraction cell size 5 mL to 200 mL
- Real time plot of temperature and pressure
- Reduced Solvent Consumption
- Lower Energy Consumption
- In Cell Sample Cleanup





Extraction Cell

Heater



Economical Extraction Cells



Easy to Use End Caps




Modular and Expandable

Expandable from 1 to 8 Modules

Parallel Extraction



Method Documentation



Pressurized Liquid Extraction for PAHs in PUF

- Large PUF plug (ca 6 cm OD x 7 cm length) spiked with native PAH standards
- Sample placed in 40 mL extraction cell and void space filled with Ottawa Sand
- Capped with disposable Teflon end caps
- 50% Dichloromethane/50% Hexane for 15 min at 100 °C and 1500 psi
- 15 min cool down
- Nitrogen flush to transfer analytes and extract to 50 mL collection tubes

Concentration of PLE extract



Extract after Cleanup

- System pre-heated to 40 °C.
- Samples evaporated at stable T under 5-6 psi nitrogen.
- Solvent exchange to hexane.
- Reduce volume to 1 mL.



Automated EZprep/+ for cleanup



EzPrep/ + Features (1)

Features:

- Programmable Flow rate and Volume
- Pressure indicator and over pressure alarm
- Real time read-out for dispensed volume and pressure
- Ability to select from 1 to six samples
- Can accommodate up to 4 solvents
- Economical & less expensive automation





EZPrep/+ Features (2)


- Closed loop system, eliminates background contaminants & exposure to chemicals
- Optimized for solvent reduction while obtaining highest possible recoveries
- Certified disposable Columns with guaranteed low contaminants background and Excellent Recoveries
- Quick connect SNAP columns simplifies system set up
- Multi pump Solvent Delivery system brings convenient automated solvent selection & dispense with controllable flow & volume
- EzPrep/+ designed with Minimum number of electronics and Electromechanical valve to lower cost and simplify the maintenance



Columns/ Oxidation Capacity



EzPrep /+ Control Panel



PAHs Cleanup

- Assemble 6 g neutral silica columns with EZPrep +
- Syringe vial at top is used for sample loading.
- **Stage 1:**
- Dilute sample extract with 9 mL hexane.
- Load sample extract onto silica column
- Elute column with 10 mL hexane (waste)
- **Stage 2:**
- Elute columns with 35 mL dichloromethane, collecting aromatic fraction (PAHs).



GC vial

Final Concentration



- Reduce volume to 1 mL.
- 24 positions in parallel.
- $T = 25\text{ }^{\circ}\text{C}$ for most volume reductions.
- Nitrogen flow 1-2 psi.

Impinger w/PAHs: use SPE (1)

LLE/CLE

Open to laboratory background

Uses >360mls solvent

Shaking / Continuous process

Forms emulsions requiring centrifuging

Little Selectivity

Requires water removal

Automated SPE

Closed system

Uses <60mls solvent

Filtration process

No emulsions formed

Wide Selectivity (adsorbent)

In-line water removal



Comparison of LLE/CLE vs SPE Methods (2)

LLE/CLE

No Separation of waste

More volume to evaporate

Massive solvent emission

CLE uses a lot of solvent

Requires lots of solvent for cleaning

Semi-Automated SPE

Separates Aqueous and Organic Waste

<60mls solvent to evaporate

6 times less solvent emission

Easily Capture Solvent

Lower solvent costs

Lower Disposal Costs



Automated SPE system



TurboTrace® SPE

System Features

- **Simultaneous Processing:** Modular, processes up to six samples in parallel.
- **Versatile Sample Handling:** Compatible with beverages, blood, milk, urine.
- **Sample Size:** Handles sample volumes from 2 mL to unlimited.
- **Programmable Operations:** Built-in computer.
- **Precise Delivery:** Positive pressure pumping for accurate solvent delivery.
- **Loading:** Incorporates vacuum or positive pressure pumps for efficient sample loading.
- **On-line evaporation:** Concentrates up to 250 mL directly into a GC vial for analysis.



SPE for PAHs (1)

- Acidify 10% Diethylene glycol/water samples with 6 N Hydrochloric Acid to pH < 2
- Add surrogates if so desired
- Load samples on automated system
- Install C-18 cartridges
- Fill rinse bottles with 35 mL dichloromethane
- Condition cartridges: 10 mL dichloromethane, 20 mL methanol and 35 mL water (waste)
- Load samples under vacuum across cartridges (waste)



SPE for PAHs (2)

- Dry with nitrogen for 3 min
- Rinse sample bottles with 35 mL dichloromethane and load across cartridges for elution and collect
- Elute 5 mL, another 5 mL, and another 2 mL of dichloromethane for soaking and elute and collect
- Blow out fraction lines with nitrogen and collect
- Dry samples over sodium sulfate (can be in-line cartridges)
- Concentration in 12-position evaporator as described.




Analysis




Agilent 7010 B Triple
Quad GC/MS/MS

PAHs recoveries from PUF

Average PAHs recoveries from PUF



PAHs recoveries from impinger



Conclusions (1)

- Work illustrates use of Sample Prep automation in monitoring stack emissions
- Solid matrices used for collection can be easily processed with Pressurized Liquid Extraction, automated cleanup with column chromatography and automated sample concentration
- PLE is green technique compared to Soxhlet because of fast process ~ 1 h, less power use
- Automated Solid Phase Extraction can provide quick processing of emissions collected in liquids using impingers. No additional cleanup is necessary.
- Automation reduces errors and produces consistent, reproducible results.
- Increased productivity.

Conclusions (2)

- Faster and easier operator training
- Automatic documentation of extraction and cleanup and concentration conditions
- Reduced errors due to mistakes eliminating manual steps and conditions.
- Reduced solvent usage and disposal costs.



Come see us at our booth
Questions?

raddink@fms-inc.com

thall@fms-inc.com