

Comparative analysis of air sampling strategies for VOC monitoring

Dr Caroline Widdowson MBA MRSC

cwiddowson@markes.com



Caroline Widdowson 
Specialist in trace chemical characterisation - Chair of BSI committee
EH2/5 Indoor Air Quality - Market Development Professional

An Indian Case Study....

India's pollution crisis linked to pre-term births and low birth weight, study warns
Just 10mg increase in PM 2.5 per cubic metre of air linked to 5 per cent rise in rate of low birth weight

New Delhi air pollution: Schools closed and construction stopped as smog worsens to levels far above WHO safety limit

A thick, toxic smog has enveloped India's capital, home to more than 33 million people. One technology firm which monitors air pollution ranked New Delhi as the most polluted city in the world.

Air pollution: Low ozone layer found over Brahmaputra river valley

Tropospheric, or ground level ozone, is created by chemical reactions between oxides of nitrogen (NOx) and volatile organic compounds (VOC).

Quarter of all deaths worldwide linked to environmental damage and pollution, UN says

"Do we continue on our current path, which will lead to a bleak future for humankind, or do we pivot to a more sustainable development pathway? That is the choice our political leaders must make, now"

Josh Gabbatore Science Correspondent | @jgabbatore |
Wednesday 13 March 2019 12:30 |



Air pollution responsible for more deaths than smoking, study says
"EU lagging a long way behind" in tackling toxic air, say researchers as study shows pollution deaths double previous estimates

Alex MacLennan-King Health Correspondent |
Tuesday 12 March 2019 10:45 |



We Made Plastic. We Depend on It. Now We're Drowning in It.

The miracle material has made modern life possible. But more than 40 percent of it is used just once, and it's choking our waterways.


Tuesday, 9 April 2019

WORLD
ECONOMIC
FORUM

In India, air pollution is cutting lives short by up to 11 years

Source: University of Chicago's Air Quality Life Index (AQLI)

Where should we be monitoring?



Where start?:

- **Establish baseline VOC levels** using advanced analytical tools
- **Deploy real-time ambient air quality system** near industrial clusters
- **Prioritize regulation enforcement** in pollution hotspots
- **Use data for policy advocacy** and public health protection

CLEAN AIR FOR HEALTH

#AirPollution



Case Study

New Bureau of Indian Standards

- BIS 5182-27
 - Air Pollution - Methods for Measurement Part 27 Vapour-Phase Organic Chemicals Vinyl Chloride to nC22 Hydrocarbons in Air and Gaseous Emissions by Diffusive (Passive) Sampling onto Sorbent Tubes or Cartridges Followed by Thermal Desorption (TD) and Capillary Gas Chromatography (GC) Analysis.
- BIS 5182-28
 - Air Pollution - Methods for Measurement Part 28 Vapour-Phase Organic Chemicals C3 to nC22 Hydrocarbons in Air and Gaseous Emissions by Pumped Sampling onto Sorbent Tubes or Cartridges Followed by Thermal Desorption (TD) and Capillary Gas Chromatography (GC) Analysis.



Analytical Instrumentation: TD-GC-MS system

How do we monitor VOCs in air?



Markes – Thermal desorber (TD)



Agilent – GCMS

What is thermal desorption (TD)?

Sample collection on sorbent tube

- Sample (e.g. air) is collected



- Compounds of interest are adsorbed on the sorbent surface




- Lighter gases such as nitrogen, argon and carbon dioxide pass through



What is thermal desorption (TD)?

Sample desorption


- Tube heated in a reversed flow of clean carrier gas ('backflushed')




What is thermal desorption (TD)?

Sample desorption

- Tube heated in a reversed flow of clean carrier gas ('backflushed')




- Compounds are released from the sorbent into the flow of carrier gas



PROBLEM:


Compounds are released
SLOWLY from the sorbent tube
Would lead to very wide
chromatographic peaks and low
sensitivity



Two stage thermal desorption

Stage 1: Sample desorption

SOLUTION: Use a narrow secondary trap




Electrically cooled
narrow bore cold trap


Two stage thermal desorption

Stage 2: Trap desorption


- Cold trap heated rapidly ($100^{\circ}\text{C}/\text{sec}$) for sharp chromatographic peaks
- **Backflush** of cold trap for greater volatility range



What is thermal desorption (TD)?



TD sampling options




Sampling methods

Passive (Diffusive) sampling

- Simple and low-cost for large scale air monitoring studies
- Unobtrusive for personal monitoring
- Single bed sorbent
- Tube inlet must have 15 mm air gap with 5 mm I.D.

Key methods:



- Fenceline monitoring (US EPA Method 325)
- Ambient air quality / Industrial air (EN 14662-4)
- Indoor, ambient and workplace air (ISO 16017-2, BIS 5182-27)



Passive sampling

Axial and Radial

- **Axial (tube based) sampling**
 - Sorbent bed exposed through diffusion cap fitted to an industry standard sorbent tube
 - The industry standard sampling route
 - Allows for short and long-term sampling (8 hours to 4 weeks)
- **Radial sampling**
 - Sorbent cartridge exposed through cylindrical (coaxial) outer surface of a diffusive body
 - Allows for a much larger exposed surface area, and a shorter diffusive path length
 - Allows faster uptake rates (typically 1 hour to 3 days)



Passive sampling

Axial

- Simple and low-cost for large scale air monitoring studies
- Unobtrusive for personal monitoring (doesn't impact human behavior)
- The mass of each VOC collected is proportional to the ambient concentration

To calculate ambient concentration we need to know:

Constant


— **U: Uptake rate** —————→

Verified, pre-published data.
Dependant on air gap and surface area.

Experimentally determined


— t: Time

— M: Mass on sorbent tube



Sampling methods

Pumped (Active) sampling



- Pump air through sorbent tube (typically 10-200 mL/min)
- Multiple sorbents \Rightarrow wide volatility range
- Targeted and untargeted compounds can be sampled and analysed in a single experiment

Key methods:

- PFAS (ASTM D8591-24)
- Ambient air (US EPA TO-17 / Chinese EPA HJ 644 / EN 14662-1)
- Stack emissions (CEN/TS 13649 / Chinese HJ 734)
- Indoor, ambient and workplace air (ISO 16017-1, BIS 5182-28)
- Polycyclic Aromatic Hydrocarbons (PAHs)
- Breath – Active sampling




Active Sampling Experiment

- Pump air through sorbent tube
 - Flow rate: **20 mL/min**
 - Time: 1 hr
 - Volume: 1.2 L
 - 4x locations:
 - Industrial Area
 - Residential Area
 - Traffic Area
 - Petrol station



Active Sampling – Chromatogram overlays



Common compounds found in all locations




Compound	RT	Mass	Formula
Isobutane	4.27	43.08	C ₄ H ₁₀
Butane	4.35	43.07	C ₄ H ₁₀
Butane, 2-methyl-	4.64	43.07	C ₅ H ₁₂
Pentane, 2-methyl-	5.44	43.07	C ₆ H ₁₄
Ethyl Acetate	6.12	43.04	C ₄ H ₈ O ₂
Benzene	7.01	78.04	C ₆ H ₆
Cyclohexane, methyl-	8.25	83.07	C ₇ H ₁₄
Toluene	9.15	91.07	C ₇ H ₈
Octane	9.75	43.07	C ₈ H ₁₈
Acetic acid, butyl ester	10.05	43.04	C ₆ H ₁₂ O ₂
Ethylbenzene	10.97	91.06	C ₈ H ₁₀
Benzene, 1,3-dimethyl-	11.11	91.05	C ₈ H ₁₀
m,p-Xylene	11.52	91.05	C ₈ H ₁₀
Benzene, propyl-	12.46	91.04	C ₉ H ₁₂
Benzene, 1,2,3-trimethyl-	13.24	105.05	C ₉ H ₁₂
Indane	13.64	117.06	C ₉ H ₁₀
Undecane	14.23	57.08	C ₁₁ H ₂₄
Dodecane	15.31	57.07	C ₁₂ H ₂₆
Tridecane	16.30	57.06	C ₁₃ H ₂₈
Tetradecane	17.22	57.07	C ₁₄ H ₃₀



[Industrial_area]
35 entities



[Traffic_area]
49 entities



[Petrol_pump]
43 entities



[Residential]
30 entities

Unique compounds for each location



Active sampling

- **Industrial area**



- Heptane, 2,2,4,6,6-pentamethyl-
- Hexadecane, 2,6,10,14-tetramethyl-
- 2-Ethoxyethyl acetate
- Tridecane, 6-methyl-

- **Petrol station area**



- Cyclohexane, 1,4-dimethyl-
- 2-Naphthalenol, decahydro-
- 2,4-Dimethyl-1-heptene
- Propane, 1,2-dichloro-
- Ethanol, 2-butoxy-
- Ethane, 1,2-dichloro-

- **Residential area**



- Trichloroethylene
- Cyclotetrasiloxane, octamethyl-
- Octanal

- **Traffic area**







- Tricyclo[5.2.1.0(2,6)]dec-4-ene, 4-methyl-
- 2-Heptene, (E)-
- Heptane, 3-methyl-
- 3a,4,5,6,7,7a-Hexahydro-4,7-methanoindene
- Linalool
- Furan, 3-methyl-
- Butane, 2,2-dimethyl-
- Cyclohexane, 1,3-dimethyl-, cis-
- Cyclohexanol, 5-methyl-2-(1-methylethyl)-
- Cyclohexane, ethyl-
- 1H-Indene, 2,3-dihydro-1,6-dimethyl-

Quantification of BTEX



Active – 1hr sampling (20ml/min; 1.2L sample)

				
Compound (ppb)	Industrial	Residential	Traffic	Petrol station
Benzene	4.94	4.20	18.60	10.37
Toluene	48.04	12.62	52.84	37.79
Ethylbenzene	45.29	3.82	18.36	13.55
m,p-Xylene	119.64	14.26	142.70	75.08


Passive Sampling Experiment




- Sorbent tubes exposed for 24-hours, 1 week and 2 weeks respectively to the air at:
 - Industrial Area
 - Residential Area



Passive Sampling – Industrial area



Passive Sampling – Residential area



Quantification of BTEX – Comparison

Passive (1-wk) v Active sampling



Compound (ppb)	Industrial		Residential	
	Passive	Active	Passive	Active
Benzene	4.03	4.94	6.17	4.20
Toluene	20.97	48.04	19.13	12.62
Ethylbenzene	29.73	45.29	8.57	3.82
m,p-Xylene	34.53	55.27	12.33	5.61


Summary

- Excellent linearity for calibrations $R^2 > 0.997$
- Great correlation between passive and active sampling
 - Some variability with passive sampling, but still giving comparable results to active
- Passive and active sampling complimentary
 - Passive sampling appropriate technique when sampling over a long period of time
 - Weight time average reduces the impact of a one-off pollution event (Average conc over 1-week period)
 - Pumped sample accurate for sampling time
 - Would need to take multiple pumped samples to generate true picture

Thank you



- Praveen Arya - Agilent India
- Dr. SK Tyagi (retd) - **Former Director & Divisional Head** (Air Toxics Lab, CAAQM, Quality Assurance, Environmental Training) & Quality Manager at **Central Pollution Control Board**
- Kiran Piduru & Hannah Calder - Markes International



Contact Markes



enquiries@markes.com



UK: +44 (0)1443 230935

USA: +1 866-483-5684 (toll-free)

Germany: +49 (0)69 6681089-10

P.R. China: +86 21 5465 1216



www.markes.com



[@MarkesInt](https://twitter.com/MarkesInt)



<https://uk.linkedin.com/company/markes-international>